

Travel in comfort and style



info@continentglobe.com



OM | P25014 | 1006060 | R01



WARNING

Failure to follow any or all safety instructions may result in serious injury or death

DO NOT ADJUST, MAINTAIN, OR OPERATE THIS TRANSPORT CHAIR WITHOUT READING, UNDERSTANDING, AND FOLLOWING THE PROPER INSTRUCTIONS INCLUDED WITH THIS PRODUCT

WARNINGS

- ▶ **Risk of tipping over** - This transport chair is not equipped with an anti-tip device. It must be used only under the supervision and with the assistance of an attendant. This product is not designed to be used as a self-propelled manual wheelchair. The user must not attempt to move by pushing against a table, wall, desk, or any other surface, as this could cause the chair to tip backward and result in serious or fatal injury.
- ▶ This transport chair has not been approved as a seat surface to be used within a vehicle.
- ▶ Do not operate this transport chair on streets or roadways.
- ▶ Do not attempt inclines unassisted.
- ▶ Engage the wheel locks upon entering or leaving the transport chair.
- ▶ Transport chair must be on a level surface when transferring user to or from seat.
- ▶ Footrests should be no less than 2 1/2 inches from the floor
- ▶ Do not stand on or put pressure on footrests when entering or exiting the transport chair.
- ▶ Do not lean over the top of the transport chair back. This can cause the transport chair to tip over.
- ▶ Do not use the wheel locks to slow your descent down and incline.
- ▶ When going down an inclined, changing directions can cause instability.

Usage Instructions

To Open

1. Pull armrests apart to unfold chair.
2. Swing folded-down backup and lock upright into place on back supports.
3. Attach footrests on each side by inserting the holes onto the vertical pegs of the frame. Swing footrests inward to lock in place.

To Fold

1. Push button release on both back uprights to fold back down behind chair.
2. Pull folding straps on outer seat edges towards each other to fold chair.
3. To remove footrests, pull silver release lever to unlock. Lift the footrests off the hooks.

To Transfer Patient

1. Push forward on brake handle to lock rear wheels and secure chair from rolling.
2. Swing footplates up to outer sides or turn both footplates upward to a vertical position.

To Adjust Length of Footrests

1. Locate tool wrench included with the footrests.
2. Turn wheels so they are in the "wheeled" position (as pictured).
3. Use the tool wrench to slightly loosen the silver hex screw located on the bottom of the footrests by turning the screw left.
4. Pull on the footplate to desired length. Footrests should be a minimum of 2 1/2" from the ground to ensure proper clearance.
5. Re-tighten the hex screw.
6. Repeat on other side.

Care and Maintenance

General Care

Wipe off all metal parts at least once a week using a soft cloth. Dry chair immediately if it is exposed to moisture. Abrasive cleaners or materials will scratch the surface of the chair.

Tires

Minimal maintenance is required for rubber tires. Wipe with a damp cloth occasionally. If tires become cracked or worn, replace immediately.

Front Casters

Inspect casters periodically and lubricate them when necessary to lubricate, remove the fork and grease the caster bearings. Remove the wheel from the fork and grease the wheel bearings. The wheel and fork can be reassembled by tightening the nuts until there is no play remaining, but the wheel rolls freely. The fork will revolve freely when adjusted correctly. The caster will utter if the nut is too loose or be difficult to steer if the nut is too tight.

Wheels

Inspect monthly for correct tightness. If the wheel seems to be loose or wobbly, release the lock nut and tighten the axle bolt. Lean chair to one side and spin adjusted wheel to check for correct adjustments.

Wheel Locks

Wheel locks should be checked for correct locking action before actual use. Pushing the handle completely forward activates the wheel locks. If the wheel locks are working correctly, they will be embedded into the tire at least 1/8". To adjust wheel locks, loosen nut on carriage bolt. Slide clamp toward rear wheel until lock shoe is embedded at least 1/8 "into rear tire when lock handle is engaged. Tighten nut to secure into position.

Care and Maintenance

1. Do not attempt to repair the transport chair yourself.
2. Call the store where you purchased the transport chair. Most dealers do repairs and will re-order parts for you.
3. If your dealer is unable to help you, please contact Continent Globe at:

info@continentglobe.com

Lifetime Limited Warranty

*One year on parts * Five years on brakes * Limited lifetime on Frame*

Continent Globe Corp warrants this product to be free from defects in material or workmanship for a period of lifetime from date of purchase to the first retail purchaser if within the specified warranty period, this product is proven to be defective, it will be repaired or replaced at the manufacturer, Continent Globes option. This warranty does not cover shipping costs incurred in replacement or repair of this product.

Continent Globe stands on our commitment to providing the highest quality products, as well as exceptional service. This product was built to precise standards and thoroughly inspected prior to shipment. This Lifetime Limited Warranty represents our confidence in the materials and workmanship of our products.

For warranty service, please contact the dealer from whom you purchased your product. In the event that you do not receive satisfactory service, please call us directly at the number below. Be prepared to indicate the nature of the defect, the product serial number, and the name, and location of your dealer. You will also be required to provide a receipt of your purchase. DO NOT return any products to our office without prior consent.

This warrant does not cover problems due to user negligence or misuse of the product. This warranty shall not extend to non-durable components, including but not limited to rubber tips, hand grips, casters, and other parts subject to normal "wear and tear."



Record this information for future reference:

Serial # _____

Date of purchase _____

Purchased from _____