
# luditta tilt-in-space wheelchair



JSA- Cod. 220169 26.05.25



**USE AND MAINTENANCE HANDBOOK** 

# ORMESA

MADE IN ITALY SINCE 1980



ORMESA thanks you for choosing JUDITTA tilting seat system.

JUDITTA is a medical device designed and manufactured by ORMESA srl.

**ORMESA s.r.l. recommends that you read this manual carefully to fully understand its contents**. In addition to help you become familiar with the device quickly, it also contains practical tips for making the best use of it safely and keeping it in good working order.

If, after reading it, you still have questions, contact your retailer, who will be happy to help you, or call AMYLIOR directly at T 1 450 424-0288 | 1 888 453-0311 or send an email to info@amylior.com.

#### Ormesa Team



### TABLE OF CONTENTS

WARNINGS	4
GUARANTEE	6
HOW TO STORE AND TRANSPORT THE MEDICAL DEVICE	6
WHAT NOT TO DO WITH JUDITTA	7
HOW TO TILT THE WHEELCHAIR AND CLIMB UP PAVEMENTS	
HOW TO FACE UP TO SLOPES IN GETTING ON AND OFF (B60 version)	
HOW TO USE JUDITTA FOR TRANSPORT ON A MOVING VEHICLE (private cars bus, etc)	s, 12
LABELS AND PLATES	16
LABELS AND PLATES	17
ABOUT JUDITTA WHEELCHAIR. INTENDED USE	18
ENVIRONMENT OF USE AND TYPE OF USERS	19
CONTRAINDICATIONS	
TECHNICAL FEATURES and MODELS	20
INITIAL PREPARATION	24
ASSEMBLY	26
HOW TO FIX SELF-PROPELLING WHEELS (B60 version)	
ACTIVATION of the ANTI-TIP DEVICE (B60 version)	
USING JUDITTA WHEELCHAIR	34
BRAKING THE WHEELCHAIR ( B30 MODEL)	34
BRAKING THE WHEELCHAIR ( B12 MODELS)	
BRAKING THE WHEELCHAIR ( B60 MODEL)	
ADJUSTING THE TILT OF THE BACKREST	
ADJUSTING THE INCLINATION OF THE SEAT BOTTOM (TILT- IN- SPACE)	
CLOSING THE FOOTRESTSADJUSTING THE FOOTRESTS IN HEIGHT	
FLEX. EXTENSION ADJUSTMENT OF THE FOOTRESTS	
TURNING THE LEGRESTS	
REMOVING THE LEGRESTS	
ADJUSTING THE LEGREST INCLINATION	
TURNING THE CALF RESTS	
DEPTH ADJUSTMENT of the CALFRESTS	



FLESSO-EXTENSION of the CALFRESTS	
ADJUSTING THE HEIGHT OF THE ARMRESTS	
OPENING THE SIDE SUPPORTS	
ADJUSTING THE HEADREST	
ADJUSTING THE DEPTH OF THE SEAT	
HEIGHT ADJUSTMENT OF THE BACKREST	
TILT ADJUSTMENT OF THE PUSH HANDLE	
REMOVING THE UPHOLSTERY FOR THE WASHING	
ASSEMBLING THE ADDITIONAL COMPONENTS	56
RECLINABLE TABLE 824-R	56
ABDUCTION BLOCK 834	63
FOOTREST LOCK 914	64
CALF REST 916	66
REMOVABLE FOOTPLATE COVER 892	68
REMOVABLE TABLE COVER 896	
HEADREST REMOVABLE COVER 936	70
TIEDOWN HOOKS (4 RINGS) 891	70
IV POLE 933 (B30 and B12 version)	71
SET OF HAND BRAKES 905 (B30 version)	73
SET OF HAND BRAKES 905 (B60 version)	77
ANTI-TIP UP DEVICE (B30 version) 932	
45° PELVIC HARNESS 894	85
MAINTENANCE	86
1. ORDINARY MAINTENANCE (monthly)	86
2. PREVENTIVE MAINTENANCE (biennial)	
3. EXTRAORDINARY MAINTENANCE (or corrective)	
CLEANING, DISINFECTION, SANITIZATION	
INFORMATION	
WARNINGS	
PROCEDURE	
ASSISTANCE AND SPARE PARTS	
DISPOSAL	91
CONDITIONS FOR DURARILITY RELISE AND REASSIGNMENT TO A NEW	IISER 92



### MARNINGS

- Carefully read the instructions in this manual before using the device; they will help you use the product safely and keep it in good working order. It is an integral part of the product and should be kept for future reference
- The manual is intended for technically qualified retailers, those who use the product and those who assist them; it provides instructions for using the product correctly and must be supplemented by adequate knowledge of the rehabilitation program for which it has been prescribed
- The manual reflects the technical state of the product at the time it was sold. ORMESA reserves the right to make any changes to the product or manual suggested by experience, technical considerations or regulatory developments without any obligation to update its previous production and relevant manual
- JUDITTA tilting seat system may be only used by persons who have read the rules and instructions in this manual
- The doctor must check that the use of JUDITTA tilting seat system and its accessories conforms to current laws, regulations and standards
- The seat system must be used by persons of the BODY SIZE AND WEIGHT specified in the paragraph "TECHNICAL FEATURES" on pages 20, 21, 22 and 23
- The seat system must be guided by persons of suitable height and weight to control its movements
- Before each use, always check that the brakes are working, the wheels are not excessively worn and that the headrest, back and footrest are locked
- When pausing, always lock both brakes, even on level ground
- Always lock the brakes before the user gets in or out of the seat system and before making any adjustments!
- The upholstery is fireproof in conformity with the UNI 1021-1/2 standard.
   Nevertheless, keep it away from cigarettes, matches, fireplaces and open flame in general
- Keeping the cushion and patient clean will prevent bedsores



#### JUDITTA IS A MECHANICAL MEDICAL DEVICE:



Perform REGULAR MAINTENANCE and

CLEANING following the instructions shown in the "Maintenance" chapter on page 86-87-88-89-90-91 and HAVE THE DEVICE INSPECTED AT THE INTERVALS INDICATED to VERIFY that it is WORKING PROPERLY and in GOOD CONDITION, otherwise warranty will fail and €€ marking will lapse



PRODUCT REPAIRS other than the regular maintenance shown on page 86-87-88-89-90-91 of the manual must be ONLY MADE by a SPECIALIZED SERVICE CENTER in the maintenance of mechanical aids for disables, otherwise warranty and C marking will be voided



Any CHANGES in the product and the REPLACEMENT WITHOUT ORIGINAL PARTS OR COMPONENTS ARE NOT AUTHORIZED.

REFER ONLY TO ORMESA, otherwise warranty and CE marking will be voided



In case of **DOUBT** about the **SAFETY** of the product or **DAMAGE** to parts or components, you are urged to **IMMEDIATELY DISCONTINUE USE** and **CONTACT** a **SPECIALIZED SERVICE CENTER** in the maintenance of mechanical aids for disables, or directly ORMESA.





**ORMESA guarantees the product for 2 years ;** in case of problems, contact the retailer where you purchased it. Always ask for original spare parts, otherwise the guarantee will decline

- ORMESA will not be liable for damage in the following cases:
  - use by an unsuitable person;
  - incorrect assembly of parts or accessories;
  - unauthorized modifications or service:
  - use of other than original replacement parts and parts subject to wear (upholstery, wheels, etc):
  - improper use (such as, transporting objects or loads larger or heavier than those shown in the instruction handbook t);
  - damage caused by incorrect use and lack of regular maintenance, as shown in the instruction handbook;
  - exceptional events:
  - failure to follow the instructions in this manual.

THE GUARANTEE DOES NOT COVER WEAR PARTS, which are subject to wear and tear, such as the upholsteries and the wheels.



#### HOW TO STORE AND TRANSPORT THE MEDICAL DEVICE

THE DEVICE MUST BE STORED AND PACKED using Ormesa original packaging materials, unless the guarantee will be voided

- -Once unpacked, THE TRANSPORT OF THE DEVICE must be done by ANCHORING ADEQUATELY it to the vehicle
- When travelling by plane, or in the car, DO NOT SUBJECT THE FOLDED FRAME TO LOADS THAT, especially with road bumps, COULD DAMAGE ITS STRUCTURE

#### THE DEVICE MUST BE PARKED / STORED IN CLOSED AND DRY PLACES

- THE OPERATING ENVIRONMENT HAS NO PARTICULAR INFLUENCE ON THE DEVICE UNLESS IT IS USED INCORRECTLY, such as by leaving it parked for a long time in direct sunlight or exposed to bad weather such as rain, or in marine environments, where salt could corrode the paint and sliding parts. IN THIS CASE, WE RECOMMEND CAREFUL CLEANING AND DRYING THE FRAME FOLLOWING THE INSTRUCTIONS SHOWN IN THE "MAINTENANCE, CLEANING AND DISINFECTION" CHAPTER ON PAGES 86-87-88-89-90-91 AND THE WARNINGS \$\frac{\Delta}{2}\$ ON PAGE 86





#### WHAT NOT TO DO WITH JUDITTA

- DO NOT USE THE DEVICE WITH SUBJECTS OF LARGER DIMENSIONS OR WEIGHT than shown on pages 21, 22 and 23
- DO NOT ALLOW CHILDREN TO USE THE DEVICE OR EVEN PLAY WITH IT
- DO NOT PLACE EXCESSIVELY HOT LIQUID CONTAINERS OR OBJECTS ON THE TABLE SURFACE that could cause damage or burns if overturned
- DO NOT ALLOW ANYONE TO STAND ON THE FOOTREST OR THE BACK OF THE WHEELCHAIR
- DO NOT LIFT THE WHEELCHAIR BY THE FRAME OF THE LEGREST OR THE FOOTREST to climb over obstacles
- NEVER LEAVE THE USER UNATTENDED, ESPECIALLY ON THE BEACH OR NEAR A RAVINE, SWIMMING POOL, SLOPE, OPEN FLAME OR STAIRS
- DO NOT PARK THE DEVICE FOR LONG LONG WAY IN THE SUN OR NEAR HEAT SOURCES: this will prevent the device from overheating and discoloration of the upholstery
- DO NOT CARRY THE USER OVER ROUGH OR STEEP SURFACES OR STAIRS
- NEVER LEAVE THE AID PARKED ON SLOPING GROUND
- DO NOT USE THE WHEELCHAIR TO CARRY THE USER UP AND DOWN STAIRS: the structure of the wheelchair was not designed for that purpose (see fig. on pag. 10)

- DO NOT CLIMB CURBS OR STEPS WITH JUST THE SIDE WHEELS (front or back) because it could tip over (see fig. on page 10)


When the SEAT is COMPLETELY TILTED BACKWARDS and the BACK is COMPLETELY RECLINED, DO NOT CLIMB OVER OBSTACLES OR STEEP SLOPES that, SINCE THEY COULD EXCEED 12°, COULD CAUSE THE WHEELCHAIR TO TIP OVER (SEE FIG. PAGE 10)

- DO NOT USE JUDITTA WHEELCHAIR IN THE SHOWER: its steel structure is not suitable for that purpose. Contact ORMESA S.r.I. or your retailer for a specific painted aluminium aid made for that purpose
- DO NOT USE THE PRODUCT IF PARTS ARE MISSING OR DAMAGED. Always insist on original replacement parts or the warranty will be void

In case of HYPERKINETIC patients, the Anti-tip up device, art. 932, is recommended when JUDITTA B30 is used tilted and reclined. (see fig. on page 85)




#### HOW TO TILT THE WHEELCHAIR AND CLIMB UP PAVEMENTS







#### CLIMB FACING FORWARD



#### DESCEND FACING BACKWARD








#### HOW TO FACE UP TO SLOPES IN GETTING ON AND OFF (B60 version)

Help is needed to get on and off ramps and slopes

The user can still get on a slight climb by himself, as JUDITTA B60 has an anti-tip device.







## HOW TO USE JUDITTA FOR TRANSPORT ON A MOVING VEHICLE (private cars, bus, etc...)

JUDITTA WHEELCHAIR was tested and conforms to the ISO 7176-19 and ANSI RESNA WC/19 standard, which defines the safety requirements for transport, facing in the direction of travel on moving vehicles (private cars, buses, etc..) following these warnings:

- THE WHEELCHAIR-SEATED OCCUPANT MUST FACE THE FRONT OF THE VEHICLE (see fig. 2, page 14). According to the ISO 7176-19 point 6.3.4, compliance with ISO 7176-19 does not preclude using the wheelchair facing rearwards in large accessible vehicles equipped and approved with rear-facing wheelchair passenger stations
- The User should weigh no more than 136 kilos because the CRASH TEST provided by the ISO and ANSI/RESNA standards only refers to wheelchairs/pushchairs for users whose weight is 136 kilos or less
- THE VEHICLE MUST BE APPROVED FOR TRANPORTING ORTHOPAEDIC WHEELCHAIRS/PUSHCHAIRS IN CONFORMITY WITH STANDARD ISO 7176-19 (see fig.2, page 14)
- THE WHEELCHAIR MUST BE ANCHORED TO THE VEHICLE WITH THE BACKREST AND THE LEGRESTS IN THE UPRIGHT POSITION using: 1) a four-point strap-type tie down system that complies with ISO 10542 supplied with the vehicle (see fig. 2, page 14) and 2) our accessory 891 SET OF 4 TIE-DOWN HOOKS (4 rings) (see fig. 1, page 13, fig.2 page 14 and fig.3 page 15)
- THE USER MUST WEAR BOTH A LAP BELT AND A DIAGONAL SHOULDER BELT ISO 7176-19 or ISO 10542 APPROVED (i.e. belts with a label of conformity to these ISO standards). THE USE OF A LAP BELT BY ITSELF IS NOT RECOMMENDED. Postural belts, such as model 894 are not type-approved for use as restraint systems and can only be used in addition.

#### **⚠** WARNINGS!

ORMESA RECOMMENDS THE USE OF A COMPLETE ISO 10542 COMPLIANT "WTORS" SYSTEM, which consists of restraint systems for both the wheelchair and its occupant

- Whenever the user is transferred on a seat of the vehicle, the unoccupied wheelchair should be stored in a cargo area or secured in the vehicle during the travel
- In order to preserve their effectiveness, **belt restraints should not be held away from the body by wheelchair components or pats,** such as the wheelchair armrests or wheels (see fig. 3 and 4, page 15)
- shoulder-belt restraints should fit over the shoulders (see fig. 3 page 15)
- belt restraints should be adjusted as tightly as possible, consistent with user comfort and belt webbing should not be twisted when in use
- in order to reduce the potential of injury to vehicle occupants, during transport the tray and all accessories should be removed and secured separately in the vehicle

THE WHEELCHAIR MUST BE INSPECTED BY THE MANUFACTURER OR ITS AUTHORIZED RETAILER OR DISTRIBUTOR BEFORE REUSE FOLLOWING INVOLVEMENT IN ANY TYPE OF VEHICLE COLLISION


ANY ALTERATION OR SUBSTITUTION OF THE 891 TIEDOWN HOOKS OR WHEELCHAIR PARTS OR COMPONENTS, WITHOUT THE WRITTEN AUTHORIZATION OF ORMESA SRL, IS FORBIDDEN AND WILL VOID THE "C" MARK AND THE WARRANTY

### (9

## 1) HOW TO FIX THE 891 TIEDOWN (4 RINGS) ON THE FRAME only for B30 and B60 version

#### WARNINGS:


- Suitable for ISO 10542-approved four-point straps provided in the vehicle. To anchor the wheelchair to the transport vehicle see the instructions on page 14
- The User has to wear both a lap belt and a diagonal shoulder belt ISO 7176-19- or ISO 10542-approved (i.e. belts with a label of conformity to these ISO standards). Postural belts such as model 894 are not approved for use as restraint systems and can be only used in addition
- The vehicle must be suitable for transporting orthopedic wheelchairs and pushchairs in conformity with current laws
- Read the notices and instructions provided on page 4 of this use and maintenance manual.
   The 891 tiedown hooks are safety devices and, as such, they must be fixed to the 869 frame by skilled personnel, such as the manufacturer or reseller that supplied the wheelchair.





### 2) HOW THE WHEELCHAIR IS TO BE SECURED IN A VEHICLE (private cars, bus, etc..)

HOOK THE APPROVED 4-POINT BELTS OF THE VEHCLE TO THE 891 4 TIEDOWN HOOKS (4 RINGS) AS SHOWN IN THE FIGURE AND IN THE DIRECTION OF TRAVEL OF







### 3) CORRECT POSITIONING OF THE OCCUPANT BELT RESTRAINTS ON THE USER

**ATTENTION!** The pelvic belt restraint should be worn low across the front of the pelvis, so that the angle of the pelvic belt restraint is within the preferred zone of 30° to 75° to the horizontal (see fig. 3). A steeper (greater) angle within the preferred zone is desirable

ATTENTION! Postural belts such as 894 are not intended to be as an occupant restraint in a moving, but they can be used in addition






#### THE INCORRECT PLACEMENT OF BELT RESTRAINTS


fig. 4



#### LABELS AND PLATES



THE **CC** MARKING CERTIFIES THAT JUDITTA CONFORMS TO THE SAFETY REQUIREMENTS defined with the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 applicable for MEDICAL DEVICE




This Medical Device is identified by a UDI (Unique Device Identifier), in compliance with Regulation (EU) 2017/745 and Code of Federal Regulations (CFR, US) Title 21. The UDI is a standardized code that allows the traceability and unique identification of the device throughout its life cycle.




#### **LABELS AND PLATES**



SECUREMENT POINT TO FOR THE ANCHORAGE OF THE ISO 10542-1 APPROVED 4-POINT STRAP-TYPE TIEDOWNS SUPPLIED BY THOSE VEHICLES WHICH ARE SUITABLE FOR TRANSPORTING ORTHOPEDIC WHEELCHAIRS IN COMPLIANCE WITH THE ISO 7176-19 STANDARD











#### ABOUT JUDITTA WHEELCHAIR. INTENDED USE

JUDITTA is a tilting posture system that allows the user to be positioned in a sitting posture up to a semi-lying position.

Intended for users between 140 cm (55,1") and 190 cm (74,8") height, suffering from disabilities resulting from neurological diseases of various kinds, systemic diseases, outcomes of cranio-encephalic trauma, brain neoplasms, non-responsiveness, geriatric diseases.

It comfortably and appropriately welcomes people who have postural control deficits, and who need prolonged use during the day.

It helps prevent the formation of decubitus, makes bed-wheelchair transfers easy and vice versa, and allows the care giver to frequently alternate the postural position and the distribution of the user's loads without effort.

The conformation and reclining of the backrest, the seat tilting and the reclining of the leg rests allow the user to be positioned lying down, semi-lying down or seated, while maintaining constant trunk-pelvis containment.

THE INCLINATION OF THE SEAT, THE CONFORMATION OF THE BACKREST (TRIANGULAR) AND THE LOAD-BEARING PADDING PREVENT SLIDING FORWARD AND FAVOUR TRANSPIRATION.

THIS ARRANGEMENT ALLOWS **THE TRUNK TO REST ON SEVERAL POINTS**, predominantly lateral with respect to the **spinal column which remains excluded from rubbing and compression** 

THE TRUNK IS KEPT IN POSITION without particular mechanical lateral supports or containment belts.

The **HEADREST** completes the user's posture.

The product must be prescribed by a specialist doctor, configured and adjusted by a healthcare professional authorized by the NHS as in the following paragraph.

To find the best postural structure, the positioning of the user must be done by a healthcare professional (Physiotherapist, Occupational Therapist, Orthopedic Technician), in order to ensure maximum comfort and maximum support for the user's needs.



Once the optimal assessment of the device has been established, the various transfers can also be carried out by non-professional personnel (the family member, for example).

#### **ENVIRONMENT OF USE AND TYPE OF USERS**

The B12 version is intended for indoor use only; the B30 version can be used both in internal and external contexts, as well as the B60 version which, however, can both allow a minimum autonomous movement by the user and allow smoother movement, especially outdoors, by the caregiver.

To find the best postural assessment, the positioning of the user must be done by a health professional (Physiotherapist, Occupational Therapist, Orthopedic Technician), in order to ensure maximum comfort and maximum support for the user's needs.

Once the optimal assessment of the device has been established, the various transfers can also be carried out by non-professional personnel (the family member, for example).

#### CONTRAINDICATIONS

The use of Juditta is not indicated if the user has a particularly significant residual ability to use the upper limbs and has control abilities, for these patients we may suggest the use of an ultra-light wheelchair (the patient with low spinal cord injury and consequent paraplegia).

If the user reports a risk of onset of skin lesions from high decubitus, identified according to the well-known international protocols (Braden, Norton, Waterflow scales, etc.), it is recommended to use an appropriate anti-decubitus device (bubble air pillow for example).





#### **TECHNICAL FEATURES and MODELS**

#### JUDITTA is available in 3 models:


- **JUDITTA B30 (tilt-in-space):** with 20 cm (8") front solid swiveling wheels and 30 cm (12") rear wheels; drum brakes are operated by a single control pedal.
- JUDITTA B12 (tilt-in-space): with four 12 cm (5") solid swiveling wheels, the rear ones with brakes.
- **JUDITTA B60 (tilt-in-space**): with 20 cm (8") front solid wheels and 56 cm (22") rear pneumatic wheels with hand brake and anti-tip device.
- FRAME made of chromed steel, painted with non-toxic epoxy powder paint. Narrow base to facilitate entering elevators.
- **SEAT**: adjustable inclination from 10° to 35°; available in four widths (36, 40, 45 and 50 cm 15,7"-16,9"-18,1"-19,2"); inclined by 10° to prevent the user sliding forward.
- **HEADREST** wide and wrap-around, padded, height-adjustable, with vertical and front-back rotation, and lateral translation.
- BACKREST wide and retentive, can be adjusted in four positions to vary the depth of the seat bottom
- SEAT low, to facilitate the user getting on and off
- GAS SPRING to adjust the inclination of the seat and backrest. The gas spring is a mechanism that supports the patient's weight during the adjustments of the backrest and of the seat tilting, with the following advantages: a) IT REDUCES THE ASSISTANT'S EFFORT IN RAISING AND LOWERING SIGNIFICANTLY; b) GUARANTEES UNIFORM MOVEMENT WITHOUT JERKING AND AVOIDS UNBALANCING THE PATIENT DURING ADJUSTMENT; c) ALLOWS STOPPING THE INCLINATION WHEREVER DESIRED
- UPHOLSTERY fireproof, breathable and removable for machine washing
- SIDES, closed, they can be opened to facilitate the user getting on and off
- ARMRESTS padded and height-adjustable
- SIDE PROTECTIONS padded for the knees
- **FOOTRESTS** foldable and adjustable in height and flexion extension.
- LEG SUPPORTS REMOVABLE, adjustable in inclination. They can be turned sideways to help the operator get close, the user get on and off. It is provided with padded calf rests.







<sup>\*</sup> footrest adjusted to the maximum length



# DIMENSIONS: in cm and Inches


# WEIGHT: in kilos and pounds


**User height:** from 140 to 190 cm (from 55.1" to 74.8")

JUDITTA B12 (Tilt in space)	size 36	size 40	size 45	size 50			
A Seat bottom width	36 cm (14,1")	40 cm (15,7")	45 cm (17,7")	50 cm (19,6")			
<b>B</b> Seat bottom depth	40	40-43-46-49 cm (15,7"-16,9"-18,1"-19,2")					
C Backrest width	47 cm (18,5")	47 cm (18,5")	52 cm (20,4")	57 cm (22,4")			
<b>D</b> Overall width	60 cm (23,6")	60 cm (23,6")	65 cm (25,5")	70 cm (27,5")			
H Seat bottom inclination (tilt)	from 10° to 35°	from 10° to 35°	from 10° to 35°	from 10° to 35°			
I Backrest inclination	from 10° to 38°	from 10° to 38°	from 10° to 38°	from 10° to 38°			
Aid weight	37 kg (81,5 lb)	37 kg (81,5 lb)	39 kg (85,9 lb)	40 kg (88,1 lb)			
Max. capacity	120 kg (264,5 lb)	120 kg (264,5 lb)	130 kg (286,6 lb)	140 kg (308,6 lb)			

Hyphen-separated numbers indicate the minimum and maximum range of adjustment







\* footrest adjusted to the maximum length


# DIMENSIONS: in cm and Inches

# WEIGHT: in kilos and pounds

**User height:** from 140 to 190 cm (from 55.1" to 74.8")

JUDITTA B30 (Tilt in space)	size 36	size 40	size 45	size 50			
A Seat bottom width	36 cm (14,1")	40 cm (15,7")	45 cm (17,7")	50 cm (19,6")			
<b>B</b> Seat bottom depth	40	40-43-46-49 cm (15,7"-16,9"-18,1"-19,2")					
C Backrest width	47 cm (18,5")	47 cm (18,5")	52 cm (20,4")	57 cm (22,4")			
<b>D</b> Overall width	62 cm (24,4")	62 cm (24,4")	67 cm (26,3")	72 cm (28,3")			
H Seat bottom inclination (tilt)	from 10° to 35°	from 10° to 35°	from 10° to 35°	from 10° to 35°			
I Backrest inclination	from 10° to 38°	from 10° to 38°	from 10° to 38°	from 10° to 38°			
Aid weight	38 kg (83,7 lb)	38 kg (83,7 lb)	41 kg (90,3 lb)	42 kg (92,5 lb)			
Max. capacity	120 kg (264,5 lb)	120 kg (264,5 lb)	130 kg (286,6 lb)	140 kg (308,6 lb)			

Hyphen-separated numbers indicate the minimum and maximum range of adjustment




size 36 size 40 size 45 size 50 JUDITTA B60 (Tilt in space) A Seat bottom width 36 cm (14,1") 40 cm (15,7") 45 cm (17,7") 50 cm (19,6") **B** Seat bottom depth 40-43-46-49 cm (15,7"-16,9"-18,1"-19,2") 47 cm (18,5") C Backrest width 47 cm (18,5") 52 cm (20,4") 57 cm (22,4") 74 cm (29,1") D Overall width 69 cm (27,1") 69 cm (27,1") 79 cm (31,1") from 10° to 35° from 10° to 35° from 10° to 35° from 10° to 35° **H** Seat bottom inclination (tilt) I Backrest inclination from 10° to 38° from 10° to 38° from 10° to 38° from 10° to 38° 40 kg (88,1 lb) 40 kg (88,1 lb) 42 kg (92,5 lb) 43 kg (94,7 lb) Aid weight 120 kg (264,5 lb) 120 kg (264,5 lb) 130 kg (286,6 lb) 140 kg (308,6 lb) Max. capacity

Hyphen-separated numbers indicate the minimum and maximum range of adjustment



### INITIAL PREPARATION



#### **PACKING INFORMATION**

### cm and kg:

JUDITTA B30	A cm	B cm	C cm	VOLUME m³	WEIGHT kg
Size 36	80	71	83	0.47	44
Size 40	80	71	83	0.47	44
Size 45	80	71	83	0.47	45
Size 50	80	71	83	0.47	46

JUDITTA B12	A cm	B cm	C cm	VOLUME m³	WEIGHT kg
Size 36	80	71	83	0.47	44
Size 40	80	71	83	0.47	44
Size 45	80	71	83	0.47	45
Size 50	80	71	83	0.47	46

JUDITTA B60	A cm	B cm	C cm	VOLUME m³	WEIGHT kg
Size 36	80	71	83	0.47	44
Size 40	80	71	83	0.47	44
Size 45	80	71	83	0.47	45
Size 50	80	71	83	0.47	46

### inches and pounds:


JUDITTA B30	A Inches	B Inches	C Inches	VOLUME m3	WEIGHT Lb
Size 36	31,4	27,9	32,6	32,6	97
Size 40	31,4	27,9	32,6	32,6	97
Size 45	31,4	27,9	32,6	32,6	99,2
Size 50	31,4	27,9	32,6	32,6	101,4

JUDITTA B12		B Inches		VOLUME m3	WEIGHT Lb
Size 36	31,4	27,9	32,6	32,6	97
Size 40	31,4	27,9	32,6	32,6	97
Size 45	31,4	27,9	32,6	32,6	99,2
Size 50	31,4	27,9	32,6	32,6	101,4

JUDITTA B60	A Inches	B Inches	C Inches	VOLUME m3	WEIGHT Lb
Size 36	31,4	27,9	32,6	32,6	97
Size 40	31,4	27,9	32,6	32,6	97
Size 45	31,4	27,9	32,6	32,6	99,2
Size 50	31,4	27,9	32,6	32,6	101,4





#### HOW TO REMOVE JUDITTA WHEELCHAIR FROM THE BOX






#### ASSEMBLY
























#### **HOW TO FIX SELF-PROPELLING WHEELS (B60 version)**




to remove the self propelling wheels, do the inverse procedure



#### **ACTIVATION of the ANTI-TIP DEVICE (B60 version)**




To deactivate the anti-tip device reverse the procedure. See the following Figure.






#### **USING JUDITTA WHEELCHAIR**


#### **BRAKING THE WHEELCHAIR (B30 MODEL)**









#### **BRAKING THE WHEELCHAIR (B12 MODELS)**









#### ADJUSTING THE TILT OF THE BACKREST








#### ADJUSTING THE INCLINATION OF THE SEAT BOTTOM (TILT-IN-SPACE)






#### **CLOSING THE FOOTRESTS**



#### **ADJUSTING THE FOOTRESTS IN HEIGHT**


- 1. UNSCREW THE NUTS AND REMOVE THE SCREWS USING THE SUPPLIED WRENCH
- 2. POSITION THE FOOTRESTS AT THE DESIRED HEIGHT
- 3. INSERT THE SCREWS AND TIGHTEN THE NUTS






# FLEX. EXTENSION ADJUSTMENT OF THE FOOTRESTS

- 1. LOOSEN THE SCREW USING THE WRENCH
- 2. ADJUST THE FOOTRESTS IN FLEX. EXTENSION
- 3. TIGHTEN THE SCREW






#### **TURNING THE LEGRESTS**






# **REMOVING THE LEGRESTS**






#### **ADJUSTING THE LEGREST INCLINATION**






# TURNING THE CALF RESTS






## **DEPTH ADJUSTMENT of the CALFRESTS**




## **FLESSO-EXTENSION of the CALFRESTS**






#### ADJUSTING THE HEIGHT OF THE ARMRESTS





#### **OPENING THE SIDE SUPPORTS**

**∕!**∖ WARNING! !




- CLOSE THE SIDE JUST ROTATING UP TO THE HORIZONTAL POSITION. A CLICK WILL CONFIRM THE RIGHT LOCK
- FOR SAFETY REASONS, THE SIDE HAS FRICTION (BRAKING) OPENING; CHECK ITS CORRECT FUNCTIONING PERIODICALLY (AS SPECIFIED ON PAGE 86)
- BRING THE SIDE TO THE OPEN POSITION AND MAKE SURE IT DOES NOT RETURN TO THE CLOSED POSITION IF NOT PUSHED WITH THE HAND




#### **ADJUSTING THE HEADREST**

# MWARNINGS!

- FIRST ADJUST THE HEADREST IN SIDE TRANSLATION, THEN GO ON WITH THE OTHER ADJUSTMENTS
- USE BOTH HANDS TO ADJUST THE HEADREST










#### ADJUSTING THE DEPTH OF THE SEAT

You can adjust the seat depth by adjusting:

I the two screws on the sides


I

- II the central bush under the seat
  For a right positioning see the chart on page 49
  - 1. LOOSEN THE TWO SCREWS ON THE RIGHT AND LEFT SIDES
  - 2. MOVE THE BACKREST IN THE DESIRED POSITION
  - 3. TIGHTEN THE TWO SCREWS






- II
- 1 OPEN THE SPRING OF THE PIN OF THE CENTRAL BUSH UNDER THE SEAT
- 2 PULL THE PIN
- 3-4 MOVE THE CENTRAL BUSH (THEN THE BACKREST) AS SPECIFIED ON THE CHART ON PAGE 49








# FOLLOW THE INSTRUCTIONS IN THE CHART A DIFFERENT POSITIONING FROM THE INDICATED ONE AFFECTS THE POSTURE AND THE STABILITY OF THE WHEELCHAIR






#### HEIGHT ADJUSTMENT OF THE BACKREST

- 1. LOOSEN THE FOUR SCREWS BEHIND THE BACKREST
- 2. ADJUST THE BACKREST HEIGHT
- 3. TIGHTEN THE FOUR SCREWS






## TILT ADJUSTMENT OF THE PUSH HANDLE






#### REMOVING THE UPHOLSTERY FOR THE WASHING





#### FOR THE WASHING, FOLLOW THE INSTRUCTIONS GIVEN ON PAGE 84





#### REAR VIEW OF THE BACK









## **ASSEMBLING THE ADDITIONAL COMPONENTS**

## RECLINABLE TABLE 824-R




4. insert the cap and the try lock in the tube of the armrest using a rubber hammer












# **B** INSERTION AND FIXING



- A. insert the tray so that the top is horizontal
  B. lock the indicated lever
- C. Before changing the inclination of the seat bottom (tilt-in-space) MAKE SURE THAT THERE ARE NO HOT OBJECTS or LIQUID CONTAINERS on the tray, that could cause







#### **C** ADJUSTMENT

## ADJUSTING THE DEPTH OF THE TABLE

- 1) press the safety lever
- 2) open the locking lever pressing the safety lever down (repeat the operation on the other side)
- 3) adjust the tray in depth
- 4) lock the right and left locking levers










## **D** SETTING OF CONTROL LOCKING LEVERS


A

should you have noticed an INEFFICIENCY, even partially, on the LOCKING of the GREEN LEVERS of the TRAY adjustments, please contact your ORTHOPAEDIC WORKSHOP to ADJUST IT WITH LITTLE INCREASES, as shown here below, up to the **ACHIEVEMENT OF THE PROPER FASTENING** 

#### depth adjustment lever




## tilt adjustment lever






# ABDUCTION BLOCK 834






# FOOTREST LOCK 914






# **B** UNLOCKING the footrest lock




# **C** - OPERATING the footrest lock






# CALF REST 916










#### REMOVABLE FOOTPLATE COVER 892






# REMOVABLE TABLE COVER 896



#### VIEW FROM BELOW OF THE TABLE






#### HEADREST REMOVABLE COVER 936



#### REAR VIEW OF THE HEADREST




# TIEDOWN HOOKS (4 RINGS) 891






JUDITTA WHEELCHAIR conforms to the ISO 7176-19 and ANSI RESNA WC/19 standard, which defines the safety requirements for transport, facing in the direction of travel on moving vehicles (private cars, buses, etc.)




FOR THE ASSEMBLY OF THE 891 ACCESSORY AND ANCHORING THE WHEELCHAIR TO THE VEHICLE, CARE-FULLY READ THE INSTRUCTIONS AND NOTICES ON PAGES 13, 14 AND 15 OF THIS MANUAL

**REAR VIEW** 














### SET OF HAND BRAKES 905 (B30 version)

### ASSEMBLE THE BRAKE LEVER AND THE DOUBLER














## ASSEMBLE THE RIGHT AND LEFT BRAKES, THE CABLE-LOCK PLATES AND THE WHEELS







### SET OF HAND BRAKES 905 (B60 version)


















## 932 ANTI-TIP UP DEVICE (B30 version)











### F ACTIVATION of the ANTI-TIP DEVICE



To deactivate the anti-tip device reverse the procedure. See the following Figure.







#### **BEWARE!**

In case of HYPERKINETIC patients, the Anti-tip up device, art. 932, is recommended when JUDITTA is used tilted and reclined.

 $\underline{\text{Outside this case,}} \text{ it is recommended to } \underline{\text{deactivate}} \text{ the device, which would make the handling of the wheelchair more difficult.}$ 



### 45° PELVIC HARNESS 894





### **MAINTENANCE**

The execution of all maintenance operations is necessary to maintain the correct functionality and safety of the device.

#### 1. ORDINARY MAINTENANCE (monthly)

JUDITTA is a mechanical device and as such must be kept in efficiency to maintain a high level of performance. This is possible with regular maintenance, including cleaning, lubrication and inspection. Routine periodic checks and maintenance should be carried out by a person with good technical competence; otherwise contact the health professional who provided it, or company specializing in maintenance.

<u> </u>	Maintenance and replacement of parts or additional components must be carried out without the user sitting in the device.
<u> </u>	Intensify all checks in marine environments, clean more often, anoint the parts exposed to salt corrosion (such as chrome, bolts and screws).
<u> </u>	Regular lubrication reduces friction between moving parts, preventing breakage or seizing and delays rust formation.

PART	DESCRIPTION/INTERVENTION	MODE
WHEELS	Clean the WHEELS from dust and dirt, to keep them smooth Check that the wheels turn freely In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component.	functional test
BRAKES	Check the effectiveness of the wheel locking/unlocking mechanism according to the instructions of the user manual.  In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component.	functional and visual test
BOLTS/SCREWS	Check the tightening of all bolts and screws Tighten them if loose.	functional test/tools
ADJUSTMENT AND LOCKING MECHANISM	Check the smoothness of the MOVING PARTS and keep them clean from dust and dirt, to avoid friction that could compromise the correct functioning In case of failure, lubricate them with commercial dry silicone oil following the instructions on the container label. After lubrication, dry the treated parts well with a soft cloth to remove excess residual grease  Check that the adjustment mechanisms are held in place.  In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component.	functional test
TILTING AND RECLINING SYSTEM	Check the functioning of gas springs and their cables. In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component.	functional test
SEAT'S VISCOELASTIC CORDS	Check the state of viscoelastic cords. In case of failure, discontinue use and contact the Healthcare Professional who supplied the device or the	visual test



	distributor for the replacement with original component.	
UPHOLSTERIES	Check the state of wear of the upholsteries. In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component.	visual test
PRODUCT INTEGRITY	Check that all components described in the user manual are present.     Check that there are no oxidized parts, and the uniformity of coating on the support elements or frame fastenings.  In case of failure, please contact the Healthcare Professional who supplied the device or the distributor for the replacement with original component	visual test/tools

#### 2. PREVENTIVE MAINTENANCE (biennial)

The Manufacturer recommends a COMPLETE REVIEW of the device, in order to verify it according to the specific use and to maintain the initial performances for all its lifetime.

This activity can be carried out by the Dealer who supplied the product or by an SERVICE CENTER, specialized in the maintenance of mechanical aids for people with disabilities and includes:

- 1) The general verification of the product, the integrity of the components and the locking of the mechanical parts, the smoothness of the moving parts
- 2) Performance verifications
- 3) The possible interventions aimed to restore the correct functionality
- 4) Mechanical control (in case of mechanical repairs)
- 5) Sanitization

The person who performed the maintenance is required to issue A REPORT WITH THE EVIDENCE OF THE INTERVENTIONS (VERIFICATIONS, REPAIRS, CONTROLS, SANITIZATIONS)

#### 3. EXTRAORDINARY MAINTENANCE (or corrective)

EXTRAORDINARY MAINTENANCE are all operations done on the device other than regular or preventive maintenance (mentioned above).

Extraordinary maintenance must be carried out by the Dealer who supplied the product or by another subject indicated by the Manufacturer or the Distributor of the ORMESA products in the country of destination.

The interventions must be performed with ORIGINAL REPLACEMENTS PARTS of the manufacturer. Modifications of the product ARE NOT ALLOWED, except for those of possible configurations provided on the brochure

For each maintenance work, the following must be performed:

- 1) The general verification of the product, the integrity of the components and the locking of the mechanical parts, the smoothness of the moving parts
- 2) Performance verifications
- 3) Mechanical control (in case of mechanical repairs)
- 4) Sanitization

THE MANUFACTURER or the AUTHORIZED PARTY must issue A REPORT WITH EVIDENCE OF THE INTERVENTIONS (CHECKS, REPAIRS, TESTS, SANITATION) CARRIED OUT.



### CLEANING, DISINFECTION, SANITIZATION

#### **INFORMATION**

**SANIFICATION** is the complex of cleaning and / or disinfection procedures and operations The **CLEANING** is a physical and mechanical process (i.e., rubbing) with which a large part of potential pathogenic microorganisms (bacteria, fungus, or virus), further to the visible dirt, is removed from the surface.

The combination of mechanical action with other factors such as the use of detergents (chemical action), temperature and duration can efficiently and sufficiently reduce the microbial load of the device.

**DISINFECTION**, after cleaning and cleansing, further reduces the number of microorganisms on a surface and it eliminates pathogenic microorganisms, ie bacteria causing disease and viruses. Products that on the label bear the authorization / registration Number of the Ministry of Health or other competent authority of an EU member State. are "disinfectants". Each disinfection action must always be preceded by a cleaning and cleansing operation, since dirt reduces the activity of the disinfectant. DRYING is essential because microbial growth can occur in the residual aqueous film

**CHEMICAL PRODUCTS WITH VIRUCIDAL-GERMICIDAL-FUNGICIDAL EFFICACY** are available on the market for hospital equipment. These products are effective in cold conditions, able to perform cleaning (elimination of dirt) and disinfection (elimination of pathogenic microorganisms, i.e. bacteria causing disease) in a single operation. and viruses).

Removable **UPHOLSTERIES** are washable respecting the symbols on the label. Washing is an extraordinary sanitation measure Replace the **REMOVABLE UPHOLSTERY** and PADDING when worn / difficult to sanitize.



#### WARNING!

I componenti 938 (WATERPROOF UPHOLSTERY) 916 (CALF REST) They are NOT WASHABLE! Only perform cleaning and disinfection with a sponge, as indicated in points 1) and 2)



#### WARNING!

The upholstery of the SIDES is HAND WASHABLE, according to the instructions on the washing label.

#### WARNINGS

**Read the TECHNICAL SHEET of the chemical product** to verify that **it is suited** to be used on chrome plated/varnished surfaces and on plastic components in PVC, PA, PP and to test on a small surface to ensure that it does not damage the aid.

For an effective operation, it is important to RESPECT THE POSOLOGY and TIMES OFACTION indicated on the product LABEL.

Use the PPE (gloves, FFP mask, visor, etc...) required on the product LABELP



The operations of sanitization must be performed without the user inside the device. Do not use compressed air, which can cause aerosol and contamination of possible virus and bacteria in the environment and on his own person.



During the Covid-19 epidemic, the contaminated PPE (gloves, gowns, glasses, masks, caps, etc...) must be **thrown away in the general waste bin**, unless otherwise indications by the belonging municipality



### **PROCEDURE**

Activities	Cadence	Description
4.1 CLEANING and CLEANSING	daily or weekly based on intensity of use and biological risk (patient with particular sweating, salivation; pandemic or endemic emergency period eg Covid-19)  Before disinfection	Soak a sponge or a clean disposable cloth with neutral detergent previously diluted in water. Rub upholstery, frame, accessories and finally wheels. Remove any traces of product by wiping with a clean damp sponge or cloth. Dry with a clean soft cloth. Proceed with any disinfection
4.2 DISINFECTION	Before Re-using the product with a new user	Spray a virucidal / germicidal / fungicidal chemical product for hospital equipment that is effective cold on a clean disposable cloth.  Rub the upholstery, the frame and the clean accessories (4.2) finally the wheels, until complete evaporation.



4.3 WASHING	Based on use	It is possib	le to wash the REMOVABLE
WASHING	24004 011 400	UPHOLSTERY of the SEAT and the	
		BACKRES	T, of the components
			rest), 892 (footplate cover), 834
			block), 896 (table cover)
		Respecting the following indications and the washing labels:	
		uic wasilli	g 140013.
			Machine wash with
		140	neutral soap .Max.
		_	temperature. 40°. Delicate cycle
			•
		X	Do not bleach
		$\overline{\Rightarrow}$	Do not iron
		$\sim$	Do not non
		$\boxtimes$	Do not tumble dry
			Dry horizontally.
			PADDED
			UPHOLSTERIES require
			EFFECTIVE DRYING IN A VENTILATED or
			HEATED
			ENVIRONMENT for at
			least 6-10 hours
		For on ANI	TIVIDAL DISINIFFCTANT action
		For an ANTIVIRAL DISINFECTANT action, a SPECIFIC CHEMICAL PRODUCT can	
			o the normal washing cycle;
			a high temperature (60 ° C) is
			s long as occasionally, as the may wear out.
		uprioistery	may wear out.
		WARNING	-
			removable padded cover of s of the leg rest inside out.
			onents <b>938</b> (waterproof
		upholstery)	916 (calf rest) They are NOT
		WASHABI	_E! Only perform cleaning and
		disinfection in points 1	on with a sponge, as indicated
		in points	i j aliu 2 j
4.4	Based on use		rbonate is a natural product of
SANITIZATION			igin that is commercially
WITH   PERCARBONATE			when dissolved in water, it
FERCARDONATE			ective oxygen already at 30 ° ectant, antibiotic and
		antibacteri	
		For a san	itizing action during the
		cleaning	of the frame and the washing
			ness and removable



upholsteries, you can add 1 teaspoon of percarbonate > 30% to the detergent:

- In the washing machine: add 1 teaspoon of percarbonate in the drum together with the detergent.

- By hand/ for the frame cleaning: dissolve 1 teaspoon of percarbonate in the bowl together with the detergent, proceed with the washing, and cleansing.

#### WARNING!

When washing / cleaning the fabric it is recommended not to mix sodium percarbonate with acids (for example: Vinegar, Lemon), as it could create chemical reactions that could damage it.

#### **ASSISTANCE AND SPARE PARTS**

For service and replacement parts, only contact the retailer where you purchased the product.

#### **DISPOSAL**



For the device disposal, contact a body authorized to recycle aluminum according to the laws in force.



# CONDITIONS FOR DURABILITY, REUSE AND REASSIGNMENT TO A NEW USER

Granted that Ormesa products should be selected, evaluated and ordered for the needs of an individual user, reuse is however possible with the respect of the following conditions.

Based on the experience of other similar sold models, on technological progress, on the guarantees of the Quality Management System certificated in 1998 according to ISO 13485, there is adequate confidence that the average lifespan of JUDITTA is about 5 years, on condition that it is used according to the directions given in the user manual.

Only when the conditions for storage and transport of chapter HOW TO STORE AND TRANSPORT THE MEDICAL DEVICE" are followed, the periods in which the product is stored at the dealer should not be considered in that time period

Factors unrelated to the product such as the development of the user, its diseases, the use and the surrounding environment can make significantly lower the duration of life of the product; on the contrary, if the indications on the use and maintenance are properly observed, the reliability of the product can extend well beyond the lifetime average above.

Prior to recycling or reassignment an already used Ormesa product, it is required that:

- a doctor or therapist verifies that the aid is appropriate and adequate to meet the dimensional, functional and postural needs of the new user, and if all its components are suitable / appropriate for him. You should also consider that the C€-marking and the manufacturer's responsibility for safety requirements for the product remain only if the original product still has not changed and only original accessories or spare parts have been applied
- 2. qualified technical personnel of a company specialized in the maintenance of technical aids for disabled people performs a detailed technical inspection to verify its condition and wear, the absence of any damage and failure of all components / adjustments, the presence of the user's manual, of the label with the date and serial number. A copy of the manual and maintenance may be always requested to the retailer that supplied the product or directly to Ormesa
- 3. The product has been thoroughly cleaned and disinfected following the directions given in the "MAINTENANCE, CLEANING AND DISINFECTION" Chapter

We recommend to keep written records on all inspections performed on the product before any assignment to the new user.



In case of doubt about the safety of the product or damage to parts or components, you are urged to immediately discontinue use and contact the Dealer who supplied you with the product; the Distributor or the Manufacturer are at your disposal for any further doubt or assistance.





Manufactured by Ormesa srl Via delle Industrie, 6/8 - Z.I. Sant'Eraclio - 06034 FOLIGNO (PG) ITALY
P.I. IT 00574020541 CCIAA Perugia 119215 Iscr. Trib. Perugia 11907
Tel. + 39 0742 22927 Fax +39 0742 22637 info@ormesa.com

#### www.ormesa.com



## JUDITTA MOD. B30 AND B.60 PASSED CRASH TEST

ACCORDING TO THE ISO 7176-19 AND ANSI RESNA WC/19

FDA REGISTERED with Listing number A980195

#### Imported and distributed by:



3190 F.-X.-Tessier, Vaudreuil-Dorion, QC J7V 5V5, Canada 1360 Boltonfield St., Columbus, OH 43228, USA T 1450 424-0288 | 1888 453-0311 F 1450 424-7211 | 1877 501-8458 www.amylior.com | info@amylior.com



OM | P25001 | 1005590 | R01

Authorized	Dealer

The technical Data and information given in this handbook are not binding. ORMESA reserves the right to carry out any necessary change without prior notice or replacements