POWER WHEELCHAIR TRANSPORTATION SAFETY GUIDANCE DOCUMENT

TABLE OF CONTENTS

1. Introduction	3
2. Transportation of a wheelchair in a motor vehicle	4-6
3. AMYLIOR wheelchairs approved for transportation in accordance with ISO 7176-19.....	7
4. History of applicable standards	8
5. Crash Test photos	9
6. Securement points for wheelchair tie-down systems on AMYLIOR products	10
7. Checklist: Transporting a person in a wheelchair	10-11
8. Manufacturer addresses of WTORS (Wheelchair Tie-down and Occupant Restraint Systems)	12

1. Introduction

This document contains information on **AMYLIOR** crash tested power wheelchairs and describes the applicable standards to which they were tested. It also includes information on transit securement of an occupied wheelchair in a motor vehicle, on tie-down systems used during crash tests, as well as images of wheelchair securement points.

Based on tests performed and results achieved, approved for transit **AMYLIOR** products withstand strong forces that occur during crash tests (in accordance with ISO 7176-19) and are therefore suitable to be used as a forward-facing seat for the transportation of a wheelchair and its user in a motor vehicle.

A simulated crash test of a frontal impact performed in a test laboratory can in no way cover all accidental situations that occur on the road. Therefore, it is recommended that the wheelchair user transfers into the vehicle seat which offers, without a doubt, the highest level of safety while the vehicle is in motion; and that the wheelchair be securely stowed in a designated area.

Because of the high number of systems offered on the market, it is not possible for **AMYLIOR** to crash test all Wheelchair Tie-down and Occupant Restraint Systems (WTORS). However, all **AMYLIOR** crash tested wheelchairs, can be transported with any ISO 10542-certified WTORS.

In order to determine whether the system you have meets these requirements, contact the wheelchair manufacturer and/or the occupant restraint system supplier (see Page 10).

This document is intended to provide transit compliance information and guidance for the safe use of tie-down systems and securement points. If you have questions on using your wheelchair for seating in a motor vehicle, please contact **AMYLIOR** Technical Support at **1 888 453-0311**.

2. Transportation of a wheelchair in a motor vehicle

A wheelchair secured in a vehicle does not offer the same level of safety as the bolted-in seating system of a vehicle. Therefore, whenever feasible, **AMYLIOR** recommends that the wheelchair user transfers or is transferred into the vehicle seat and uses the vehicle-installed restraint systems (seat belt) and that the unoccupied wheelchair is stored in a cargo area or secured in the vehicle during travel. If this is not possible and the wheelchair occupant must remain in the wheelchair while traveling in a vehicle, then these instructions must be followed:

1. Before entering a vehicle, all power options must be back to their original position. This implies that the wheelchair is without elevation and/or tilt, fully retracted center mount and backrest in an upright position.


Ease of access to, and maneuverability in, motor vehicles can be significantly affected by wheelchair size and turning radius, and that smaller wheelchairs and/or wheelchairs with a shorter turning radius will generally provide greater ease of vehicle access and maneuverability to a forward-facing position.

2. The wheelchair and occupant must be forward-facing in the direction of travel and must be secured using the wheelchair **tie-down system**. The user's **occupant restraint system** must meet the ISO 10542 or SAE J2249 requirements. These systems must be fitted in accordance with their manufacturer's instructions as well as the wheelchair manufacturer's instructions (See checklist on page 9).

Assuming the wheelchair conforms with the requirements of this document and, as such, has been designed and tested for use only as a forward-facing seat in a motor vehicle.

The test standard in accordance with ISO 7176-19 only covers the transportation of wheelchairs in a forward-facing direction. This means that the wheelchair must never be transported in a side facing direction. (Fig. 1).

Fig. 1

POWER WHEELCHAIR TRANSPORTATION SAFETY GUIDANCE DOCUMENT


3. All add-on components or accessories should be removed from the wheelchair and stowed safely. These include but are not limited to:

Loose cushion – Tray Table – Vent Tray – Crutches – Hip and Thigh Supports – Bags

4. Alterations or substitutions must not be made to the wheelchair securement points or to the chassis/frame components without consulting the manufacturer. If modifications are made without the consent of the manufacturer, the wheelchair is no longer adequate to be transported in a vehicle and the manufacturer warranty is void.
5. When transported in a vehicle, a power wheelchair must be equipped with leak-proof, sealed batteries (e.g. gel batteries).
6. Should there be an accident or impact, the wheelchair must be inspected by an authorised wheelchair dealer before it is used again.
7. A combination of a pelvic-belt and a shoulder-belt must be used to restrain the wheelchair occupant (Fig. 2). The shoulder-belt restraint must not lie across the neck or underneath the armpit but across the chest. The pelvic-belt restraint should be positioned just above the hipbone.

The use of only a pelvic belt may compromise the performance of the wheelchair tie-down and occupant restraint system (WTORS), and increase the risk of serious or fatal injuries to the wheelchair occupant.

Fig. 2

8. A headrest suitable for transportation should be properly fitted and positioned at all times during transportation in a vehicle.
9. Standard wheelchair lap belts are not suitable.


Postural supports (lap straps, lap belts) should not be used or relied upon as occupant restraints in a moving vehicle unless they are labelled as meeting the requirements specified in RESNA WC-4, Section 18 (with successful performance in a 48 km/h (30 mph) frontal impact test).

OCCUPANT RESTRAINT INSTRUCTIONS

Before fitting occupant restraint, make sure that all power options are back to their original position i.e. the wheelchair is without elevation and/or tilt, fully retracted center mount or legrests and backrest in an upright position. This position should also be used when the wheelchair is in a motor vehicle.




1. The wheelchair tie-down and occupant restraint systems (WTORS) must be fitted on the lower front side of the hipbone (Fig. 5) in a way that the angle of the pelvic-belt restraint is within the range of 30° to 75° to the horizontal (Fig. 3).

A steeper (greater) angle of the belt is desirable but must never exceed 75°.

NOTE: Only the designated securement points on the wheelchair should be used to attach it to the vehicle. The wheelchair has been tested with and must be secured by manufacturer-designated wheelchair-mounted belt restraints.

2. Restraint belts must not be placed over wheelchair parts or obstructed by wheelchair components such as armrests or wheels so that they are being held away from the body (Fig. 4).
3. The shoulder-belt restraint must be fitted diagonally across the chest over the collar bone (Fig. 6).
4. The occupant restraint belts must be as tight as possible, without restricting the user.
5. Restraint belts must not be twisted when in use.
6. Exercise caution when securing the occupant restraint, ensuring that the seatbelt buckle is positioned in a way that, in the event of a crash, a wheelchair component will not contact the release button.
7. Suitable headrests (strongly recommended) should be positioned correctly (Fig. 5).

3. AMYLIOR power wheelchairs approved for transportation in accordance with ISO 7176-19

AMYLIOR wheelchairs are tested in accordance with ISO 7176-19 with a hybrid III test dummy of 100 kg (220 lb), forward-facing in the direction of travel with a frontal impact and a wheelchair tie-down system consisting of a 4-point tie-down system with a 3-point vehicle anchored belt occupant restraint system.

Some of the variants of the models tested were assessed and deemed that they fulfill the dynamic test requirements because they have the same construction (i.e. the same frame resistance cover materials, rigidity, structural integrity of the components and connections) as well as geometrical similarity with the models which have been approved for transportation in a vehicle.

The products shown in the following table were tested in their standard configuration and with standard settings. It would be impossible to test **AMYLIOR's** entire range of adjustable configurations and options.

The **AMYLIOR** power wheelchairs tested with a modified configuration, must be **selected at the initial purchase** to comply with ISO 7176-19 standards.

The following list of crash tested **AMYLIOR** products expands constantly. Please refer to the most recent version of this table on the *WC Transportation Safety* website and click the tab "Crash Tested Product Lists".

AMYLIOR POWER WHEELCHAIRS ARE FULLY COMPLIANT WITH ISO 7176-19 STANDARDS AS OF JANUARY 2017

The following link features the latest ISO 7176-19 compliance information:

https://wc-transportation-safety.umtri.umich.edu/wp-content/uploads/sites/517/2024/03/UMTRI-Wheelchair-Compliance_March2024.pdf

Name and Model	Configuration and Options	User Weight Capacity Range	Wheelchair Weight
<u>Alltrack Signature Series</u> Alltrack M3 Series Alltrack R3 Series Alltrack P3 Series	<u>Signature Seat</u> • Static seat (no options) • With Power Tilt • With Power Recline • With Power Elevate • With Power Tilt & Recline • With Power Tilt & Recline & Elevate	51 lb – 300 lb (23 kg – 136 kg)	Up to 400 lb (181 kg)
<u>Alltrack Signature Series</u> Alltrack M3-HD Series Alltrack R3-HD Series	<u>Signature Seat</u> • Static seat (no options) • With Power Tilt • With Power Recline • With Power Elevate • With Power Tilt & Recline • With Power Tilt & Recline & Elevate	51 lb – 450 lb (23 kg – 205 kg)	Up to 430 lb (195 kg)
<u>Alltrack Signature Series</u> Alltrack HD3 Series	<u>Signature Seat</u> • Static seat (no options) • With Power Tilt • With Power Recline • With Power Tilt & Recline	51 lb – 550 lb (23 kg – 250 kg)	Up to 480 lb (218 kg)

4. History of applicable standards

1. **ISO 7176-19 Crash Test Standard** – The wheelchair crash test was derived from the crash test already used in the car industry. Wheelchairs are crash tested at a speed of 48 kph (30 mph) and an impact deceleration of 20 G. The test dummy is limited to a maximum weight of 100 lb (220 lb). Only frontal impact has so far been simulated for wheelchairs.
2. **Wheelchair Tie-down and Occupant Restraint System (WTORS) in accordance with ISO 10542-1** – There are various versions of WTORS that meet ISO 10542 requirements. The wheelchair is attached with loops, hooks, carabiner clips, buckle tongues or other means of securement. The WTORS itself is also tested by a dynamic crash test. To do this, a "surrogate wheelchair" is used. With this WTORS, particular attention is given to the weight that is being secured (wheelchair weight + user weight). For this reason, power wheelchairs (considered very heavy) must be secured in a vehicle using a 4-point tie-down system.

20 G FRONTAL CRASH

- ISO 10542-1, 3-point tie-down occupant restraint system;
- Person (100 lb/220 lb hybrid dummy) ISO 7176-19 Frontal Crash Test for wheelchairs as a seat in a vehicle (48 kph/30 mph, 20 G deceleration);
- ISO 10542-2 wheelchair tie-down system.


Occupant restraint system

Power wheelchair tie-down system

5. Crash Test Photos





These photos show you what happens during a wheelchair crash test.

6. Securement points for wheelchair tie-down systems on AMYLOR products

This section shows location of securement points for approved **AMYLOR** products.

In accordance with ISO 7176-19, all securement points on the wheelchair are labeled with the relevant standard hook symbol. These labels show where to attach the wheelchair tie-down system hardware (hooks, carabiner clips, loops, etc.)

SIGNATURE ALLTRACK MWD SERIES	SIGNATURE ALLTRACK HWD SERIES
4 securement points on the base (2 on each side)	
	
	

7. Checklist: Transportation of a person in a wheelchair

For those who transport and/or accompany the wheelchair user, we advise going through this checklist:

1. Whenever feasible, it is recommended that the wheelchair user transfers or is transferred into the vehicle seat and uses the vehicle-installed restraint systems (seat belt), and that the wheelchair be securely stowed in a designated area.
2. If the occupant must remain seated in the power wheelchair during transportation, consider the following on the next page.

POWER WHEELCHAIR TRANSPORTATION SAFETY GUIDANCE DOCUMENT

- a) The wheelchair must be tested in accordance with ISO 7176-19 requirements. All seating components added to the wheelchair, must be compliant with ISO 16840-4. **AMYLOR** wheelchairs that have been tested carry the relevant labels.

- b) Use only a suitable wheelchair tie-down and occupant restraint system (WTORS) in accordance with ISO 10542. This requires a 4-point vehicle tie-down system with a 3-point vehicle-anchored belt occupant restraint system.
- c) Make sure that all power options are back to their original position i.e. the wheelchair is without elevation and/or tilt, fully retracted center mount and backrest in an upright position.
- d) The wheelchair and user must be forward-facing in the direction of travel and symmetrically aligned over the vehicle's securement rails.
- e) When possible, all removable and loose parts (water bottles, crutches, etc.) or other auxiliary wheelchair equipment should be either secured to the wheelchair or removed from the wheelchair and stowed away in the vehicle during travel so that it does not break free and cause injury to vehicle occupants in the event of a collision. For a tray that is permanently installed on the wheelchair, there must be energy-absorbing padding placed between the tray and the occupant.
- f) If the wheelchair is equipped with freewheel or brake release levers, they need to be engaged.
- g) The wheelchair tie-down system is to be fitted to the wheelchair in accordance with the manufacturer's recommendations. The hook labels indicate the position of all securement points on the wheelchair. Attach the front straps on the wheelchair first. Release the levers. Pull the wheelchair back toward the rear to create tension and tighten the straps. Then attach the back straps making sure that they are tight by the wheelchair's traction. Re-engage the levers. The rear straps set the system under tension.
- h) After the wheelchair is firmly attached to the vehicle floor, the occupant restraint system should be fitted in accordance with the manufacturer's recommendations. When doing this, make sure that the shoulder-belt and pelvic-belt restraints are fitted securely. The belts must not be twisted nor pass over wheelchair components such as the sideguards, armrests or other edges, they must closely hug the occupant. The shoulder-belt should be positioned diagonally across the chest and over the collar bone away from the neck. The pelvic-belt must always be positioned just above the hipbone, not over the abdomen.
- i) A headrest, tested according to ISO 16840-4 standard, should be fitted and positioned correctly.
- j) Finally, belts and straps should be checked one more time to make sure that they are correctly positioned:
- Are the wheelchair straps firmly tensioned and fitted in the right place?
 - Is the pelvic-belt restraint positioned just above the hipbone?
 - Is the shoulder-belt restraint positioned diagonally across the chest and over the collar bone away from the neck?
 - Do the belts closely hug the body?

8. Manufacturers of WTORS (Wheelchair Tie-down and Occupant Restraint Systems)

Q'Straint Europe

72-76 John Wilson Business Park,
Whitstable Kent,
CT5 3QT, UK

Tel.: 44 (0) 1227 773035
qstraint.com

Q'Straint America

4031, NE 12th Terrace,
Oakland Park,
FL 33334, USA

Tel.: 1 954 986-6665
qstraint.com

Q'Straint Australia

Tramanco Pty Ltd.
21 Shoebury St.,
Rocklea, QLD 4106 Australia

Tel.: 61 7 3892 2311
qstraint.com

EZ Lock Inc.

2001 Wooddale Blvd.,
Baton Rouge,
LA 70806, USA

Tel.: 1 225 214-4620
ezlock.net

This document is available in a printable PDF on our website at amylior.com

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