

INSTRUCTION GUIDE

YOMPER+ EASYGO SETTINGS FOR PROFESSIONALS



THIS DOCUMENT IS STRICTLY RESERVED FOR HEALTHCARE PROFESSIONALS OR AUTHORIZED PROVIDERS



1. To access the Settings Menu, start the Yomper+ and the EasyGo control box, then wait a few moments for the Bluetooth pairing to complete.
2. Press & hold the main button and the right (-) button simultaneously for a few seconds (see P01).



P01

The menu for PROFESSIONALS will appear, giving access to the settings for **maximum speed** (2.1), **"Magnet"** (2.2), and **"Stalling"** (2.3).

2.1 Maximum speed setting: Select the maximum speed using the arrows to limit the Yomper (and the max. speed value in the USER menu) to the desired speed under no load (see P02).

A "2" must be entered at the end of the speed for it to be taken into account.

Example: 92 to limit to 9 km/h (5 ½ miles), 122 to limit to 12 km/h (7 ½ miles).



P02

Press the main button to go to the next setting (see P03).



P03

YOMPER 2.6 – EasyGo Settings for Professionals

2.2 “Magnet” Setting: Default number of magnets is 1, 12 for external magnets (see P04).



Press the main button to go to the next setting (see P05).



2.3 “Stalling” Setting: Internal setting is set to 6 by default. Do not modify unless explicitly instructed by the manufacturer (see P06).



Press the main button to go to the next setting (see P07).



3. Press the main button once again to confirm the changes and the main screen will appear (see P08).

