

INFORMATION BULLETIN

NEW PARTNERSHIP WITH ORMESA IN THE UNITED STATES

Pediatric/Adult Mobility and Seating Solutions



New Partnership with ORMESA Mobility Solutions

Expanding our mobility product offering!

INTRODUCTION

Amylior is proud to announce a new partnership with Ormesa, a globally recognized manufacturer of premium pediatric and adult mobility and seating equipment. This collaboration brings three of Ormesa's flagship products — the Bug, Grillo Gait Trainer, and Juditta — to the U.S. market. For more than 45 years, this Italian company has made children the focus of everything they do.

BRIEF PRODUCT PRESENTATION

Bug

The Bug is a safe, modular, lightweight, compact, and versatile stroller for children with special needs, featuring an integrated tilt-in-space seat. Highly engineered and adaptable, it offers a wide range of adjustments, making it a long-lasting mobility aid that grows with the child. The Bug can be mounted on four different bases: one for outdoor use and three for indoor use or placed securely on any flat surface.



Grillo Gait Trainer

The Grillo Gait Trainer is a versatile support walker for children and adults with reduced mobility, designed to encourage movement and promote independence. Available in both front- and rear-directic configurations, it offers exceptional customization and ease of use with a wide range of adjustments and support options. Its ergonomic, independently adjustable pelvic and thoracic supports provide optimal positioning and stability at every step.



New Partnership with ORMESA Mobility Solutions

Juditta

The Juditta is a tilting and reclining comfort wheelchair designed for adults and seniors with reduced mobility. Its backrest promotes proper posture, providing stability while helping to prevent contractures and the development of deformities. The seat's viscoelastic cords conform to the user's anatomy, reducing friction and lowering the risk of pressure injuries. Juditta's unique design also ensures breathability, hygiene, and lasting comfort.




RESOURCES

Our website contains all the information you need. By clicking on *Products*, then **Ormesa**, you will be on the landing page. On the page you can click on the product to see:

- Order Forms
- Owner's Manual(s)
- Brochure
- Picture Gallery

Videos are integrated with the descriptions



ADDITIONAL INFORMATION

Availability and Pricing

Both products are available immediately in the United States. Contact Customer Service for pricing.

Corresponding updates

For
Parts:

PARTS MANUAL	
CHOOSE YOUR SECTION	
SECTIONS	1 MODULES 1.1 - POWER TEST MODULES 1.2 - POWER BLUETOOTH MODULES 1.4 - ADJUSTABLE POWER STANDING MODULE 1.5 - SHOULDER 1.6 - LIGHT MODULES
	2 SEAT 2.1 - BACKS 2.2 - SEAT PANS 2.3 - ARMRESTS 2.4 - CENTER MOUNT FOOTBOARD 2.5 - LEG RESTS 2.6 - HEADRESTS 2.7 - POSITIONING 2.8 - BACKS AND CUSHIONS 2.9 - ACCESSORIES
	3 POWER BASE 3.1 - ATRACK SERIES
	4 INTERFACE BRACKETS 4.1 - AMTOSYSTEMS 4.2 - QUADROCK 4.3 - QUANTUM REHAB 4.4 - PIVOTABLE 4.5 - INNOVATION IN MOTION 4.6 - OTHER MANUFACTURERS
	5 USER INTERFACE BRACKET 5.1 - STANDARD JOYSTICK MOUNTING 5.2 - MODULAR MOUNTING 5.3 - SIMULTANEOUS DEVICE MOUNTING KIT 5.4 - LIGHT MOUNTS 5.5 - SWITCH MOUNTS
	6 POWER BASE ELECTRONICS 6.1 - VBE 6.2 - BARET 6.3 - SIGETS 6.4 - PROGRAMMERS 6.5 - BATTERIES 6.6 - JOYSTICK ACCESSORIES
	7 POWER SEATING ELECTRONICS 7.1 - ELECTRONICS
	9 MOBILITY DEVICES 9.1 - ROLLATORS
	10 ORMESA'S PRODUCTS For any parts requests, please contact info@amylior.com

10

ORMESA'S PRODUCTS

For any parts requests, please contact
info@amylior.com