

INFORMATION BULLETIN

Update Notice: Alltrack M3/P3 (MWD) Narrow N2221 EZ-Lock & Q'Straint

INTRODUCTION

We are pleased to announce that vehicle docking solutions for the Alltrack M3/P3 Narrow N2221 are now available for both Q'Straint and EZ Lock.

This compatibility improves the accessibility of our product, facilitates integration into different types of adapted vehicles, and strengthens our position in the market.

It also officially provides users with safe, validated securement options from these two companies, which were not previously available for this model.

WHAT'S NEW

EZ-LOCK

The bracket for the Alltrack M3/P3 (MWD) Narrow N2221 is available on [EZ-Lock](#) under the reference number: **BKT-AS-MWD-22-N**

Q'STRAIT

The brackets for the Alltrack M3/P3 (MWD) Narrow N2221 are also available through [Q'Straint](#) in two configurations: front row and second row.

- Front row = driver or front passenger, forward-facing dock – **Q04BR156**
- Second row = rear section, side-facing dock – **Q04BR155**

ADDITIONAL INFORMATION

The companies also used this project as an opportunity to add more information to their websites, making it easier to identify the correct bracket based on the wheelchair model and its base configuration.