

INFORMATION BULLETIN

June 2026

NEW SIGNATURE HARDWARE KIT

INTRODUCTION

We are excited to introduce a new dealer-level hardware kit (1005533) for the Signature Power Positioning System. Designed to complement the hardware already included with every Signature system (1005532), this new offering gives dealers even more flexibility and convenience when mounting aftermarket backs and accessories to Signature.

These kits make it easier to keep the right fasteners on hand for a variety of installation needs. Please note that the fasteners in these kits fit only in the extruded slots of the Signature back canes. T-nuts for our seat pan are a different size.

WHAT'S NEW

The new **Signature Hardware Kit 1005533** is packaged in a clearly labeled, transparent storage box for easy organization and quick access. It includes:

- M6 screws in 10 mm, 12 mm and 14 mm lengths (x 24 ea.)
- Single t-nuts (x 12)
- Short Double t-nuts (x 12)
- M6 washers (x 36)

Kit 1005533



The new kit (1005533) is designed to supplement **Hardware Kit 1005532**, which is included with every Signature power wheelchair and is also available through Customer Service. Packaged in a clear zip-lock bag, Kit 1005532 includes:

- M6 screws in 10mm, 12mm and 14mm lengths (x4 ea.)
- Single t-nuts (x4)
- Short Double t-nuts (x4)
- M6 washers (x4)

Kit 1005532



ADDITIONAL INFORMATION

To order your kit(s) or learn more, please contact our Customer Service department at info@amylor.com or our Technical Service department at techsupport@amylor.com.