

## Folding Aluminum Walker

This **ultra lightweight aluminum folding walker** has been designed to help people with reduced mobility in their daily activities.

**4081 Series - 7/8" frame**

The walker's sturdy folding frame is chrome plated to resist scratches and features two different wheel position options: exterior for increased stability or interior for limited motion space.

### Included features:

- Adjustable floor/wrist height
- Hypoallergenic, comfortable rubber grips
- Scratch resistant chrome plated folding frame
- Bearing wheels (smooth and quiet motion)
- Wheel stem inlet on both sides  
(to alternate interior/exterior wheel position without removing the wheels)



**For cases of ataxia, it is possible to add 6 lb (2.7 kg) to the walker's weight.**

	#ITEM	DESCRIPTION	HEIGHT	USER HEIGHT (approx.)
WHEELS / SKIS	4081AW5	ADULT with wheels / skis	33" – 37" (840 mm – 940 mm)	5'8" – 6" (1.73 m – 1.83 m)
	4081JW5	JUNIOR with wheels / skis	29" – 33" (735 mm – 840 mm)	5'4" – 5'8" (1.63 m – 1.73 m)
	4081PW5	JUNIOR PETITE with wheels / skis	26" – 30" (660 mm – 760 mm)	5'1" – 5'5" (1.55 m – 1.65 m)
RUBBER TIPS	4081A	ADULT rubber tips	33" – 37" (840 mm – 940 mm)	5'8" – 6" (1.73 m – 1.83 m)
	4081J	JUNIOR rubber tips	29" – 33" (735 mm – 840 mm)	5'4" – 5'8" (1.63 m – 1.73 m)
	4081P	JUNIOR PETITE rubber tips	25" – 29" (635 mm – 735 mm)	5' – 5'4" (1.52 m – 1.63 m)

### LIMITED WARRANTY

**5 years:** frame

**1 year:** every other component  
(except for wheels / skis)

SPECIFICATIONS	
Overall width (exterior wheel position)	26" (660 mm)
Overall width (interior wheel position)	22.75" (580 mm)
Overall width (rear-skis)	22.5" (570 mm)
Inner width (between handles)	17" (430 mm)
Weight	6.5 lb (3 kg)
Weight capacity	330 lb (149.7 kg)