

Bariatric Folding Walker

This **bariatric folding walker** has been designed to help people with reduced mobility in their daily activities.

The walker's reliably sturdy folding frame is chrome plated to resist scratches and features two different wheel position options: exterior for increased stability or interior for limited motion space.

Included features:

- Adjustable floor/wrist height
- Hypoallergenic, comfortable rubber grips
- Sturdy 5" x 1.75" wheels
- Reinforced, double pin adjustment tubes
- Chrome plated, scratch resistant folding frame (400 lb capacity)



4082 Series - 1" frame

For cases of ataxia, it is possible to add 9 lb (4.1 kg) to the walker's weight.

| | #ITEM | DESCRIPTION | HEIGHT | USER HEIGHT (approx.) |
|---------------|-----------|--------------------|------------------------------|-------------------------------|
| WHEELS / SKIS | 4082AOW5 | ADULT with wheels | 33" - 37" (840 mm - 940 mm) | 5'8" - 6" (1.73 m - 1.83 m) |
| | 4082A+OW5 | ADULT+ with wheels | 37" - 41" (940 mm - 1040 mm) | 6' - 6'4" (1.83 m - 1.93 m) |
| | 4082JOW5 | JUNIOR with wheels | 29" - 33" (735 mm - 840 mm) | 5'4" - 5'8" (1.63 m - 1.73 m) |
| RUBBER TIPS | 4082AO | ADULT rubber tips | 33" - 37" (840 mm - 940 mm) | 5'8" - 6' (1.73 m - 1.83 m) |
| | 4082A+OW5 | ADULT+ rubber tips | 37" - 41" (940 mm - 1040 mm) | 6' - 6'4" (1.83 m - 1.93 m) |
| | 4082JO | JUNIOR rubber tips | 29" - 33" (735 mm - 840 mm) | 5'4" - 5'8" (1.63 m - 1.73 m) |

LIMITED WARRANTY

5 years: frame

1 year: every other component
(except for wheels / skis)

| SPECIFICATIONS | |
|---|-------------------|
| Overall width (exterior wheel position) | 30.5" (775 mm) |
| Overall width (interior wheel position) | 26.5" (675 mm) |
| Overall width (rear-skis) | 25" (635 mm) |
| Inner width (between handles) | 22" (560 mm) |
| Weight | 12.5 lb (5.7 kg) |
| Weight capacity | 400 lb (181.4 kg) |