

Folding Aluminum Walker

This **ultra lightweight aluminum folding walker** has been designed to help people with reduced mobility in their daily activities.

4082 Series - 1" frame

The walker's sturdy folding frame is chrome plated to resist scratches and features two different wheel position options: exterior for increased stability or interior for limited motion space.

Included features:

- Adjustable floor/wrist height
- Hypoallergenic, comfortable rubber grips
- Scratch resistant chrome plated folding frame
- Bearing wheels (smooth and quiet motion)
- Wheel stem inlet on both sides
(to alternate interior/exterior wheel position without removing the wheels)



For cases of ataxia, it is possible to add 9 lb (4.1 kg) to the walker's weight

	#ITEM	DESCRIPTION	HEIGHT	USER HEIGHT (approx.)
WHEELS / SKIS	4082AW5	ADULT with wheels / skis	32" – 38" (810 mm – 965 mm)	5' 7" – 6' 1" (1.70 m – 1.85 m)
	4082A+W5	ADULT+ with wheels / skis	35" – 41" (890 mm – 1040 mm)	5' 10" – 6' 4" (1.78 m – 1.93 m)
	4082JW5	JUNIOR with wheels / skis	29" – 33" (735 mm – 840 mm)	5' 4" – 5' 8" (1.63 m – 1.73 m)
RUBBER TIPS	4082A	ADULT rubber tips	32" – 38" (810 mm – 965 mm)	5' 7" – 6' 1" (1.70 m – 1.85 m)
	4082A+	ADULT+ rubber tips	35" – 41" (890 mm – 1040 mm)	5' 10" – 6' 4" (1.78 m – 1.93 m)
	4082J	JUNIOR rubber tips	29" – 33" (735 mm – 840 mm)	5' 4" – 5' 8" (1.63 m – 1.73 m)

LIMITED WARRANTY

5 years: frame

1 year: every other component
(except for wheels / skis)

SPECIFICATIONS	
Overall width (exterior wheel position)	25.75" (655 mm)
Overall width (interior wheel position)	23" (585 mm)
Overall width (rear-skis)	23.25" (590 mm)
Inner width (between handles)	17" (430 mm)
Weight	7.5 lb (3.4 kg)
Weight capacity	330 lb (149.7 kg)