



1. ACCOUNT INFORMATION

P.O. # or Quote:	_____		
Dealer Contact Name:	_____		
Email Address:	_____		
Phone Number:	_____	Fax Number:	_____
Bill to (Name):	_____		
Bill to (Location):	_____		
Ship to (Name):	_____		
Address:	_____		
City:	_____		
Province:	_____	Postal Code:	_____
User Name or Number:	_____		
User Weight:	_____		
User Height:	_____		
Seat Cushion Thickness:	_____		
Knee to Heel Length: Left:	_____	Right:	_____



THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:

- A user name or number, the client's measurements and weight - Amylior is not responsible for configuration or size discrepancies resulting from customer errors **or omissions** on the order form.
- Order confirmations will be provided by Amylior summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

***IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES**

2. STANDARD OPTIONS

Alltrack HD3 Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.4mph motors, Weight Capacity = 550 lb

3. ALLTRACK HD3 POWER WHEELCHAIR

ATHD3	ALLTRACK HD3 Power Wheelchair	\$11,295
--------------	-------------------------------	----------

3.5 OWNER'S MANUAL

SOOMATENG	English / Anglais	N/C	SOOMATFR	Français / French	N/C
------------------	-------------------	-----	-----------------	-------------------	-----

4. USER WEIGHT PACKAGES

WCA55	550 lb Maximum user weight capacity	STD
--------------	-------------------------------------	-----

5. STANDARD COLORS

SCMOBM	Onyx Black Matte	N/C	SCMCAB	Cabernet	N/C
SCMCTB	Cobalt Blue	N/C	SCMBGM	Bubble Gum Metallic	N/C
SCMIOB	Indigo Blue	N/C	SCMHAB	Habanero Metallic	N/C
SCMCRE	Candy Red	N/C	SCMLIM	SubLime	N/C

▼ OPTIONAL COLORS

FCHDD-M	Desert Digital Camo HydroGraphics - (Military Tan base)	\$795
FCHCF-M	Carbon fiber HydroGraphics - (Matte black base)	\$795
FCHWL-M	Woodland Camo HydroGraphics - (Military Green base)	\$795



5.5 WHEEL INSERT COLORS

SCCTB-WI	Cobalt Blue	N/C	SCHAB-WI	Habanero Metallic	N/C
SCIOB-WI	Indigo Blue	N/C	SCLIM-WI	SubLime	N/C
SCCRE-WI	Candy Red	N/C	SCBGM-WI	Bubble Gum Metallic	N/C
SCCAB-WI	Cabernet	N/C			
▲ <i>Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)</i>					
▲ <i>Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)</i>					
▲ <i>Wheel inserts are not available in Onyx Black Matte</i>					

7. TRANSPORT OPTION

TTD4	4-Point wheelchair transport brackets	STD
▲ <i>4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.</i>		

8. MOTOR PACKAGE

MP544B	5.4 MPH 4-Pole Motor Package / with Brake Release	N/C
▲ <i>Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)</i>		
MP544-OPT	5.4 MPH 4-Pole Motor Package / with COG Release	CLIENT COST \$395

9. ASSISTIVE DRIVING PACKAGE

SMARTT	Smart-Track Assistive Driving Package	\$1,420
▲ <i>Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.</i>		
▲ <i>Requires the R-net PM 120 Amp Expandable Controller (ECR120)</i>		
▲ <i>This option requires the MP544-OPT motor package.</i>		
▲ <i>Recommended for non-proportional Drive Controls</i>		

10. BATTERY CHARGER

BCOB	8A Offboard Battery Charger	N/C
-------------	-----------------------------	-----

11. BATTERY TRAY AND BASE WIDTH

FS24B	Battery Tray for Group 24 Batteries – (Total base width = 27 3/8")	N/C
--------------	--	-----

12. BATTERIES - REQUIRES A SELECTION

BAT24GEL	GEL Battery – Group 24 – 70 Amps (Each)	AADL CONTRACT	RFQ
BATO	Omission of the batteries		N/C

14. DRIVE WHEEL

DW14PT	14" Cast Wheel Pneumatic Tires	STD
DW14PI	14" Cast Wheel Pneumatic with Foam Filled Tires	N/C
DW144PT	14" x 4" Extra Wide Cast Wheels, Pneumatic Tires	CLIENT COST \$200
DW144PI	14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires	CLIENT COST \$595
▲ <i>Increases the overall width of the base by 2 3/4"</i>		

15. FRONT CASTERS

FC8PTM	8" Pneumatic Casters	STD
FC8PIM	8" Pneumatic Casters with Foam Filled Tires	N/C



16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± 1/4")

SEAT TO FLOOR HEIGHT TABLE SEAT ANGLE	NO POWER OPTIONS			TILT OR TILT/RECLINE
	0°	3°	6°	
STFH1700 (17.00" STFH)				
STFH1800 (18.00" STFH)				
STFH1900 (19.00" STFH)				
STFH2000 (20.00" STFH)				

17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)

SOMHSS	Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges.	RFQ
▲ <i>Required at: _____" (± 1/4" Range)</i>		

18. PRE-SET SEAT FRAME ANGLE

SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

19. LIGHTING PACKAGE

EOLP	LED Light Package – (Includes front and rear lights with integrated flashers)	CLIENT COST	\$1,425
▲ <i>Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL, ECRJCL, ECRJC35B or EDIB).</i>			

21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

ECEL90	R-net EL 90 Amp Expandable Controller	\$175
ECR120	R-net PM 120 Amp Expandable Controller	\$275

22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

▶ SELECTION OF ONE JOYSTICK OPTION BELOW IS REQUIRED WITH THE EXPANDABLE CONTROLLER ECEL90 OR ECR120			
ECO	Omit Joystick – (For use with an Omni 2 Display)		N/C
ECELJ	R-net - (LED) – Joystick		STD
▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>			
ECRJC	R-net Joystick with Color screen		\$155
▶ JOYSTICKS WITH LIGHTING OPERATION			
ECELJL	R-net – (LED) – Joystick with Lighting Operation	CLIENT COST	\$125
▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>			
ECRJCL	R-net Joystick with Color screen and Lighting Operation		\$180
▶ ADVANCED JOYSTICK			
ECRJC35B	R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	CLIENT COST	\$1,105

23. JOYSTICK MOUNT

JSMF	Fixed Joystick Mount	N/C
-------------	----------------------	-----



23. JOYSTICK MOUNT (CONTINUED)

JSMSA	Swing Away Joystick Mount – Height adjustable, Angle adjustable and Multi-Axis		\$235
Select:	Set-up with standard installation to swing-away sideways	Set-up at 35° to swing-away sideways and downwards	
Select:	Mechanical Lock (STD)	Magnetic Lock	
JSMCM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount		\$275

24. JOYSTICK MOUNT SUPPORT

JSMFH	Joystick mount support Bracket	STD
--------------	--------------------------------	-----

25. JOYSTICK HAND MOUNT

HML	Install Joystick on Left hand side	HMR	Install Joystick on Right hand side
------------	------------------------------------	------------	-------------------------------------

26. JOYSTICK HANDLE

JSH	Standard Joystick Gimble	STD	JBAT	Stick Handle 3" in Length	\$130
JBH	Ball Handle	\$85	JFBH	Foam Ball Handle	\$70
JTH	T-Handle	\$155	JMUSH	Bodypoint Mushroom Handle	\$210
JGP	Bodypoint 4" U-shaped Handle	\$210			

29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)

EDIB	Omni 2 – (Enhanced Display & Digital Interface Box with Infrared)		\$1,915
	▲	<i>Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120)</i>	
	▲	<i>The Omni 2 display mount must be selected from the option below.</i>	

30. OMNI 2 ENHANCED DISPLAY MOUNT

EDSAM	User Swing Away Mount	CODE W100	N/C
	EDML – (Install Display on LH side)	EDMR – (Install Display on RH side)	

31. ENVIRONMENTAL CONTROL MODULE

ECU	Environment control module – (Output Module)	CLIENT COST	\$995
	▲	<i>Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).</i>	

33. AMYSEAT REHAB SEAT / POWER SEATING SYSTEMS AND MODULES

SAABAHD-550	Bariatric Amyseat Rehab Seat – (Seat Frame only)		\$1,045
SAPTXVHD	Bariatric CG Power Tilt	CODE W995	\$8,995
	▲	<i>This option includes the Amyseat Rehab Seat – (SAABAHD)</i>	
SAPRXVHD	Bariatric Power Recline with Mechanical Shear Reduction (MSR)	CODE W995TR	\$9,845
	▲	<i>This option includes the Amyseat Rehab Seat – (SAABAHD)</i>	
SAPTPRXVHD	Bariatric CG Power Tilt and Power Recline with Mechanical Shear Reduction (MSR)	CODE W995TR	\$14,495
	▲	<i>This option includes the Amyseat Rehab Seat – (SAABAHD)</i>	



34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER

▶ SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (For both expandable and non-expandable systems)										
SOLO	Single function controller (Includes an external toggle switch)								STD	
▼ OPTIONAL SWITCHES										
SWTOGGLE	Single Toggle Switch						Attendant Activation	\$100		
SWPUSH1	Single Push Button Switch – (Single Mode: SAC)				User Activation	Attendant Activation	\$90			
SWPUSH2	Dual Push Button Switch – (Dual Mode: SAC)				User Activation	Attendant Activation	\$180			
SSBS	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$230		
SSES	Egg Switch-Black				User Activation	Attendant Activation	\$116			
SSML	Ultra Light Switch				User Activation	Attendant Activation	\$295			
▶ SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL										
EPSEL	Single Function Through Drive Control with R-net 90-amp controller								N/C	
	▲ <i>Allows for activation of single power option through R-net 90-amp controller</i>									
	▲ <i>Compatible with all R-net joysticks</i>									
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>									
	▲ <i>Can not operate more than one power option</i>									
	▲ <i>If selecting memory seating, choose SAC below</i>									
	▲ <i>If selecting 120-amp R-net controller, choose SAC below</i>									
SAC	Single Function Through Drive Control with R-net 90-amp or 120-amp controller								N/C	
	▲ <i>Allows for activation of single power option through R-net drive controls</i>									
	▲ <i>Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls</i>									
▶ SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE)										
SACMEMORY	Memory Seat Program						CLIENT COST	\$1,495		
	▲ <i>This option is not compatible with LED joysticks (ECELJ and ECELJL)</i>									
	▲ <i>Memory seat programming is not compatible with the Very Heavy Duty Tilt option (SAPT XVHD and SAPTPRXVHD)</i>									
▶ SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH										
SACTOGGLE	Single Toggle Switch						Attendant Activation	CODE W100	\$410	
SACPUSH1	Single Push Button Switch – (Single Mode: SAC)				User Activation	Attendant Activation	CODE W100	\$400		
SACPUSH2	Dual Push Button Switch – (Dual Mode: SAC)				User Activation	Attendant Activation	CODE W100	\$490		
SACSSBS	Buddy Button Switch				User Activation	Attendant Activation	CODE W100	\$540		
	▲	Black	Red	Green	Yellow					
SACSSSES	Egg Switch-Black				User Activation	Attendant Activation	CODE W100	\$426		
SACSSML	Ultra Light Switch				User Activation	Attendant Activation	CODE W100	\$605		



35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE

▶ MULTIPLE SEAT FUNCTIONS THROUGH R-NET					
MAC	E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators				\$2,295
▶ MAC OPTIONS					
MACMEMORY	Memory Seat Program				\$175
▶ SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH					
MACKEYPAD-P	Keypad w/ 5 Push Button Switches	User Activation	Attendant Activation	CODE W100	\$445
MACKEYPAD-T	Keypad w/ 5 Toggle Switches	User Activation	Attendant Activation	CODE W100	\$445
MACMULTI	4-Way Toggle Switch – (Up to 4 Activations)	User Activation	Attendant Activation	CODE W100	\$180
▲ By selecting this switch, memory recall is only available through the Joystick.					
MACMULTIT	Triple Toggle Switch – (Up to 5 Activations)	User Activation	Attendant Activation	CODE W100	\$135
▲ By selecting this switch, memory recall is only available through the Joystick.					
MACPUSH4	4-Way Push Button Switch	User Activation	Attendant Activation	CODE W100	\$225
▲ By selecting this switch, memory recall is only available through the Joystick.					
MACPUSH5	5-Way Push Button Switch	User Activation	Attendant Activation	CODE W100	\$240
▲ By selecting this switch, memory recall is only available through the Joystick.					

36. SWITCH ACTIVATION HAND MOUNT

Switch Installation:	TSML – (Install Switch on LH side)	TSMR – (Install Switch on RH side)	N/C
----------------------	---	---	-----

37. PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)

SOBAA	Select one:	86°	90°	94°	98°	102°	106°	110°	N/C
--------------	-------------	-----	-----	-----	-----	------	------	------	-----

39. LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19” AND 28”

SOWN19-24-22	Width adjustment from 19” to 24” and Depth range of 18” to 22”								N/C
SOWN24-28-24	Width adjustment from 24” to 28” and Depth range of 18” to 24”								\$475
19” (SOWN19L)	20” (SOWN20L)	21” (SOWN21)	22” (SOWN22)	23” (SOWN23)	24” (SOWN24)				
24” (SOWN24HDL)	25” (SOWN25HDL)	26” (SOWN26HDL)	27” (SOWN27HDL)	28” (SOWN28HDL)					
▶ CUSTOM WIDTH ABOVE 28”									
SOWCUS	Custom Seat Width (May require the signature of a waiver) ◀ Specify: _____								RFQ

40. ADJUSTABLE SEAT DEPTH RANGE BETWEEN 18” AND 24”

18” (SODN18)	19” (SODN19)	20” (SODN20)	21” (SODN21)					N/C	
22” (SODN22)	23” (SODN23)	24” (SODN24)					N/C		
▶ CUSTOM DEPTH ABOVE 24”									
SODCUS	Custom Seat Depth (May require the signature of a waiver) ◀ Specify: _____								RFQ

41. BACK CANE HEIGHT (THE BACK CANE HEIGHT MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE PUSH HANDLES)

▼ 22” is the minimum back height required with a power recline.

16” (SOBCH16)	18” (SOBCH18)	20” (SOBCH20)	22” (SOBCH22)	24” (SOBCH24)	26” (SOBCH26)
----------------------	----------------------	----------------------	----------------------	----------------------	----------------------



42. BACK SUPPORT

DSTR-wwhh	Adjustable Tension Back Upholstery, WW x HH	N/C
	▲ <i>Only available with Dual Post armrests (SOADP)</i>	
STDSB-wwhh	Contour back support up to 21" wide, with a maximum shell height of 22".	SEATING BENEFIT \$605
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
STDSB-wwhh	Contour back support up to 21" wide for non-standard sizes with shell height ≥ 23".	SEATING BENEFIT \$755
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
STDSB-wwhh	Contour back support for Heavy-Duty from 22" wide and higher, with any height.	SEATING BENEFIT \$1,075
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
▶ MATRX BACKS		
PBE-wwhh-hd	MATRX HEAVY DUTY ELITE BACK – (21" to 24" Wide - Back height of 16" and 20")	SEATING BENEFIT \$1,150
	▲ <i>Back shell height is available in 16" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>Weight capacity up to 500 lb</i>	
EDB-wwhh-hd	MATRX ELITE DEEP HD BACK – (21" to 24" Wide - Back height of 16", 18" and 20")	SEATING BENEFIT \$1,195
	▲ <i>Back shell height is available in 16", 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Weight capacity up to 500 lb</i>	
▶ NXT BACKS		
NOACB-wwhh	NXT OPTIMA FH AC CONTOUR BACK – (15" to 21" Wide)	SEATING BENEFIT \$690
	▲ <i>Weight capacity up to 300 lb</i>	
NOACB-wwhh	NXT OPTIMA FH AC CONTOUR BACK – (22" to 28" Wide)	SEATING BENEFIT \$720
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
NOPCB-wwhh	NXT OPTIMA FH PC DEEP CONTOUR BACK – (15" to 21" Wide)	SEATING BENEFIT \$790
	▲ <i>Weight capacity up to 300 lb</i>	
NOPCB-wwhh	NXT OPTIMA FH PC DEEP CONTOUR BACK – (22" to 28" Wide)	SEATING BENEFIT \$820
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
▶ OMISSION OF THE BACK SUPPORT		
STDSBO	Omission of the Contour Back Support ◀ Specify: _____	N/C

43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

▼ <i>22" is the minimum back height required with a power recline.</i>						
16" (SOBH16)	18" (SOBH18)	20" (SOBH20)	22" (SOBH22)	24" (SOBH24)	26" (SOBH26)	28" (SOBH28)
This height is not available with the NXT Optima AC and PC options ▲						



44. SEAT CUSHIONS (SEATING BENEFIT)

▶ **FOR CUSTOM SIZES AND CONFIGURATIONS, PLEASE REFER TO THE PERFORMA CG AIR ORDER FORM.**

▶ **PERFORMA CG AIR - THERAPEUTIC AIR CUSHION**

Custom Cushion (Please use the Performa CG Air order form)

Cell Height (Choose 1): 2" 3" 4"

▶ **SINGLE (1) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING**

▶ **PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER**

PCGA1VSCxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$812
----------------------	------------------------------------	------------------------	-------

PCGA1VSCxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$914
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER**

PCGA1V3Dxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$857
----------------------	------------------------------------	------------------------	-------

PCGA1V3Dxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$963
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER**

PCGA1VICxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$915
----------------------	------------------------------------	------------------------	-------

PCGA1VICxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$969
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **DUAL (2) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING**

▶ **PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER**

PCGA2VSCxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$812
----------------------	------------------------------------	------------------------	-------

PCGA2VSCxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$914
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER**

PCGA2V3Dxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$857
----------------------	------------------------------------	------------------------	-------

PCGA2V3Dxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$963
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER**

PCGA2VICxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$915
----------------------	------------------------------------	------------------------	-------

PCGA2VICxxx-x	22" + (560 mm +) wide	SEATING BENEFIT	\$969
----------------------	-----------------------	------------------------	-------

▲ **Cushion sizes only come in even numbers for width and depth measurements.**

▶ **DUAL VALVE CUSHION SEGMENTATION (PLEASE SELECT 1 BOX)**

PCGA-2VLR	Dual Valves with Left / Right Location	SEATING BENEFIT	N/C
------------------	--	------------------------	-----

PCGA-2VF6	Dual Valves with Front / Back Location, Front distance 6 in (150 mm) from the front	SEATING BENEFIT	N/C
------------------	---	------------------------	-----

PCGA-2VF8	Dual Valves with Front / Back Location, Front distance 8 in (205 mm) from the front	SEATING BENEFIT	N/C
------------------	---	------------------------	-----

PCGA-2VF10	Dual Valves with Front / Back Location, Front distance 10 in (255 mm) from the front	SEATING BENEFIT	N/C
-------------------	--	------------------------	-----

▶ **PERFORMA CG AIR OPTIONS**

PCGA-GZ	GripZone Anti-Slip Technology	SEATING BENEFIT	\$86
----------------	-------------------------------	------------------------	------

PCGA-TT-WW	Transfer Trellis	SEATING BENEFIT	\$100
-------------------	------------------	------------------------	-------

PCGA-AC	Anatomic Contour Seat Insert	SEATING BENEFIT	\$100
----------------	------------------------------	------------------------	-------

▲ **Only available with 2" and 3" cell heights**

PCGA-SP	Seat Pouch	SEATING BENEFIT	\$55
----------------	------------	------------------------	------



44. SEAT CUSHIONS (CONTINUED)

▶ **MATRX CUSHION**

PSP-wydd	MATRX PSP Seat Cushion – (11" to 20" Wide)	SEATING BENEFIT	\$645
	▲ The maximum user weight capacity is 300 lb		
	▲ Cushion sizes only come in even numbers for width and depth measurements.		
PSP-wydd	MATRX PSP Seat Cushion – (22" to 26" Wide)	SEATING BENEFIT	\$995
	▲ The maximum user weight capacity is 550 lb		
	▲ Cushion sizes only come in even numbers for width and depth measurements.		

▶ **NXT BIO-FIT CUSHION**

NXTBIOFIT-wydd	NXT BIO FIT Cushion – (15" to 21" Wide)	SEATING BENEFIT	\$595
	▲ The maximum user weight capacity is 300 lb		
NXTBIOFIT-wydd	NXT BIO FIT Cushion – (22" to 26" Wide)	SEATING BENEFIT	\$625
	▲ The maximum user weight capacity is 400 lb		
	▲ Cushion sizes only come in even numbers for width and depth measurements.		

45. CANTILEVER – HEIGHT ADJUSTABLE ARMRESTS – OPENS WITH A 5° ANGLE (CHOOSE 2)

SOAFBDLH	Flip-Back Cantilever – Desk Length Left Side – (Height Adjustable from 7" to 15" with 1" increments)	STD Length 12"	N/C
SOAFBDRH	Flip-Back Cantilever – Desk Length Right Side – (Height Adjustable from 7" to 15" with 1" increments)	STD Length 12"	N/C
SOAFBFLH	Flip-Back Cantilever – Full Length Left Side – (Height Adjustable from 7" to 15" with 1" increments)	STD Length 16"	STD
SOAFBFRH	Flip-Back Cantilever – Full Length Right Side – (Height Adjustable from 7" to 15" with 1" increments)	STD Length 16"	STD

46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)

SOADPDLH	Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
SOADPDRH	Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
SOADPFLH	Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
SOADPFRH	Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168

47. RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (CHOOSE 2)

SOANIDLH	Reclining armrest – Desk Length Left Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments)	STD Length 12"	N/C		
SOANIDRH	Reclining armrest – Desk Length Right Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments)	STD Length 12"	N/C		
SOANIFLH	Reclining armrest – Full Length Left Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments)	STD Length 16"	N/C		
SOANIFRH	Reclining armrest – Full Length Right Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments)	STD Length 16"	N/C		
CUST-48247	Custom Armrest Length – Select Length for the Left side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C
CUST-48247	Custom Armrest Length – Select Length for the Right side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C



48. ARMREST HEIGHT ADJUSTMENT AT FACTORY

▶	Set the Armrest Height at: _____ “
▲	The Armrest Height is taken from the seat pan to the top of the armpad.
▲	Not required for Dual Post.

49. ARMPADS

▶ FOAM TYPE SELECTION

SOWAPD	Desk Length Waterfall – (8 ½” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	STD
▲	<i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>				
SOWAPF	Full Length Waterfall – (13” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	STD
▲	<i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>				

▶ HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)

SOGELD	Desk Length Hybrid Gel and Foam Flat armpad – (9 ¾” x 2”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$40
SOGELF	Full Length Hybrid Gel and Foam Flat armpad – (13 ½” x 2”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$40
SOWGELD	Desk Length Hybrid Gel & Foam Waterfall armpad – (8 ½” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$95
▲	<i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>					
SOWGELF	Full Length Hybrid Gel and Foam Waterfall armpad – (13” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$95
▲	<i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>					

▶ SPECIALTY TYPE SELECTION

SOCAR	Arm Trough with Hand Support (18” x 4 ½”)	Select side:	Left Hand Side	Right Hand Side	\$290
▲	<i>Arm Trough Armpads may interfere with Armrests Flipping back.</i>				

50. FOOTRESTS

SOAAF	Signature Heavy Duty Swing-Away Footrests				\$175
	Set at:	60°		70°	
	Select Range:	SOFLHD12 (12 ½” – 17 ½”)		SOFLHD14 (14 ½” – 19 ½”)	
ADDML	Signature Manual Elevating & Articulating Legrests - Pair				\$595
	Select Range:	SOFLME12 (12 ½” – 17 ½”)		SOFLME14 (14 ½” – 19 ½”)	
▲	<i>Does not include calf pads. Please select calf pads in section 51.</i>				
ADDPLHD	Signature Heavy Duty Power Elevating & Articulating Legrests			CLIENT COST	\$2,745
	Select Function:	Individual		Combined	
	Select Range:	SOFLPW12 (12 ½” – 17 ½”)		SOFLPW14 (14 ½” – 19 ½”)	
▲	<i>This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.</i>				

51. FOOTPLATES & ACCESSORIES

SOSFAF	Flip-up Footplates	STD
SOAFP	Angle Adjustable & Flip-up Footplates – (Each)	\$65
SOAHL	Heel Loops – (Each)	\$20



51. FOOTPLATES & ACCESSORIES (CONTINUED)

SOACPF	Calf Pads for Adult Heavy-Duty Swing-Away footrests – (SOAAF) – with adjustable hardware – (Each)										\$93
SOAMA	Swing-away residual limb Support	Select:	Left leg		Right leg						\$455

52. FOOTREST LENGTH ADJUSTMENT AT FACTORY

Footrest / Legrest Length	Set at:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"
		16.5"	17"	17.5"	18"	18.5"	19"	19.5"	

▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

53. CENTER MOUNTS

SOAAF90	Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°)										\$200
	Select Foot Platform Range:	SOFBA12 (12 ½" – 16 ½")		SOFBA16 (16 ½" – 19 ½")							
	▲	<i>When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 18.75"</i>									
	▲	<i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>									
	▲	<i>With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.</i>									
ADDPFVHD	90° Heavy Duty Power Center Mount Elevating & Articulating Foot Platform										\$2,745
	Select Footrest Range:	SOFBPW12 (10 ½" - 16")		SOFBPW16 (16" - 20")							
	▲	<i>When selecting SOFBPW16 (16" - 20") range, the lowest possible seat to floor height is 18.75"</i>									
	▲	<i>This option includes calf pads. The power center mount counts as one power function.</i>									
	▲	<i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>									
	▲	<i>For leg lengths between 10 ½" and 12 ½", there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.</i>									

54. FOOT PLATFORMS & ACCESSORIES

SOFUFP	Flip-up one piece Foot Platform										\$140
SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount – (Each)										\$65
SOACPCM	Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each)										\$93

55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY

Fixed Center Mount Foot Platform Length	Set at:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"
		17.5"	18"	18.5"	19"	19.5"					

Power Center Mount Foot Platform Length	Set at:	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"
		15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"

▲ The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates

56. SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19-ISO OPTION)

SOAUTO80-STD	80" Seat Positioning Strap										STD
SOBPHP2C	Bodypoint padded Hip Belt – (Two Point – Center Pull)										\$170
SOBPHP2C-L82	XL 82" Bodypoint padded Hip Belt – (Two Point – Center Pull)										\$190
SOPHB	Padded Hip Belt 60", 2-PT Magnetic Center Pull										\$194
SOPCS	Padded Chest Strap 60", 2-PT Magnetic Center Pull										\$262



57. HEADREST

▼ **SELECT ONE OF THE HEADREST OPTIONS WITH THE MOUNTING HARDWARE AND SPECIFY WHERE THE HEADREST WILL BE INSTALLED.**

SONUAHR	Using headrest mounting hardware other than SOMAH located below. ▲ Specify: _____	N/C
SOFHRS	Small Curved Headrest Pad – 4" x 9"	\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"	\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"	\$155
SOMAH	Multi-axis Removable Mounting Hardware	\$220
SOHRBC	Bracket to Mount on the Contour Back Support – (Without Recline)	N/C
SOHRPR	Bracket to Mount on the Seat Frame – (For Reclining Backs)	N/C

58. LATERAL TRUNK SUPPORTS

▼ **SELECT ONE OF THE LATERAL PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

SOSLS3	Adjustable Laterals with Small Curved Pads – (Each) – 3" x 5"	SEATING BENEFIT	\$100
SOSLM4	Adjustable Laterals with Medium Curved Pads – (Each) – 4" x 6"	SEATING BENEFIT	\$100
SOSLL5	Adjustable Laterals with Large Curved Pads – (Each) – 5" x 7"	SEATING BENEFIT	\$100
SOSLHM	Adjustable Swing-Away Hardware Mounts – (Each)	SEATING BENEFIT	\$220
▲ Waterfall armpads (SOWAPx & SOWGELx) may interfere with flipping back the armrests when combined with Lateral Trunk supports			

59. HIP & THIGH SUPPORTS

▼ **SELECT ONE OF THE HIP & THIGH SUPPORT PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

SOHTSS10	Adjustable Hip & Thigh Support Pads – (Each) – Small 10" x 4"	SEATING BENEFIT	\$130
SOHTSM12	Adjustable Hip & Thigh Support Pads – (Each) – Medium 12" x 5"	SEATING BENEFIT	\$130
SOHTSL14	Adjustable Hip & Thigh Support Pads – (Each) – Large 14" x 6"	SEATING BENEFIT	\$130
SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware – (Each)	SEATING BENEFIT	\$136

60. ACCESSORIES

SOCH	Cup Holder – (Self-Leveling)	Install on:	Left Side	Right Side	\$115
SOCPH	Ram Mounts X-Grip Cell Phone Holder	Install on:	Left Side	Right Side	\$190
SOTHS	Ram Mounts X-Grip Small Tablet Holder	Install on:	Left Side	Right Side	\$205
SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	Install on:	Left Side	Right Side	\$360
▲ Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.					
▲ With multiple accessories installed on the same side, a custom upgrade will be required.					

61. REFLECTOR KIT

REFLECT	Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner's Manual)	N/C
----------------	---	-----

62. BACK PACK AND ARMREST SIDE POUCHES

SOAP	Armrest Side Pouch (Comes with durable straps that can be hooked around the armrest tubes)	Order:	Left Side	Right Side	CLIENT COST	\$53
▲ If only choosing one (1) armrest pouch, please show what side it should be installed						
SOBP	Back Pack – (Comes with durable straps that can be hooked around the back canes)				CLIENT COST	\$90



63. SWITCHES AND OPTIONS

ASCPWIFI	PGDT - Wireless programmer, WIFI dongle	NOT COVERED BY AADL	\$844
ASCPCP	PGDT R-NET Dongle and PC programmer kit	NOT COVERED BY AADL	\$600
SSBS	Buddy Button Switch		\$230
SSES	Egg Switch-Black		\$116
SSML	Ultra Light Switch		\$295
ESS-RED	Emergency STOP Switch - Red		\$314

SPECIALTY DRIVER CONTROL OPTIONS

► Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)

64. THERAFIN SIP-N-PUFF KIT (CODE W100)

SCSNPK	Therafin Sip-N-Puff Kit				\$1,975
SCSNPBTK	Therafin Sip-N-Puff Breath Tube Kit				N/C
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRMSN	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	SWL (Swings on Left)	SWR (Swings on Right)		\$300
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Therafin Sip-N-Puff kit.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Therafin Sip-N-Puff kit.				
► CHOOSE ONE (1) HEADREST PAD					
SOFHRS	Small Curved Headrest Pad – 4" x 9"				\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"				\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"				\$155

65. PROXIMITY HEAD ARRAY (CODE W100)

CUSTOM	Adult – 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse				RFQ
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Proximity Head Array.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Proximity Head Array.				



66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL (CODE W100)

CUSTOM	Light touch mini Proportional Joystick			RFQ
SCSRMMP	Stealth Rock-n-Lock Swingaway Arm for Chin MP	SWL (Swings on Left)	SWR (Swings on Right)	\$870
	▲ Required with the Light touch mini Proportional Joystick.			
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section			
JCC	Chin Cup			\$55
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Light touch mini Proportional Joystick.			
▶	CHOOSE ONE (1) HEADREST PAD			
SOFHRS	Small Curved Headrest Pad – 4" x 9"			\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"			\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"			\$155

67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT (CODE W100)

SCPGMJS	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$2,369
SCGMJSM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for JC200	MLH (Mount on Left)	MRH (Mount on Right)	\$595
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
CUSTOM	Light touch mini Proportional Joystick			RFQ
SCGMPPM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for Light Touch Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$825
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$850
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT (CODE W100)

SCPGMJS-M	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$2,369
SCARMJSM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
CUSTOM	Light touch mini Proportional Joystick			RFQ
SCARMPM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$375
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$850
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

69. PROPORTIONAL HEAD CONTROL PACKAGE (CODE W100)

SCPHC	Proportional Head Control Package			\$4,495
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510



69. PROPORTIONAL HEAD CONTROL PACKAGE (CONTINUED) (CODE W100)

SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section			
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Proportional RIM Head Control.			

CUSTOM SECTION

- ▶ **Lead times will increase depending on the custom order pipeline at the time the purchase order is received.**
- ▶ **Please call for Quote and specifications for these options.**

70. O² HOLDER

SUPOXY	Custom Non-Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$295
	▲ Not available with a power recline			
SUPOXYA	Custom Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$1,195
	▲ Required with a power recline			
	D-Size	Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve		
	E-Size	Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve		

71. VENT TRAY (PLEASE COMPLETE IN FULL)

▼ **AMYLIOR DOES NOT SUPPLY THE BATTERY, CHARGER OR THE WIRING FOR THE VENT. THESE ITEMS ARE NORMALLY PROVIDED BY THE VENT MANUFACTURER.**

SOVENT	Custom Vent Tray			
	Articulating Vent Box (Required with Power Tilt and/or Recline)			\$1,360
	Fixed Vent Box			\$400
	SPECIFY MODEL			
	Trilogy 100/200			
	LTV			
	LTV with Sprintpack			
	Other, specify: _____			
SOVTBE	Extended Battery Compartment			\$495
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries			
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries			

72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)

▼ **PLEASE SPECIFY DETAILS AND SPECIFICATIONS FOR THESE OPTIONS.**

CUSTOM	USB Charger Port ◀ Specify the location for the installation: _____	CLIENT COST	RFQ			
CUSTOM	Accessory Hooks located at the rear of the seat structure (2 hooks)	CLIENT COST	RFQ			
CUSTOM	Heavy-Duty Wheelchair Cover, Large	CLIENT COST	RFQ			
CUSTOM	Transfer Handles – (Pair) Select a height range: 5" to 8" 6" to 9"	CLIENT COST	RFQ			
CUSTOM	Folding Backrest for Additional Storage possibilities and Travel accessibility	CLIENT COST	RFQ			
CUSTOM	Power Anterior Tilt	CLIENT COST	RFQ			
CUSTOM	Power Lateral Tilt 15°	Raise the Left Side only	Raise the Right Side only	Raise the Left & Right Side	CLIENT COST	RFQ



72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE) (CONTINUED)

CUSTOM	Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				CLIENT COST	RFQ
CUSTOM	Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				CLIENT COST	RFQ
CUSTOM	Rearview mirror	Left Side	Right Side		CLIENT COST	RFQ
CUSTOM	Rearview Camera with 5" Monitor				CLIENT COST	RFQ
CUSTOM	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan				CLIENT COST	RFQ
CUSTOM	Endomorphic Back Canes (Pair)	1" per side towards the inside = 2"	2" per side towards the inside = 4"		CLIENT COST	RFQ
CUSTOM	Foot Box	Specify	Length: _____	Width: _____	Please confirm with CSR on sizing possibilities	CLIENT COST RFQ
CUSTOM #1	Other custom options. Specify in notes section below				CLIENT COST	RFQ
CUSTOM #2	Other custom options. Specify in notes section below				CLIENT COST	RFQ

73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS