



1. ACCOUNT INFORMATION

P.O. # or Quote: _____

Dealer Contact Name: _____

Email Address: _____

Phone Number: _____ Fax Number: _____

Bill to (Name): _____

Bill to (Location): _____

Ship to (Name): _____

Address: _____

City: _____

Province: _____ Postal Code: _____

User Name or Number: _____

User Weight: _____

User Height: _____

Seat Cushion Thickness: _____

Knee to Heel Length: Left: _____ " Right: _____ "



THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:

- A user name or number, the client's measurements and weight - Amylor is not responsible for configuration or size discrepancies resulting from customer errors *or omissions* on the order form
 - Order confirmations will be provided by Amylor summarizing the items selected. Please review this carefully.
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.
- *IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES*

2. STANDARD OPTIONS

Alltrack M3 Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 6" front casters, and standard auto style buckle seat positioning strap, 6.5 mph motors, Weight Capacity = 300 lb

3. ALLTRACK M3 POWER WHEELCHAIR

ATM3	ALLTRACK M3 Power Wheelchair	\$8,020
-------------	------------------------------	---------

3.5 OWNER'S MANUAL

SOOMATENG	English / Anglais	N/C	SOOMATFR	Français / French	N/C
------------------	-------------------	-----	-----------------	-------------------	-----

4. USER WEIGHT PACKAGES

WCA25	250 lb Maximum user weight capacity	STD	WCA3	300 lb Maximum user weight capacity	STD
--------------	-------------------------------------	-----	-------------	-------------------------------------	-----

5. STANDARD COLORS

SCMOBM	Onyx Black Matte	N/C	SCMCAB	Cabernet	N/C
SCMCTB	Cobalt Blue	N/C	SCMBGM	Bubble Gum Metallic	N/C
SCMIOB	Indigo Blue	N/C	SCMHAB	Habanero Metallic	N/C
SCMCRE	Candy Red	N/C	SCMLIM	SubLime	N/C

▼ OPTIONAL COLORS

FCHDD-M	Desert Digital Camo HydroGraphics	\$795
FCHCF-M	Carbon fiber HydroGraphics	\$795
FCHWL-M	Woodland Camo HydroGraphics	\$795



5.5 WHEEL INSERT COLORS

SCCTB-WI	Cobalt Blue	N/C	SCHAB-WI	Habanero Metallic	N/C
SCIOB-WI	Indigo Blue	N/C	SCLIM-WI	SubLime	N/C
SCCRE-WI	Candy Red	N/C	SCBGM-WI	Bubble Gum Metallic	N/C
SCCAB-WI	Cabernet	N/C			
▲ <i>Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)</i>					
▲ <i>Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)</i>					
▲ <i>Wheel inserts are not available in Onyx Black Matte</i>					

7. TRANSPORT OPTION

TTD4	Tie-down attachment points compliant with ISO 7176-19. See owner's manual for additional details	STD
-------------	--	-----

8. MOTOR PACKAGE

MP65B	6.5 MPH 4-Pole Motor Package / with Brake Release	STD
▲ <i>Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)</i>		
MP65-OPT	6.5 MPH 4-Pole Motor MP65-OPT Package / with COG Release	\$395
▲ <i>Not compatible with the narrow base configuration (FS2221B)</i>		

9. ASSISTIVE DRIVING PACKAGE

SMARTT	Smart-Track Assistive Driving Package	\$1,420
▲ <i>Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.</i>		
▲ <i>Requires the R-net PM 120 Amp Expandable Controller (ECR120).</i>		
▲ <i>This option requires the MP65-OPT motor package.</i>		
▲ <i>Not compatible with the narrow base configuration (FS2221B)</i>		
▲ <i>Recommended for non-proportional Drive Controls</i>		

10. BATTERY CHARGER

BCOB	8A Offboard Battery Charger	N/C
-------------	-----------------------------	-----

11. BATTERY TRAY AND BASE WIDTH

FS2221B	N2221 - Battery Tray for Group 22 Batteries with a Narrow 21" base width (Total base width = 21")	N/C
▶ NOT AVAILABLE WITH THE FOLLOWING OPTIONS:		
▲ <i>6.5 MPH 4-Pole Motor Package / with COG Release option (MP65-OPT)</i>		
▲ <i>Smart-Track option (SMARTT)</i>		
▲ <i>8" front casters (FC8PTM and FC8PIM)</i>		
▲ <i>14" x 4" Extra Wide Cast Wheels, Pneumatic Tires (DW144PT and DW144PI)</i>		
▲ <i>Not compatible with 23" and 24" seat widths selections (SOWN23 and SOWN24)</i>		
FS22B	Battery Tray for Group 22NF Batteries – (Total base width = 23 3/8")	STD
FS2424B	N2424 - Battery Tray for Group 24 Batteries with a Narrow 24" base width (Total base width = 24")	N/C
FS24B	Battery Tray for Group 24 Batteries – (Total base width = 25 3/8")	N/C



12. BATTERIES - REQUIRES A SELECTION

BAT22NFGEL	Premium CG Power - GEL Battery – Group 22NF – 50 Amps (Each)	\$335
BAT24GEL	Premium CG Power - GEL Battery – Group 24 – 70 Amps (Each)	\$465
BATO	Omission of the batteries	N/C

14. DRIVE WHEEL

DW14PT	14" Cast Wheel Pneumatic Tires	STD
DW14PI	14" Cast Wheel Pneumatic with Foam Filled Tires	N/C
DW144PT	14" x 4" Extra Wide Cast Wheels, Pneumatic Tires	\$200
DW144PI	14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires	\$595
	▲ Increases the overall width of the base by 2 ¾"	
	▲ Not compatible with the narrow base configuration (FS2221B)	
DW12PT	12 ½" Cast Wheel Pneumatic Tires	N/C
DW12PI	12 ½" Cast Wheel Pneumatic with Foam Filled Tires	N/C
	▲ 12 ½" Drive wheels lower ground clearance and speed, but increase motor torque.	

15. FRONT CASTERS

FC6S	6" Solid Casters – (Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI)	STD
FC8PTM	8" Pneumatic Casters	N/C
FC8PIM	8" Pneumatic Casters with Foam Filled Tires	N/C
	▲ The 8" casters are NOT compatible with the 12 ½" Drive wheel configuration (DW12PT and DW12PI)	
	▲ Not compatible with the narrow base configuration (FS2221B)	

16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± ¼")

▶ YELLOW BOXES REPRESENT OPTIONS THAT CAN ONLY BE SELECTED WITH 12 ½" DRIVE WHEELS (SECTION 14)					
SEAT TO FLOOR HEIGHT TABLE SEAT ANGLE	NO POWER OPTIONS			TILT <u>OR</u> ELEVATE (WITH OR WITHOUT RECLINE)	TILT <u>AND</u> ELEVATE (WITH OR WITHOUT RECLINE)
	0°	3°	6°		
STFH1550 (15.50" STFH)					
STFH1650 (16.50" STFH)					
STFH1750 (17.50" STFH)					
STFH1850 (18.50" STFH)					
STFH1950 (19.50" STFH)					

17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)

SOMHSS	Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges.	RFQ
	▲ Required at: _____" (± ¼" Range)	

18. PRE-SET SEAT FRAME ANGLE

SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C



19. LIGHTING PACKAGE

EOLP	LED Light Package – (Includes front and rear lights with integrated flashers)	\$1,425
▲ <i>Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL, ECELJL, ECRJCL, ECRJC35B or EDIB).</i>		

20. NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

▶ <i>VR2 Electronics can only accommodate specific power seat functions and combinations such as Tilt / Elevate / Recline / Power Center Mount / Combined Power Legs / Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs / Tilt + Elevate (Custom)</i>		
ECV902A	VR2 90 Amp Drive Controller – (Requires the selection of VR2 joystick below)	STD
ECVJ4K	VR2 Standard 4-Key Joystick	STD
ECVJ6K	VR2 6-Key Joystick (Allows for single or select dual power options to operate through this joystick)	N/C
ECVJ9KL	VR2 9-Key Joystick with Lighting Operation (Allows for single or select dual power options to operate through this joystick)	\$150

21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

ECEL90	R-net EL 90 Amp Expandable Controller	\$175
ECR120	R-net PM 120 Amp Expandable Controller	\$275

22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

▶ <i>SELECTION OF ONE JOYSTICK OPTION BELOW IS REQUIRED WITH THE EXPANDABLE CONTROLLER ECEL90 OR ECR120</i>		
ECO	Omit Joystick – (For use with an Omni 2 Display)	N/C
ECELJ	R-net - (LED) – Joystick	STD
▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>		
ECRJC	R-net Joystick with Color screen	\$155
▶ <i>JOYSTICKS WITH LIGHTING OPERATION</i>		
ECELJL	R-net – (LED) – Joystick with Lighting Operation	\$125
▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>		
ECRJCL	R-net Joystick with Color screen and Lighting Operation	\$180
▶ <i>ADVANCED JOYSTICK</i>		
ECRJC35B	R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,105

23. JOYSTICK MOUNT

JSMF	Fixed Joystick Mount	N/C
JMSA	Swing Away Joystick Mount – Height adjustable, Angle adjustable and Multi-Axis	\$235
Select:	<i>Set-up with standard installation to swing-away sideways</i>	
	<i>Set-up at 35° to swing-away sideways and downwards</i>	
▲ <i>Set-up at 35° angle may cause interference with the user's limbs when the new Signature cantilever armrest (SOAFBxxx) is internally rotated.</i>		
JSMCM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount	\$275
▲ <i>This option is not compliant with ISO 7176-19</i>		

24. JOYSTICK MOUNT SUPPORT

JSMFH	Joystick mount support Bracket	STD
--------------	--------------------------------	-----

25. JOYSTICK HAND MOUNT

HML	Install Joystick on Left hand side	HMR	Install Joystick on Right hand side
------------	------------------------------------	------------	-------------------------------------



26. JOYSTICK HANDLE

JSH	Standard Joystick Gimble	STD	JBAT	Stick Handle 3" in Length	\$130
JBH	Ball Handle	\$85	JFBH	Foam Ball Handle	\$70
JTH	T-Handle	\$155	JMUSH	Bodypoint Mushroom Handle	\$210
JGP	Bodypoint 4" U-shaped Handle	\$210			

27. ATTENDANT DRIVE CONTROL OPTION

EADC	Attendant Hand Control	\$1,150
-------------	------------------------	---------

28. ATTENDANT JOYSTICK HAND MOUNT

EADCL	Install Attendant Joystick on Rear Left hand side	EADCR	Install Attendant Joystick on Rear Right hand side
--------------	---	--------------	--

29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)

EDIB	Omni 2 – (Enhanced Display & Digital Interface Box with Bluetooth and Infrared)	\$1,915
	▲ <i>Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120)</i>	
	▲ <i>The Omni 2 display mount must be selected from the option below.</i>	

30. OMNI 2 ENHANCED DISPLAY MOUNT

EDSAM	User Swing Away Mount	N/C
	EDML – (Install Display on LH side)	EDMR – (Install Display on RH side)

31. ENVIRONMENTAL CONTROL MODULE

ECU	Environment control module – (Output Module)	\$995
	▲ <i>Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).</i>	

33. SIGNATURE TAILORED POWER POSITIONING SOLUTIONS / REHAB SEAT

SAABA	Signature Rehab Seat – (Seat Frame only)	\$880
SAPTX	Signature CG Power Tilt	\$6,195
	▲ <i>This option includes the Signature Rehab Seat – (SAABA)</i>	
SAPRX	Signature 165° Power Recline with Mechanical Shear Reduction (MSR)	\$7,845
	▲ <i>This option includes the Signature Rehab Seat – (SAABA)</i>	
SAPTPRX	Signature CG Power Tilt and 165° Power Recline with Mechanical Shear Reduction (MSR)	\$11,495
	▲ <i>This option includes the Signature Rehab Seat – (SAABA)</i>	
SAPEX	ALTA 11" Power Seat Elevation	\$4,495
	▲ <i>This option includes the Signature Rehab Seat – (SAABA)</i>	
ADDPEX	Module – Add ALTA 11" Power Seat Elevation to either SAPTX, SAPRX or SAPTPRX	\$2,495
ADDMR	Module – Manual Recline with Mechanical Shear Reduction – (MSR)	\$1,125
	▲ <i>The maximum user weight capacity with this option is 250 lb.</i>	
	▲ <i>This option is not compliant with ISO 7176-19</i>	



34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER

▶ SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (For both expandable and non-expandable systems)									
SOLO	Single function controller (Includes an external toggle switch)								STD
▼ OPTIONAL SWITCHES									
SWTOGGLE	Single Toggle Switch							Attendant Activation	\$100
SWPUSH1	Single Push Button Switch – (Single Mode: SAC)					User Activation	Attendant Activation	\$90	
SWPUSH2	Dual Push Button Switch – (Dual Mode: SAC)					User Activation	Attendant Activation	\$180	
SSBS	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$230	
SSES	Egg Switch-Black					User Activation	Attendant Activation	\$116	
SSML	Ultra Light Switch					User Activation	Attendant Activation	\$295	
▶ SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL									
EPSV	Single Function Through VR2 Non-Expandable Drive Control								N/C
	▲ <i>Allows for activation of a single power option through VR2 controller</i>								
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>								
	▲ <i>Control of power seat elevation with tilt is available as a custom. Check appropriate box in custom section 72 in addition to EPSV.</i>								
	▲ <i>VR2 Electronics can only accommodate specific power seat functions such as Tilt / Recline / Power Center Mount / Combined Power Legs / Tilt + Elevate (Custom)</i>								
	▲ <i>Choose EPSV2 to operate select power combinations</i>								
EPSEL	Single Function Through Drive Control with R-net 90-amp controller								N/C
	▲ <i>Allows for activation of single power option through R-net 90-amp controller</i>								
	▲ <i>Compatible with all R-net joysticks</i>								
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>								
	▲ <i>Can not operate more than one power option</i>								
	▲ <i>If selecting memory seating, choose SAC below</i>								
	▲ <i>If selecting 120-amp R-net controller, choose SAC below</i>								
SAC	Single Function Through Drive Control with R-net 90-amp or 120-amp controller								N/C
	▲ <i>Allows for activation of single power option through R-net drive controls</i>								
	▲ <i>Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls</i>								
	▲ <i>Capable of operating tilt with power seat elevation & compatible with memory seating</i>								
▶ SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE)									
SACMEMORY	Memory Seat Program								\$1,495
	▲ <i>This option is not compatible with any VR2 or LED joysticks (ECELJ and ECELJL)</i>								
▶ SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH									
SACTOGGLE	Single Toggle Switch							Attendant Activation	\$410
SACPUSH1	Single Push Button Switch – (Single Mode: SAC)					User Activation	Attendant Activation	\$400	
SACPUSH2	Dual Push Button Switch – (Dual Mode: SAC)					User Activation	Attendant Activation	\$490	
SACSSBS	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$540	
SACSES	Egg Switch-Black					User Activation	Attendant Activation	\$426	
SACSSML	Ultra Light Switch					User Activation	Attendant Activation	\$605	



35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE

► DUAL SEAT FUNCTIONS THROUGH DRIVE CONTROL

EPSV2	Dual Function Through Drive Control with VR2 Non-Expandable Electronics	N/C
	▲ <i>Allows for the activation of 2 power options through the VR2 Controller. Combination of TILT + RECLINE or ELEVATE + RECLINE are not possible with EPSV2</i>	
	▲ <i>Does not allow for control of power options through drive control and separate switches</i>	
	▲ <i>VR2 Electronics can only accommodate specific power seat combinations such as Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs</i>	

► MULTIPLE SEAT FUNCTIONS THROUGH R-NET

MAC	E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators	\$2,295
	▲ <i>Required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2221B)</i>	

► MAC OPTIONS

MACMEMORY	Memory Seat Program	\$175
------------------	---------------------	-------

► SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH

MACKEYPAD-P	Keypad w/ 5 Push Button Switches	User Activation	Attendant Activation	\$445
MACKEYPAD-T	Keypad w/ 5 Toggle Switches	User Activation	Attendant Activation	\$445
MACMULTI	4-Way Toggle Switch – (Up to 4 Activations)	User Activation	Attendant Activation	\$180
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
MACMULTIT	Triple Toggle Switch – (Up to 5 Activations)	User Activation	Attendant Activation	\$135
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
MACPUSH4	4-Way Push Button Switch	User Activation	Attendant Activation	\$225
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
MACPUSH5	5-Way Push Button Switch	User Activation	Attendant Activation	\$240
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			

36. SWITCH ACTIVATION HAND MOUNT

Switch Installation:	TSML – (Install Switch on LH side)	TSMR – (Install Switch on RH side)	N/C
----------------------	---	---	-----

37. SIGNATURE TAILORED SOLUTIONS - PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)

SOBAA	Select one:	86°	90°	94°	98°	102°	106°	110°	N/C
--------------	-------------	-----	-----	-----	-----	------	------	------	-----

38. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT WIDTH RANGE BETWEEN 14" AND 22" (355 - 560 MM)

STS14-19-19	Width adjustment from 14" to 19" (355-480 MM) and Depth range of 15" to 19" (380 - 480 MM)								N/C
STS14-19-22	Width adjustment from 14" to 19" (355-480 MM) and Depth range of 18" to 22" (455 - 560 MM)								N/C
	▲ <i>The minimum seat depth available with this option is 18" (455 MM)</i>								
STS17-22-19	Width adjustment from 17" to 22" (430-560 MM) and Depth range of 15" to 19" (380 - 480 MM)								N/C
STS17-22-22	Width adjustment from 17" to 22" (430-560 MM) and Depth range of 18" to 22" (455 - 560 MM)								N/C
	▲ <i>The minimum seat depth available with this option is 18" (455 MM)</i>								
14" (SOSW14)	15" (SOSW15)	16" (SOSW16)	17" (SOSW17)	18" (SOSW18)	19" (SOSW19)				
20" (SOSW20)	21" (SOSW21)	22" (SOSW22)							



39. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT WIDTH RANGE BETWEEN 20" AND 24" (510 - 610 MM)

STS20-25-22	Width adjustment from 20" to 25" (510-635 MM) and Depth range of 18" to 22" (455 - 560 MM)	N/C
	▲ Seat width of 25" is a custom dimension and require an (RFQ)	
	▲ 23" and 24" Seat widths are not compatible with the N2221 - Narrow 21" base width option (FS2221B)	
	▲ The minimum seat depth available with this option is 18" (455 MM)	
20" (SOSW20)	21" (SOSW21)	22" (SOSW22)
	23" (SOSW23)	24" (SOSW24)

40. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22" (380 - 560 MM)

15" (SODN15)	16" (SODN16)	17" (SODN17)	18" (SODN18)	19" (SODN19)	20" (SODN20)	N/C
21" (SODN21)	22" (SODN22)		▲ 18" is the minimum possible seat depth when choosing a 22" long Seat Pan (STS14-19-22, STS17-22-22 and STS20-25-22)			N/C

41. SIGNATURE TAILORED SOLUTIONS - EXTRUDED BACK CANE HEIGHT (THE MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE BACK CANE)

▶ **Not compatible with a power or manual recline**

16" (SOBCH16)	17" (SOBCH17)	18" (SOBCH18)	19" (SOBCH19)			
▼ 20" is the minimum back height required with a power or manual recline						
20" (SOBCH20)	21" (SOBCH21)	22" (SOBCH22)	23" (SOBCH23)	24" (SOBCH24)	25" (SOBCH25)	
26" (SOBCH26)	27" (SOBCH27)	28" (SOBCH28)				

41.5 SIGNATURE TAILORED SOLUTIONS - BACK CANE OPTIONS

▼ PUSH HANDLES

SOPH	Push Handles (pair)	N/C
-------------	---------------------	-----

▼ ENDOMORPH BACK CANES FOR NON-RECLINING BACKS

STSEND0-F1	Endomorph Back Canes for non-reclining backs, 1 IN (25 MM) total	\$120
STSEND0-F2	Endomorph Back Canes for non-reclining backs, 2 IN (50 MM) total	\$120
STSEND0-F3	Endomorph Back Canes for non-reclining backs, 3 IN (75 MM) total	\$120
STSEND0-F4	Endomorph Back Canes for non-reclining backs, 4 IN (100 MM) total	\$120

▼ ENDOMORPH BACK CANES FOR RECLINING BACKS

STSEND0-R1	Endomorph Back Canes for reclining backs, 1 IN (25 MM) total	N/C
STSEND0-R2	Endomorph Back Canes for reclining backs, 2 IN (50 MM) total	N/C

▼ FOLDING BACKREST (FOR NON-RECLINING BACKS)

STSF0	Folding Backrest for Storage/Travel Access for non-reclining backs	\$240
--------------	--	-------

42. BACK SUPPORT

DSTR-WWHH	Adjustable Tension Back Upholstery, WW x HH	N/C
	▲ Only available with Dual Post armrests (SOADP)	
STDSB-WWHH	Contour back support up to 21" wide, with a maximum shell height of 22".	\$605
	▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back	
STDSB-WWHH	Contour back support up to 21" wide for non-standard sizes with shell height ≥ 23".	\$755
	▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back	
STDSB-WWHH	Contour back support for Heavy-Duty from 22" wide and higher, with any height.	\$1,075
	▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back	



42. BACK SUPPORT (CONTINUED)

▶ MATRX BACKS

E2S-WWHH	MATRX ELITE E2S BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	\$995
	▲ <i>Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>Weight capacity up to 300 lb</i>	
E2D-WWHH	MATRX DEEP CONTOUR ELITE E2D BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	\$1,155
	▲ <i>Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Weight capacity up to 300 lb</i>	
PBE-WWHH-TR	MATRX ELITE TR BACK – (15" to 20" Wide - Back height of 14", 16", 18" and 20")	\$910
	▲ <i>Back shell height is available in 14", 16" 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>Weight capacity up to 300 lb</i>	
PBE-WWHH-HD	MATRX HEAVY DUTY ELITE BACK – (21" to 24" Wide - Back height of 16" and 20")	\$1,150
	▲ <i>Back shell height is available in 16" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>Weight capacity up to 500 lb</i>	
EDB-WWHH-HD	MATRX ELITE DEEP HD BACK – (21" to 24" Wide - Back height of 16", 18" and 20")	\$1,195
	▲ <i>Back shell height is available in 16", 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Weight capacity up to 500 lb</i>	

▶ NXT BACKS

NOACB-WWHH	NXT OPTIMA FH AC CONTOUR BACK – (15" to 21" Wide)	\$690
	▲ <i>Weight capacity up to 300 lb</i>	
NOACB-WWHH	NXT OPTIMA FH AC CONTOUR BACK – (22" to 28" Wide)	\$720
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
NOPCB-WWHH	NXT OPTIMA FH PC DEEP CONTOUR BACK – (15" to 21" Wide)	\$790
	▲ <i>Weight capacity up to 300 lb</i>	
NOPCB-WWHH	NXT OPTIMA FH PC DEEP CONTOUR BACK – (22" to 28" Wide)	\$820
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	

▶ OMISSION OF THE BACK SUPPORT

STDSBO	Omission of the Back Support - (Using an Aftermarket Back)	N/C
	▲ <i>Back support MUST be omitted before selecting aftermarket back hardware</i>	
	▲ <i>Upcharges apply when both a back support AND aftermarket back hardware are selected</i>	
CUST-BHW	Aftermarket Back Hardware - Information on the back must be specified:	\$185
	<i>Back Manufacturer:</i> _____	
	<i>Product Name:</i> _____	
	<i>Hardware being used:</i> _____	
	▲ <i>Required information</i>	
	▲ <i>Upcharges apply when both a back support AND aftermarket back hardware are selected</i>	
	▲ <i>Hardware will not be provided IF the required information for the aftermarket back support is omitted, but a standard T-Nut kit will be supplied</i>	



43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

						▼ 20" is the minimum back height required with a reclining back.
16" (SOBH16)	18" (SOBH18)	20" (SOBH20)	22" (SOBH22)	24" (SOBH24)	26" (SOBH26)	28" (SOBH28)
This height is not available with the NXT Optima AC and PC options						▲

44. SEAT CUSHIONS

► FOR CUSTOM SIZES AND CONFIGURATIONS, PLEASE REFER TO THE PERFORMA CG AIR ORDER FORM.

► PERFORMA CG AIR - THERAPEUTIC AIR CUSHION

Custom Cushion (Please use the Performa CG Air order form)

Cell Height (Choose 1): 2" 3" 4"

► SINGLE (1) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING

► PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER

PCGA1VSCxxx-x 14" to 20" (355 mm to 510 mm) wide \$812

PCGA1VSCxxx-x 22" + (560 mm +) wide \$914

▲ Cushion sizes only come in even numbers for width and depth measurements.

► PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER

PCGA1V3Dxxx-x 14" to 20" (355 mm to 510 mm) wide \$857

PCGA1V3Dxxx-x 22" + (560 mm +) wide \$963

▲ Cushion sizes only come in even numbers for width and depth measurements.

► PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER

PCGA1VICxxx-x 14" to 20" (355 mm to 510 mm) wide \$915

PCGA1VICxxx-x 22" + (560 mm +) wide \$969

▲ Cushion sizes only come in even numbers for width and depth measurements.

► DUAL (2) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING

► PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER

PCGA2VSCxxx-x 14" to 20" (355 mm to 510 mm) wide \$812

PCGA2VSCxxx-x 22" + (560 mm +) wide \$914

▲ Cushion sizes only come in even numbers for width and depth measurements.

► PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER

PCGA2V3Dxxx-x 14" to 20" (355 mm to 510 mm) wide \$857

PCGA2V3Dxxx-x 22" + (560 mm +) wide \$963

▲ Cushion sizes only come in even numbers for width and depth measurements.

► PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER

PCGA2VICxxx-x 14" to 20" (355 mm to 510 mm) wide \$915

PCGA2VICxxx-x 22" + (560 mm +) wide \$969

▲ Cushion sizes only come in even numbers for width and depth measurements.

► DUAL VALVE CUSHION SEGMENTATION (PLEASE SELECT 1 BOX)

PCGA-2VLR Dual Valves with Left / Right Location N/C

PCGA-2VF6 Dual Valves with Front / Back Location, Front distance 6 in (150 mm) from the front N/C

PCGA-2VF8 Dual Valves with Front / Back Location, Front distance 8 in (205 mm) from the front N/C

PCGA-2VF10 Dual Valves with Front / Back Location, Front distance 10 in (255 mm) from the front N/C

► PERFORMA CG AIR OPTIONS

PCGA-GZ GripZone Anti-Slip Technology \$86

PCGA-TT-WW Transfer Trellis \$100

PCGA-AC Anatomic Contour Seat Insert \$100

▲ Only available with 2" and 3" cell heights

PCGA-SP Seat Pouch \$55



44. SEAT CUSHIONS (CONTINUED)

▶ CONTOUR SEAT CUSHION

SCPSC-WWDD	Contour Seat Cushion up to 19" Wide	\$620
SCPSC-WWDD	Contour Seat Cushion for Non-Standard Sizes of 20" to 24" Wide	\$735

▶ MATRX CUSHION

PSP-WWHH	MATRX PSP Seat Cushion – (11" to 20" Wide)	\$645
	▲ The maximum user weight capacity is 300 lb	
	▲ Cushion sizes only come in even numbers for width and depth measurements.	
PSP-WWHH	MATRX PSP Seat Cushion – (22" to 26" Wide)	\$995
	▲ The maximum user weight capacity is 550 lb	
	▲ Cushion sizes only come in even numbers for width and depth measurements.	

▶ NXT BIO-FIT CUSHION

NXTBIOFIT-WWHH	NXT BIO FIT Cushion – (15" to 21" Wide)	\$595
	▲ The maximum user weight capacity is 300 lb	
NXTBIOFIT-WWHH	NXT BIO FIT Cushion – (22" to 26" Wide)	\$625
	▲ The maximum user weight capacity is 400 lb	
	▲ Cushion sizes only come in even numbers for width and depth measurements.	

45. SIGNATURE TAILORED SOLUTIONS - CANTILEVER ARMRESTS - ADJUSTABLE HEIGHT, INWARD/OUTWARD ROTATION AND VERTICAL ANGLE. ALLOWS ARMPAD TO MOVE 1" FORWARD (CHOOSE 2)

FOR NON-RECLINING BACKS

- ▼ **Height Adjustable from 7.25" (185 mm) to 15.25" (385 mm) in 1" (25 mm) increments**
- ▼ **Not compatible with Waterfall Armpads (SOWAPx et SOWGELx)**

SOAFBDLH	Signature Cantilever Armrest – Desk Length 12" (305 mm) - Left Side	N/C
SOAFBDRH	Signature Cantilever Armrest – Desk Length 12" (305 mm) - Right Side	N/C
SOAFBFLH	Signature Cantilever Armrest – Full Length 16" (405 mm) - Left Side	N/C
SOAFBFRH	Signature Cantilever Armrest – Full Length 16" (405 mm) - Right Side	N/C
SOAFBELH	Signature Cantilever Armrest – Extended Length 18" (455 mm) - Left Side	N/C
SOAFBERH	Signature Cantilever Armrest – Extended Length 18" (455 mm) - Right Side	N/C

▼ CUSTOM LENGTH

SOAFBCLH	Signature Cantilever Armrest – Custom Length 20" (510 mm) - Left Side	\$95
SOAFBCRH	Signature Cantilever Armrest – Custom Length 20" (510 mm) - Right Side	\$95
SOAFBCLH	Signature Cantilever Armrest – Custom Length 22" (560 mm) - Left Side	\$95
SOAFBCRH	Signature Cantilever Armrest – Custom Length 22" (560 mm) - Right Side	\$95

FOR RECLINING BACKS

- ▼ **Height Adjustable from 10.25" (260 mm) to 16.25" (415 mm) in 1" (25 mm) increments**
- ▼ **Choose a Height Range :** **Low Range 10.25" to 13.25" (260 - 335 mm)** **High Range 14.25" to 16.25" (360 - 415 mm)**
- ▼ **Not compatible with Waterfall Armpads (SOWAPx and SOWGELx)**

SOAFBDLH-R	Signature Cantilever Armrest for Reclining Back – Desk Length 12" (305 mm) - Left Side	N/C
SOAFBDRH-R	Signature Cantilever Armrest for Reclining Back – Desk Length 12" (305 mm) - Right Side	N/C
SOAFBFLH-R	Signature Cantilever Armrest for Reclining Back – Full Length 16" (405 mm) - Left Side	N/C
SOAFBFRH-R	Signature Cantilever Armrest for Reclining Back – Full Length 16" (405 mm) - Right Side	N/C
SOAFBELH-R	Signature Cantilever Armrest for Reclining Back – Extended Length 18" (455 mm) - Left Side	N/C
SOAFBERH-R	Signature Cantilever Armrest for Reclining Back – Extended Length 18" (455 mm) - Right Side	N/C

▼ CUSTOM LENGTH

SOAFBCLH-R	Signature Cantilever Armrest for Reclining Back – Custom Length 20" (510 mm) - Left Side	\$95
SOAFBCRH-R	Signature Cantilever Armrest for Reclining Back – Custom Length 20" (510 mm) - Right Side	\$95
SOAFBCLH-R	Signature Cantilever Armrest for Reclining Back – Custom Length 22" (560 mm) - Left Side	\$95
SOAFBCRH-R	Signature Cantilever Armrest for Reclining Back – Custom Length 22" (560 mm) - Right Side	\$95



46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)

▼ *NOT compatible with the extra-long armpads (16" x 3 1/2") (SOXWFAP16 or SOXWGAP16) and the Arm Trough with Hand Support (SOCAR)*

SOADPDLH	Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPDRH	Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPFLH	Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPFRH	Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168

48. ARMREST HEIGHT ADJUSTMENT AT FACTORY

Cantilever Armrest Range (SOAFB for Non-Recl. backs)					Set at:					
7.25"	8.25"	9.25"	10.25"	11.25"	12.25"	13.25"	14.25"	15.25"	16.25"	

Cantilever Armrest Range (SOAFB for Reclining backs)					Set at:					
			10.25"	11.25"	12.25"	13.25"	14.25"	15.25"	16.25"	

- ▲ The Armrest Height is taken from the seat pan to the top of the armpad.
- ▲ Not required for Dual Post.

49. SIGNATURE TAILORED SOLUTIONS - ARMPADS

FOAM TYPE SELECTION

SOFAPD	Desk Length Flat Armpad – (10 1/2" x 2 1/2")	Select side:	Left Hand Side	Right Hand Side	STD
SOFAPF	Full Length Flat Armpad – (14" x 2 1/2")	Select side:	Left Hand Side	Right Hand Side	STD
SOXWFAP10	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 10" (255 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C
SOXWFAP14	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 14" (355 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C
SOXWFAP16	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 16" (405 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C

- ▲ *The X-Wide 16" x 3 1/2" Armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more*
- ▲ *The X-Wide 16" x 3 1/2" Armpads are NOT compatible with Dual Post Armrests (SOADPxx)*

SOWAPD	Desk Length Waterfall – (8 1/2" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	N/C
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i>				
SOWAPF	Full Length Waterfall – (13" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	N/C
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i>				

HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)

SOGELD	Desk Length Hybrid Gel and Foam Flat armpad – (9 3/4" x 2")	Select side:	Left Hand Side	Right Hand Side	\$40
SOGELF	Full Length Hybrid Gel and Foam Flat armpad – (13 1/2" x 2")	Select side:	Left Hand Side	Right Hand Side	\$40
SOXWGAP10	X-Wide 3.5" (90 mm) Flat GEL Armpad, 10" (255 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$95
SOXWGAP14	X-Wide 3.5" (90 mm) Flat GEL Armpad, 14" (355 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$95
SOXWGAP16	X-Wide 3.5" (90 mm) Flat GEL Armpad, 16" (405 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$105

- ▲ *The X-Wide 16" x 3 1/2" Armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more*
- ▲ *The X-Wide 16" x 3 1/2" Armpads are NOT compatible with Dual Post Armrests (SOADPxx)*



49. SIGNATURE TAILORED SOLUTIONS - ARMPADS (CONTINUED)

SOWGELD	Desk Length Hybrid Gel & Foam Waterfall armpad – (8 ½" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i>				
SOWGELF	Full Length Hybrid Gel and Foam Waterfall armpad – (13" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i>				
▶ SPECIALTY TYPE SELECTION					
SOCAR	Arm Trough with Hand Support (18" x 4 ½")	Select side:	Left Hand Side	Right Hand Side	\$290
	▲ <i>Arm Trough Armpads may interfere with Armrests Flipping back.</i>				
▲ <i>SOCAR armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more</i>					

50. SIGNATURE TAILORED SOLUTIONS - FOOTRESTS

SOSAAF	Signature Swing-Away Footrests				STD
		Set at:	60°	70°	
		Select Range:	SOFLA12 (12 ½" - 17 ½")	SOFLA14 (14 ½" - 19 ½")	
	▲ <i>This option is not compliant with ISO 7176-19</i>				
SOAAF	Signature Heavy Duty Swing-Away Footrests				\$175
		Set at:	60°	70°	
		Select Range:	SOFLHD12 (12 ½" - 17 ½")	SOFLHD14 (14 ½" - 19 ½")	
ADDML	Signature Manual Elevating & Articulating Legrests - Pair				\$595
		Select Range:	SOFLME12 (12 ½" - 17 ½")	SOFLME14 (14 ½" - 19 ½")	
	▲ <i>Does not include calf pads. Please select calf pads in section 51.</i>				
ADDPL	Signature Power Elevating & Articulating Legrests				\$2,495
		Select Function:	Individual	Combined	
		Select Range:	SOFLPW12 (12 ½" - 17 ½")	SOFLPW14 (14 ½" - 19 ½")	
	▲ <i>This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.</i>				

51. SIGNATURE TAILORED SOLUTIONS - FOOTPLATES & ACCESSORIES

SOSFAF	Flip-up Footplates				STD
SOAFP	Angle Adjustable & Flip-up Footplates – (Each)				\$65
SOAHL	Heel Loops – (Each)				\$20
SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests – (SOAAF) – with adjustable hardware – (Each)				\$93
SOAMA	Swing-away residual limb Support	Select:	Left leg	Right leg	\$455

52. FOOTREST LENGTH ADJUSTMENT AT FACTORY

Footrest / Legrest Length	Set at:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"
		16.5"	17"	17.5"	18"	18.5"	19"	19.5"	

▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

53. SIGNATURE TAILORED SOLUTIONS - CENTER MOUNTS

SOAAF90	Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°)				\$200
		Select Foot Platform Range:	SOFBA12 (12 ½" - 16 ½")	SOFBA16 (16 ½" - 19 ½")	
	▲ <i>When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 18.75"</i>				
	▲ <i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>				
	▲ <i>With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.</i>				



53. SIGNATURE TAILORED SOLUTIONS - CENTER MOUNTS (CONTINUED)

ADDPF	90° Power Center Mount Elevating & Articulating Foot Platform	\$2,495
	Select Footrest Range: SOFBPW12 (10 ½" - 16") SOFBPW16 (16" - 20")	
▲	<i>When selecting SOFBPW16 (16" - 20") range, the lowest possible seat to floor height is 18.75"</i>	
▲	<i>This option includes calf pads. The power center mount counts as one power function.</i>	
▲	<i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>	
▲	<i>For leg lengths between 10 ½" and 12 ½", there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.</i>	
▲	<i>The starting angle is reduced from 90° to 85° for some seat configurations when used with the narrow base configuration (FS2424B)</i>	
▲	<i>The E-BX (MAC) option is required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2424B)</i>	

54. SIGNATURE TAILORED SOLUTIONS - FOOT PLATFORMS & ACCESSORIES

SOFUFP	Flip-up one piece Foot Platform	\$140
SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount – (Each)	\$65
SOACPCM	Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each)	\$93

55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY

Fixed Center Mount Foot Platform Length		Set at:									
12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"		
17.5"	18"	18.5"	19"	19.5"							
Power Center Mount Foot Platform Length		Set at:									
10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"		
15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		
▲	The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates										

56. SEAT POSITIONING

SOAUTO	55" Seat Positioning Strap	STD
SOBPHP2C	Bodypoint padded Hip Belt – (Two Point – Center Pull)	\$170
SOPHB	Padded Hip Belt, 2-PT Magnetic Center Pull	\$194
SOBPCS	Bodypoint padded Chest Strap	\$250
SOPCS	Padded Chest Strap, 2-PT Magnetic Center Pull	\$262

57. HEADREST

▼	SELECT ONE OF THE HEADREST OPTIONS WITH THE MOUNTING HARDWARE AND SPECIFY WHERE THE HEADREST WILL BE INSTALLED.	
SONUAHR-12	Using headrest mounting hardware other than SOMAH below. Available with both Fixed and Reclining backrests. 1/2 IN (12.7 MM) Mounts	N/C
SONUAHR-15	Using headrest mounting hardware other than SOMAH below. Available with Reclining backrests only. 19/32 IN (15 MM) Mounts	N/C
SOFHRS	Small Curved Headrest Pad – 4" x 9"	\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"	\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"	\$155
SOMAH	Multi-axis Removable Mounting Hardware	\$220
SOHRBC	Bracket to Mount on the Contour Back Support – (Without Recline)	N/C
SOHRPR	Bracket to Mount on the Seat Frame – (For Reclining Backs)	N/C



58. LATERAL TRUNK SUPPORTS

▼ SELECT ONE OF THE LATERAL PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.

SOSLS3	Adjustable Laterals with Small Curved Pads – (Each) – 3" x 5"	\$100
SOSLM4	Adjustable Laterals with Medium Curved Pads – (Each) – 4" x 6"	\$100
SOSLL5	Adjustable Laterals with Large Curved Pads – (Each) – 5" x 7"	\$100
SOSLHM	Adjustable Swing-Away Hardware Mounts – (Each)	\$220
▲	<i>The Lateral Trunk Supports are not compatible with the NXT Optima FH PC 6" Deep Contour Fixed Thoracic Back Support (NOPCB).</i>	
▲	<i>Waterfall armpads (SOWAPx & SOWGELx) may interfere with flipping back the armrests when combined with Lateral Trunk supports</i>	

59. HIP & THIGH SUPPORTS

▼ SELECT ONE OF THE HIP & THIGH SUPPORT PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.

SOHTSS10	Adjustable Hip & Thigh Support Pads – (Each) – Small 10" x 4"	\$130
SOHTSM12	Adjustable Hip & Thigh Support Pads – (Each) – Medium 12" x 5"	\$130
SOHTSL14	Adjustable Hip & Thigh Support Pads – (Each) – Large 14" x 6"	\$130
SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware – (Each)	\$136

60. ACCESSORIES

SOCH	Cup Holder – (Self-Leveling)	Install on:	Left Side	Right Side	\$115
SOCPH	Ram Mounts X-Grip Cell Phone Holder	Install on:	Left Side	Right Side	\$190
SOTHS	Ram Mounts X-Grip Small Tablet Holder	Install on:	Left Side	Right Side	\$205
SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	Install on:	Left Side	Right Side	\$360
▲	<i>Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.</i>				
▲	<i>With multiple accessories installed on the same side, a custom upgrade will be required.</i>				

61. REFLECTOR KIT

REFLECT	Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner's Manual)	N/C
----------------	---	-----

62. BACK PACK AND ARMREST SIDE POUCHES

SOBP	Back Pack – (Comes with durable straps that can be hooked around the back canes)				\$90
SOAP	Armrest Side Pouch – (Comes with durable straps that can be hooked around the armrest tubes)	Order:	Left Side	Right Side	\$53
▲	<i>If only choosing one (1) armrest pouch, please show what side it should be installed</i>				

63. SWITCHES AND OPTIONS

ASCPWIFI	PGDT - Wireless programmer, WIFI dongle	\$844
ASCPCP	PGDT R-NET Dongle and PC programmer kit	\$600
ASCPVR2	PGDT VR2 with PPI Dongle and PC programmer kit	\$1,450
SSBS	Buddy Button Switch	\$230
SSES	Egg Switch-Black	\$116
SSML	Ultra Light Switch	\$295
ESS-RED	Emergency STOP Switch - Red	\$314



SPECIALTY DRIVER CONTROL OPTIONS

- ▶ Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)
- ▶ Specialty Control options in this section are not compliant with ISO 7176-19

64. THERAFIN SIP-N-PUFF KIT

▼ *This option is not compliant with ISO 7176-19*

SCSNPK	Therafin Sip-N-Puff Kit				\$1,975
SCSNPBTK	Therafin Sip-N-Puff Breath Tube Kit				N/C
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRMSNP	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	SWL (Swings on Left)	SWR (Swings on Right)		\$300
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Therafin Sip-N-Puff kit.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Therafin Sip-N-Puff kit.				
▶ CHOOSE ONE (1) HEADREST PAD					
SOFHRS	Small Curved Headrest Pad – 4" x 9"				\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"				\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"				\$155

65. PROXIMITY HEAD ARRAY

▼ *This option is not compliant with ISO 7176-19*

CUSTOM	Adult – 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse				RFQ
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Proximity Head Array.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Proximity Head Array.				

66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL

▼ *This option is not compliant with ISO 7176-19*

CUSTOM	Light touch mini Proportional Joystick				RFQ
SCSRMMP	Stealth Rock-n-Lock Swingaway Arm for Chin MP	SWL (Swings on Left)	SWR (Swings on Right)		\$870
	▲ Required with the Light touch mini Proportional Joystick.				
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section				
JCC	Chin Cup				\$55



66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL (CONTINUED)

SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ <i>Required with the Light touch mini Proportional Joystick.</i>				
▶ CHOOSE ONE (1) HEADREST PAD					
SOFHRS	Small Curved Headrest Pad – 4" x 9"				\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"				\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"				\$155

67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT

▼ *This option is not compliant with ISO 7176-19*

SCPGMJS	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount				\$2,369
SCGMJSM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for JC200	MLH (Mount on Left)	MRH (Mount on Right)		\$595
	▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i>				
CUSTOM	Light touch mini Proportional Joystick				RFQ
SCGMPPM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for Light Touch Joystick	MLH (Mount on Left)	MRH (Mount on Right)		\$825
	▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i>				
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)		\$850
	▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i>				

68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT

▼ *This option is not compliant with ISO 7176-19*

SCPGMJS-M	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount				\$2,369
SCARMJSM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)		\$325
	▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i>				
CUSTOM	Light touch mini Proportional Joystick				RFQ
SCARMPM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)		\$375
	▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i>				
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)		\$850
	▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i>				

69. PROPORTIONAL HEAD CONTROL PACKAGE

▼ *This option is not compliant with ISO 7176-19*

SCPHC	Proportional Head Control Package				\$4,495
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i>				
	▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section</i>				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ <i>Required with the Proportional RIM Head Control.</i>				



CUSTOM SECTION

▶ Lead times will increase depending on the custom order pipeline at the time the purchase order is received.
 ▶ Please call for Quote and specifications for these options.

70. O² HOLDER

▼ *This option is not compliant with ISO 7176-19*

SUPOXY	Custom Non-Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$295
	▲ <i>Not available with a power recline</i>			
SUPOXYA	Custom Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$1,195
	▲ <i>Required with a power recline</i>			
	D-Size	Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve		
	E-Size	Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve		

71. VENT TRAY (PLEASE COMPLETE IN FULL)

▼ **AMYLOR DOES NOT SUPPLY THE BATTERY, CHARGER OR THE WIRING FOR THE VENT. THESE ITEMS ARE NORMALLY PROVIDED BY THE VENT MANUFACTURER.**

▼ *This option is not compliant with ISO 7176-19*

SOVENT	Custom Vent Tray		
	Articulating Vent Box (Required with Power Tilt and/or Recline)		\$1,360
	Fixed Vent Box		\$400
	SPECIFY MODEL		
	Trilogy 100/200		
	LTV		
	LTV with Sprintpack		
	Other, specify: _____		
SOVTBE	Extended Battery Compartment		\$495
	▲ <i>The Vent Tray option is not available on narrow options N2424 & N2221 (FS2424B & FS2221B)</i>		
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries		
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries		
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries		
	The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries		
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries		



72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)

▼ **PLEASE SPECIFY DETAILS AND SPECIFICATIONS FOR THESE OPTIONS.**

▼ **These following custom options and accessories are ISO 7176-19 compliant**

CUSTOM	USB Charger Port ◀ Specify the location for the installation: _____	RFQ
CUSTOM	Tilt + Elevate Function Through VR2 Non-Expandable Drive Control	RFQ
CUSTOM	Accessory Hooks located at the rear of the seat structure (2 hooks)	RFQ
CUSTOM	Heavy-Duty Wheelchair Cover, Large	RFQ
CUSTOM	Transfer Handles – (Pair) Select a height range: 5" to 8" 6" to 9"	RFQ

▼ **The following custom options and accessories below are not compliant with ISO 7176-19**

CUSTOM	CUSTOM Lithium ion - 54 Amps (represents one 24V Battery pack, equivalent to two Group 22 Batteries)	RFQ
CUSTOM	CUSTOM Lithium ion - 78 Amps (represents one 24V Battery pack, equivalent to two Group 24 Batteries)	RFQ
CUSTOM	CUSTOM 8A Offboard Battery Charger for Lithium ion batteries	N/C
	▲ This option is required when choosing custom Lithium ion batteries	
CUSTOM	CUSTOM Lithium ion Battery capacity LCD display gauge	N/C
	▲ This option is required when choosing custom Lithium ion batteries	
CUSTOM	Folding Backrest for Storage/Travel Access for Reclining backs	RFQ
CUSTOM	Power Anterior Tilt	RFQ
CUSTOM	Power Lateral Tilt 15° Raise the Left Side only Raise the Right Side only Raise the Left & Right Side	RFQ
CUSTOM	Power Latitude System (Features 14" of forward travel and a 5" seat to floor height from the pan to the ground) (200 Lb max. user weight capacity)	RFQ
CUSTOM	Power Elevating Stretcher – (Please call CSR for availability with configuration and user weight capacity)	RFQ
CUSTOM	Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity)	RFQ
CUSTOM	Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity)	RFQ
CUSTOM	Rearview Camera with 5" Monitor	RFQ
CUSTOM	Rearview mirror Left Side Right Side	RFQ
CUSTOM	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan	RFQ
CUSTOM	Foot Box ◀ Specify Length: _____ Width: _____ Please confirm with CSR on sizing possibilities	RFQ
CUSTOM #1	Other custom options. Specify in notes section below	RFQ
CUSTOM #2	Other custom options. Specify in notes section below	RFQ



73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS

A large, empty rectangular area intended for handwritten or typed notes, special requests, or custom specifications.