



1. ACCOUNT INFORMATION

| | | | |
|----------------------------|-------|--------------|-------|
| P.O. # or Quote: | _____ | | |
| Dealer Contact Name: | _____ | | |
| Email Address: | _____ | | |
| Phone Number: | _____ | Fax Number: | _____ |
| Bill to (Name): | _____ | | |
| Bill to (Location): | _____ | | |
| Ship to (Name): | _____ | | |
| Address: | _____ | | |
| City: | _____ | | |
| Province: | _____ | Postal Code: | _____ |
| User Name or Number: | _____ | | |
| User Weight: | _____ | | |
| User Height: | _____ | | |
| Seat Cushion Thickness: | _____ | | |
| Knee to Heel Length: Left: | _____ | Right: | _____ |



THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:

- A user name or number, the client's measurements and weight - Amylor is not responsible for configuration or size discrepancies resulting from customer errors **or omissions** on the order form
 - Order confirmations will be provided by Amylor summarizing the items selected. Please review this carefully.
 - Once the order is shipped, any changes are subject to the returns policy and may be prohibited.
- *IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES**

2. STANDARD OPTIONS

Alltrack M3 (HD) Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 6" front casters, and standard auto style buckle seat positioning strap, 6.5 mph motors, Weight Capacity = 450 lb

3. ALLTRACK M3HD POWER WHEELCHAIR

| | | |
|-----------------|-----------------------------------|---------|
| ATM3-3HD | ALLTRACK M3 (HD) Power Wheelchair | \$8,545 |
|-----------------|-----------------------------------|---------|

3.5 OWNER'S MANUAL

| | | | | | |
|------------------|-------------------|-----|-----------------|-------------------|-----|
| SOOMATENG | English / Anglais | N/C | SOOMATFR | Français / French | N/C |
|------------------|-------------------|-----|-----------------|-------------------|-----|

4. USER WEIGHT PACKAGES

| | | | | | |
|-------------|-------------------------------------|-----|--------------|-------------------------------------|-----|
| WCA4 | 400 lb Maximum user weight capacity | STD | WCA45 | 450 lb Maximum user weight capacity | STD |
|-------------|-------------------------------------|-----|--------------|-------------------------------------|-----|

5. STANDARD COLORS

| | | | | | |
|---------------|------------------|-----|---------------|---------------------|-----|
| SCMOBM | Onyx Black Matte | N/C | SCMCAB | Cabernet | N/C |
| SCMCTB | Cobalt Blue | N/C | SCMBGM | Bubble Gum Metallic | N/C |
| SCMIOB | Indigo Blue | N/C | SCMHAB | Habanero Metallic | N/C |
| SCMCRE | Candy Red | N/C | SCMLIM | SubLime | N/C |

▼ OPTIONAL COLORS

| | | |
|----------------|-----------------------------------|-------|
| FCHDD-M | Desert Digital Camo HydroGraphics | \$795 |
| FCHCF-M | Carbon fiber HydroGraphics | \$795 |
| FCHWL-M | Woodland Camo HydroGraphics | \$795 |



5.5 WHEEL INSERT COLORS

| | | | | | |
|---|-------------|-----|-----------------|---------------------|-----|
| SCCTB-WI | Cobalt Blue | N/C | SCHAB-WI | Habanero Metallic | N/C |
| SCIOB-WI | Indigo Blue | N/C | SCLIM-WI | SubLime | N/C |
| SCCRE-WI | Candy Red | N/C | SCBGM-WI | Bubble Gum Metallic | N/C |
| SCCAB-WI | Cabernet | N/C | | | |
| ▲ <i>Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)</i> | | | | | |
| ▲ <i>Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)</i> | | | | | |
| ▲ <i>Wheel inserts are not available in Onyx Black Matte</i> | | | | | |

7. TRANSPORT OPTION

| | | |
|-------------|--|-----|
| TTD4 | Tie-down attachment points compliant with ISO 7176-19. See owner's manual for additional details | STD |
|-------------|--|-----|

8. MOTOR PACKAGE

| | | |
|--|--|-------|
| MP65B | 6.5 MPH 4-Pole Motor Package / with Brake Release | STD |
| ▲ <i>Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)</i> | | |
| ▲ <i>The Maximum user weight with this option is 400lb</i> | | |
| MP65-OPT | 6.5 MPH 4-Pole Motor MP65-OPT Package / with COG Release | \$395 |
| ▲ <i>The Maximum user weight with this option is 400lb</i> | | |
| MP544B | 5.4 MPH 4-Pole Motor Package / with Brake Release | N/C |
| ▲ <i>Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)</i> | | |
| MP544-OPT | 5.4 MPH 4-Pole Motor Package / with COG Release | \$395 |

9. ASSISTIVE DRIVING PACKAGE

| | | |
|--|---------------------------------------|---------|
| SMARTT | Smart-Track Assistive Driving Package | \$1,420 |
| ▲ <i>Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.</i> | | |
| ▲ <i>Requires the R-net PM 120 Amp Expandable Controller (ECR120).</i> | | |
| ▲ <i>This option requires the MP65-OPT or MP544-OPT motor package.</i> | | |
| ▲ <i>Recommended for non-proportional Drive Controls</i> | | |

10. BATTERY CHARGER

| | | |
|-------------|-----------------------------|-----|
| BCOB | 8A Offboard Battery Charger | N/C |
|-------------|-----------------------------|-----|

11. BATTERY TRAY AND BASE WIDTH

| | | |
|--------------|--|-----|
| FS24B | Battery Tray for Group 24 Batteries – (Total base width = 25 3/8") | N/C |
|--------------|--|-----|

12. BATTERIES - REQUIRES A SELECTION

| | | |
|-----------------|--|-------|
| BAT24GEL | Premium CG Power - GEL Battery – Group 24 – 70 Amps (Each) | \$465 |
| BATO | Omission of the batteries | N/C |



14. DRIVE WHEEL

| | | |
|--|---|-------|
| DW14PT | 14" Cast Wheel Pneumatic Tires | STD |
| DW14PI | 14" Cast Wheel Pneumatic with Foam Filled Tires | N/C |
| DW144PT | 14" x 4" Extra Wide Cast Wheels, Pneumatic Tires | \$200 |
| DW144PI | 14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires | \$595 |
| ▲ <i>Increases the overall width of the base by 2 3/4"</i> | | |
| DW12PT | 12 1/2" Cast Wheel Pneumatic Tires | N/C |
| DW12PI | 12 1/2" Cast Wheel Pneumatic with Foam Filled Tires | N/C |
| ▲ <i>12 1/2" Drive wheels lower ground clearance and speed, but increase motor torque.</i> | | |

15. FRONT CASTERS

| | | |
|---|---|-----|
| FC6S | 6" Solid Casters – (Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI) | STD |
| FC8PTM | 8" Pneumatic Casters | N/C |
| FC8PIM | 8" Pneumatic Casters with Foam Filled Tires | N/C |
| ▲ <i>The 8" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)</i> | | |

16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± 1/4")

▶ **YELLOW BOXES REPRESENT OPTIONS THAT CAN ONLY BE SELECTED WITH 12 1/2" DRIVE WHEELS (SECTION 14)**

| SEAT TO FLOOR HEIGHT TABLE SEAT ANGLE | NO POWER OPTIONS | | | TILT (WITH OR WITHOUT RECLINE) | HD ELEVATE (WITH OR WITHOUT RECLINE) | TILT AND HD ELEVATE (WITH OR WITHOUT RECLINE) |
|--|------------------|----|----|--------------------------------|--------------------------------------|---|
| | 0° | 3° | 6° | | | |
| STFH1550 (15.50" STFH) | | | | | | |
| STFH1650 (16.50" STFH) | | | | | | |
| STFH1700 (17" STFH) | | | | | | |
| STFH1750 (17.50" STFH) | | | | | | |
| STFH1800 (18" STFH) | | | | | | |
| STFH1850 (18.50" STFH) | | | | | | |
| STFH1900 (19" STFH) | | | | | | |
| STFH1950 (19.50" STFH) | | | | | | |
| STFH2000 (20" STFH) | | | | | | |

17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)

| | | |
|---|---|-----|
| SOMHSS | Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges. | RFQ |
| ▲ <i>Required at: _____" (± 1/4" Range)</i> | | |

18. PRE-SET SEAT FRAME ANGLE

| | | |
|-------|---|-----|
| SFE-0 | 0° Pre-set seat frame seat angle (Must be selected when using power options) | STD |
| SFE-3 | 3° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |
| SFE-6 | 6° Pre-set seat frame seat angle (Only available with a seat configuration without power options) | N/C |



19. LIGHTING PACKAGE

| | | |
|---|---|---------|
| EOLP | LED Light Package – (Includes front and rear lights with integrated flashers) | \$1,425 |
| ▲ <i>Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL, ECRJCL, ECRJC35B or EDIB).</i> | | |

20. NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

| | | |
|--|--|-------|
| ▶ <i>VR2 Electronics can only accommodate specific power seat functions and combinations such as Tilt / Elevate / Recline / Power Center Mount / Combined Power Legs / Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs / Tilt + Elevate (Custom)</i> | | |
| ECV902A | VR2 90 Amp Drive Controller – (Requires the selection of VR2 joystick below) | STD |
| ECVJ4K | VR2 Standard 4-Key Joystick | STD |
| ECVJ6K | VR2 6-Key Joystick (Allows for single or select dual power options to operate through this joystick) | N/C |
| ECVJ9KL | VR2 9-Key Joystick with Lighting Operation (Allows for single or select dual power options to operate through this joystick) | \$150 |

21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

| | | |
|---------------|--|-------|
| ECEL90 | R-net EL 90 Amp Expandable Controller | \$175 |
| ECR120 | R-net PM 120 Amp Expandable Controller | \$275 |

22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

| | | |
|---|--|---------|
| ▶ <i>SELECTION OF ONE JOYSTICK OPTION BELOW IS REQUIRED WITH THE EXPANDABLE CONTROLLER ECEL90 OR ECR120</i> | | |
| ECO | Omit Joystick – (For use with an Omni 2 Display) | N/C |
| ECELJ | R-net - (LED) – Joystick | STD |
| ▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i> | | |
| ECRJC | R-net Joystick with Color screen | \$155 |
| ▶ <i>JOYSTICKS WITH LIGHTING OPERATION</i> | | |
| ECELJL | R-net – (LED) – Joystick with Lighting Operation | \$125 |
| ▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i> | | |
| ECRJCL | R-net Joystick with Color screen and Lighting Operation | \$180 |
| ▶ <i>ADVANCED JOYSTICK</i> | | |
| ECRJC35B | R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks | \$1,105 |

23. JOYSTICK MOUNT

| | | |
|----------------|---|-------|
| JSMF | Fixed Joystick Mount | N/C |
| JMSA | Swing Away Joystick Mount – Height adjustable, Angle adjustable and Multi-Axis | \$235 |
| Select: | <i>Set-up with standard installation to swing-away sideways</i> | |
| | <i>Set-up at 35° to swing-away sideways and downwards</i> | |
| | ▲ <i>Set-up at 35° angle may cause interference with the user's limbs when the new Signature cantilever armrest (SOAFBxxx) is internally rotated.</i> | |
| JSMCM | Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount | \$275 |
| | ▲ <i>This option is not compliant with ISO 7176-19</i> | |

24. JOYSTICK MOUNT SUPPORT

| | | |
|--------------|--------------------------------|-----|
| JSMFH | Joystick mount support Bracket | STD |
|--------------|--------------------------------|-----|

25. JOYSTICK HAND MOUNT

| | | | |
|------------|------------------------------------|------------|-------------------------------------|
| HML | Install Joystick on Left hand side | HMR | Install Joystick on Right hand side |
|------------|------------------------------------|------------|-------------------------------------|



26. JOYSTICK HANDLE

| | | | | | |
|------------|------------------------------|-------|--------------|---------------------------|-------|
| JSH | Standard Joystick Gimble | STD | JBAT | Stick Handle 3" in Length | \$130 |
| JBH | Ball Handle | \$85 | JFBH | Foam Ball Handle | \$70 |
| JTH | T-Handle | \$155 | JMUSH | Bodypoint Mushroom Handle | \$210 |
| JGP | Bodypoint 4" U-shaped Handle | \$210 | | | |

27. ATTENDANT DRIVE CONTROL OPTION

| | | |
|-------------|------------------------|---------|
| EADC | Attendant Hand Control | \$1,150 |
|-------------|------------------------|---------|

28. ATTENDANT JOYSTICK HAND MOUNT

| | | | |
|--------------|---|--------------|--|
| EADCL | Install Attendant Joystick on Rear Left hand side | EADCR | Install Attendant Joystick on Rear Right hand side |
|--------------|---|--------------|--|

29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)

| | | |
|-------------|--|---------|
| EDIB | Omni 2 – (Enhanced Display & Digital Interface Box with Infrared) | \$1,915 |
| | ▲ <i>Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120)</i> | |
| | ▲ <i>The Omni 2 display mount must be selected from the option below.</i> | |

30. OMNI 2 ENHANCED DISPLAY MOUNT

| | | |
|--------------|--|--|
| EDSAM | User Swing Away Mount | N/C |
| | EDML – (Install Display on LH side) | EDMR – (Install Display on RH side) |

31. ENVIRONMENTAL CONTROL MODULE

| | | |
|------------|---|-------|
| ECU | Environment control module – (Output Module) | \$995 |
| | ▲ <i>Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).</i> | |

33. SIGNATURE TAILORED POWER POSITIONING SOLUTIONS / REHAB SEAT

| | | |
|--------------------|---|----------|
| SAABAHD-450 | Signature HD Rehab Seat – (Seat Frame only) | \$1,045 |
| SAPTXHRD | Signature Heavy-Duty CG Power Tilt | \$6,195 |
| | ▲ <i>This option includes the Signature Rehab Seat – (SAABAHD)</i> | |
| SAPRXHD | Signature Heavy-Duty 165° Power Recline with Mechanical Shear Reduction (MSR) | \$7,845 |
| | ▲ <i>This option includes the Signature Rehab Seat – (SAABAHD)</i> | |
| SAPTPRXHRD | Signature Heavy-Duty CG Power Tilt and 165° Power Recline with Mechanical Shear Reduction (MSR) | \$12,995 |
| | ▲ <i>This option includes the Signature Rehab Seat – (SAABAHD)</i> | |
| SAPESXHRD | Heavy-Duty ALTA 11" Power Seat Elevation | \$5,495 |
| | ▲ <i>This option includes the Signature Rehab Seat – (SAABAHD)</i> | |
| | ▲ <i>The maximum user weight with this option is 400 lb</i> | |
| ADDPESXHRD | Module – Add Heavy-Duty ALTA 11" Power Seat Elevation to either SAPTXHD, SAPRXHD or SAPTPRXHD | \$3,995 |
| | ▲ <i>The maximum user weight with this option is 400 lb</i> | |



34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER

| | | | | | | | | | |
|---|---|-------|-----|-------|--------|-----------------|----------------------|----------------------|---------|
| ▶ SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (For both expandable and non-expandable systems) | | | | | | | | | |
| SOLO | Single function controller (Includes an external toggle switch) | | | | | | | | STD |
| ▼ OPTIONAL SWITCHES | | | | | | | | | |
| SWTOGGLE | Single Toggle Switch | | | | | | | Attendant Activation | \$100 |
| SWPUSH1 | Single Push Button Switch – (Single Mode: SAC) | | | | | User Activation | Attendant Activation | \$90 | |
| SWPUSH2 | Dual Push Button Switch – (Dual Mode: SAC) | | | | | User Activation | Attendant Activation | \$180 | |
| SSBS | Buddy Button Switch | Black | Red | Green | Yellow | User Activation | Attendant Activation | \$230 | |
| SSES | Egg Switch-Black | | | | | User Activation | Attendant Activation | \$116 | |
| SSML | Ultra Light Switch | | | | | User Activation | Attendant Activation | \$295 | |
| ▶ SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL | | | | | | | | | |
| EPSV | Single Function Through VR2 Non-Expandable Drive Control | | | | | | | | N/C |
| | ▲ <i>Allows for activation of a single power option through VR2 controller</i> | | | | | | | | |
| | ▲ <i>Does not allow for control of power option through drive control and separate switches</i> | | | | | | | | |
| | ▲ <i>Control of power seat elevation with tilt is available as a custom. Check appropriate box in custom section 72 in addition to EPSV.</i> | | | | | | | | |
| | ▲ <i>VR2 Electronics can only accommodate specific power seat functions such as Tilt / Recline / Power Center Mount / Combined Power Legs /Tilt +Elevate (Custom)</i> | | | | | | | | |
| | ▲ <i>Choose EPSV2 to operate select power combinations</i> | | | | | | | | |
| EPSEL | Single Function Through Drive Control with R-net 90-amp controller | | | | | | | | N/C |
| | ▲ <i>Allows for activation of single power option through R-net 90-amp controller</i> | | | | | | | | |
| | ▲ <i>Compatible with all R-net joysticks</i> | | | | | | | | |
| | ▲ <i>Does not allow for control of power option through drive control and separate switches</i> | | | | | | | | |
| | ▲ <i>Can not operate more than one power option</i> | | | | | | | | |
| | ▲ <i>If selecting memory seating, choose SAC below</i> | | | | | | | | |
| | ▲ <i>If selecting 120-amp R-net controller, choose SAC below</i> | | | | | | | | |
| SAC | Single Function Through Drive Control with R-net 90-amp or 120-amp controller | | | | | | | | N/C |
| | ▲ <i>Allows for activation of single power option through R-net drive controls</i> | | | | | | | | |
| | ▲ <i>Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls</i> | | | | | | | | |
| | ▲ <i>Capable of operating tilt with power seat elevation & compatible with memory seating</i> | | | | | | | | |
| ▶ SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE) | | | | | | | | | |
| SACMEMORY | Memory Seat Program | | | | | | | | \$1,495 |
| | ▲ <i>This option is not compatible with any VR2 or LED joysticks (ECELJ and ECELJL)</i> | | | | | | | | |
| | ▲ <i>Memory seat programming is not compatible with the Heavy Duty Seat Elevate option (SAPESXHD and ADDPESXHD)</i> | | | | | | | | |
| ▶ SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH | | | | | | | | | |
| SACTOGGLE | Single Toggle Switch | | | | | | | Attendant Activation | \$410 |
| SACPUSH1 | Single Push Button Switch – (Single Mode: SAC) | | | | | User Activation | Attendant Activation | \$400 | |
| SACPUSH2 | Dual Push Button Switch – (Dual Mode: SAC) | | | | | User Activation | Attendant Activation | \$490 | |
| SACSSBS | Buddy Button Switch | Black | Red | Green | Yellow | User Activation | Attendant Activation | \$540 | |
| SACSES | Egg Switch-Black | | | | | User Activation | Attendant Activation | \$426 | |
| SACSSML | Ultra Light Switch | | | | | User Activation | Attendant Activation | \$605 | |



35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE

▶ **DUAL SEAT FUNCTIONS THROUGH DRIVE CONTROL**

| | | |
|--------------|--|-----|
| EPSV2 | Dual Function Through Drive Control with VR2 Non-Expandable Electronics | N/C |
| | ▲ <i>Allows for the activation of 2 power options through the VR2 Controller. Combination of TILT + RECLINE or ELEVATE + RECLINE are not possible with EPSV2</i> | |
| | ▲ <i>Does not allow for control of power options through drive control and separate switches</i> | |
| | ▲ <i>VR2 Electronics can only accommodate specific power seat combinations such as Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs</i> | |

▶ **MULTIPLE SEAT FUNCTIONS THROUGH R-NET**

| | | |
|------------|--|---------|
| MAC | E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators | \$2,295 |
|------------|--|---------|

▶ **MAC OPTIONS**

| | | |
|------------------|---|-------|
| MACMEMORY | Memory Seat Program | \$175 |
| | ▲ <i>Memory seat programming is not compatible with the Heavy Duty Seat Elevate option (SAPESXHD and ADDPESXHD)</i> | |

▶ **SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH**

| | | | | |
|--------------------|--|-----------------|----------------------|-------|
| MACKEYPAD-P | Keypad w/ 5 Push Button Switches | User Activation | Attendant Activation | \$445 |
| MACKEYPAD-T | Keypad w/ 5 Toggle Switches | User Activation | Attendant Activation | \$445 |
| MACMULTI | 4-Way Toggle Switch – (Up to 4 Activations) | User Activation | Attendant Activation | \$180 |
| | ▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i> | | | |
| MACMULTIT | Triple Toggle Switch – (Up to 5 Activations) | User Activation | Attendant Activation | \$135 |
| | ▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i> | | | |
| MACPUSH4 | 4-Way Push Button Switch | User Activation | Attendant Activation | \$225 |
| | ▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i> | | | |
| MACPUSH5 | 5-Way Push Button Switch | User Activation | Attendant Activation | \$240 |
| | ▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i> | | | |

36. SWITCH ACTIVATION HAND MOUNT

| | | | |
|----------------------|---|---|-----|
| Switch Installation: | TSML – (Install Switch on LH side) | TSMR – (Install Switch on RH side) | N/C |
|----------------------|---|---|-----|

37. SIGNATURE TAILORED SOLUTIONS - PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)

| | | | | | | | | | |
|--------------|-------------|-----|-----|-----|-----|------|------|------|-----|
| SOBAA | Select one: | 86° | 90° | 94° | 98° | 102° | 106° | 110° | N/C |
|--------------|-------------|-----|-----|-----|-----|------|------|------|-----|

38. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT WIDTH RANGE BETWEEN 17” AND 22” (430 - 560 MM)

| | | |
|---------------------|--|---------------------|
| STS17-22-19 | Width adjustment from 17” to 22” (430-560 MM) and Depth range of 15” to 19” (380 - 480 MM) | N/C |
| STS17-22-22 | Width adjustment from 17” to 22” (430-560 MM) and Depth range of 18” to 22” (455 - 560 MM) | N/C |
| | ▲ <i>The minimum seat depth available with this option is 18” (455 MM)</i> | |
| 17” (SOSW17) | 18” (SOSW18) | 19” (SOSW19) |
| 20” (SOSW20) | 21” (SOSW21) | 22” (SOSW22) |

39. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT WIDTH RANGE BETWEEN 20” AND 26” (510 - 660 MM)

| | | |
|---------------------|--|---------------------|
| STS20-25-22 | Width adjustment from 20” to 25” (510-635 MM) and Depth range of 18” to 22” (455 - 560 MM) | N/C |
| STS23-28-22 | Width adjustment from 23” to 28” (585-710 MM) and Depth range of 18” to 22” (455 - 560 MM) | \$475 |
| | ▲ <i>Seat widths of 27” and 28” are custom dimensions and require an (RFQ)</i> | |
| | ▲ <i>The minimum seat depth available with this option is 18” (455 MM)</i> | |
| 20” (SOSW20) | 21” (SOSW21) | 22” (SOSW22) |
| 23” (SOSW23) | 24” (SOSW24) | 25” (SOSW25) |
| 26” (SOSW26) | | |



40. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22" (380 - 560 MM)

| | | | | | | |
|--------------|--------------|--------------|---|--------------|--------------|-----|
| 15" (SODN15) | 16" (SODN16) | 17" (SODN17) | 18" (SODN18) | 19" (SODN19) | 20" (SODN20) | N/C |
| 21" (SODN21) | 22" (SODN22) | | ▲ 18" is the minimum possible seat depth when choosing a 22" long Seat Pan (STS14-19-22, STS14-22-22 and STS20-25-22) | | | N/C |

41. SIGNATURE TAILORED SOLUTIONS - EXTRUDED BACK CANE HEIGHT (THE MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE BACK CANE)

► **Not compatible with a power or manual recline**

| | | | | | |
|--|---------------|---------------|---------------|---------------|---------------|
| 16" (SOBCH16) | 17" (SOBCH17) | 18" (SOBCH18) | 19" (SOBCH19) | | |
| ▼ 20" is the minimum back height required with a power or manual recline | | | | | |
| 20" (SOBCH20) | 21" (SOBCH21) | 22" (SOBCH22) | 23" (SOBCH23) | 24" (SOBCH24) | 25" (SOBCH25) |
| 26" (SOBCH26) | 27" (SOBCH27) | 28" (SOBCH28) | | | |

41.5 SIGNATURE TAILORED SOLUTIONS - BACK CANE OPTIONS

▼ **PUSH HANDLES**

| | | |
|------|---------------------|-----|
| SOPH | Push Handles (pair) | N/C |
|------|---------------------|-----|

▼ **ENDOMORPH BACK CANES FOR NON-RECLINING BACKS**

| | | |
|------------|---|-------|
| STSEND0-F1 | Endomorph Back Canes for non-reclining backs, 1 IN (25 MM) total | \$120 |
| STSEND0-F2 | Endomorph Back Canes for non-reclining backs, 2 IN (50 MM) total | \$120 |
| STSEND0-F3 | Endomorph Back Canes for non-reclining backs, 3 IN (75 MM) total | \$120 |
| STSEND0-F4 | Endomorph Back Canes for non-reclining backs, 4 IN (100 MM) total | \$120 |

▼ **ENDOMORPH BACK CANES FOR RECLINING BACKS**

| | | |
|------------|--|-----|
| STSEND0-R1 | Endomorph Back Canes for reclining backs, 1 IN (25 MM) total | N/C |
| STSEND0-R2 | Endomorph Back Canes for reclining backs, 2 IN (50 MM) total | N/C |

▼ **FOLDING BACKREST (FOR NON-RECLINING BACKS)**

| | | |
|-------|--|-------|
| STSFb | Folding Backrest for Storage/Travel Access for non-reclining backs | \$240 |
|-------|--|-------|

42. BACK SUPPORT

| | | |
|-----------|---|-----|
| DSTR-WWHH | Adjustable Tension Back Upholstery, WW x HH | N/C |
|-----------|---|-----|

▲ **Only available with Dual Post armrests (SOADP)**

| | | |
|------------|--|-------|
| STDSB-WWHH | Contour back support up to 21" wide, with a maximum shell height of 22". | \$605 |
|------------|--|-------|

▲ **The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back**

| | | |
|------------|---|-------|
| STDSB-WWHH | Contour back support up to 21" wide for non-standard sizes with shell height ≥ 23". | \$755 |
|------------|---|-------|

▲ **The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back**

| | | |
|------------|--|---------|
| STDSB-WWHH | Contour back support for Heavy-Duty from 22" wide and higher, with any height. | \$1,075 |
|------------|--|---------|

▲ **The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back**

► **MATRX BACKS**

| | | |
|-------------|--|---------|
| PBE-WWHH-HD | MATRX HEAVY DUTY ELITE BACK – (21" to 24" Wide - Back height of 16" and 20") | \$1,150 |
|-------------|--|---------|

▲ **Back shell height is available in 16" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back**

▲ **Weight capacity up to 500 lb**

| | | |
|-------------|--|---------|
| EDB-WWHH-HD | MATRX ELITE DEEP HD BACK – (21" to 24" Wide - Back height of 16", 18" and 20") | \$1,195 |
|-------------|--|---------|

▲ **Back shell height is available in 16", 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back**

▲ **This product is not compatible with Lateral Trunk Supports**

▲ **Weight capacity up to 500 lb**



42. BACK SUPPORT (CONTINUED)

| | |
|--------------------|---|
| ▶ NXT BACKS | |
| NOACB-WWHH | NXT OPTIMA FH AC CONTOUR BACK – (15” to 21” Wide) \$690 |
| | ▲ <i>Weight capacity up to 300 lb</i> |
| NOACB-WWHH | NXT OPTIMA FH AC CONTOUR BACK – (22” to 28” Wide) \$720 |
| | ▲ <i>Weight capacity up to 550 lb</i> |
| | ▲ <i>Back shell heights only come in 13”, 16” or 19” high. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i> |
| NOPCB-WWHH | NXT OPTIMA FH PC DEEP CONTOUR BACK – (15” to 21” Wide) \$790 |
| | ▲ <i>Weight capacity up to 300 lb</i> |
| NOPCB-WWHH | NXT OPTIMA FH PC DEEP CONTOUR BACK – (22” to 28” Wide) \$820 |
| | ▲ <i>Weight capacity up to 550 lb</i> |
| | ▲ <i>This product is not compatible with Lateral Trunk Supports</i> |
| | ▲ <i>Back shell heights only come in 13”, 16” or 19” high. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i> |

| | |
|---------------------------------------|---|
| ▶ OMISSION OF THE BACK SUPPORT | |
| STDSBO | Omission of the Back Support - (Using an Aftermarket Back) N/C |
| | ▲ <i>Back support MUST be omitted before selecting aftermarket back hardware</i> |
| | ▲ <i>Upcharges apply when both a back support AND aftermarket back hardware are selected</i> |
| CUST-BHW | Aftermarket Back Hardware - Information on the back must be specified: \$185 |
| | <i>Back Manufacturer:</i> _____ |
| | <i>Product Name:</i> _____ |
| | <i>Hardware being used:</i> _____ |
| | ▲ <i>Required information</i> |
| | ▲ <i>Upcharges apply when both a back support AND aftermarket back hardware are selected</i> |
| | ▲ <i>Hardware will not be provided IF the required information for the aftermarket back support is omitted, but a standard T-Nut kit will be supplied</i> |

43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

| | | | | | | |
|---|---|---------------------|---------------------|---------------------|---------------------|---------------------|
| | ▼ <i>20” is the minimum back height required with a reclining back.</i> | | | | | |
| 16” (SOBH16) | 18” (SOBH18) | 20” (SOBH20) | 22” (SOBH22) | 24” (SOBH24) | 26” (SOBH26) | 28” (SOBH28) |
| <i>This height is not available with the NXT Optima AC and PC options ▲</i> | | | | | | |

44. SEAT CUSHIONS

| | |
|---|--|
| ▶ FOR CUSTOM SIZES AND CONFIGURATIONS, PLEASE REFER TO THE PERFORMA CG AIR ORDER FORM. | |
| ▶ PERFORMA CG AIR - THERAPEUTIC AIR CUSHION | |
| Custom Cushion (Please use the Performa CG Air order form) | |
| Cell Height (Choose 1): | 2” 3” 4” |
| ▶ SINGLE (1) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING | |
| ▶ PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER | |
| PCGA1VSCxxxx-x | 14” to 20” (355 mm to 510 mm) wide \$812 |
| PCGA1VSCxxxx-x | 22” + (560 mm +) wide \$914 |
| | ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> |



44. SEAT CUSHIONS

| | | |
|--|---|--------------|
| ▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER | | |
| PCGA1V3Dxxxx-x | 14" to 20" (355 mm to 510 mm) wide | \$857 |
| PCGA1V3Dxxxx-x | 22" + (560 mm +) wide | \$963 |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| ▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER | | |
| PCGA1VICxxxx-x | 14" to 20" (355 mm to 510 mm) wide | \$915 |
| PCGA1VICxxxx-x | 22" + (560 mm +) wide | \$969 |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| ▶ DUAL (2) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING | | |
| ▶ PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER | | |
| PCGA2VSCxxxx-x | 14" to 20" (355 mm to 510 mm) wide | \$812 |
| PCGA2VSCxxxx-x | 22" + (560 mm +) wide | \$914 |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| ▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER | | |
| PCGA2V3Dxxxx-x | 14" to 20" (355 mm to 510 mm) wide | \$857 |
| PCGA2V3Dxxxx-x | 22" + (560 mm +) wide | \$963 |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| ▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER | | |
| PCGA2VICxxxx-x | 14" to 20" (355 mm to 510 mm) wide | \$915 |
| PCGA2VICxxxx-x | 22" + (560 mm +) wide | \$969 |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| ▶ DUAL VALVE CUSHION SEGMENTATION (PLEASE SELECT 1 BOX) | | |
| PCGA-2VLR | Dual Valves with Left / Right Location | N/C |
| PCGA-2VF6 | Dual Valves with Front / Back Location, Front distance 6 in (150 mm) from the front | N/C |
| PCGA-2VF8 | Dual Valves with Front / Back Location, Front distance 8 in (205 mm) from the front | N/C |
| PCGA-2VF10 | Dual Valves with Front / Back Location, Front distance 10 in (255 mm) from the front | N/C |
| ▶ PERFORMA CG AIR OPTIONS | | |
| PCGA-GZ | GripZone Anti-Slip Technology | \$86 |
| PCGA-TT-WW | Transfer Trellis | \$100 |
| PCGA-AC | Anatomic Contour Seat Insert | \$100 |
| ▲ <i>Only available with 2" and 3" cell heights</i> | | |
| PCGA-SP | Seat Pouch | \$55 |
| ▶ CONTOUR SEAT CUSHION | | |
| SCPSC-WWDD | Contour Seat Cushion up to 19" Wide | \$620 |
| SCPSC-WWDD | Contour Seat Cushion for Non-Standard Sizes of 20" to 24" Wide | \$735 |
| SCPSC-WWDD | Contour Seat Cushion for Heavy-Duty of ≥ 25" Wide ▲ Specify: Width: _____" Depth: _____" | \$865 |
| ▶ MATRX CUSHION | | |
| PSP-WWDD | MATRX PSP Seat Cushion – (11" to 20" Wide) | \$645 |
| ▲ <i>The maximum user weight capacity is 300 lb</i> | | |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |
| PSP-WWDD | MATRX PSP Seat Cushion – (22" to 26" Wide) | \$995 |
| ▲ <i>The maximum user weight capacity is 550 lb</i> | | |
| ▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i> | | |



44. SEAT CUSHIONS (CONTINUED)

▶ **NXT BIO-FIT CUSHION**

| | | |
|-----------------------|--|-------|
| NXTBIOFIT-WWDD | NXT BIO FIT Cushion – (15” to 21” Wide) | \$595 |
| | ▲ The maximum user weight capacity is 300 lb | |
| NXTBIOFIT-WWDD | NXT BIO FIT Cushion – (22” to 26” Wide) | \$625 |
| | ▲ The maximum user weight capacity is 400 lb | |
| | ▲ Cushion sizes only come in even numbers for width and depth measurements. | |

45. SIGNATURE TAILORED SOLUTIONS - CANTILEVER ARMRESTS - ADJUSTABLE HEIGHT, INWARD/ OUTWARD ROTATION AND VERTICAL ANGLE. ALLOWS ARMPAD TO MOVE 1” FORWARD (CHOOSE 2)

FOR NON-RECLINING BACKS

▼ **Height Adjustable from 7.25” (185 mm) to 15.25” (385 mm) in 1” (25 mm) increments**

▼ **Not compatible with Waterfall Armpads (SOWAPx et SOWGELx)**

| | | |
|------------------------|--|------|
| SOAFBDLH | Signature Cantilever Armrest – Desk Length 12” (305 mm) - Left Side | N/C |
| SOAFBDRH | Signature Cantilever Armrest – Desk Length 12” (305 mm) - Right Side | N/C |
| SOAFBFLH | Signature Cantilever Armrest – Full Length 16” (405 mm) - Left Side | N/C |
| SOAFBFRH | Signature Cantilever Armrest – Full Length 16” (405 mm) - Right Side | N/C |
| SOAFBELH | Signature Cantilever Armrest – Extended Length 18” (455 mm) - Left Side | N/C |
| SOAFBERH | Signature Cantilever Armrest – Extended Length 18” (455 mm) - Right Side | N/C |
| ▼ CUSTOM LENGTH | | |
| SOAFBCLH | Signature Cantilever Armrest – Custom Length 20” (510 mm) - Left Side | \$95 |
| SOAFBCRH | Signature Cantilever Armrest – Custom Length 20” (510 mm) - Right Side | \$95 |
| SOAFBCLH | Signature Cantilever Armrest – Custom Length 22” (560 mm) - Left Side | \$95 |
| SOAFBCRH | Signature Cantilever Armrest – Custom Length 22” (560 mm) - Right Side | \$95 |

FOR RECLINING BACKS

▼ **Height Adjustable from 10.25” (260 mm) to 16.25” (415 mm) in 1” (25 mm) increments**

▼ **Choose a Height Range :** **Low Range 10.25” to 13.25” (260 - 335 mm)** **High Range 14.25” to 16.25” (360 - 415 mm)**

▼ **Not compatible with Waterfall Armpads (SOWAPx and SOWGELx)**

| | | |
|------------------------|---|------|
| SOAFBDLH-R | Signature Cantilever Armrest for Reclining Back – Desk Length 12” (305 mm) - Left Side | N/C |
| SOAFBDRH-R | Signature Cantilever Armrest for Reclining Back – Desk Length 12” (305 mm) - Right Side | N/C |
| SOAFBFLH-R | Signature Cantilever Armrest for Reclining Back – Full Length 16” (405 mm) - Left Side | N/C |
| SOAFBFRH-R | Signature Cantilever Armrest for Reclining Back – Full Length 16” (405 mm) - Right Side | N/C |
| SOAFBELH-R | Signature Cantilever Armrest for Reclining Back – Extended Length 18” (455 mm) - Left Side | N/C |
| SOAFBERH-R | Signature Cantilever Armrest for Reclining Back – Extended Length 18” (455 mm) - Right Side | N/C |
| ▼ CUSTOM LENGTH | | |
| SOAFBCLH-R | Signature Cantilever Armrest for Reclining Back – Custom Length 20” (510 mm) - Left Side | \$95 |
| SOAFBCRH-R | Signature Cantilever Armrest for Reclining Back – Custom Length 20” (510 mm) - Right Side | \$95 |
| SOAFBCLH-R | Signature Cantilever Armrest for Reclining Back – Custom Length 22” (560 mm) - Left Side | \$95 |
| SOAFBCRH-R | Signature Cantilever Armrest for Reclining Back – Custom Length 22” (560 mm) - Right Side | \$95 |



46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)

▼ **NOT compatible with the extra-long armpads (16" x 3 1/2") (SOXWFAP16 or SOXWGAP16) and the Arm Trough with Hand Support (SOCAR)**

| | | |
|-----------------|---|-------|
| SOADPDLH | Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments) | \$168 |
| SOADPDRH | Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments) | \$168 |
| SOADPFLH | Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments) | \$168 |
| SOADPFRH | Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments) | \$168 |

48. ARMREST HEIGHT ADJUSTMENT AT FACTORY

| | | | | | | | | | | |
|---|-------|-------|--------|--------|---------|--------|--------|--------|--------|--|
| Cantilever Armrest Range (SOAFB for Non-Recl. backs) | | | | | Set at: | | | | | |
| 7.25" | 8.25" | 9.25" | 10.25" | 11.25" | 12.25" | 13.25" | 14.25" | 15.25" | 16.25" | |
| Cantilever Armrest Range (SOAFB for Reclining backs) | | | | | Set at: | | | | | |
| | | | 10.25" | 11.25" | 12.25" | 13.25" | 14.25" | 15.25" | 16.25" | |

- ▲ The Armrest Height is taken from the seat pan to the top of the armpad.
- ▲ Not required for Dual Post.

49. SIGNATURE TAILORED SOLUTIONS - ARMPADS

FOAM TYPE SELECTION

| | | | | | |
|---|---|--------------|----------------|-----------------|-----|
| SOFAPD | Desk Length Flat Armpad – (10 1/2" x 2 1/2") | Select side: | Left Hand Side | Right Hand Side | STD |
| SOFAPF | Full Length Flat Armpad – (14" x 2 1/2") | Select side: | Left Hand Side | Right Hand Side | STD |
| SOXWFAP10 | X-Wide 3.5" (90 mm) Flat FOAM Armpad, 10" (255 mm) Long | Select side: | Left Hand Side | Right Hand Side | N/C |
| SOXWFAP14 | X-Wide 3.5" (90 mm) Flat FOAM Armpad, 14" (355 mm) Long | Select side: | Left Hand Side | Right Hand Side | N/C |
| SOXWFAP16 | X-Wide 3.5" (90 mm) Flat FOAM Armpad, 16" (405 mm) Long | Select side: | Left Hand Side | Right Hand Side | N/C |
| ▲ The X-Wide 16" x 3 1/2" Armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more | | | | | |
| ▲ The X-Wide 16" x 3 1/2" Armpads are NOT compatible with Dual Post Armrests (SOADPxx) | | | | | |
| SOWAPD | Desk Length Waterfall – (8 1/2" x 2 3/4") | Select side: | Left Hand Side | Right Hand Side | N/C |
| ▲ May interfere with flipping back the armrests when combined with Lateral Trunk supports. | | | | | |
| ▲ NOT compatible with Cantilever (SOAFBxxx) | | | | | |
| SOWAPF | Full Length Waterfall – (13" x 2 3/4") | Select side: | Left Hand Side | Right Hand Side | N/C |
| ▲ May interfere with flipping back the armrests when combined with Lateral Trunk supports. | | | | | |
| ▲ NOT compatible with Cantilever (SOAFBxxx) | | | | | |

HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)

| | | | | | |
|---|--|--------------|----------------|-----------------|-------|
| SOGELD | Desk Length Hybrid Gel and Foam Flat armpad – (9 3/4" x 2") | Select side: | Left Hand Side | Right Hand Side | \$40 |
| SOGELF | Full Length Hybrid Gel and Foam Flat armpad – (13 1/2" x 2") | Select side: | Left Hand Side | Right Hand Side | \$40 |
| SOXWGAP10 | X-Wide 3.5" (90 mm) Flat GEL Armpad, 10" (255 mm) Long | Select side: | Left Hand Side | Right Hand Side | \$95 |
| SOXWGAP14 | X-Wide 3.5" (90 mm) Flat GEL Armpad, 14" (355 mm) Long | Select side: | Left Hand Side | Right Hand Side | \$95 |
| SOXWGAP16 | X-Wide 3.5" (90 mm) Flat GEL Armpad, 16" (405 mm) Long | Select side: | Left Hand Side | Right Hand Side | \$105 |
| ▲ The X-Wide 16" x 3 1/2" Armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more | | | | | |
| ▲ The X-Wide 16" x 3 1/2" Armpads are NOT compatible with Dual Post Armrests (SOADPxx) | | | | | |



49. SIGNATURE TAILORED SOLUTIONS - ARMPADS (CONTINUED)

| | | | | | |
|---|---|--------------|----------------|-----------------|-------|
| SOWGELD | Desk Length Hybrid Gel & Foam Waterfall armpad – (8 ½" x 2 ¾") | Select side: | Left Hand Side | Right Hand Side | \$95 |
| | ▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i> | | | | |
| | ▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i> | | | | |
| SOWGELF | Full Length Hybrid Gel and Foam Waterfall armpad – (13" x 2 ¾") | Select side: | Left Hand Side | Right Hand Side | \$95 |
| | ▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i> | | | | |
| | ▲ <i>NOT compatible with Cantilever (SOAFBxxx)</i> | | | | |
| ▶ SPECIALTY TYPE SELECTION | | | | | |
| SOCAR | Arm Trough with Hand Support (18" x 4 ½") | Select side: | Left Hand Side | Right Hand Side | \$290 |
| | ▲ <i>Arm Trough Armpads may interfere with Armrests Flipping back.</i> | | | | |
| ▲ <i>SOCAR armpads require Cantilever (SOAFBExx) armrests with an extended length of 18" (455 mm) or more</i> | | | | | |

50. SIGNATURE TAILORED SOLUTIONS - FOOTRESTS

| | | | | | |
|----------------|---|--------------------------|--|--------------------------|---------|
| SOAAF | Signature Heavy Duty Swing-Away Footrests | | | | \$175 |
| | Set at: | 60° | | 70° | |
| | Select Footrest Range: | SOFLHD12 (12 ½" – 17 ½") | | SOFLHD14 (14 ½" – 19 ½") | |
| ADDML | Signature Manual Elevating & Articulating Legrests - Pair | | | | \$595 |
| | Select Footrest Range: | SOFLME12 (12 ½" – 17 ½") | | SOFLME14 (14 ½" – 19 ½") | |
| | ▲ <i>Does not include calf pads. Please select calf pads in section 51.</i> | | | | |
| ADDPLHD | Signature Heavy Duty Adult Power Elevating & Articulating Legrests | | | | \$2,745 |
| | Select Function: | Individual | | Combined | |
| | Select Footrest Range: | SOFLPW12 (12 ½" – 17 ½") | | SOFLPW14 (14 ½" – 19 ½") | |
| | ▲ <i>This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.</i> | | | | |

51. SIGNATURE TAILORED SOLUTIONS - FOOTPLATES & ACCESSORIES

| | | | | | |
|----------------|---|---------|----------|-----------|-------|
| SOSFAF | Flip-up Footplates | | | | STD |
| SOAFP | Angle Adjustable & Flip-up Footplates – (Each) | | | | \$65 |
| SOAHL | Heel Loops – (Each) | | | | \$20 |
| SOACPFR | Calf Pads for Adult Heavy-Duty Swing-Away footrests – (SOAAF) – with adjustable hardware – (Each) | | | | \$93 |
| SOAMA | Swing-away residual limb Support | Select: | Left leg | Right leg | \$455 |

52. FOOTREST LENGTH ADJUSTMENT AT FACTORY

| | | | | | | | | | |
|----------------------------------|---------|-------|-----|-------|-----|-------|-----|-------|-----|
| Footrest / Legrest Length | Set at: | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" |
| | | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | |

▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

53. SIGNATURE TAILORED SOLUTIONS - CENTER MOUNTS

| | | | | | |
|----------------|---|-------------------------|--|-------------------------|-------|
| SOAAF90 | Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°) | | | | \$200 |
| | Select Foot Platform Range: | SOFBA12 (12 ½" – 16 ½") | | SOFBA16 (16 ½" – 19 ½") | |
| | ▲ <i>When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 18.75"</i> | | | | |
| | ▲ <i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i> | | | | |
| | ▲ <i>With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.</i> | | | | |



53. SIGNATURE TAILORED SOLUTIONS - CENTER MOUNTS (CONTINUED)

| | | |
|----------------|---|---------|
| ADDPFHD | 90° Heavy Duty Power Center Mount Elevating & Articulating Foot Platform | \$2,745 |
| | Select Footrest Range: SOFBPW12 (10 ½"- 16") SOFBPW16 (16"- 20") | |
| ▲ | <i>When selecting SOFBPW16 (16"- 20") range, the lowest possible seat to floor height is 18.75"</i> | |
| ▲ | <i>This option includes calf pads. The power center mount counts as one power function.</i> | |
| ▲ | <i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i> | |
| ▲ | <i>For leg lengths between 10 ½" and 12 ½", there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.</i> | |

54. SIGNATURE TAILORED SOLUTIONS - FOOT PLATFORMS & ACCESSORIES

| | | |
|----------------|--|-------|
| SOFUFP | Flip-up one piece Foot Platform | \$140 |
| SOAFPFB | Angle Adjustable & Flip-up Split Footplates for center mount – (Each) | \$65 |
| SOACPCM | Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each) | \$93 |

55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY

| | | | | | | | | | | | |
|--|---------|-------|-----|-------|-----|-------|-----|-------|-----|--|--|
| Fixed Center Mount Foot Platform Length | Set at: | | | | | | | | | | |
| 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | | |
| 17.5" | 18" | 18.5" | 19" | 19.5" | | | | | | | |
| Power Center Mount Foot Platform Length | Set at: | | | | | | | | | | |
| 10.5" | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | | |
| 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | | |

▲ The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates

56. SEAT POSITIONING

| | | |
|---------------------|--|-------|
| SOAUTO | 55" Seat Positioning Strap | STD |
| SOAUTO80 | 80" Seat Positioning Strap | \$75 |
| SOBPHP2C | Bodypoint padded Hip Belt – (Two Point – Center Pull) | \$170 |
| SOBPHP2C-L82 | XL 82" Bodypoint padded Hip Belt – (Two Point – Center Pull) | \$190 |
| SOPHB | Padded Hip Belt, 2-PT Magnetic Center Pull | \$194 |
| SOBPCS | Bodypoint padded Chest Strap | \$250 |
| SOPCS | Padded Chest Strap, 2-PT Magnetic Center Pull | \$262 |

57. HEADREST

▼ **SELECT ONE OF THE HEADREST OPTIONS WITH THE MOUNTING HARDWARE AND SPECIFY WHERE THE HEADREST WILL BE INSTALLED.**

| | | |
|-------------------|---|-------|
| SONUAHR-12 | Using headrest mounting hardware other than SOMAH below. Available with both Fixed and Reclining backrests. 1/2 IN (12.7 MM) Mounts | N/C |
| SONUAHR-15 | Using headrest mounting hardware other than SOMAH below. Available with Reclining backrests only. 19/32 IN (15 MM) Mounts | N/C |
| SOFHRS | Small Curved Headrest Pad – 4" x 9" | \$155 |
| SOFHRM | Medium Curved Headrest Pad – 5" x 10" | \$155 |
| SOFHRL | Large Curved Headrest Pad – 6" x 11" | \$155 |
| SOMAH | Multi-axis Removable Mounting Hardware | \$220 |
| SOHRBC | Bracket to Mount on the Contour Back Support – (Without Recline) | N/C |
| SOHRPR | Bracket to Mount on the Seat Frame – (For Reclining Backs) | N/C |



58. LATERAL TRUNK SUPPORTS

▼ **SELECT ONE OF THE LATERAL PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

| | | |
|---------------|---|-------|
| SOSLS3 | Adjustable Laterals with Small Curved Pads – (Each) – 3” x 5” | \$100 |
| SOSLM4 | Adjustable Laterals with Medium Curved Pads – (Each) – 4” x 6” | \$100 |
| SOSLL5 | Adjustable Laterals with Large Curved Pads – (Each) – 5” x 7” | \$100 |
| SOSLHM | Adjustable Swing-Away Hardware Mounts – (Each) | \$220 |
| ▲ | <i>The Lateral Trunk Supports are not compatible with the NXT Optima FH PC 6” Deep Contour Fixed Thoracic Back Support (NOPCB).</i> | |
| ▲ | <i>Waterfall armpads (SOWAPx & SOWGELx) may interfere with flipping back the armrests when combined with Lateral Trunk supports</i> | |

59. HIP & THIGH SUPPORTS

▼ **SELECT ONE OF THE HIP & THIGH SUPPORT PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

| | | |
|-----------------|--|-------|
| SOHTSS10 | Adjustable Hip & Thigh Support Pads – (Each) – Small 10” x 4” | \$130 |
| SOHTSM12 | Adjustable Hip & Thigh Support Pads – (Each) – Medium 12” x 5” | \$130 |
| SOHTSL14 | Adjustable Hip & Thigh Support Pads – (Each) – Large 14” x 6” | \$130 |
| SOHTSARH | Adjustable and Removable Hip & Thigh Support Hardware – (Each) | \$136 |

60. ACCESSORIES

| | | | | | |
|--------------|---|-------------|-----------|------------|-------|
| SOCH | Cup Holder – (Self-Leveling) | Install on: | Left Side | Right Side | \$115 |
| SOCPH | Ram Mounts X-Grip Cell Phone Holder | Install on: | Left Side | Right Side | \$190 |
| SOTHS | Ram Mounts X-Grip Small Tablet Holder | Install on: | Left Side | Right Side | \$205 |
| SOTHL | Ram Mounts X-Grip Large Tablet Holder with Swivel | Install on: | Left Side | Right Side | \$360 |
| ▲ | <i>Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.</i> | | | | |
| ▲ | <i>With multiple accessories installed on the same side, a custom upgrade will be required.</i> | | | | |

61. REFLECTOR KIT

| | | |
|----------------|---|-----|
| REFLECT | Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner’s Manual) | N/C |
|----------------|---|-----|

62. BACK PACK AND ARMREST SIDE POUCHES

| | | | | | |
|-------------|--|--------|-----------|------------|------|
| SOBP | Back Pack – (Comes with durable straps that can be hooked around the back canes) | | | | \$90 |
| SOAP | Armrest Side Pouch – (Comes with durable straps that can be hooked around the armrest tubes) | Order: | Left Side | Right Side | \$53 |
| ▲ | <i>If only choosing one (1) armrest pouch, please show what side it should be installed</i> | | | | |

63. SWITCHES AND OPTIONS (CONTINUED)

| | | |
|-----------------|--|---------|
| ASCPWIFI | PGDT - Wireless programmer, WIFI dongle | \$844 |
| ASCPCP | PGDT R-NET Dongle and PC programmer kit | \$600 |
| ASCPVR2 | PGDT VR2 with PPI Dongle and PC programmer kit | \$1,450 |
| SSBS | Buddy Button Switch | \$230 |
| SSES | Egg Switch-Black | \$116 |
| SSML | Ultra Light Switch | \$295 |
| ESS-RED | Emergency STOP Switch - Red | \$314 |



SPECIALTY DRIVER CONTROL OPTIONS

▶ Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)
 ▶ Specialty Control options in this section are not compliant with ISO 7176-19

64. THERAFIN SIP-N-PUFF KIT

▼ *This option is not compliant with ISO 7176-19*

| | | | | | |
|--------------------------------------|--|----------------------|-----------------------|--|---------|
| SCSNPK | Therafin Sip-N-Puff Kit | | | | \$1,975 |
| SCSNPBTK | Therafin Sip-N-Puff Breath Tube Kit | | | | N/C |
| | ▲ <i>Required with the Therafin Sip-N-Puff kit.</i> | | | | |
| SCSRMSN | Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff | SWL (Swings on Left) | SWR (Swings on Right) | | \$300 |
| | ▲ <i>Required with the Therafin Sip-N-Puff kit.</i> | | | | |
| SCSRM | Stealth Rock-n-Lock Swingaway Mount for Mode Switch | SWL (Swings on Left) | SWR (Swings on Right) | | \$510 |
| SSES | Egg Switch for Changing Modes – Black | MLH (Mount on Left) | MRH (Mount on Right) | | \$116 |
| | ▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i> | | | | |
| | ▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section</i> | | | | |
| ESS-RED | Emergency STOP Switch - Red | User Activation | Attendant Activation | | \$314 |
| | ▲ <i>Required with the Therafin Sip-N-Puff kit.</i> | | | | |
| SOMAH | Multi-axis Headrest Removable Mounting Hardware | | | | \$220 |
| | ▲ <i>Required with the Therafin Sip-N-Puff kit.</i> | | | | |
| ▶ CHOOSE ONE (1) HEADREST PAD | | | | | |
| SOFHRS | Small Curved Headrest Pad – 4" x 9" | | | | \$155 |
| SOFHRM | Medium Curved Headrest Pad – 5" x 10" | | | | \$155 |
| SOFHRL | Large Curved Headrest Pad – 6" x 11" | | | | \$155 |

65. PROXIMITY HEAD ARRAY

▼ *This option is not compliant with ISO 7176-19*

| | | | | | |
|----------------|---|----------------------|-----------------------|--|-------|
| CUSTOM | Adult – 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse | | | | RFQ |
| SCSRM | Stealth Rock-n-Lock Swingaway Mount for Mode Switch | SWL (Swings on Left) | SWR (Swings on Right) | | \$510 |
| SSES | Egg Switch for Changing Modes – Black | MLH (Mount on Left) | MRH (Mount on Right) | | \$116 |
| | ▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i> | | | | |
| | ▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.</i> | | | | |
| ESS-RED | Emergency STOP Switch - Red | User Activation | Attendant Activation | | \$314 |
| | ▲ <i>Required with the Proximity Head Array.</i> | | | | |
| SOMAH | Multi-axis Headrest Removable Mounting Hardware | | | | \$220 |
| | ▲ <i>Required with the Proximity Head Array.</i> | | | | |



66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|--|---|----------------------|-----------------------|-------|
| CUSTOM | Light touch mini Proportional Joystick | | | RFQ |
| SCSRMMP | Stealth Rock-n-Lock Swingaway Arm for Chin MP | SWL (Swings on Left) | SWR (Swings on Right) | \$870 |
| ▲ <i>Required with the Light touch mini Proportional Joystick.</i> | | | | |
| SCSRM | Stealth Rock-n-Lock Swingaway Mount for Mode Switch | SWL (Swings on Left) | SWR (Swings on Right) | \$510 |
| SSES | Egg Switch for Changing Modes – Black | MLH (Mount on Left) | MRH (Mount on Right) | \$116 |
| ▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i> | | | | |
| ▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section</i> | | | | |
| JCC | Chin Cup | | | \$55 |
| SOMAH | Multi-axis Headrest Removable Mounting Hardware | | | \$220 |
| ▲ <i>Required with the Light touch mini Proportional Joystick.</i> | | | | |
| ▶ CHOOSE ONE (1) HEADREST PAD | | | | |
| SOFHRS | Small Curved Headrest Pad – 4" x 9" | | | \$155 |
| SOFHRM | Medium Curved Headrest Pad – 5" x 10" | | | \$155 |
| SOFHRL | Large Curved Headrest Pad – 6" x 11" | | | \$155 |

67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|--|---|---------------------|----------------------|---------|
| SCPGMJS | PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount | | | \$2,369 |
| SCGMJSM | Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for JC200 | MLH (Mount on Left) | MRH (Mount on Right) | \$595 |
| ▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i> | | | | |
| CUSTOM | Light touch mini Proportional Joystick | | | RFQ |
| SCGMPPM | Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for Light Touch Joystick | MLH (Mount on Left) | MRH (Mount on Right) | \$825 |
| ▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i> | | | | |
| SCBTM | Bullet Hand Tray for Light Touch mini Proportional Joystick | MLH (Mount on Left) | MRH (Mount on Right) | \$850 |
| ▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i> | | | | |

68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|--|--|---------------------|----------------------|---------|
| SCPGMJS-M | PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount | | | \$2,369 |
| SCARMJSM | Swing-away and Multi-Axis armrest mount | MLH (Mount on Left) | MRH (Mount on Right) | \$325 |
| ▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i> | | | | |
| CUSTOM | Light touch mini Proportional Joystick | | | RFQ |
| SCARMPM | Swing-away and Multi-Axis armrest mount | MLH (Mount on Left) | MRH (Mount on Right) | \$375 |
| ▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i> | | | | |
| SCBTM | Bullet Hand Tray for Light Touch mini Proportional Joystick | MLH (Mount on Left) | MRH (Mount on Right) | \$850 |
| ▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i> | | | | |



69. PROPORTIONAL HEAD CONTROL PACKAGE

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|--------------|--|----------------------|-----------------------|---------|
| SCPHC | Proportional Head Control Package | | | \$4,495 |
| SCSRM | Stealth Rock-n-Lock Swingaway Mount for Mode Switch | SWL (Swings on Left) | SWR (Swings on Right) | \$510 |
| SSES | Egg Switch for Changing Modes – Black | MLH (Mount on Left) | MRH (Mount on Right) | \$116 |
| | ▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i> | | | |
| | ▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section</i> | | | |
| SOMAH | Multi-axis Headrest Removable Mounting Hardware | | | \$220 |
| | ▲ <i>Required with the Proportional RIM Head Control.</i> | | | |

CUSTOM SECTION

- ▶ Lead times will increase depending on the custom order pipeline at the time the purchase order is received.
- ▶ Please call for Quote and specifications for these options.

70. O² HOLDER

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|----------------|---|--|----------------------|---------|
| SUPOXY | Custom Non-Articulating Oxygen Tank Holder | MLH (Mount on Left) | MRH (Mount on Right) | \$295 |
| | ▲ <i>Not available with a power recline</i> | | | |
| SUPOXYA | Custom Articulating Oxygen Tank Holder | MLH (Mount on Left) | MRH (Mount on Right) | \$1,195 |
| | ▲ <i>Required with a power recline</i> | | | |
| | D-Size | Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve | | |
| | E-Size | Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve | | |

71. VENT TRAY (PLEASE COMPLETE IN FULL)

▼ **AMYLOR DOES NOT SUPPLY THE BATTERY, CHARGER OR THE WIRING FOR THE VENT. THESE ITEMS ARE NORMALLY PROVIDED BY THE VENT MANUFACTURER.**

▼ *This option is not compliant with ISO 7176-19*

| | | | | |
|---------------|---|--|--|---------|
| SOVENT | Custom Vent Tray | | | |
| | Articulating Vent Box (Required with Power Tilt and/or Recline) | | | \$1,360 |
| | Fixed Vent Box | | | \$400 |
| | SPECIFY MODEL | | | |
| | Trilogy 100/200 | | | |
| | LTV | | | |
| | LTV with Sprintpack | | | |
| | Other, specify: _____ | | | |

71. VENT TRAY (PLEASE COMPLETE IN FULL) (CONTINUED)

| | | | | |
|---------------|---|--|--|-------|
| SOVTBE | Extended Battery Compartment | | | \$495 |
| | The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries | | | |
| | The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries | | | |



72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)

▼ **PLEASE SPECIFY DETAILS AND SPECIFICATIONS FOR THESE OPTIONS.**

▼ **These following custom options and accessories are ISO 7176-19 compliant**

| | | |
|---------------|--|-----|
| CUSTOM | USB Charger Port ◀ Specify the location for the installation: _____ | RFQ |
| CUSTOM | Tilt + Elevate Function Through VR2 Non-Expandable Drive Control | RFQ |
| CUSTOM | Accessory Hooks located at the rear of the seat structure (2 hooks) | RFQ |
| CUSTOM | Heavy-Duty Wheelchair Cover, Large | RFQ |
| CUSTOM | Transfer Handles – (Pair) Select a height range: 5" to 8" 6" to 9" | RFQ |

▼ **The following custom options and accessories below are not compliant with ISO 7176-19**

| | | |
|------------------|---|-----|
| CUSTOM | Folding Backrest for Storage/Travel Access for Reclining backs | RFQ |
| CUSTOM | Power Anterior Tilt | RFQ |
| CUSTOM | Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity) | RFQ |
| CUSTOM | Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity) | RFQ |
| CUSTOM | Rearview Camera with 5" Monitor | RFQ |
| CUSTOM | Rearview mirror Left Side Right Side | RFQ |
| CUSTOM | Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan | RFQ |
| CUSTOM | Foot Box ◀ Specify Length: Width: Please confirm with CSR on sizing possibilities | RFQ |
| CUSTOM #1 | Other custom options. Specify in notes section below | RFQ |
| CUSTOM #2 | Other custom options. Specify in notes section below | RFQ |

73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS