

Please check off the interested products: ☒ _____

☐ size **Small**

☐ size **Medium**

Frame colour: **black**

Upholstery colour: ☐ **dark grey/green**

☐ **dark grey/pink**

☐ **light grey/green**

☐ **light grey/pink**



CRASH TESTED
according to the ISO 7176/19
and ANSI RESNA WC/19

Product registered **FDA (USA)**

Product marked **CE**

SEAT UNIT, the only one with a tilt-in-space system incorporated in the seat

- HEADREST incorporated in the backrest: your child's head is always supported during the tilt adjustments of the seat and of the backrest, also without accessories.
- Rigid BACKREST, reclinable and adjustable in height. The hinge coincides with the child's pelvis physiological fulcrum.
- Tilt-in-space SEAT, adjustable in width and depth. Reversible (size Small, base 869).
- Fireproof, recyclable FABRIC. Colours durable and highly resistant to sweat, light, dry cleaning and washing. Maximum resistance to rubbing. Oeko-Tex Certificate - for textile products intended for all fabrics that are in direct contact with the skin, a guarantee of quality and safety for health.
- Tip-up FOOTREST, flexion-extension adjustable, with legrest adjustable in height and tilt; the hinge coincides with the child's knee physiological fulcrum.

869 OUTDOOR four wheels base - CRASH TESTED ISO 7176/19 and ANSI RESNA WC/19

- Foldable aluminium FRAME. Push-handle adjustable in tilt: it adapts to the mother's height and reduces the encumbrance of the pushchair.
- Puncture proof WHEELS, light and easy to push on every surface: front 17 cm with directional locks; rear 25 cm (size Small) or 30 cm (size Medium), with shock absorbers and independent drum brakes.

856 - 851 indoor high-low bases, 850 indoor compact floor base

- Foldable FRAME in steel. Adjustable in height. It is easy and safe to use, thanks to a gas spring mechanism, hand lever (base 851) and foot pedal (base 856)
- 10 cm swivelling WHEELS with brakes (for base 856)
- 7.5 cm swivelling WHEELS with brakes (for base 851)



Attention!
Read carefully the user and
maintenance manual



MEASURES



Stroller closed with seat removed



Stroller closed with seat attached



	A	B	C	D	E	F	G	J	H	I	L	M	N	O	P	Q	R	S	max. load	frame weight	seat weight	total weight
Small	20-30 cm	52-62 cm	37 cm	52 cm	20-30 cm	14-30 cm	0°-35°	35°	100-110 cm	57 cm	67-96 cm	35 cm	88 cm	33 cm	82 cm	50 cm	81 cm	90-120 cm	35 kg	9,5 kg	8,5 kg	18 kg
	7,8-11,8 in	20,4-24,4 in	14,5 in	20,4 in	7,8-11,8 in	5,5-11,8 in			39,3-43,3 in	22,4 in	26,3-37,7 in	13,7 in	34,6 in	12,9 in	32,2 in	19,6 in	31,8 in	35,4-47,2 in	77,1 lb	20,9 lb	18,7 lb	39,6 in
Medium	25-36 cm	60-75 cm	46 cm	55 cm	30-40 cm	20-38 cm	5°-30°	40°	112-127 cm	62 cm	80-107 cm	37 cm	112 cm	39 cm	94 cm	52 cm	94 cm	115-140 cm	45 kg	11 kg	12,7 kg	23,7 kg
	9,8-14,1 in	23,6-29,5 in	18,1 in	21,6 in	11,8-15,7 in	7,8-14,9 in			44-50 in	24,4 in	31,4-42,1 in	14,5 in	44 in	15,3 in	37 in	20,4 in	37 in	45,2-55,1 in	99,2 lb	24,2 lb	27,9 lb	52,2 lb

MKT.F108 | P25001 | USA | July 2025 | R01

Hyphen-separated numbers indicate the minimum and maximum range of adjustment

Distributed by:



Dealer name

Date of order

Order ref./name of the User

ADDITIONAL COMPONENTS

☐ 863

Multadjustable Headrest with parietal supports



☐ 852

Headrest with parietal supports



942

Ergonomic headrest for Bug Small

☐ S

☐ M



942

Ergonomic headrest for Bug Medium

☐ M

☐ L



☐ 835

Shaped headrest



☐ 934

Suspension headrest attachment. Only available for 863 headrest



☐ 838

Adjustable lateral supports



☐ 868

Wrappable and flexible trunk supports



☐ 834

Padded abduction block




☐ 834R

Adjustable padded abduction block



☐ 853

Vest harness




☐ 903

Five point vest harness



☐ 906

Five point harness



☐ 894

45° pelvic belt



☐ 920

4 point belt



☐ 827

Foot straps




☐ 839

Adjustable front handle



☐ 824

Adjustable transparent tray



☐ 824B

Wide tray only for 850-851-856 indoor bases.



☐ 905

Hand brake for 869 base



☐ 935

Oxygen bottle holder

Maximum dimensions of the oxygen bottle: height 600 mm (valve included) - diameter 140 mm.



819 Canopy

☐ Black/green ☐ Black/pink



☐ 825

Rain cover (to fix at the canopy 819) for 869 base



☐ 818

Thermal cover



☐ 858

Net basket for 869 base



☐ 911

Vent tray



☐ 892

Padded covering for footrest



☐ 891

Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private/cars, buses...). Only for 869 base



CONFORME ISO 7176/19 e ANSI RESNA WC/19

☐ 922

Bottle holder for 869 base



☐ 851

Indoor high-low compact base



☐ 856

Indoor high-low feeding base



MKT.F108 | P25001 | USA | July 2025 | R01

Amylior - T 1 888 453-0311 - F 1 877 501-8458 www.amylior.com - info@amylior.com