

Distributed by:



It stimulates and helps the step pattern. It trains the user to maintain the balance of the trunk and pelvis. It organizes a wider and structured vision of space, enhances the perception of distances and environment that surrounds the user. It develops physical self-sufficiency and psychological independence by experimenting the possibility of acting independently without the help of other people.

Grillo compensates the user lack of balance thanks to its ergonomic support and wide base, it allows the user to concentrate his/her residual abilities on walking.

The wide range of components and supports offers maximum versatility. All accessories are interchangeable quickly and without the use of tools to suit different clinical pictures.



Product registered FDA (USA)

Product marked CE



- **LIGHT ALUMINUM FRAME**, foldable for an easy transportation with PATENTED HEIGHT ADJUSTMENT SYSTEM.
- **PELVIC SUPPORT and THORACIC SUPPORT**, padded, adjustable in height, width and circumference. They allow an easy access. They hold the user just like two HANDS giving him/her a sense of security and stability allowing the user to concentrate on walking without any fear. The Thoracic support is adjustable in tilt and backward-forward.
- **ERGONOMIC HARNESS**, comfortable, adjustable and easily removable for washing. It supports the user in case of necessity without generating painful pressures.
- **HANDLEBAR**, removable and adjustable in height, it is used as a reference and as an armrest.



Attention!
Read carefully the user and maintenance manual

FRONT DRIVE

Inner circumference
pelvic and thoracic support
from 74 to 120 cm
from 29,1-47,2 in



WEIGHT OF THE AID 22,5 kg / 49,6 lb
MAX LOAD 90 kg / 198,4 lb

REAR DRIVE


Inner circumference
pelvic and thoracic support
from 74 to 120 cm
from 29,1-47,2 in



WEIGHT OF THE AID 21,7 kg / 47,8 lb
MAX LOAD 90 kg / 198,4 lb

USER DIMENSIONS:

USER HEIGHT (C):
110-165 cm / 43,3-64,9 in



IMPORTANT:
Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the supports.



Dealer name
Date of order
Order ref./name of the User

ORDER FORM



Front Drive Rear Drive

CHOOSE THE CONFIGURATION >		PT	PA	P	PT	PA	P
ADDITIONAL COMPONENTS▼		Including:					
925, 924, 930: you can choose other size		> thoracic support > pelvic support > ergonomic harness > handlebar					
865 Multiadjustable headrest		NA	NA	NA	 available in combination with thoracic support 925	 available in combination with thoracic support 925	 available in combination with thoracic support 925
925 Thoracic support		●	 ●	 ●	 ●	 ●	 ●
924 Pelvic support		●	●	●	●	●	●
930 Ergonomic harness		●	●	●	●	●	●
926 Adjustable handlebar		●	 ○	●	 NA	 NA	 NA
809 Arm supports		 ○	●	 ○	 ○ with joints	●	 ○ with joints
943 Arm straps		 ○ available in combination with arm supports 809	 ○	 ○ available in combination with arm supports 809	 ○ available in combination with arm supports 809	 ○	 ○ available in combination with arm supports 809
939 Ergonomic handles		 ○	●	 ○	 ○ with joints	●	 ○ with joints
927 Removable knobs		NA	NA	NA	●	 ○	●
919 Lumbar thrust support		NA	NA	NA	 ○	 ○	 ○
918 Folding seat		NA	NA	NA	 ○	 ○	 ○
815 Proximal abductor		 ○	 ○	 ○	 ○	 ○	 ○
946 Ergonomic saddle		 ○	 ○	 ○	 ○	 ○	 ○
890-C Leg divider with thigh loops		 ○	 ○	 ○	 ○	 ○	 ○
890-SC Distal abductor		 ○	 ○	 ○	 ○	 ○	 ○
810 Weighting bars		 ○	 ○	 ○	 ○	 ○	 ○
923 Additional assistant push handle		 ○	 ○	 ○	 ○	 ○	 ○
928 Rear narrow base wheels with parking brake pedals. ONLY FOR INDOOR USE		 ○	 ○	 ○	 ○	 ○	 ○