

Distributed by:



ORDER FORM

GRILLO **SIZE mini**

FRONT DRIVE and
REAR DRIVE gait trainer

GRILLO
mini



Product registered FDA (USA)

Product marked **CE**




- **ALUMINUM FRAME**, foldable for easy transportation and convertible, from **front configuration to rear configuration and vice versa**.
- **PADDED PELVIC and THORACIC SUPPORTS**, independently adjustable in height, width and circumference; allow an easy access thanks to the back opening. They hold the user just like two hands and this gives him/her a sense of security and stability allowing the user to concentrate only on walking without any fear.
- **GRADUATED BARS** as a reference on height adjustments.
- **ERGONOMIC HARNESS**, ergonomic, comfortable, adjustable and easily removable for washing. It supports the users in case of necessity without generating painful pressures.
- **WHEELS**, light and puncture-proof: 2 front swiveling, Ø15 and 2 rear solid, Ø cm 20, for indoor and outdoor use.
- **DIRECTIONAL LOCKS** on front wheels.
- **ANTI-ROLLBACK** device, **ADJUSTABLE BRAKES**, parking brake on rear wheels, easy to operate by the red pedal.
- **SEMI-ELASTIC BUMPERS** on front wheels to protect against bumps.



Attention!
Read carefully the user and
maintenance manual

FRONT DRIVE

Inner **circumference pelvic and thoracic support**
from 53 to 79 cm
from 20,8 to 31,1 in




WEIGHT kg 15/33 lb

MAX. LOAD kg 35/77,1 lb

REAR DRIVE

Inner **circumference pelvic and thoracic support**
from 53 to 79 cm
from 20,8 to 31,1 in




WEIGHT kg 15/33 lb

MAX. LOAD kg 35/77,1 lb

USER DIMENSIONS:

USER HEIGHT (C):
min 75 cm/29,5 in
max 95 cm/37,4 in



IMPORTANT:
Make sure that measures A and B
indicated in the figure do not exceed the
suggested range indicated for the supports.

28.07.25

Dealer name

Date of order

Order ref./name of the User

ORDER FORM

GRILLO



Front Drive

Rear Drive

CHOOSE THE CONFIGURATION >	PT	PA	P	PTM	PT	PA	P
ADDITIONAL COMPONENTS >	 Required: > thoracic support > pelvic support > ergonomic harness > handlebar	 Required: > pelvic support > ergonomic harness > arm supports	 Required: > pelvic support > ergonomic harness > handlebar	 Required: > transparent tray > handlebar > pelvic support > ergonomic harness	 Required: > thoracic support > pelvic support > ergonomic harness > handles	 Required: > pelvic support > ergonomic harness > arm supports	 Required: > pelvic support > ergonomic harness > handles
865 Multiadjustable headrest	NA	NA	NA	NA	☐	available in combination with thoracic support 925	available in combination with thoracic support 925
925 Thoracic support	●	☐	☐	☐ when transparent tray is removed	●	☐	☐
924 Pelvic support	●	●	●	●	●	●	●
930 Ergonomic harness	●	●	●	●	●	●	●
894 Rear table belt for grillo mini PTM	NA	NA	NA	☐	NA	NA	NA
926 Adjustable handlebar	●	☐	●	NA	NA	NA	NA
809 mini Arm supports	☐	●	☐	NA	☐	●	☐
809 mini Arm supports with curved bars	NA	●	NA	☐ when handlebar and tray are removed	NA	●	NA
943 mini Arm straps for size mini	☐ available in combination with arm supports 809	☐	☐ available in combination with arm supports 809	☐ available in combination with arm supports 809	☐ available in combination with arm supports 809	☐	☐ available in combination with arm supports 809
937 Ergonomic handles for size mini	NA	NA	NA	☐ when handlebar is removed	NA	NA	NA
927 Removable knobs	NA	NA	NA	NA	●	☐	●
946 Ergonomic saddle	☐	☐	☐	☐	☐	☐	☐
890-SC Distal abductor	☐	☐	☐	☐	☐	☐	☐
923 Additional assistant push handle	☐	☐	☐	☐	☐	☐	☐

MKT.F109 | P25001 | USA | January 2026 | R02 ● > REQUIRED IN THE CHOSEN CONFIGURATION

☐ > ACCESSORY

NA > NOT AVAILABLE