

Distributed by:



## ORDER FORM

TILT-IN-SPACE WHEELCHAIR

# JUDITTA B30

The alternative to being bedridden



Attention!  
Read carefully the user  
and maintenance  
manual



- FEATURES**
- **FRAME** made of chromed steel, painted with non-toxic epoxy powder paint. Narrow base to facilitate entering elevators.
  - **The seat and backrest are covered with visco-elastic fire-proof material. It is resistant, breathable and anti-slip. It is made of a special and exclusive ORMESA formula:**
    - 1) It adapts to the user's body shapes and it prevents its sliding, with the result of an **EXCELLENT COMFORT and KEEPING POSTURE**, even when changing the backrest recline and tilt (it adapts to the user's body shapes maintaining its support function).
    - 2) Easily washable to keep the **HYGIENE** of the user; it lets **AIR PASS** keeping the **OPTIMAL MICROCLIMATE** that **HELPS PREVENTING BED SORES**.
    - 3) Durable and fireproof, **GUARANTEE OF DURABILITY and SAFETY**.
    - 4) Ergonomic, breathable, non-toxic, fireproof **UPHOLSTERY**; it does not interfere with the benefits of the visco-elastic cover. The upholstery is fixed with buttons, easy to remove for machine washing.
  - Ergonomic **BACKREST**, it holds the shapes of the user. **The special structure of the backrest gives a proper alignment to the trunk and it doesn't generate any pressure on the spine.** Adjustable in tilt by levers with safety pin. The pin allows blocking the levers to maintain the proper adjustment for the user, height adjustable, folding for transportation.
  - Low **SEAT** and slightly inclined. Adjustable in depth on four positions.
  - Wide and wrap **HEADREST**, padded, multi adjustable (in height, tilt, depth, lateral traslation); thus assuring the **proper alignment between head and trunk**, without interfering with the user's view. Anti-bacterial and Hygienic.
  - **ARMRESTS**, padded, easily and stably adjustable in height thanks to the lever.
  - **PUSH HANDLE** adjustable in tilt/ height, it helps the care.
  - **CONTINUOUS AND WITHOUT EFFORT ADJUSTMENTS**: tilting, backrest recline and legrest tilt are operated by independent gas springs.
  - **SIDES** with scissor opening. Thanks to the special joint, the armrests are positioned behind the backrest, eliminating completely all obstacles to the transfer.
  - **LEG RESTS**: ergonomic, removable, turnable and adjustable in tilt without effort, activated with gas springs.
  - **CALFRETS** adjustable in depth and free tilt to accommodate the movements of the leg.
  - **FOOT RESTS**, adjustable in flexion-extension and height.
  - **8" front solid swiveling wheels and 12" rear wheels**; drum brakes are operated by a single control pedal.



Tie down hooks (4 rings) for the front facing transport of the user in moving vehicles (private cars, conformed buses etc.) **AS STANDARD**

Product registered **FDA (USA)**

Product marked **CE**

Please check off the interested products:



SIZE




UPHOLSTERY COLOUR:



FRAME COLOUR:



MEASURES



\*100 cm / \*39,3 in

\*footrest adjusted to the maximum length

	A seat bottom width	B seat bottom depth	C backrest width	D overall width	H seat bottom inclination-tilt	I backrest inclination	Total weight	Max. Capacity	User height
size 36	36 cm 14,1 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	47 cm 18,5 in	62 cm 24,4 in	from 10° to 35°	from 10° to 38°	38 kg 83,7 lb	120 kg 264,5 lb	from 140 to 190 cm from 55,1 to 74,8 in
size 40	40 cm 15,7 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	47 cm 18,5 in	62 cm 24,4 in	from 10° to 35°	from 10° to 38°	38 kg 83,7 lb	120 kg 264,5 lb	from 140 to 190 cm from 55,1 to 74,8 in
size 45	45 cm 17,7 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	52 cm 20,4 in	67 cm 26,3 in	from 10° to 35°	from 10° to 38°	41 kg 90,3 lb	130 kg 286,6 lb	from 140 to 190 cm from 55,1 to 74,8 in
size 50	50 cm 19,6 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	57 cm 22,4 in	72 cm 28,3 in	from 10° to 35°	from 10° to 38°	42 kg 92,5 lb	140 kg 308,6 lb	from 140 to 190 cm from 55,1 to 74,8 in

MKT.F105 | P25001 | USA | July 2025 | R01

Hyphen separated numbers indicate a minimum and maximum range of adjustment



AMYLIOR

Date of order

Order ref./name of the User

**824-R**  
Reclinable tray,  
adjustable in tilt,  
backwards -  
forwards, with  
security block

It is made by resistant and hygienic material, design with no edges, it doesn't keep residual food.

**igienico - hygienic**  
**resistente e leggero**  
durable and lightweight  
**ergonomico - ergonomic**  
**regolazioni semplici e sicure**  
simple and safe adjustments

**New design!**

**894** 45° Pelvic belt

### 834 Padded abduction block

**892** Footplate covering

**905** Seat of hand  
brakes with single lever

 **914** Footrest lock

 **896** Tray covering

**916** Padded leg panel

**932 Anti-tip device**

933 IV Pole

**936** Headrest  
cover

 **938**  
Water-proof cover

Notes:

MKT.F105 | P25001 | USA | July 2025 | R01

**Amylior** - T 1 888 453-0311 - F 1 877 501-8458 [www.amylior.com](http://www.amylior.com) - [info@amylior.com](mailto:info@amylior.com)