

ORDERING INFORMATION: PHONE: (888) 453-0311, FAX: (877) 501-8458, EMAIL: info@amylior.com

**ACCOUNT INFORMATION**

|                              |       |                        |       |
|------------------------------|-------|------------------------|-------|
| <b>P.O. or Quote Number:</b> | _____ | <b>Bill to (Name):</b> | _____ |
| <b>Dealer Contact Name:</b>  | _____ | <b>Address:</b>        | _____ |
| <b>Email Address:</b>        | _____ | <b>City:</b>           | _____ |
| <b>Phone Number:</b>         | _____ | <b>Province:</b>       | _____ |
| <b>Fax Number:</b>           | _____ | <b>Postal Code:</b>    | _____ |
| <b>User Name or Number:</b>  | _____ |                        |       |

**YOMPER+ / ADD-ON POWER-ASSIST DEVICE FOR MANUAL WHEELCHAIRS**
**YOMPER+**

- AKYOMPEREG-** **Yomper+ with Easy-Go Bluetooth control** \$7,990
- ▲ Included: Carrying case, Yomper charger, mounting bracket and strap, and USB cable for Easy-Go.
  - ▲ Maximum user weight capacity is 330 lb (150 kg)


**BATTERY SELECTION (the battery will be installed in the Yomper+ at the factory)**

- AKYPB160** **Standard Battery 24V 6 Ah / 160Wh lithium-ion** \$955
- AKYPB240** **High-Capacity Battery 24V 10 Ah / 240Wh lithium-ion** \$995
- ▲ This option is available by special order only and may result in additional delays.

**FRAME SUPPORT AND SPECIAL ADAPTATION**

▼ Depending on the configuration, excessive depth may require the addition of a support bar under the seat.

▼ Select a support and adaptation kit for a Motion Composites wheelchair

|                                     |            |                 | Compatible seat width selection / Central bar              |  |  |   |       |
|-------------------------------------|------------|-----------------|--|--|--|---|-------|
| Wheelchair                          | Code       | Wheel Size (in) | Central Bar <b>XS</b>                                      | Central Bar <b>L</b>                                     | Central Bar <b>XL</b>                                      | Yomper external Info.   | PRICE |
| <input type="checkbox"/> Apex A     | AKYPMCAXA  | 24,25,26        | N/A  | N/A  | N/A  | N/A   | N/C   |
| <input type="checkbox"/> Apex C     | AKYPMCAXC  | 24,25,26        | N/A  | N/A  | N/A  | N/A   | N/C   |
| <input type="checkbox"/> Helio A6   | AKYPMCHA6- | 22,24,25,26     | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | <b>CAUTION:</b> Yomper is only compatible with the rear-most wheel position       | \$679 |
| <input type="checkbox"/> Helio A7   | AKYPMCHA7- | 22,24,25,26     | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | Specific support.<br>Compatible if: Camber = 0° and <b>see photo (Appendix A)</b> | \$679 |
| <input type="checkbox"/> Helio C2   | AKYPMCHC2- | 22,24,25,26     | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | Specific support.<br>Compatible if: Camber = 0° and <b>see photo (Appendix A)</b> | \$679 |
| <input type="checkbox"/> Helio Kids | AKYPMCHKD- | 22,24           | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | Specific support.<br>Compatible if: Camber = 0° and <b>see photo (Appendix A)</b> | \$679 |
| <input type="checkbox"/> Move       | AKYPMCMOV- | 22,24,25,26     | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | <b>CAUTION:</b> Yomper is only compatible with the rear-most wheel position       | \$679 |
| <input type="checkbox"/> Veloce     | AKYPMCVLC- | 22,24,25,26     | <input type="checkbox"/> 13.9" - 15.3"<br>(35,4 - 38,9 cm) | <input type="checkbox"/> 15.3" - 18.2"<br>(39 - 46,2 cm) | <input type="checkbox"/> 18.2" - 21.1"<br>(46,3 - 53,6 cm) | Specific support.<br>Compatible if: Camber = 0° and <b>see photo (Appendix A)</b> | \$679 |

▼ Select a support and adaptation kit for a Permobil TiLite wheelchair

|  |            |                 | Compatible seat width selection / Central bar            |   |   |   |       |
|--|------------|-----------------|--|---|---|---|-------|
| Wheelchair                             | Code       | Wheel Size (in) | Central Bar <b>XS</b>                                    | Central Bar <b>L</b>  | Central Bar <b>XL</b>                                     | Yomper external Info.   | PRICE |
| <input type="checkbox"/> TiLite 2GX    | AKYPTL2GX- | 22,24,26        | <input type="checkbox"/> 13" - 14.8"<br>(33 - 37,5 cm)   | <input type="checkbox"/> 14.8" - 17.6"<br>(37,6 - 44,8 cm)  | <input type="checkbox"/> 17.6" to 20.5"<br>(44,9 - 52 cm) |   | \$679 |
| <input type="checkbox"/> TiLite Aero T | AKYPTLAOT  | 24,25,26        | N/A  | N/A   | N/A   |   | N/C   |
| <input type="checkbox"/> TiLite Aero X | AKYPTLAOX- | 22,24,26        | <input type="checkbox"/> 13" - 14.4"<br>(33,1 - 36,7 cm) | <input type="checkbox"/> 14.5" to 17.2"<br>(36,8 - 43,9 cm) | <input type="checkbox"/> 17.3" to 21.8"<br>(44 - 55,4 cm) |   | \$679 |
| <input type="checkbox"/> TiLite Aero Z | AKYPTLAOZ  | 24              | N/A  | N/A   | N/A   |   | N/C   |
| <input type="checkbox"/> TiLite TR     | AKYPTLTR   | 24,25,26        | N/A  | N/A   | N/A   | <b>CAUTION:</b> available with carbon frame bar; if selected, add 2 x 1005164 | N/C   |
| <input type="checkbox"/> TiLite ZRA    | AKYPTLZRA  | 24,25,26        | N/A  | N/A   | N/A   |   | N/C   |

▼ Select a support and adaptation kit for a Ki Mobility (Etac) wheelchair

| Wheelchair                                 | Code        | Wheel Size (in) | Compatible seat width selection / Central bar          |  |   | Yomper external Info. | PRICE |
|--|-------------|-----------------|--|--|---|-----------------------|-------|
|  |             |                 | Central Bar <b>XS</b>                                  | Central Bar <b>L</b>                                       | Central Bar <b>XL</b>                                     |                       |       |
| <input type="checkbox"/> CATALYST 5        | AKYPKMC5-   | 22,24,25,26     | <input type="checkbox"/> 13" - 14.8"<br>(33 - 37,5 cm) | <input type="checkbox"/> 14.8" - 17.6"<br>(37,6 - 44,8 cm) | <input type="checkbox"/> 17.6" to 20.5"<br>(44,9 - 52 cm) |                       | \$679 |
| <input type="checkbox"/> CATALYST 5Vx      | AKYPKMC5VX- | 22,24,25,26     | <input type="checkbox"/> 13" - 14.8"<br>(33 - 37,5 cm) | <input type="checkbox"/> 14.8" - 17.6"<br>(37,6 - 44,8 cm) | <input type="checkbox"/> 17.6" to 20.5"<br>(44,9 - 52 cm) |                       | \$679 |
| <input type="checkbox"/> CATALYSTE 4       | AKYPKMC4-   | 22,24           | <input type="checkbox"/> 13" - 14.8"<br>(33 - 37,5 cm) | <input type="checkbox"/> 14.8" - 17.6"<br>(37,6 - 44,8 cm) | <input type="checkbox"/> 17.6" to 20.5"<br>(44,9 - 52 cm) |                       | \$679 |
| <input type="checkbox"/> CATALYSTE 4C      | AKYPKMC4C-  | 22,24,25,26     | <input type="checkbox"/> 13" - 14.8"<br>(33 - 37,5 cm) | <input type="checkbox"/> 14.8" - 17.6"<br>(37,6 - 44,8 cm) | <input type="checkbox"/> 17.6" to 20.5"<br>(44,9 - 52 cm) |                       | \$679 |
| <input type="checkbox"/> LITTLE WAVE CLICK | AKYPKMLWC   | 24              | N/A  | N/A  | N/A   |                       | \$679 |
| <input type="checkbox"/> ROGUE 2           | AKYPKMR2    | 24,25,26        | N/A  | N/A  | N/A   |                       | N/C   |
| <input type="checkbox"/> ROGUE XP          | AKYPKMRXP   | 24,25           | N/A  | N/A  | N/A   |                       | N/C   |

▼ Select a support and adaptation kit for a Sunrise Medical (Quickie-RGK) wheelchair

| Wheelchair                              | Code       | Wheel Size (in) | Compatible seat width selection / Central bar |                      |                       | Yomper external Info. | PRICE |
|---|------------|-----------------|---|----------------------|-----------------------|-----------------------|-------|
|   |            |                 | Central Bar <b>XS</b>                         | Central Bar <b>L</b> | Central Bar <b>XL</b> |                       |       |
| <input type="checkbox"/> RGK Tiga       | AKYPSUNT   | 24,25,26        | N/A   | N/A                  | N/A                   |                       | N/C   |
| <input type="checkbox"/> RGK Tiga F2    | AKYPSUNTF2 | 24,25           | N/A   | N/A                  | N/A                   |                       | N/C   |
| <input type="checkbox"/> RGK Tiga Fx    | AKYPSUNTFX | 24,25           | N/A   | N/A                  | N/A                   |                       | N/C   |
| <input type="checkbox"/> RGK Tiga Lite  | AKYPSUNTLT | 24,25,26        | N/A   | N/A                  | N/A                   |                       | N/C   |
| <input type="checkbox"/> RGK Tiga Sub 4 | AKYPSUNTS4 | 24,25           | N/A   | N/A                  | N/A                   |                       | N/C   |

**ATTENDANT CONTROL**

- AKYPCOMASSIST **Yomper for attendant control** \$360
- ▲ Can be used with or without Bluetooth control box
- ▲ Includes the mounting of the control unit on the backrest (clip) mounting (1005161).
- ▲ Includes the extension of the attendant control – Length 39.4" (100 cm) (1005162) – Extension available upon request.

**APPENDIX A**

▼ COMPATIBILITY WITH MOTION COMPOSITES HELIO A7, C2, KIDS, AND VELOCE

