

The light weight, easy to use **STELLA** rollator includes features to better help people with reduced mobility in their daily activities.

Utilizing an aluminum built frame, the **STELLA** weighs only 14 lb (6.4 kg) in addition to all of its features as well as the oversized soft material backrest which offers enhanced comfort when seated.

Included features :

- Reinforced seat
- Four 8" (205 mm) caster with smooth urethane Tires
- Locking system when folded
- 20% thicker brake cable
- *Cane and cup holder included
- Only 14 lb (6.4 kg)
- Oversized height adjustable backrest strap
- Adjustable handles with quick release system
- Removable bag
- Cell phone pouch



* Cane not included

Options :

Candy red, blue or black

MODEL	SEAT PAN HEIGHT	HANDLE HEIGHT	DISTANCE BTW HANDLES	OVERALL WIDTH	WHEELBASE	CAPACITY
STELLA 18	18" (455 mm)	29.5" - 33.5" (750 mm - 850 mm)	17.5" (445 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 20	20" (510 mm)	30.5" - 34.5" (775 mm - 875 mm)	17.5" (445 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 23	23" (585 mm)	33.5" - 37.5" (850 mm - 950 mm)	17.5" (445 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA HD	20" (510 mm)	30.5" - 34.25" (775 mm - 870 mm)	19.5" (495 mm)	26" (660 mm)	28.5" (725 mm)	400 lb (182 kg)

LIMITED WARRANTY

Life: frame and cross-braces

5 years: handles

1 year: every other component