

## DELUXE CGT Transport Chair

To improve user experience the **DELUXE CGT transport chair** is equipped with retractable and removable armrests which feature two settings; full length or desk option (to get near tables). The footrests are both retractable and adjustable in length.

The 12" (305 mm) back wheels improve the chair's rolling motion and obstacle handling while the fold-down backrest makes storage for transport so easy.

Dynamic brakes allow the attendant to slow down prior to coming to a complete stop and the wheel locking mechanism makes user transfers safer.



### Included features:

- Dynamic brakes
- 12" (305 mm) back wheels
- Quick release, fold-down backrest
- Flip-up, removable and dual-function armrests (full length and desk length)
- Wheel lock mechanism
- Swing-away/adjustable length footrest

### Available options :

- Width 17", 19" or 22" (430 mm, 485 mm or 560 mm)
- Burgundy or blue

MODEL	WIDTH	WEIGHT	WEIGHT CAPACITY
TC1217DLX	17" (430 mm)	28 lb (13 kg)	285 lb (129 kg)
TC1219DLX	19" (485 mm)	30 lb (14 kg)	285 lb (129 kg)
TC1222DLX	22" (560 mm)	32 lb (15 kg)	400 lb (181 kg)

### LIMITED WARRANTY

**5 years:** frame

**1 year:** every other component