

## STANDARD **CGT** Transport Chair

Find comfort, peace of mind and safety in the standard model of **CGT transport chair**.

Always aiming for the user's comfort, **CGT transport chair** features full length armrests, retractable and length adjustable footrests as well as a quality design and materials.

Aluminum built with a quick-release foldable backrest, this transport chair is ready to travel.

Also, the chair's wheel locking mechanism makes transferring easy and safe.

### Included features:

- Quick-release foldable backrest
- 4 x 8" (205 mm) solid wheels
- Fixed full-length armrests
- Swing-away and length adjustable footrests
- Wheel lock mechanism

### Available options :

- 17" or 19" width
- Burgundy or blue

### LIMITED WARRANTY

**5 years:** frame

**1 year:** every other component



MODEL	WIDTH	WEIGHT	WEIGHT CAPACITY
TC817S	17" (430 mm)	20 lb (9 kg)	285 lb (129 kg)
TC819S	19" (485 mm)	22 lb (10 kg)	285 lb (129 kg)