

Stella



Stella

Meet Stella, there for you every day, enabling easy, confident, real-life mobility.

The light weight, easy to use **STELLA** rollator includes features to help people with reduced mobility in their daily activities.

Built on an aluminum frame, the **STELLA** weighs as little as 14.75 lb (6,7 kg) in its standard configuration.* The oversized soft material backrest offers enhanced comfort and stability when seated.

*Weight does not include optional equipment.

Choice of colors:



LIGHTWEIGHT ROLLATOR

Standard features:

- Reinforced seat
- Four 8" (205 mm) casters with smooth urethane tires
- Locking system when folded
- Stronger brake cables
- Starting at 14.75 lb (6,7 kg)
- Oversized height adjustable backrest strap
- Height adjustable handles with quick release system
- Left and right curb climbers
- Cup holder
- Cane holder** (not compatible with armrests)
- Removable storage bag



**Cane not included

Optional equipment:



- Armrests (Pair: 1.9 lb/0,86 kg)



- Footrests (Pair: 2.8 lb/1,27 kg)



- Slow-down brakes (Pair: 1.1 lb/0,5 kg)



- Transport bag

MODEL	SEAT HEIGHT	HANDLE HEIGHT	SEAT WIDTH	FOOTREST LENGTH	OVERALL WIDTH	OVERALL LENGTH	WEIGHT CAPACITY
STELLA 18	18" (455 mm)	29.5" - 33.5" (750 mm - 850 mm)	17.5" (445 mm)	15" - 17.5" (380 mm - 445 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 20	20" (510 mm)	30.5" - 34.5" (775 mm - 875 mm)	17.5" (445 mm)	17" - 19.5" (430 mm - 495 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 23	23" (585 mm)	33.5" - 37.5" (850 mm - 950 mm)	17.5" (445 mm)	20" - 22.5" (510 mm - 570 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 20 HD	20" (510 mm)	30.5" - 34.5" (775 mm - 875 mm)	19.5" (495 mm)	17" - 19.5" (430 mm - 495 mm)	26" (660 mm)	28.5" (725 mm)	400 lb (182 kg)
STELLA 23 HD	23" (585 mm)	33.5" - 37.5" (850 mm - 950 mm)	19.5" (495 mm)	20" - 22.5" (510 mm - 570 mm)	26" (660 mm)	28.5" (725 mm)	400 lb (182 kg)

NOTE - Some metric measurements are rounded to the nearest 5 mm

LIMITED WARRANTY

Life: Frame and cross-braces | 5 years: Brake handles | 1 year: All other components
Exclusions: Parts subject to wear or misuse such as upholstery, hand-grips, and wheels.



For more information:



1.450.424.0288

1.888.453.0311

info@amylor.com

www.amylor.com



BR | TMK260302 | 1000099 | R09