

CG4216 Bariatric Rollator

Like all of Continent Globe products, the **CG4216 bariatric rollator** was designed in collaboration with therapists, users and medical equipment specialists whose main concern is to alleviate daily challenges of people with reduced mobility.

This sturdy rollator walker is equipped with static and dynamic brakes, 8" urethane tires, hand brakes as well as ergonomic, adjustable handles.

Include features:

- Reinforced rollator
- 8" urethane tires
- Static and dynamic brakes
- Ergonomic and height adjustable handles
- Padded seat
- Basket

Available options :

- Burgundy or blue

LIMITED WARRANTY

5 years: Frame

1 year: All other components



TECHNICAL SPECIFICATIONS	
Overall width	28" (710 mm)
Distance between handles	23" (585 mm)
Handle adjustable height	31" - 37" (785 mm - 940 mm)
Floor / seat height	22" (560 mm)
Weight	28 lb (13 kg)
Weight capacity	500 lb (227 kg)