

TANDEM Wheelchair

Proudly presenting the **TANDEM** aluminum wheelchair by Continent Globe.

The many adjustment features of this wheelchair manage to meet the needs of most users and **TANDEM** is also a great asset for the attendant.

The brake system helps the attendant by slowing the wheelchair to a stop with ease and the adjustable height on the push handles provides comfort and contributes to avoid injury for the attendant.



Included features:

- Push handles, adjustable in height
- Drum brakes for attendant
- Fast adjustment fork angle
- Floor/ seat height adjustable from 17" to 20"
- Tension adjustable backrest
- Flip-up / height adjustable armrest
- Seat width 14", 16", 18" and 20"
- Seat depth 16" and 18"
- Height adjustable legrest
- Retractable footrest (inward and outward)
- Adjustable angle footrest platform
- 24" reinforced wheels with urethane tires
- 7" front wheels with urethane tires
- Quick release axles

SPECIFICATIONS	WC14DB247	WC16DB247	WC18DB247	WC20DB247
Seat width	14" (355 mm)	16" (405 mm)	18" (455 mm)	20" (50.8 cm)
Seat depth	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)
Floor/ seat height	17" - 18" - 19" - 20" (430 - 510 mm)	17" - 18" - 19" - 20" (430 - 510 mm)	17" - 18" - 19" - 20" (430 - 510 mm)	17" - 18" - 19" - 20" (430 - 510 mm)
Backrest height	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)	16" - 18" (405 - 455 mm)
Armrest height	9" - 10" - 11" (230 - 280 mm)	9" - 10" - 11" (230 - 280 mm)	9" - 10" - 11" (230 - 280 mm)	9" - 10" - 11" (230 - 280 mm)
Overall width	23" (585 mm)	25" (635 mm)	27" (685 mm)	29" (735 mm)
Width (folded)	12" (305 mm)	12" (305 mm)	12" (305 mm)	12" (305 mm)
Weight	36.5 lb (16.5 kg)	37 lb (16.9 kg)	38 lb (17.2 kg)	40 lb (18.1 kg)
Weight capacity	330 lb (150 kg)	330 lb (150 kg)	330 lb (150 kg)	330 lb (150 kg)

Available options:

- Plastic coated hand rims
- Elevating legrests
- Burgundy or blue

LIMITED WARRANTY

5 years: frame

1 year: every other component