

INFORMATION BULLETIN

HIP & THIGH SUPPORT – HARDWARE KITS & PADS

INTRODUCTION

There are two changes to the HIP & THIGH SUPPORT:

1. Support Hardware Kit
2. Thigh & Hip Pads

1. NEW HARDWARE KITS

Effective December 2022 or when the current inventory runs out, the Imperial-dimensioned hip & thigh guard support hardware kit, part number C-9612, will be discontinued and replaced with **two (2) Metric-dimensioned kits**:

- a) CURVED Hip & Thigh support hardware kit for Reclining and Cantilever Armrests
 - Serviceable part # 1004020
 - On-chair part # SOHTSARH-C
- b) STRAIGHT Hip & Thigh support hardware kit for Dual Post Armrests
 - Serviceable part # 1004021
 - On-chair part # SOHTSARH-S

Metric curved
hip & thigh support
#1004020
#SOHTSARH-C
Armrest type:
Cantilever: SOAFB
Reclining: SOANI



Adding
Metric straight
hip & thigh support
#1004021
#SOHTSARH-S
Armrest type:
Single Post: SOSP
Dual Post: SOADP



2. THIGH & HIP PADS

The foam density on the pads has changed to a more cushioned feel. We are keeping the same part numbers, but under a new revision. Both supports above are compatible with all old & new thigh pads.

H00097
(Small)
10" x 4"
(255 mm x
100 mm)



H00098
(Medium)
12" x 5"
(305 mm x
125 mm)



H00099
(Large)
14" x 6"
(355 mm x
150 mm)



Corresponding updates

There is no change to the Order Forms. The configurator will choose the appropriate Hip & Thigh Support according to the type of armrests selected. The Parts Manual has been modified and is online.