

INFORMATION BULLETIN

MATRIX SEATING COMPONENTS NOW AVAILABLE ON-CHAIR

INTRODUCTION

Please note that, effective June 5, 2023, select MATRIX seating components will be included on our Alltrack power wheelchairs Order Forms as an on-chair option. Along with our current back supports and seat cushions (CG Air, NXT and Standard), MATRIX seating components below are NOW AVAILABLE.

Note: HCPCS codes are for US use only



SELECT MODELS

ORDER FORM UPDATES (BACKS)

NEW MATRIX BACKS	
<input type="checkbox"/> E2S-WWHH <small>NEW</small>	<p>MATRIX ELITE E2S BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments) – (E2615)</p> <ul style="list-style-type: none"> ▲ The actual back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back ▲ Weight capacity up to 300 lb
<input type="checkbox"/> E2D-WWHH <small>NEW</small>	<p>MATRIX DEEP CONTOUR ELITE E2D BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments) – (E2620)</p> <ul style="list-style-type: none"> ▲ The actual back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back ▲ This product is not compatible with Lateral Trunk Supports ▲ Weight capacity up to 300 lb
<input type="checkbox"/> PBE-WWHH-TR <small>NEW</small>	<p>MATRIX PB ELITE TR BACK – (15" to 20" Wide - Back height of 14", 16", 18" and 20") – (E2620)</p> <ul style="list-style-type: none"> ▲ The actual back shell height is available in 14", 16" 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back ▲ Weight capacity up to 300 lb
<input type="checkbox"/> PBE-WWHH-HD <small>NEW</small>	<p>MATRIX HEAVY DUTY ELITE BACK – (21" to 24" Wide - Back height of 16" and 20") – (E2620 / E2621)</p> <ul style="list-style-type: none"> ▲ The actual back shell height is available in 16" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back ▲ Weight capacity up to 500 lb
<input type="checkbox"/> EDB-WWHH-HD <small>NEW</small>	<p>MATRIX ELITE DEEP HD BACK – (21" to 24" Wide - Back height of 16", 18" and 20") – (E2620 / E2621)</p> <ul style="list-style-type: none"> ▲ The actual back shell height is available in 16", 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back ▲ This product is not compatible with Lateral Trunk Supports ▲ Weight capacity up to 500 lb

New MATRX Seating Components (cont'd)

BACKS

MATRX Elite E2S (E2615)



- Width: 14" to 20", in 1" increments
- Height: 6" to 20", in 2" increments

MATRX Deep Contour Elite E2D (E2620)



- Width: 14" to 20", in 1" increments
- Height: 6" to 20", in 2" increments

MATRX PB Elite TR (E2620)



- Width: 15" to 20", in 1" increments
- Height: 14", 16", 18", and 20"

MATRX Heavy Duty Elite HD (E2620/E2621)



- Width: 21" to 24", in 1" increments
- Height: 16" and 20"

MATRX Elite Deep HD (E2620/E2621)



- Width: 21" to 24"
- Height: 16", 18", and 20"

New MATRX Seating Components (cont'd)

ORDER FORM UPDATES (CUSHIONS)

▶ NEW MATRX CUSHION	
<input type="checkbox"/> PSP-WWDD <small>NEW</small>	MATRX PSP Seat Cushion – (11" to 20" Wide) – (E2607)
	▲ The maximum user weight capacity is 300 lb
	▲ The actual cushion sizes only come in even numbers for width and depth measurements.
<input type="checkbox"/> PSP-WWDD <small>NEW</small>	MATRX PSP Seat Cushion – (22" to 26" Wide) – (E2608)
	▲ The maximum user weight capacity is 550 lb
	▲ The actual cushion sizes only come in even numbers for width and depth measurements.

MATRX PSP Seat Cushion



E2607

Max user weight: 300 lb

Width: 11" to 20"

Depth: 13" to 21"

E2608

Max user weight: 550 lb

Width: 22" to 26"

Depth: 18" to 23"

OTHER IMPLICATIONS

Availability and Pricing

On-chair pricing according to Order Forms. Additional sizes and models available through Customer Service.

Installation

Alltrack power wheelchairs ordered from Amylior will be shipped with backs fully installed and ready for delivery.

Corresponding updates

Order Forms and Parts Manual have been modified and will be online June 5th, 2023.