

INFORMATION BULLETIN

NEW STELLA 23 HD ROLLATOR

INTRODUCTION

As of August 28th, we are adding a new Stella to our line up. The Stella 23 HD is now in stock and available to order.



DETAILS ON THE NEW ADDITION

The **NEW Stella 23 HD** comes with a 23" seat height, and is offered in addition to our existing Stella 20 HD with 20" seat height. Both the Stella 20 HD and Stella 23 HD have a 19.5" seat width, and 400 lb (182 kg) weight capacity. Please see the table below for specifications on our full line of Stella rollators

MODEL	SEAT HEIGHT	HANDLE HEIGHT	SEAT WIDTH	FOOTREST LENGTH	OVERALL WIDTH	OVERALL LENGTH	WEIGHT CAPACITY
STELLA 18	18" (455 mm)	29.5" - 33.5" (750 mm - 850 mm)	17.5" (445 mm)	15" - 17.5" (380 mm - 445 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 20	20" (510 mm)	30.5" - 34.5" (775 mm - 875 mm)	17.5" (445 mm)	17" - 19.5" (430 mm - 495 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 23	23" (585 mm)	33.5" - 37.5" (850 mm - 950 mm)	17.5" (445 mm)	20" - 22.5" (510 mm - 570 mm)	24" (610 mm)	28.5" (725 mm)	300 lb (136 kg)
STELLA 20 HD	20" (510 mm)	30.5" - 34.25" (775 mm - 870 mm)	19.5" (495 mm)	17" - 19.5" (430 mm - 495 mm)	26" (660 mm)	28.5" (725 mm)	400 lb (182 kg)
STELLA 23 HD	23" (585 mm)	33.5" - 38.5" (850 mm - 980 mm)	19.5" (495 mm)	20" - 22.5" (510 mm - 570 mm)	26" (660 mm)	28.5" (725 mm)	400 lb (182 kg)

New Stella 23 HD**ADDITIONAL INFORMATION****Availability and Pricing**

The new Stella 23 HD Rollator carries the same pricing as the Stella 20 HD. Please consult Amylior.com for the most current order form for your location.

Color availability

Black only.

Corresponding updates

- Order Forms and website are updated.
- Flyer is now online and available to order with part ID 1000099-R07.
https://amylior.com/wp-content/uploads/2023/08/FL_P23027_1000099_Stella-HD-23_EN_2023-06_R07-web.pdf
- The Parts Manual is not affected by the additional size.