

INFORMATION BULLETIN

FRONT LIGHTS NOW BRIGHTER


INTRODUCTION

Owing to recent feedback regarding the brightness of the front lights in our optional lighting package, we've taken action to significantly enhance their intensity. We have upgraded the front lights to provide a brighter and more efficient illumination.



BRIGHTER FRONT LIGHTS

After conducting internal tests with previous front lights and today's new front lights, results showed that the luminous intensity (95-100 cd) of the new lights is 8 times higher than the previous lights (11-12 cd). These upgraded lights not only enhance the user's visibility during low-light conditions but also significantly improve how well others can see the user.

Lighting the wall with FORMER front lights



Lighting the wall with NEW front lights



New Light Package (cont'd)**OTHER CONSIDERATIONS****Availability and Pricing**

The new lights will be included in the optional LED Lighting Package, for Alltrack power wheelchairs in early January 2024. The same kit numbers (hardware included) will be used: #1004183 left and #1004184 right. On-chair pricing remains unchanged. Parts pricing requires an RFQ.