

INFORMATION BULLETIN

November 2023

REPLACEMENT SHROUDS WITH OLD COLOR SCHEME

INTRODUCTION

As of October 1st, 2023, replacement shrouds for Alltrack power wheelchairs with painted frames and black/gray shrouds (built before December 2020), will be completely Onyx Black Matte instead of Black and Gray. Please see below images.

VISUAL IMPACT

Alltrack power wheelchair built prior to December 2020 with old color scheme.

ALLTRACK M3 MID-WHEEL DRIVE

ALLTRACK R3 HYBRID-WHEEL DRIVE

Former shroud color in Black and Gray



Replacement shrouds are now completely Onyx Black Matte




Replacement shroud with old color scheme (cont'd)



Former Battery Shroud in Gray

Replacement Battery Shrouds in Onyx Black Matte

Front



Rear



ADDITIONAL CONSIDERATIONS

Availability and Pricing

Original Black and Gray shrouds are no longer available, Onyx Black Matte shrouds are currently in stock. Replacement part numbers are available in the Parts Manual and pricing can be obtained by RFQ.

Corresponding updates

The Parts Manual has been modified and is online.