

# INFORMATION BULLETIN

**AMYLIOR'S 11" POWER SEAT ELEVATION OPTION IS NOW**



## INTRODUCTION

Please note that, effective IMMEDIATELY, our 11" power seat elevation options will be marketed under a new name.



May 2023, United States CMS decision:

*"The Centers for Medicare & Medicaid Services (CMS) finds that power seat elevation equipment on Medicare-covered power wheelchairs (PWCs) falls within the benefit category for durable medical equipment (DME). CMS finds in this national coverage analysis that the evidence is sufficient to determine that **power seat elevation equipment is reasonable and necessary for individuals using complex rehabilitative power-driven wheelchairs** when certain conditions are met."* We couldn't agree more and are thrilled that more beneficiaries will be able to realize the functional and clinical benefits provided by power seat elevation. Amylior continues to work with industry groups as coding and reimbursement guidelines are finalized for this technology.

## PROVEN FOR OVER 15 YEARS

Power Seat Elevation has been available on Alltrack chairs since 2008. With 11" of vertical elevation, this extremely solid, stable and durable design has been helping people reach new heights for 15 years. Available in both standard and heavy-duty configurations, Alta Power Seat Elevation can support user weights up to 400 pounds and allows users to drive safely at up to 25% of the maximum programmed forward speed. All this function is delivered with no additional mechanical lockouts or modifications to the standard suspension, ensuring each user can safely enjoy the full capabilities of their Alltrack base.

## ALTA POWER SEAT ELEVATION (cont'd)

### ADVERTISING IN PROGRESS

The "ALTA" name was chosen carefully. It means HIGH or ELEVATED in Latin and speaks to the importance of the recent decision by CMS to open access for more people to these systems.

**ALTA Power Seat Elevation** will take you to new heights!

**ALLTRACK** with **ALTA** POWER SEAT ELEVATION

**NARROW MEETS PERFORMANCE**

**M3 N2221**

**P3 N2221**

**AMYLIOR**

Like and follow us on social media

Use QR code to visit our website

**N2221 and Alta Highlights**

- 21" base width available on Alltrack M3 or P3 with 14" drive wheels and 6" casters
- Full-size performance in compact package!
- Available with 11" of power seat elevation, 50° of power tilt and 160° of power recline
- Able to drive while fully elevated
- Compatible with V82 (non-expandable) and R-net (expandable) electronics
- Patented Dual-Active Suspension – Same performance and function as full-size M3 and P3

AM | TMK230519 | R01 01

www.amylior.com | Phone: 888-453-0311 | Fax: 877-501-8458 | info@amylior.com

### OTHER IMPLICATIONS

#### Order Forms Updates

##### US version:

<b>SAPESX</b>	<b>ALTA 11" Power Seat Elevation – (E2300)</b> ▲ This option includes the Amyseat Rehab Seat – (SAABA)
<b>ADDPEX</b>	Module – Add <b>ALTA 11" Power Seat Elevation</b> to either <b>SAPT</b> X, <b>SAPR</b> X or <b>SAPTPR</b> X – (E2300)

**Note:** HCPCS codes are for US use only and are reflective of current, established HCPCS codes

##### Canada and International version:

<b>SAPESX</b>	<b>ALTA 11" Power Seat Elevation</b> ▲ This option includes the Amyseat Rehab Seat – (SAABA)
<b>ADDPEX</b>	Module – Add <b>ALTA 11" Power Seat Elevation</b> to either <b>SAPT</b> X, <b>SAPR</b> X or <b>SAPTPR</b> X