


# INFORMATION BULLETIN

November 2023

## *NEW CUSTOM OPTION FOR CG AIR CUSHION*

### INTRODUCTION

Amylior now offers an optional 1" (25mm) thick EVA foam "trellis" that is easily installed around the front rows of cells on the CG Air cushion. The trellis helps to protect the base of these cells, which are often subjected to the most stress and movement during sliding transfers. The EVA material also provides a firmer surface for hand and/or sliding board placement.



### ADDITIONAL DETAILS

The trellis is available as an option when ordering a CG Air cushion. It is defined by the cushion's width and is available in either a 2-cell depth (4.5" or 115mm) or 3 cell depth (6.5" or 165mm). The trellis is only offered in a 1" (25 mm) thickness.

### ADDITIONAL CONSIDERATIONS

#### Availability and Pricing

The trellis is a line-item choice on the CG Air Order form. Refer to it for pricing.

#### Corresponding updates

All Order Forms for CG Air cushions have been updated with this new option and are published on our website. This item will be in the Parts Manual in its next release. An Instruction Guide detailing its installation will be published on our website in the Support/For Pros/Instruction Guides section.