



1. ACCOUNT INFORMATION

P.O. # or Quote: _____

Email Address: _____

Phone Number: _____ Fax Number: _____

Bill to (Name): _____

Bill to (Location): _____

Ship to (Name): _____

Address: _____

City: _____

Province: _____ Postal Code: _____

User Name or Number: _____

User Weight: _____

User Height: _____

Seat Cushion Thickness: _____

Knee to Heel Length: Left: _____ " Right: _____ "



THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:

- A user name or number, the client's measurements and weight - Amylor is not responsible for configuration or size discrepancies resulting from customer errors or omissions on the order form.
- Order confirmations will be provided by Amylor summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

***IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES**

2. STANDARD OPTIONS

Alltrack P3 Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 6" front casters on MWD or 8" on HWD, and standard auto style buckle seat positioning strap, 6.5mph motors, Weight Capacity = 300 lb

3. ALLTRACK P3 POWER WHEELCHAIR

ATPM3	ALLTRACK P3 Power Wheelchair (MWD - MID-WHEEL DRIVE CONFIGURATION)	\$8,020
ATPR3	ALLTRACK P3 Power Wheelchair (HWD - HYBRID-WHEEL DRIVE CONFIGURATION)	\$8,020

3.5 OWNER'S MANUAL

SOOMATENG	English / Anglais	N/C	SOOMATFR	Français / French	N/C
------------------	-------------------	-----	-----------------	-------------------	-----

4. USER WEIGHT PACKAGES

WCA25	250 lb Maximum user weight capacity	STD	WCA3	300 lb Maximum user weight capacity	STD
--------------	-------------------------------------	-----	-------------	-------------------------------------	-----

5. STANDARD COLORS

SCMOBM / SCROBM	Onyx Black Matte	N/C	SCMCAB / SCRCAB	Cabernet	N/C
SCMCTB / SCRCTB	Cobalt Blue	N/C	SCMBGM / SCRBMG	Bubble Gum Metallic	N/C
SCMIOB / SCRIOB	Indigo Blue	N/C	SCMHAB / SCRHAB	Habanero Metallic	N/C
SCMCRE / SCRCRE	Candy Red	N/C	SCMLIM / SCRLIM	SubLime	N/C

▼ OPTIONAL COLORS

FCHDD-M / FCHDD-R	Desert Digital Camo HydroGraphics - (Military Tan base)	\$795
FCHWL-M / FCHWL-R	Carbon fiber HydroGraphics - (Matte black base)	\$795
FCHCF-M / FCHCF-R	Woodland Camo HydroGraphics - (Military Green base)	\$795



5.5 WHEEL INSERT COLORS

SCCTB-WI	Cobalt Blue	N/C	SCHAB-WI	Habanero Metallic	N/C
SCIOB-WI	Indigo Blue	N/C	SCLIM-WI	SubLime	N/C
SCCRE-WI	Candy Red	N/C	SCBGM-WI	Bubble Gum Metallic	N/C
SCCAB-WI	Cabernet	N/C			
▲ <i>Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)</i>					
▲ <i>Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)</i>					
▲ <i>Wheel inserts are not available in Onyx Black Matte</i>					

7. TRANSPORT OPTION

TTD4	4-Point wheelchair transport brackets			STD
▲ <i>4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.</i>				
WC19-ISO	WC-19 and ISO 7176-19 Occupied Transit Package (includes an approved pelvic safety belt)		CLIENT COST	\$595
▲ <i>To confirm if your overall wheelchair configuration is WC-19 and ISO 7176-19 compliant, please refer to this website: https://docs.google.com/spreadsheets/d/1KtcNIEA5OQJCS5z3OD1adWEIaFLN1XO3ZcTc6DI_kaE/edit#gid=1917752765</i>				
▲ <i>Other restrictions may apply, please contact CSR for details.</i>				
▲ <i>This option is only available on new power wheelchairs as an "on-chair" option. It is not available to install on a delivered wheelchair</i>				
▲ <i>Not compatible with the narrow base configurations (FS2424B and FS2221B)</i>				

8. MOTOR PACKAGE

MP65B	6.5 MPH 4-Pole Motor Package / with Brake Release			STD
MP65-OPT	6.5 MPH 4-Pole Motor Package / with COG Release		CLIENT COST	\$395
▲ <i>Not compatible with the narrow base configuration (FS2221B)</i>				
MP85	8.5 MPH 4-Pole High Speed Motor Package / with COG Release		CLIENT COST	\$395
▲ <i>Maximum user weight for this option is 250 lb and requires the R-net 120 Amp Expandable Controller (ECR120).</i>				
▲ <i>This motor package is only available with the HWD configuration (ATPR3)</i>				

9. ASSISTIVE DRIVING PACKAGE

SMARTT	Smart-Track Assistive Driving Package			\$1,420
▲ <i>Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.</i>				
▲ <i>Requires the R-net PM 120 Amp Expandable Controller (ECR120).</i>				
▲ <i>This option requires either the MP65-OPT or MP85 motor package.</i>				
▲ <i>Not compatible with the narrow base configuration (FS2221B)</i>				
▲ <i>Recommended for non-proportional Drive Controls</i>				

10. BATTERY CHARGER

BCOB	8A Offboard Battery Charger			N/C
-------------	-----------------------------	--	--	-----

11. BATTERY TRAY AND BASE WIDTH

FS2221B	N2221 - Battery Tray for Group 22 Batteries with a Narrow 21" base width (Total base width = 21")			N/C
▶ <i>This option is only available with the MWD configuration (ATPM3). ▶ NOT AVAILABLE WITH THE FOLLOWING OPTIONS:</i>				
▲ <i>WC-19 and ISO 7176-19 Occupied Transit Package (WC19-ISO)</i>				
▲ <i>6.5 MPH 4-Pole Motor Package / with COG Release option (MP65-OPT)</i>				
▲ <i>Smart-Track option (SMARTT)</i>				
▲ <i>8" front casters (FC8PTM and FC8PIM)</i>				
▲ <i>14" x 4" Extra Wide Cast Wheels, Pneumatic Tires (DW144PT and DW144PI)</i>				



11. BATTERY TRAY AND BASE WIDTH (CONTINUED)

FS22B	Battery Tray for Group 22NF Batteries – (Total base width = 23 3/8")	STD
▲	<i>For the HWD configuration, the 23 3/8" compact base is not compatible with most of the Center Mount Front Rigging Options due to mechanical conflicts</i>	
▲	<i>The Straight Front Arm option is not available with any of the Center Mount Front Rigging Options due to mechanical conflicts</i>	
FS2424B	N2424 - Battery Tray for Group 24 Batteries with a Narrow 24" base width (Total base width = 24")	N/C
▲	<i>This option is only available with the MWD configuration (ATPM3).</i>	
▶	SELECT ONE FRONT CASTER ARM OPTION FOR THE HWD CONFIGURATION:	Flared out (STD) Straight (Optional)

12. BATTERIES - REQUIRES A SELECTION

BAT22NFGEL	GEL Battery – Group 22NF – 50 Amps (Each)	AADL CONTRACT	RFQ
BAT24GEL	GEL Battery – Group 24 – 70 Amps (Each)	AADL CONTRACT	RFQ
BATO	Omission of the batteries		N/C

14. DRIVE WHEEL

DW14PT	14" Cast Wheel Pneumatic Tires		STD
DW14PI	14" Cast Wheel Pneumatic with Foam Filled Tires		N/C
DW144PT	14" x 4" Extra Wide Cast Wheels, Pneumatic Tires	CLIENT COST	\$200
DW144PI	14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires	CLIENT COST	\$595
▲	<i>Increases the overall width of the base by 2 3/4"</i>		
▲	<i>Not compatible with the narrow base configuration (FS2221B)</i>		
DW12PT	12 1/2" Cast Wheel Pneumatic Tires		N/C
DW12PI	12 1/2" Cast Wheel Pneumatic with Foam Filled Tires		N/C
▲	<i>12 1/2" Drive wheels lower ground clearance and speed, but increase motor torque.</i>		

15. FRONT CASTERS

▶ MWD CONFIGURATION

FC6S	6" Solid Casters – (Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI)	STD
FC8PTM	8" Pneumatic Casters	N/C
FC8PIM	8" Pneumatic Casters with Foam Filled Tires	N/C
▲	<i>The 8" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)</i>	
▲	<i>Not compatible with the narrow base configuration (FS2221B)</i>	

▶ HWD CONFIGURATION

FC8PTR	8" Pneumatic Casters	STD
FC8PIR	8" Pneumatic Casters with Foam Filled Tires	N/C
FC9PT	9" Pneumatic Casters	\$150
FC9PI	9" Pneumatic Casters with Foam Filled Tires	\$175
▲	<i>The 9" casters will increase the minimum seat to floor height by 1/2"</i>	
▲	<i>The 9" casters are NOT compatible with the 12 1/2" Drive wheel configuration (DW12PT and DW12PI)</i>	
▲	<i>The 9" casters are NOT compatible with many front rigging options. Please see details in the Front Rigging sections.</i>	



16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± ¼")

▶ **YELLOW BOXES REPRESENT HEIGHTS ONLY AVAILABLE WITH 12 ½" DRIVE WHEELS (SECTION 14)**

▶ **SEAT-TO-FLOOR HEIGHT FOR PEDIATRIC SEAT PAN (STS11-16-16 IN SECTION 37.5)**

SEAT TO FLOOR HEIGHT TABLE SEAT ANGLE	NO POWER OPTIONS			TILT <u>OR</u> ELEVATE (WITH OR WITHOUT RECLINE)	TILT <u>AND</u> ELEVATE (WITH OR WITHOUT RECLINE)
	0°	3°	6°		
STFH1675 (16.75" STFH)					
STFH1775 (17.75" STFH)					
STFH1875 (18.75" STFH)					
STFH1975 (19.75" STFH)					
STFH2075 (20.75" STFH)					

▶ **SEAT-TO-FLOOR HEIGHT FOR STANDARD SEAT PAN (STS14-19-19 IN SECTION 38)**

SEAT TO FLOOR HEIGHT TABLE SEAT ANGLE	NO POWER OPTIONS			TILT <u>OR</u> ELEVATE (WITH OR WITHOUT RECLINE)	TILT <u>AND</u> ELEVATE (WITH OR WITHOUT RECLINE)
	0°	3°	6°		
STFH1550 (15.50" STFH)					
STFH1650 (16.50" STFH)					
STFH1750 (17.50" STFH)					
STFH1850 (18.50" STFH)					
STFH1950 (19.50" STFH)					

17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)

SOMHSS	Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges.	RFQ
▲	Required at: _____ " (± ¼" Range)	

18. PRE-SET SEAT FRAME ANGLE

SFE-0	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
SFE-3	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
SFE-6	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

19. LIGHTING PACKAGE

EOLP	LED Light Package – (Includes front and rear lights with integrated flashers)	CLIENT COST	\$1,425
▲	Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL, ECELJL, ECRJCL, ECRJC35B or EDIB).		

21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS

ECEL90	R-net EL 90 Amp Expandable Controller	\$175
ECR120	R-net PM 120 Amp Expandable Controller	\$275

22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

▶	SELECTION OF ONE JOYSTICK OPTION BELOW IS REQUIRED WITH THE EXPANDABLE CONTROLLER ECEL90 OR ECR120	
ECO	Omit Joystick – (For use with an Omni 2 Display)	N/C
ECELJ	R-net - (LED) – Joystick	STD
▲	This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)	
ECRJC	R-net Joystick with Color screen	\$155



22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS (CONTINUED)

ECELJL	R-net – (LED) – Joystick with Lighting Operation	CLIENT COST	\$125
	▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>		
ECRJCL	R-net Joystick with Color screen and Lighting Operation		\$180
▶ ADVANCED JOYSTICK			
ECRJ35B	R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	CLIENT COST	\$1,105

23. JOYSTICK MOUNT

JSMF	Fixed Joystick Mount		N/C
JMSA	Swing Away Joystick Mount – Height adjustable, Angle adjustable and Multi-Axis		\$235
Select:	Set-up with standard installation to swing-away sideways	Set-up at 35° to swing-away sideways and downwards	
Select:	Mechanical Lock (STD)	Magnetic Lock	
JSMCM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount		\$275

24. JOYSTICK MOUNT SUPPORT

JSMFH	Joystick mount support Bracket		STD
--------------	--------------------------------	--	-----

25. JOYSTICK HAND MOUNT

HML	Install Joystick on Left hand side	HMR	Install Joystick on Right hand side
------------	------------------------------------	------------	-------------------------------------

26. JOYSTICK HANDLE

JSH	Standard Joystick Gimble	STD	JBAT	Stick Handle 3" in Length	\$130
JBH	Ball Handle	\$85	JFBH	Foam Ball Handle	\$70
JTH	T-Handle	\$155	JMUSH	Bodypoint Mushroom Handle	\$210
JGP	Bodypoint 4" U-shaped Handle	\$210			

29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)

EDIB	Omni 2 – (Enhanced Display & Digital Interface Box with Infrared)		\$1,915
	▲ <i>Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120)</i>		
	▲ <i>The Omni 2 display mount must be selected from the option below.</i>		

30. OMNI 2 ENHANCED DISPLAY MOUNT

EDSAM	User Swing Away Mount	CODE W100	N/C
	EDML – (Install Display on LH side)	EDMR – (Install Display on RH side)	

31. ENVIRONMENTAL CONTROL MODULE

ECU	Environment control module – (Output Module)	CLIENT COST	\$995
	▲ <i>Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).</i>		

33. AMYSEAT REHAB SEAT / POWER SEATING SYSTEMS AND MODULES

SAABA	Amyseat Rehab Seat – (Seat Frame only)		\$880
SAPTX	CG Power Tilt	CODE W995	\$6,195
	▲ <i>This option includes the Amyseat Rehab Seat – (SAABA)</i>		



33. AMYSEAT REHAB SEAT / POWER SEATING SYSTEMS AND MODULES (CONTINUED)

SAPRX	Power Recline with Mechanical Shear Reduction (MSR)	CODE W995TR	\$7,845
	▲ <i>This option includes the Amyseat Rehab Seat – (SAABA)</i>		
	▲ <i>This option is not available with the pediatric seat frame SOWN11-16</i>		
SAPTPRX	CG Power Tilt and Power Recline with Mechanical Shear Reduction (MSR)	CODE W995TR	\$11,495
	▲ <i>This option includes the Amyseat Rehab Seat – (SAABA)</i>		
	▲ <i>This option is not available with the pediatric seat frame SOWN11-16</i>		
SAPESX	ALTA 11" Power Seat Elevation	CLIENT COST	\$4,495
	▲ <i>This option includes the Amyseat Rehab Seat – (SAABA)</i>		
ADDPESX	Module – Add ALTA 11" Power Seat Elevation to either SAPTX, SAPRX or SAPTPRX	CLIENT COST	\$2,495
ADDMR	Module – Manual Recline with Mechanical Shear Reduction – (MSR)		\$1,125
	▲ <i>The maximum user weight capacity with this option is 250 lb.</i>		
	▲ <i>This option is not available with the pediatric seat frame SOWN11-16</i>		

34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER

▶ SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (For both expandable and non-expandable systems)								
SOLO	Single function controller (Includes an external toggle switch)			STD				
▼ OPTIONAL SWITCHES								
SWTOGGLE	Single Toggle Switch			Attendant Activation \$100				
SWPUSH1	Single Push Button Switch – (Single Mode: SAC)	User Activation	Attendant Activation	\$90				
SWPUSH2	Dual Push Button Switch – (Dual Mode: SAC)	User Activation	Attendant Activation	\$180				
SSBS	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$230
SSES	Egg Switch-Black			User Activation	Attendant Activation	\$116		
SSML	Ultra Light Switch			User Activation	Attendant Activation	\$295		
▶ SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL								
EPSEL	Single Function Through Drive Control with R-net 90-amp controller			N/C				
	▲ <i>Allows for activation of single power option through R-net 90-amp controller</i>							
	▲ <i>Compatible with all R-net joysticks</i>							
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>							
	▲ <i>Can not operate more than one power option</i>							
	▲ <i>If selecting memory seating, choose SAC below</i>							
	▲ <i>If selecting 120-amp R-net controller, choose SAC below</i>							
SAC	Single Function Through Drive Control with R-net 90-amp or 120-amp controller			N/C				
	▲ <i>Allows for activation of single power option through R-net drive controls</i>							
	▲ <i>Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls</i>							
	▲ <i>Capable of operating tilt with power Seat Elevate & compatible with memory seating</i>							
▶ SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE)								
SACMEMORY	Memory Seat Program			CLIENT COST	\$1,495			
	▲ <i>This option is not compatible with LED joysticks (ECELJ and ECELJL)</i>							



34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER (CONTINUED)

▶ **SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH**

SACTOGGLE	Single Toggle Switch					Attendant Activation	CODE W100	\$410
SACPUSH1	Single Push Button Switch – (Single Mode: SAC)				User Activation	Attendant Activation	CODE W100	\$400
SACPUSH2	Dual Push Button Switch – (Dual Mode: SAC)				User Activation	Attendant Activation	CODE W100	\$490
SACSSBS	Buddy Button Switch				User Activation	Attendant Activation	CODE W100	\$540
	▲	Black	Red	Green	Yellow			
SACSSSES	Egg Switch-Black				User Activation	Attendant Activation	CODE W100	\$426
SACSSML	Ultra Light Switch				User Activation	Attendant Activation	CODE W100	\$605

35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE

▶ **MULTIPLE SEAT FUNCTIONS THROUGH R-NET**

MAC	E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators							\$2,295
	▲ <i>Required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2221B)</i>							
MAC OPTIONS								
MACMEMORY	Memory Seat Program							\$175
SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH								
MACKEYPAD-P	Keypad w/ 5 Push Button Switches			User Activation		Attendant Activation	CODE W100	\$445
MACKEYPAD-T	Keypad w/ 5 Toggle Switches			User Activation		Attendant Activation	CODE W100	\$445
MACMULTI	4-Way Toggle Switch – (Up to 4 Activations)			User Activation		Attendant Activation	CODE W100	\$180
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>							
MACMULTIT	Triple Toggle Switch – (Up to 5 Activations)			User Activation		Attendant Activation	CODE W100	\$135
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>							
MACPUSH4	4-Way Push Button Switch			User Activation		Attendant Activation	CODE W100	\$225
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>							
MACPUSH5	5-Way Push Button Switch			User Activation		Attendant Activation	CODE W100	\$240
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>							

36. SWITCH ACTIVATION HAND MOUNT

Switch Installation:	TSML – (Install Switch on LH side)	TSMR – (Install Switch on RH side)	N/C
-----------------------------	---	---	-----

37. PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)

SOBAA	Select one:	86°	90°	94°	98°	102°	106°	110°	N/C
--------------	-------------	-----	-----	-----	-----	------	------	------	-----



37.5 PEDIATRIC ADJUSTABLE SEAT WIDTH RANGE BETWEEN 11" AND 16"

SOWN11-16	Width adjustment from 11" to 16" with a maximum depth of 16".					N/C
11" (SOWN11)	12" (SOWN12)	13" (SOWN13)	14" (SOWN14)	15" (SOWN15)	16" (SOWN16)	

38. STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"

SOWN15-20-19	Width adjustment from 15" to 20" with a maximum Depth of 19"					N/C
15" (SOWN15)	16" (SOWN16)	17" (SOWN17)	18" (SOWN18)	19" (SOWN19S)	20" (SOWN20S)	

40. ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 19"

▶ PEDIATRIC ADJUSTABLE SEAT DEPTH RANGE BETWEEN 11" AND 16"						
11" (SODN11)	12" (SODN12)	13" (SODN13)	14" (SODN14)	15" (SODN15)	16" (SODN16)	N/C
▶ ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 19"						
15" (SODN15)	16" (SODN16)	17" (SODN17)	18" (SODN18)	19" (SODN19)		N/C

41. BACK CANE HEIGHT (THE BACK CANE HEIGHT MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE PUSH HANDLES)

16" (SOBCH16)	18" (SOBCH18)	20" (SOBCH20)	22" (SOBCH22)	24" (SOBCH24)	26" (SOBCH26)	
----------------------	----------------------	----------------------	----------------------	----------------------	----------------------	--

42. BACK SUPPORT

DSTR-WWHH	Adjustable Tension Back Upholstery, WW x HH		N/C
	▲ <i>Only available with Dual Post armrests (SOADP)</i>		
STDSB-WWHH	Contour back support up to 21" wide, with a maximum shell height of 22".	SEATING BENEFIT	\$605
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>		
▶ MATRX BACKS			
E2S-WWHH	MATRX ELITE E2S BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	SEATING BENEFIT	\$995
	▲ <i>Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>		
	▲ <i>Weight capacity up to 300 lb</i>		
E2D-WWHH	MATRX DEEP CONTOUR ELITE E2D BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	SEATING BENEFIT	\$1,155
	▲ <i>Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>		
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>		
	▲ <i>Weight capacity up to 300 lb</i>		
PBE-WWHH-TR	MATRX ELITE TR BACK – (15" to 20" Wide - Back height of 14", 16", 18" and 20")	SEATING BENEFIT	\$910
	▲ <i>Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>		
	▲ <i>Weight capacity up to 300 lb</i>		
▶ NXT BACKS			
NOACB-WWHH	NXT OPTIMA FH AC CONTOUR BACK – (15" to 21" Wide)	SEATING BENEFIT	\$690
	▲ <i>Weight capacity up to 300 lb</i>		
NOPCB-WWHH	NXT OPTIMA FH PC DEEP CONTOUR BACK – (15" to 21" Wide)	SEATING BENEFIT	\$790
	▲ <i>Weight capacity up to 300 lb</i>		
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>		
	▲ <i>Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>		
▶ OMISSION OF THE BACK SUPPORT			
STDSBO	Omission of the Contour Back Support ◀ Specify: _____		N/C



43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

							▼ 22" is the minimum back height required with a power recline.
16" (SOBH16)	18" (SOBH18)	20" (SOBH20)	22" (SOBH22)	24" (SOBH24)	26" (SOBH26)	28" (SOBH28)	
This height is not available with the NXT Optima AC and PC options							▲

44. SEAT CUSHIONS (SEATING BENEFIT)

▶ FOR CUSTOM SIZES AND CONFIGURATIONS, PLEASE REFER TO THE PERFORMA CG AIR ORDER FORM.

▶ PERFORMA CG AIR - THERAPEUTIC AIR CUSHION

Custom Cushion (Please use the Performa CG Air order form)

Cell Height (Choose 1): 2" 3" 4"

▶ SINGLE (1) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING

▶ PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER

PCGA1VSCxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$812
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER

PCGA1V3Dxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$857
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER

PCGA1VICxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$915
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ DUAL (2) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING

▶ PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER

PCGA2VSCxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$812
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER

PCGA2V3Dxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$857
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER

PCGA2VICxxx-x	14" to 20" (355 mm to 510 mm) wide	SEATING BENEFIT	\$915
---------------	------------------------------------	------------------------	--------------

▲ Cushion sizes only come in even numbers for width and depth measurements.

▶ DUAL VALVE CUSHION SEGMENTATION (PLEASE SELECT 1 BOX)

PCGA-2VLR	Dual Valves with Left / Right Location	SEATING BENEFIT	N/C
-----------	--	------------------------	-----

PCGA-2VF6	Dual Valves with Front / Back Location, Front distance 6 in (150 mm) from the front	SEATING BENEFIT	N/C
-----------	---	------------------------	-----

PCGA-2VF8	Dual Valves with Front / Back Location, Front distance 8 in (205 mm) from the front	SEATING BENEFIT	N/C
-----------	---	------------------------	-----

PCGA-2VF10	Dual Valves with Front / Back Location, Front distance 10 in (255 mm) from the front	SEATING BENEFIT	N/C
------------	--	------------------------	-----

▶ PERFORMA CG AIR OPTIONS

PCGA-GZ	GripZone Anti-Slip Technology	SEATING BENEFIT	\$86
---------	-------------------------------	------------------------	-------------

PCGA-TT-WW	Transfer Trellis	SEATING BENEFIT	\$100
------------	------------------	------------------------	--------------

PCGA-AC	Anatomic Contour Seat Insert	SEATING BENEFIT	\$100
---------	------------------------------	------------------------	--------------

▲ Only available with 2" and 3" cell heights

PCGA-SP	Seat Pouch	SEATING BENEFIT	\$55
---------	------------	------------------------	-------------

▶ CONTOUR SEAT CUSHION

SCPSC-WWDD	Contour Seat Cushion up to 19" Wide	SEATING BENEFIT	\$620
------------	-------------------------------------	------------------------	--------------

SCPSC-WWDD	Contour Seat Cushion for Non-Standard Sizes of 20" to 24" Wide	SEATING BENEFIT	\$735
------------	--	------------------------	--------------

▶ MATRX CUSHION

PSP-WWDD	MATRX PSP Seat Cushion – (11" to 20" Wide)	SEATING BENEFIT	\$645
----------	--	------------------------	--------------

▲ The maximum user weight capacity is 300 lb

▲ Cushion sizes only come in even numbers for width and depth measurements.



44. SEAT CUSHIONS (CONTINUED) (SEATING BENEFIT)

▶ **NXT BIO-FIT CUSHION**

NXTBIOFIT-WWDD	NXT BIO FIT Cushion – (15” to 21” Wide)	SEATING BENEFIT	\$595
	▲ The maximum user weight capacity is 300 lb		
	▲ Cushion sizes only come in even numbers for width and depth measurements.		

45. CANTILEVER – HEIGHT ADJUSTABLE ARMRESTS – OPENS WITH A 5° ANGLE (CHOOSE 2)

SOAFBDLH	Flip-Back Cantilever – Desk Length Left Side – (Height Adjustable from 7” to 15” with 1” increments)	STD Length 12”	N/C		
SOAFBDRH	Flip-Back Cantilever – Desk Length Right Side – (Height Adjustable from 7” to 15” with 1” increments)	STD Length 12”	N/C		
SOAFBFLH	Flip-Back Cantilever – Full Length Left Side – (Height Adjustable from 7” to 15” with 1” increments)	STD Length 16”	STD		
SOAFBFRH	Flip-Back Cantilever – Full Length Right Side – (Height Adjustable from 7” to 15” with 1” increments)	STD Length 16”	STD		
CUST-112250	Custom Armrest Length – Select Length for the Left side:	Custom Length 18”	Custom Length 20”	Custom Length 22”	N/C
CUST-112250	Custom Armrest Length – Select Length for the Right side:	Custom Length 18”	Custom Length 20”	Custom Length 22”	N/C
CUST-SOAFB	Select this code when using Cantilever armrests with Reclining backs				N/C
	▲ The maximum user weight capacity is 300 lb				

46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)

SOADPPLH	Pediatric Dual Post Height Adjustable – Left Side – (Height Adjustable from 7 1/4” to 10 1/4” with 1/2” increments)	\$168
SOADPPRH	Pediatric Dual Post Height Adjustable – Right Side – (Height Adjustable from 7 1/4” to 10 1/4” with 1/2” increments)	\$168
SOADPDLH	Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 1/4” to 15 1/4” with 1/2” increments)	\$168
SOADPDRH	Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 1/4” to 15 1/4” with 1/2” increments)	\$168
SOADPFLH	Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 1/4” to 15 1/4” with 1/2” increments)	\$168
SOADPFRH	Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 1/4” to 15 1/4” with 1/2” increments)	\$168

47. RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (CHOOSE 2)

SOANIDLH	Reclining armrest – Desk Length Left Side – (Height Adjustable from 10 1/2” to 14” with 1/2” increments)	STD Length 12”	N/C		
SOANIDRH	Reclining armrest – Desk Length Right Side – (Height Adjustable from 10 1/2” to 14” with 1/2” increments)	STD Length 12”	N/C		
SOANIFLH	Reclining armrest – Full Length Left Side – (Height Adjustable from 10 1/2” to 14” with 1/2” increments)	STD Length 16”	N/C		
SOANIFRH	Reclining armrest – Full Length Right Side – (Height Adjustable from 10 1/2” to 14” with 1/2” increments)	STD Length 16”	N/C		
CUST-48247	Custom Armrest Length – Select Length for the Left side:	Custom Length 18”	Custom Length 20”	Custom Length 22”	N/C
CUST-48247	Custom Armrest Length – Select Length for the Right side:	Custom Length 18”	Custom Length 20”	Custom Length 22”	N/C



48. ARMREST HEIGHT ADJUSTMENT AT FACTORY

- ▶ Set the Armrest Height at: _____ “
- ▲ The Armrest Height is taken from the seat pan to the top of the armpad.
- ▲ Not required for Dual Post.

49. ARMPADS

▶ FOAM TYPE SELECTION

SOWAPD	Desk Length Waterfall – (8 ½” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	STD
	▲ <i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>				
SOWAPF	Full Length Waterfall – (13” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	STD
	▲ <i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>				

▶ HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)

SOGELD	Desk Length Hybrid Gel and Foam Flat armpad – (9 ¾” x 2”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$40
SOGELF	Full Length Hybrid Gel and Foam Flat armpad – (13 ½” x 2”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$40
SOWGELD	Desk Length Hybrid Gel & Foam Waterfall armpad – (8 ½” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$95
	▲ <i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>					
SOWGELF	Full Length Hybrid Gel and Foam Waterfall armpad – (13” x 2 ¾”)	Select side:	Left Hand Side	Right Hand Side	SEATING BENEFIT	\$95
	▲ <i>May interfere with flipping back the armrests when using RECLINING armrests (SOANlxxx) and/or combined with Lateral Trunk supports.</i>					

▶ SPECIALTY TYPE SELECTION

SOCAR	Arm Trough with Hand Support (18” x 4 ½”)	Select side:	Left Hand Side	Right Hand Side	\$290
	▲ <i>Arm Trough Armpads may interfere with Armrests Flipping back.</i>				

50. FOOTRESTS

SOPA AF	Signature Pediatric Swing-Away Footrests					STD
	Set at:	60°	70°	80°	90°	
	Select Range:	SOFLP5 (5” - 10”)		SOFLP5L (5” - 12”)		
	▲ <i>Depending on the seat configuration for both MWD and HWD configurations, the 80° and 90° angles may not be available.</i>					
	▲ <i>On the HWD version, 70°, 80° and 90° angles are not compatible with the 9” front casters (FC9PT and FC9PI), and may interfere with the casters when selecting the straight front caster arms option.</i>					
	▲ <i>When selecting SOFLP5 or SOFLP5L ranges and the Calf pads option (SOACPFR), the height range is 6 3/8” - 10” for SOFLP5, and 6 3/8” - 12” for SOFLP5L</i>					
SOSA AF	Signature Swing-Away Footrests					STD
	Set at:	60°		70°		
	Select Range:	SOFLA12 (12 ½” - 17 ½”)		SOFLA14 (14 ½” - 19 ½”)		
	▲ <i>This option is not available with the pediatric seat frame (SOWN11-16).</i>					
	▲ <i>On the HWD version, the 70° footrest angle is not compatible with the 9” front casters (FC9PT and FC9PI), and may interfere with the casters when selecting the straight front caster arm option.</i>					



50. FOOTRESTS (CONTINUED)

SOAAF	Signature Heavy Duty Swing-Away Footrests									\$175
		Set at:	60°		70°					
		Select Range:	SOFLHD12 (12 ½" – 17 ½")		SOFLHD14 (14 ½" – 19 ½")					
	▲	<i>This option is not available with the pediatric seat frame (SOWN11-16).</i>								
	▲	<i>On the HWD version, the 70° footrest angle is not compatible with the 9" front casters (FC9PT and FC9PI), and may interfere with the casters when selecting the straight front caster arm option.</i>								
ADDML	Signature Manual Elevating & Articulating Legrests - Pair									\$595
		Select Range:	SOFLME12 (12 ½" – 17 ½")		SOFLME14 (14 ½" – 19 ½")					
	▲	<i>Does not include calf pads. Please select calf pads in section 51.</i>								
ADDPL	Signature Power Elevating & Articulating Legrests								CLIENT COST	\$2,495
		Select Function:	Individual		Combined					
		Select Range:	SOFLPW12 (12 ½" – 17 ½")		SOFLPW14 (14 ½" – 19 ½")					
	▲	<i>This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.</i>								
	▲	<i>This option is not available with the pediatric seat frame (SOWN11-16)</i>								

51. FOOTPLATES & ACCESSORIES

SOSFAF	Flip-up Footplates									STD
SOAFP	Angle Adjustable & Flip-up Footplates – (Each)									\$65
SOAHL	Heel Loops – (Each)									\$20
SOACPFR	Calf Pads for Adult Heavy-Duty Swing-Away footrests – (SOAAF) – with adjustable hardware – (Each)									\$93
SOAMA	Swing-away residual limb Support	Select:	Left leg	Right leg						\$455

52. FOOTREST LENGTH ADJUSTMENT AT FACTORY

Footrest Length	Set at:	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"	
Footrest / Legrest Length	Set at:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"
		16.5"	17"	17.5"	18"	18.5"	19"	19.5"	

▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

53. CENTER MOUNTS

SOAAF90	Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°)									\$200
	Select Foot Platform Range:	SOFBP7 (8 ½" – 12 ½")	SOFBA12 (12 ½" – 16 ½")	SOFBA16 (16 ½" – 19 ½")						
	▲	<i>When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 18.75" (MWD)</i>								
	▲	<i>When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 19" (HWD)</i>								
	▲	<i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>								
	▲	<i>With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.</i>								
	▲	<i>This option is not compatible with the 9" Caster option on HWD Configuration</i>								



53. CENTER MOUNTS (CONTINUED)

ADDPF	90° Power Center Mount Elevating & Articulating Foot Platform	\$2,495
	Select Footrest Range: SOFBPW12 (10 ½" - 16") SOFBPW16 (16" - 20")	
▲	<i>When selecting SOFBPW16 (16" - 20") range, the lowest possible seat to floor height is 18.75" (MWD)</i>	
▲	<i>When selecting SOFBPW16 (16" - 20") range, the lowest possible seat to floor height is 19" (HWD)</i>	
▲	<i>This option includes calf pads. The power center mount counts as one power function.</i>	
▲	<i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB)</i>	
▲	<i>This option is not available with the pediatric seat frame (SOWN11-16).</i>	
▲	<i>The power center mount option is not available with the straight front caster arms (HWD)</i>	
▲	<i>This option is not compatible with the 9" Caster option on HWD Configuration</i>	
▲	<i>For the HWD configuration, the 23 3/8" compact frame (FS22B) is not compatible with most of the Center Mount Front Rigging options due to mechanical conflicts</i>	
▲	<i>For leg lengths between 10 ½" and 12 ½", there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.</i>	
▲	<i>The starting angle is reduced from 90° to 85° for some seat configurations when used with the narrow base configuration (FS2424B)</i>	
▲	<i>The E-BX (MAC) option is required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2424B) on the MWD configuration.</i>	

54. FOOT PLATFORMS & ACCESSORIES

SOFUFP	Flip-up one piece Foot Platform	\$140
SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount – (Each)	\$65
SOACPCM	Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each)	\$93

55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY

Fixed Center Mount Foot Platform Length		Set at:									
		8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"		
12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"		
17.5"	18"	18.5"	19"	19.5"							
Power Center Mount Foot Platform Length		Set at:									
10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"		
15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		

▲ The Foot Platform Length is taken from the top of the seat pan to the top of the foot platform or split foot plates

56. SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19-ISO OPTION)

SOAUTO	55" Seat Positioning Strap	STD
SOBPHP2C	Bodypoint padded Hip Belt – (Two Point – Center Pull)	\$170
SOPHB	Padded Hip Belt, 2-PT Magnetic Center Pull	\$194
SOBPCS	Bodypoint padded Chest Strap	\$250
SOPCS	Padded Chest Strap, 2-PT Magnetic Center Pull	\$262

57. HEADREST

▼	SELECT ONE OF THE HEADREST OPTIONS WITH THE MOUNTING HARDWARE AND SPECIFY WHERE THE HEADREST WILL BE INSTALLED.	
SONUAHR	Using headrest mounting hardware other than SOMAH located below. ▲ Specify: _____	N/C
SOFHRS	Small Curved Headrest Pad – 4" x 9"	\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"	\$155



57. HEADREST (CONTINUE)

SOFHRL	Large Curved Headrest Pad – 6" x 11"				\$155
SOMAH	Multi-axis Removable Mounting Hardware				\$220
SOHRBC	Bracket to Mount on the Contour Back Support – (Without Recline)				N/C
SOHRPR	Bracket to Mount on the Seat Frame – (For Reclining Backs)				N/C

58. LATERAL TRUNK SUPPORTS

▼ SELECT ONE OF THE LATERAL PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.

SOSLS3	Adjustable Laterals with Small Curved Pads – (Each) – 3" x 5"				SEATING BENEFIT	\$100
SOSLM4	Adjustable Laterals with Medium Curved Pads – (Each) – 4" x 6"				SEATING BENEFIT	\$100
SOSLL5	Adjustable Laterals with Large Curved Pads – (Each) – 5" x 7"				SEATING BENEFIT	\$100
SOSLHM	Adjustable Swing-Away Hardware Mounts – (Each)				SEATING BENEFIT	\$220
	▲ Waterfall arm pads (SOWAPx & SOWGELx) may interfere with flipping back the armrests when combined with Lateral Trunk supports					

59. HIP & THIGH SUPPORTS

▼ SELECT ONE OF THE HIP & THIGH SUPPORT PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.

SOHTSS10	Adjustable Hip & Thigh Support Pads – (Each) – Small 10" x 4"				SEATING BENEFIT	\$130
SOHTSM12	Adjustable Hip & Thigh Support Pads – (Each) – Medium 12" x 5"				SEATING BENEFIT	\$130
SOHTSL14	Adjustable Hip & Thigh Support Pads – (Each) – Large 14" x 6"				SEATING BENEFIT	\$130
SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware – (Each)				SEATING BENEFIT	\$136

60. ACCESSORIES

SOCH	Cup Holder – (Self-Leveling)	Install on:	Left Side	Right Side	\$115
SOCPH	Ram Mounts X-Grip Cell Phone Holder	Install on:	Left Side	Right Side	\$190
SOTHS	Ram Mounts X-Grip Small Tablet Holder	Install on:	Left Side	Right Side	\$205
SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	Install on:	Left Side	Right Side	\$360
	▲ Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.				
	▲ With multiple accessories installed on the same side, a custom upgrade will be required.				

61. REFLECTOR KIT

REFLECT	Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner's Manual)				N/C
----------------	---	--	--	--	-----

62. BACK PACK AND ARMREST SIDE POUCHES

SOBP	Back Pack (Comes with durable straps that can be hooked around the back canes)				CLIENT COST	\$90
SOAP	Armrest Side Pouch (Comes with durable straps that can be hooked around the armrest tubes)	Order:	Left Side	Right Side	CLIENT COST	\$53
	▲ If only choosing one (1) armrest pouch, please show what side it should be installed					

63. SWITCHES AND OPTIONS

ASCPWIFI	PGDT - Wireless programmer, WIFI dongle				NOT COVERED BY AADL	\$844
ASCPCP	PGDT R-NET Dongle and PC programmer kit				NOT COVERED BY AADL	\$600
SSBS	Buddy Button Switch					\$230
SSES	Egg Switch-Black					\$116
SSML	Ultra Light Switch					\$295
ESS-RED	Emergency STOP Switch - Red					\$314



SPECIALTY DRIVER CONTROL OPTIONS

► Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)

64. THERAFIN SIP-N-PUFF KIT (CODE W100)

SCSNPK	Therafin Sip-N-Puff Kit				\$1,975
SCSNPBTK	Therafin Sip-N-Puff Breath Tube Kit				N/C
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRMSN	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	SWL (Swings on Left)	SWR (Swings on Right)		\$300
	▲ Required with the Therafin Sip-N-Puff kit.				
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Therafin Sip-N-Puff kit.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Therafin Sip-N-Puff kit.				
► CHOOSE ONE (1) HEADREST PAD					
SOFHRS	Small Curved Headrest Pad – 4" x 9"				\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"				\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"				\$155

65. PROXIMITY HEAD ARRAY (CODE W100)

CUSTOM	Adult – 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse				RFQ
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.				
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.				
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ Required with the Proximity Head Array.				
SOMAH	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ Required with the Proximity Head Array.				

66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL (CODE W100)

CUSTOM	Light touch mini Proportional Joystick				RFQ
SCSRMMP	Stealth Rock-n-Lock Swingaway Arm for Chin MP	SWL (Swings on Left)	SWR (Swings on Right)		\$870
	▲ Required with the Light touch mini Proportional Joystick.				
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510



66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL (CONTINUED) (CODE W100)

SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section			
JCC	Chin Cup			\$55
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Light touch mini Proportional Joystick.			
▶ CHOOSE ONE (1) HEADREST PAD				
SOFHRS	Small Curved Headrest Pad – 4" x 9"			\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"			\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"			\$155

67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT (CODE W100)

SCPGMJS	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$2,369
SCGMJSM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for JC200	MLH (Mount on Left)	MRH (Mount on Right)	\$595
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
CUSTOM	Light touch mini Proportional Joystick			RFQ
SCGMPPM	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount for Light Touch Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$825
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$850
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT (CODE W100)

SCPGMJS-M	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$2,369
SCARMJSM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
CUSTOM	Light touch mini Proportional Joystick			RFQ
SCARMPM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$375
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$850
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

69. PROPORTIONAL HEAD CONTROL PACKAGE (CODE W100)

SCPHC	Proportional Head Control Package			\$4,495
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section			
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Proportional RIM Head Control.			



CUSTOM SECTION

- ▶ Lead times will increase depending on the custom order pipeline at the time the purchase order is received.
- ▶ Please call for Quote and specifications for these options.

70. O² HOLDER

SUPOXY	Custom Non-Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$295
	▲ <i>Not available with a power recline</i>			
SUPOXYA	Custom Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$1,195
	▲ <i>Required with a power recline</i>			
	D-Size	Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve		
	E-Size	Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve		

71. VENT TRAY (PLEASE COMPLETE IN FULL)

▼ **AMYLOR DOES NOT SUPPLY THE BATTERY, CHARGER OR THE WIRING FOR THE VENT. THESE ITEMS ARE NORMALLY PROVIDED BY THE VENT MANUFACTURER.**

SOVENT	Custom Vent Tray			
	Articulating Vent Box (Required with Power Tilt and/or Recline)			\$1,360
	Fixed Vent Box			\$400
	SPECIFY MODEL			
	Trilogy 100/200			
	LTV			
	LTV with Sprintpack			
	Other, specify: _____			
SOVTBE	Extended Battery Compartment			\$495
	▲ <i>The Vent Tray option is not available on narrow options N2424 & N2221 (FS2424B & FS2221B)</i>			
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries			
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries			
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries			
	The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries			
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries			

72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)

▼ **PLEASE SPECIFY DETAILS AND SPECIFICATIONS FOR THESE OPTIONS.**

CUSTOM	USB Charger Port ◀ Specify the location for the installation: _____	CLIENT COST	RFQ
CUSTOM	Accessory Hooks located at the rear of the seat structure (2 hooks)	CLIENT COST	RFQ
CUSTOM	Heavy-Duty Wheelchair Cover, Large	CLIENT COST	RFQ
CUSTOM	Transfer Handles – (Pair) Select a height range: 5" to 8" 6" to 9"	CLIENT COST	RFQ
CUSTOM	CUSTOM Lithium ion - 54 Amps (represents one 24V Battery pack, equivalent to two Group 22 Batteries)	CLIENT COST	RFQ
CUSTOM	CUSTOM 8A Offboard Battery Charger for Lithium ion batteries	CLIENT COST	N/C
	▲ <i>This option is required when choosing custom Lithium ion batteries</i>		
CUSTOM	CUSTOM Lithium ion Battery capacity LCD display gauge	CLIENT COST	N/C
	▲ <i>This option is required when choosing custom Lithium ion batteries</i>		



72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE) (CONTINUED)

CUSTOM	Folding Backrest for Additional Storage possibilities and Travel accessibility				CLIENT COST	RFQ
CUSTOM	Power Anterior Tilt				CLIENT COST	RFQ
CUSTOM	Power Lateral Tilt 15°	Raise the Left Side only	Raise the Right Side only	Raise the Left & Right Side	CLIENT COST	RFQ
CUSTOM	Power Latitude System (Features 14" of forward travel and a 5" seat to floor height from the pan to the ground)				CLIENT COST	RFQ
	▲ This option is only available with MWD (200 Lb max. user weight capacity)					
CUSTOM	Power Elevating Stretcher – (Please call CSR for availability with configuration and user weight capacity)				CLIENT COST	RFQ
	▲ This option is not available with the Straight Front Caster Arms on HWD					
CUSTOM	Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				CLIENT COST	RFQ
	▲ This option is not available with the Straight Front Caster Arms on HWD					
CUSTOM	Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				CLIENT COST	RFQ
	▲ This option is not available with the Straight Front Caster Arms on HWD					
CUSTOM	Rearview mirror	Left Side	Right Side		CLIENT COST	RFQ
CUSTOM	Rearview Camera with 5" Monitor				CLIENT COST	RFQ
CUSTOM	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan				CLIENT COST	RFQ
CUSTOM	Endomorphic Back Canes (Pair)	1" per side towards the inside = 2"	2" per side towards the inside = 4"		CLIENT COST	RFQ
CUSTOM	Transfer Handles – (Pair)	Select a height range:	5" to 8"	6" to 9"	CLIENT COST	RFQ
CUSTOM	Foot Box	Specify	Length: _____	Width: _____	Please confirm with CSR on sizing possibilities	CLIENT COST RFQ
CUSTOM #1	Other custom options. Specify in notes section below				CLIENT COST	RFQ
CUSTOM #2	Other custom options. Specify in notes section below				CLIENT COST	RFQ

73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS