



1. ACCOUNT INFORMATION

P.O. # or Quote: _____
 Dealer Contact Name: _____
 Email Address: _____
 Phone Number: _____ Fax Number: _____
 Bill to (Name): _____
 Bill to (Location): _____
 Ship to (Name): _____
 Address: _____
 City: _____
 Province: _____ Postal Code: _____
 User Name or Number: _____
 User Weight: _____
 User Height: _____
 Seat Cushion Thickness: _____
 Knee to Heel Length: Left: _____ " Right: _____ "



THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:

- A user name or number, the client's measurements and weight - Amylor is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amylor summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

***IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES**

2. STANDARD OPTIONS

Alltrack M3 Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A onboard charger, 6" front casters, and standard auto style buckle seat positioning strap, 6.5mph motors, **Weight Capacity = 300 lb**

3. ALLTRACK M3 POWER WHEELCHAIR

ATM3	ALLTRACK M3 Power Wheelchair	\$8,020
-------------	------------------------------	---------

3.5 OWNER'S MANUAL

SOOMATENG	English / Anglais	N/C	SOOMATFR	Français / French	N/C
------------------	-------------------	-----	-----------------	-------------------	-----

4. USER WEIGHT PACKAGES

WCA25	250 lb Maximum user weight capacity	STD
WCA3	300 lb Maximum user weight capacity	STD

5. STANDARD COLORS

SCMOBM	Onyx Black Matte	N/C	SCMCAB	Cabernet	N/C
SCMCTB	Cobalt Blue	N/C	SCMBGM	Bubble Gum Metallic	N/C
SCMIOB	Indigo Blue	N/C	SCMHAB	Habanero Metallic	N/C
SCMCRE	Candy Red	N/C	SCMLIM	Sublime	N/C

OPTIONAL COLORS

FCHDD-M	Desert Digital Camo HydroGraphics - (Military Tan base)	\$795	FCHWL-M	Woodland Camo HydroGraphics - (Military Green base)	\$795
FCHCF-M	Carbon fiber HydroGraphics - (Matte black base)	\$795			

5.5 WHEEL INSERT COLORS

SCCTB-WI	Cobalt Blue	N/C	SCHAB-WI	Habanero Metallic	N/C
SCIOB-WI	Indigo Blue	N/C	SCLIM-WI	Sublime	N/C
SCCRE-WI	Candy Red	N/C	SCBGM-WI	Bubble Gum Metallic	N/C
SCCAB-WI	Cabernet	N/C			

- ▲ Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)
- ▲ Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)
- ▲ Wheel inserts are not available in Onyx Black Matte



7. TRANSPORT OPTION

TTD4	4-Point wheelchair transport brackets ▲ 4-Point wheelchair transport brackets have not yet received WC-19 crash-testing approval.	STD
WC19-ISO	WC-19 and ISO 7176-19 Occupied Transit Package ▲ To confirm if your overall wheelchair configuration is WC-19 and ISO 7176-19 compliant, please refer to this website: https://docs.google.com/spreadsheets/d/1K1cNIEA5OQJCS5z3OD1adWElaFLN1XO3ZcTc6DL_kaE/edit#gid=1917752765 ▲ Other restrictions may apply, please contact CSR for details. ▲ This option is only available on new power wheelchairs as an "on-chair" option. It is not available to install on a delivered wheelchair ▲ Not compatible with the narrow base configurations (FS2424B and FS2221B)	\$595

8. MOTOR PACKAGE

MP65B	6.5MPH 4-Pole Motor Package / with Brake Release ▲ Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)	STD
MP65-OPT	6.5MPH 4-Pole Motor MP65-OPT Package / with COG Release ▲ Not compatible with the narrow base configuration (FS2221B)	\$395

9. ASSISTIVE DRIVING PACKAGE

SMARTT	Smart-Track Assistive Driving Package ▲ Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control. ▲ Requires the R-net PM 120 Amp Expandable Controller (ECR120). ▲ This option requires the MP65-OPT motor package. ▲ Not compatible with the narrow base configuration (FS2221B) ▲ Recommended for non-proportional Drive Controls	\$1,420
---------------	---	---------

10. BATTERY CHARGER

BCOB	8A Offboard Battery Charger	N/C
-------------	-----------------------------	-----

11. BATTERY TRAY AND BASE WIDTH

FS2221B	N2221 - Battery Tray for Group 22 Batteries with a Narrow 21" base width (Total base width = 21") ▶ NOT AVAILABLE WITH THE FOLLOWING OPTIONS: ▲ WC-19 and ISO 7176-19 Occupied Transit Package (WC19-ISO) ▲ 6.5MPH 4-Pole Motor Package / with COG Release option (MP65-OPT) ▲ Smart-Track option (SMARTT) ▲ 8" front casters (FC8PTM and FC8PIM) ▲ 14" x 4" Extra Wide Cast Wheels, Pneumatic Tires (DW144PT and DW144PI) ▲ Not compatible with 23" and 24" seat widths selections (SOWN23 and SOWN24)	N/C
FS22B	Battery Tray for Group 22NF Batteries – (Total base width = 23 3/8")	STD
FS2424B	N2424 - Battery Tray for Group 24 Batteries with a Narrow 24" base width (Total base width = 24") ▲ The starting angle is reduced from 90° to 85° for some seat configurations when used with the narrow base configuration (FS2424B).	N/C
FS24B	Battery Tray for Group 24 Batteries – (Total base width = 25 3/8")	N/C

12. BATTERIES - REQUIRES A SELECTION

BAT22NFGEL	Premium CG Power - GEL Battery – Group 22NF – 50 Amps (Each)	\$335
BAT24GEL	Premium CG Power - GEL Battery – Group 24 – 70 Amps (Each)	\$465
BATO	Omission of the batteries	N/C

14. DRIVE WHEEL

DW14PT	14" Cast Wheel Pneumatic Tires	STD
DW14PI	14" Cast Wheel Pneumatic with Foam Filled Tires	N/C
DW144PT	14" x 4" Extra Wide Cast Wheels, Pneumatic Tires	\$200
DW144PI	14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires ▲ Increases the overall width of the base by 2 3/4" ▲ Not compatible with the narrow base configuration (FS2221B)	\$595
DW12PT	12 1/2" Cast Wheel Pneumatic Tires	N/C
DW12PI	12 1/2" Cast Wheel Pneumatic with Foam Filled Tires ▲ 12 1/2" Drive wheels reduce overall seat to floor height by 3/4", lowers ground clearance and speed, but increases motor torque.	N/C



15. FRONT CASTERS

FC6S	6" Solid Casters – (Not available with the 14" x 4" Extra Wide Cast Wheel options DW144PT or DW144PI)	STD
FC8PTM	8" Pneumatic Casters	N/C
FC8PIM	8" Pneumatic Casters with Foam Filled Tires ▲ The 8" casters are NOT compatible with the 12 ½" Drive wheel configuration (DW12PT and DW12PI) ▲ Not compatible with the narrow base configuration (FS2221B)	N/C

16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± ¼")

▶ OPTIONS FOR A REHAB SEAT ONLY - NO POWER SEATING

STFH155	15 ½" – (Possible Seat angle of 0° only) ▲ Requires the 12 ½" Drive wheels (DW12PT and DW12PI)		
STFH1625	16 ¼" – (Possible Seat angle of 0° only)	STFH1925	19 ¼" – (Possible Seat angle of 0°, 3° and 6°)
STFH1725	17 ¼" – (Possible Seat angle of 0° and 3°)	STFH2025	20 ¼" – (Possible Seat angle of 0°, 3°, 6°, 9°)
STFH1825	18 ¼" – (Possible Seat angle of 0° and 3°)		

▶ OPTIONS FOR A POWER TILT, POWER TILT / RECLINE OR ELEVATE ONLY MODULE

STFH155	15 ½" – (Possible Seat angle of 0° only) ▲ Requires the 12 ½" Drive wheels (DW12PT and DW12PI)		
STFH1625	16 ¼" – (Possible Seat angle of 0° only)	STFH1825	18 ¼" – (Possible Seat angle of 0° only)
STFH1725	17 ¼" – (Possible Seat angle of 0° only)	STFH1875	18 ¾" – (Possible Seat angle of 0° only)
STFH1775	17 ¾" – (Possible Seat angle of 0° only)		

▶ OPTIONS FOR A COMBO - POWER TILT OR POWER TILT / RECLINE WITH ELEVATE MODULE

STFH175	17 ½" – (Possible Seat angle of 0° only) ▲ Requires the 12 ½" Drive wheels (DW12PT and DW12PI)		
STFH1825	18 ¼" – (Possible Seat angle of 0° only)		

17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)

SOMHSS	Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges. ▲ Required at: _____" (± ¼" Range)	RFQ
---------------	---	-----

18. PRE-SET SEAT FRAME ANGLE

SFE-0	0° Pre-set seat frame seat angle – (Must be selected when using power options)	STD
SFE-3	3° Pre-set seat frame seat angle – (Only available with a seat configuration without power options)	N/C
SFE-6	6° Pre-set seat frame seat angle – (Only available with a seat configuration without power options)	N/C
SFE-9	9° Pre-set seat frame seat angle – (Only available with a seat configuration without power options)	N/C

19. LIGHTING PACKAGE

EOLP	LED Light Package – (Includes front and rear lights with integrated flashers) ▲ Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL , ECELJL , ECRJCL , ECRJC35B or EDIB)	\$1,425
-------------	--	---------

20. NON-EXPANDABLE ELECTRONICS - PG DRIVES VR2 CONTROLLER AND JOYSTICKS

▶ VR2 Electronics can only accommodate specific power seat functions and combinations such as Tilt / Elevate / Recline / Power Center Mount / Combined Power Legs / Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs / Tilt + Elevate (Custom)		
ECV902A	VR2 90 Amp Drive Controller – (Requires the selection of VR2 joystick below.)	STD
ECVJ4K	VR2 Standard 4-Key Joystick	STD
ECVJ6K	VR2 6-Key Joystick – (Allows for single or select dual power options to operate through this joystick)	N/C
ECVJ9KL	VR2 9-Key Joystick with Lighting Operation – (Allows for single or select dual power options to operate through this joystick)	\$150

21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET POWER MODULES

ECEL90	R-net EL 90 Amp Expandable Controller	\$175
ECR120	R-net PM 120 Amp Expandable Controller	\$275

22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS

▶ Selection of one joystick option below is required with the expandable controller ECEL90 or ECR120		
ECO	Omit Joystick – (For use with an Omni 2 Display)	N/C
ECELJ	R-net - (LED) – Joystick ▲ This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY).	STD
ECRJCL	R-net Joystick with Color screen	\$155
▶ JOYSTICKS WITH LIGHTING OPERATION		
ECELJL	R-net – (LED) – Joystick with Lighting Operation ▲ This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY).	\$125
ECRJCL	R-net Joystick with Color screen and Lighting Operation	\$180
▶ ADVANCED JOYSTICK		
ECRJC35B	R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,105

MKT.F70 | OF | PD-A20602 | R06



23. JOYSTICK MOUNT

JSMF	Fixed Joystick Mount	N/C
JMSA	Swing Away Joystick Mount - Height adjustable, Angle adjustable and Multi-Axis Set-up with standard installation to swing-away sideways Set-up at 35° to swing-away sideways and downwards	\$235
JSMCM	Midline Mounted Swing Away Joystick Mount	\$275

24. JOYSTICK MOUNT SUPPORT

JSMFH	Joystick mount support Bracket	STD
--------------	--------------------------------	-----

25. JOYSTICK HAND MOUNT

HML	Install Joystick on Left hand side	HMR	Install Joystick on Right hand side
------------	------------------------------------	------------	-------------------------------------

26. JOYSTICK HANDLE

JSH	Standard Joystick Gimble	STD	JBAT	Stick Handle 3" in Length	\$130
JBH	Ball Handle	\$85	JFBH	Foam Ball Handle	\$70
JTH	T-Handle	\$155	JMUSH	Bodypoint Mushroom Handle	\$210
JGP	Bodypoint 4" U-shaped Handle	\$210			

27. ATTENDANT DRIVE CONTROL OPTION

EADC	Attendant Hand Control for R-net EL & PM Controls	\$1,150
-------------	---	---------

28. ATTENDANT JOYSTICK HAND MOUNT

EADCL	Install Attendant Joystick on Rear Left hand side	EADCR	Install Attendant Joystick on Rear Right hand side
--------------	---	--------------	--

29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)

EDIB	Omni 2 – (Enhanced Display & Digital Interface Box with Bluetooth and Infrared) ▲ Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120) ▲ The Omni 2 display mount must be selected from the option below.	\$1,915
-------------	--	---------

30. OMNI 2 ENHANCED DISPLAY MOUNT

EDSAM	User Swing Away Mount	N/C
	EDML – (Install Display on LH side) EDMR – (Install Display on RH side)	

31. ENVIRONMENTAL CONTROL MODULE

ECU	Environment control module – (Output Module) ▲ Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).	\$995
------------	--	-------

33. AMYSEAT REHAB SEAT / POWER SEATING SYSTEMS AND MODULES

SAABA	Amyseat Rehab Seat – (Seat Frame only)	\$880
SAPTX	CG Power Tilt ▲ This option includes the Amyseat Rehab Seat – (SAABA)	\$6,195
SAPRX	Power Recline with Mechanical Shear Reduction (MSR) ▲ This option includes the Amyseat Rehab Seat – (SAABA)	\$7,845
SAPTRX	CG Power Tilt and Power Recline with Mechanical Shear Reduction (MSR) ▲ This option includes the Amyseat Rehab Seat – (SAABA)	\$11,495
SAPESX	ALTA 11" Power Seat Elevation ▲ This option includes the Amyseat Rehab Seat – (SAABA)	\$4,495
ADDPESX	Module – Add ALTA 11" Power Seat Elevation to either SAPTX, SAPRX or SAPTRX	\$2,495
ADDMR	Module – Manual Recline with Mechanical Shear Reduction – (MSR) ▲ The maximum user weight capacity with this option is 250 lb.	\$1,125

34. ELECTRONICS - SINGLE SEAT FUNCTION CONTROL TYPE

► SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (FOR BOTH EXPANDABLE AND NON-EXPANDABLE SYSTEMS)

SOLO	Single function controller (Includes an external toggle switch)	STD
▼ OPTIONAL SWITCHES		
SWTOGGLE	Single Toggle Switch	Attendant Activation \$100
SWPUSH1	Single Push Button Switch – (Single Mode: SAC)	User Activation Attendant Activation \$90
SWPUSH2	Dual Push Button Switch – (Dual Mode: SAC)	User Activation Attendant Activation \$180
SSBS	Buddy Button Switch Black Red Green Yellow	User Activation Attendant Activation \$230
SSES	Egg Switch-Black	User Activation Attendant Activation \$116
SSML	Ultra Light Switch	User Activation Attendant Activation \$295



► SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL

EPSV	Single Function Through VR2 Non-Expandable Drive Control ▲ Allows for activation of a single power option through VR2 controller ▲ Does not allow for control of power option through drive control and separate switches ▲ Control of power elevating seat with tilt is available as a custom. Check appropriate box in custom section 72 in addition to EPSV. ▲ VR2 Electronics can only accommodate specific power seat functions such as Tilt / Elevate / Recline / Power Center Mount / Combined Power Legs / Tilt + Elevate (Custom) ▲ Choose EPSV2 to operate select power combinations	N/C
EPSEL	Single Function Through Drive Control with R-net 90-amp controller ▲ Allows for activation of single power option through R-net 90-amp controller ▲ Compatible with all R-net joysticks ▲ Does not allow for control of power option through drive control and separate switches ▲ Can not operate more than one power option ▲ If selecting memory seating, choose SAC below ▲ If selecting 120-amp R-net controller, choose SAC below	N/C
SAC	Single Function Through Drive Control with R-net 90-amp or 120-amp controller ▲ Allows for activation of single power option through R-net drive controls ▲ Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls ▲ Capable of operating tilt with power elevating seat ▲ Compatible with memory seating	N/C

► SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE)

SACMEMORY	Memory Seat Program ▲ This option is not compatible with any VR2 or LED joysticks (ECEJL and ECELJL)	\$1,495
------------------	---	---------

► SELECTION OF SWITCH BELOW GIVES ACCES TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH

SACTOGGLE	Single Toggle Switch	Attendant Activation	\$410	
SACPUSH1	Single Push Button Switch – (Single Mode: SAC)	User Activation	Attendant Activation	\$400
SACPUSH2	Dual Push Button Switch – (Dual Mode: SAC)	User Activation	Attendant Activation	\$490
SACSSBS	Buddy Button Switch: Black Red Green Yellow	User Activation	Attendant Activation	\$540
SACSSSE	Egg Switch-Black	User Activation	Attendant Activation	\$426
SACSSML	Ultra Light Switch	User Activation	Attendant Activation	\$605

35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE

► DUAL SEAT FUNCTIONS THROUGH DRIVE CONTROL

EPSV2	Dual Function Through Drive Control with VR2 Non-Expandable Electronics ▲ Allows for the activation of 2 power options through the VR2 Controller. Combination of TILT + RECLINE or ELEVATE + RECLINE are not possible with EPSV2 ▲ Does not allow for control of power options through drive control and separate switches ▲ VR2 Electronics can only accommodate specific power seat combinations such as Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs	N/C
--------------	--	-----

► MULTIPLE SEAT FUNCTIONS THROUGH R-NET

MAC	E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators ▲ Required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2221B)	\$2,295
------------	--	---------

► MAC OPTIONS

MACMEMORY	Memory Seat Program	\$175
------------------	---------------------	-------

► SELECTION OF SWITCH BELOW GIVES ACCES TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH

MACKEYPAD-P	Keypad w/ 5 Push Button Switches	User Activation	Attendant Activation	\$445
MACKEYPAD-T	Keypad w/ 5 Toggle Switches	User Activation	Attendant Activation	\$445
MACMULTI	4-Way Toggle Switch – (Up to 4 Activations) ▲ By selecting this switch, memory recall is only available through the Joystick.	User Activation	Attendant Activation	\$180
MACMULTIT	Triple Toggle Switch – (Up to 5 Activations) ▲ By selecting this switch, memory recall is only available through the Joystick.	User Activation	Attendant Activation	\$135
MACPUSH4	4-Way Push Button Switch ▲ By selecting this switch, memory recall is only available through the Joystick.	User Activation	Attendant Activation	\$225
MACPUSH5	5-Way Push Button Switch ▲ By selecting this switch, memory recall is only available through the Joystick.	User Activation	Attendant Activation	\$240

36. SWITCH ACTIVATION HAND MOUNT

Switch Installation:	TSML – (Install Switch on LH side)	TSMR – (Install Switch on RH side)	N/C
-----------------------------	---	---	-----



37. PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)

SOBAA	Select one:	86°	90°	94°	98°	102°	106°	110°	114°	118°	122°
▲ Back angles of 102° or greater are not compatible with the OMNI 2.											

38. STANDARD ADJUSTABLE SEAT WIDTH RANGE BETWEEN 15" AND 20"

SOWN15-20-19	– (Width adjustment from 15" to 20" with a maximum Depth of 19")	N/C			
SOWN15-20-22	– (Width adjustment from 15" to 20" with a maximum Depth of 22") – (The minimum seat depth with this seat option is 18" Deep)	N/C			
15" (SOWN15)	16" (SOWN16)	17" (SOWN17)	18" (SOWN18)	19" (SOWN19S)	20" (SOWN20S)

39. LARGE ADJUSTABLE SEAT WIDTH RANGE BETWEEN 19" AND 24"

SOWN19-24-22	– (Width adjustment from 19" to 24" and Depth range of 18" to 22") – (The minimum seat depth with this seat option is 18" Deep)	N/C			
19" (SOWN19L)	20" (SOWN20L)	21" (SOWN21)	22" (SOWN22)	23" (SOWN23)	24" (SOWN24)
▲ 23" and 24" Seat widths are not compatible with the N2221 - Narrow 21" base width option (FS2221B)					

40. ADJUSTABLE SEAT DEPTH RANGE BETWEEN 15" AND 22"

15" (SODN15)	16" (SODN16)	17" (SODN17)	18" (SODN18)	N/C
19" (SODN19)	20" (SODN20)	21" (SODN21)	22" (SODN22)	N/C
▲ With a 22" long Seat Pan SOWN15-20-22 and SOWN19-24-22, the minimum possible seat depth is 18".				

41. BACK CANE HEIGHT (THE BACK CANE HEIGHT MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE PUSH HANDLES)

▼ 22" is the minimum back height required with a power or manual recline.					
16" (SOBCH16)	18" (SOBCH18)	20" (SOBCH20)	22" (SOBCH22)	24" (SOBCH24)	26" (SOBCH26)

42. BACK SUPPORT

DSTR-WWHH	NEW Adjustable Tension Back Upholstery, WW x HH	N/C
▲ Only available with Dual Post armrests (SOADP)		
STDSB-WWHH	Contour back support up to 21" wide, with a maximum shell height of 22".	\$605
▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
STDSB-WWHH	Contour back support up to 21" wide for non-standard sizes with shell height ≥ 23".	\$755
▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
STDSB-WWHH	Contour back support for Heavy-Duty from 22" wide and higher, with any height.	\$1,075
▲ The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▶ MATRIX BACKS		
E2S-WWHH	MATRIX ELITE E2S BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	\$995
▲ Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▲ Weight capacity up to 300 lb		
E2D-WWHH	MATRIX DEEP CONTOUR ELITE E2D BACK – (14" to 20" Wide - Back height of 6" to 20", in 2" increments)	\$1,155
▲ Back shell height is available in 6" to 20", in 2" increments. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▲ This product is not compatible with Lateral Trunk Supports		
▲ Weight capacity up to 300 lb		
PBE-WWHH-TR	MATRIX PB ELITE TR BACK – (15" to 20" Wide - Back height of 14", 16", 18" and 20")	\$910
▲ Back shell height is available in 14", 16" 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▲ Weight capacity up to 300 lb		
PBE-WWHH-HD	MATRIX HEAVY DUTY ELITE BACK – (21" to 24" Wide - Back height of 16" and 20")	\$1,150
▲ Back shell height is available in 16" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▲ Weight capacity up to 500 lb		
EDB-WWHH-HD	MATRIX ELITE DEEP HD BACK – (21" to 24" Wide - Back height of 16", 18" and 20")	\$1,195
▲ Back shell height is available in 16", 18" or 20". The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back		
▲ This product is not compatible with Lateral Trunk Supports		
▲ Weight capacity up to 500 lb		



▶ **NXT BACKS**

NOACB-WWHH	NXT OPTIMA FH AC CONTOUR BACK – (15" to 21" Wide) ▲ Weight capacity up to 300 lb	\$690
NOACB-WWHH	NXT OPTIMA FH AC CONTOUR BACK – (22" to 28" Wide) ▲ Weight capacity up to 550 lb ▲ Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back	\$720
NOPCB-WWHH	NXT OPTIMA FH PC DEEP CONTOUR BACK – (15" to 21" Wide) ▲ Weight capacity up to 300 lb	\$790
NOPCB-WWHH	NXT OPTIMA FH PC DEEP CONTOUR BACK – (22" to 28" Wide) ▲ Weight capacity up to 550 lb ▲ This product is not compatible with Lateral Trunk Supports ▲ Back shell heights only come in 13", 16" or 19" high. The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back	\$820

▶ **OMISSION OF THE BACK SUPPORT**

SDSBO	Omission of the Contour Back Support ◀ Specify Model of Back Support: _____	N/C
--------------	---	-----

43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)

▼ 22" is the minimum back height required with a power or manual recline.

16" (SOBH16)	18" (SOBH18)	20" (SOBH20)	22" (SOBH22)	24" (SOBH24)	26" (SOBH26)
28" (SOBH28)					

44. SEAT CUSHIONS

▶ **CG AIR - THERAPEUTIC AIR CUSHION**

Single Valve 15" to 20" Wide Note: For custom Cushions (Please use the CG Air Order Form)					
	Cell Height:	2"	3"	4"	
CGA1VSCwwdd-h	CG Air Single Valve – (with Standard 2-way Stretch Cover)				\$600
CGA1V3Dwwdd-h	CG Air Single Valve – (with 3D Mesh Cover)				\$690
CGA1VICwwdd-h	CG Air Single Valve – (with Incontinence Cover) ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$690
Dual Valve 15" to 20" Wide Note: For custom Cushions (Please use the CG Air Order Form)					
	Cell Height:	2"	3"	4"	
CGA2VSCwwdd-h	CG Air Dual Valve – (with Standard 2-way Stretch Cover)				\$600
CGA2V3Dwwdd-h	CG Air Dual Valve – (with 3D Mesh Cover)				\$690
CGA2VICwwdd-h	CG Air Dual Valve – (with Incontinence Cover) ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$690
Single Valve 22" to 24" Wide Note: For custom Cushions (Please use the CG Air Order Form)					
	Cell Height:	2"	3"	4"	
CGA1VSCwwdd-h	CG Air Single Valve – (with Standard 2-way Stretch Cover)				\$890
CGA1V3Dwwdd-h	CG Air Single Valve – (with 3D Mesh Cover)				\$910
CGA1VICwwdd-h	CG Air Single Valve – (with Incontinence Cover) ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$910
Dual Valve 22" to 24" Wide Note: For custom Cushions (Please use the CG Air Order Form)					
	Cell Height:	2"	3"	4"	
CGA2VSCwwdd-h	CG Air Dual Valve – (with Standard 2-way Stretch Cover)				\$890
CGA2V3Dwwdd-h	CG Air Dual Valve – (with 3D Mesh Cover)				\$910
CGA2VICwwdd-h	CG Air Dual Valve – (with Incontinence Cover) ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$910
▶ CONTOUR SEAT CUSHION					
SCPSC-wwdd	Contour Seat Cushion up to 19" Wide				\$620
SCPSC-wwdd	Contour Seat Cushion for Non-Standard Sizes of 20" to 24" Wide				\$735
▶ MATRX CUSHION					
PSP-WWDD	MATRX PSP Seat Cushion – (11" to 20" Wide) ▲ The maximum user weight capacity is 300 lb ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$645
PSP-WWDD	MATRX PSP Seat Cushion – (22" to 26" Wide) ▲ The maximum user weight capacity is 550 lb ▲ Cushion sizes only come in even numbers for width and depth measurements.				\$995



▶ **NXT BIO-FIT CUSHION**

NXTBIOFIT -WWDD	NXT BIO FIT Cushion – (15" to 21" Wide) ▲ The maximum user weight capacity is 300 lb	\$595
NXTBIOFIT -WWDD	NXT BIO FIT Cushion – (22" to 26" Wide) ▲ The maximum user weight capacity is 400 lb ▲ Cushion sizes only come in even numbers for width and depth measurements.	\$625

45. CANTILEVER – HEIGHT ADJUSTABLE ARMRESTS – OPENS WITH A 5° ANGLE (CHOOSE 2)

SOAFBDLH	Flip-Back Cantilever – Desk Length Left Side – (Height Adjustable from 7" to 15" with 1" increments) –	STD Length 12"	N/C		
SOAFBDRH	Flip-Back Cantilever – Desk Length Right Side – (Height Adjustable from 7" to 15" with 1" increments) –	STD Length 12"	N/C		
SOAFBFLH	Flip-Back Cantilever – Full Length Left Side – (Height Adjustable from 7" to 15" with 1" increments) –	STD Length 16"	STD		
SOAFBFRH	Flip-Back Cantilever – Full Length Right Side – (Height Adjustable from 7" to 15" with 1" increments) –	STD Length 16"	STD		
CUST-112250	Custom Armrest Length – Select Length for the Left side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C
CUST-112250	Custom Armrest Length – Select Length for the Right side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C
CUST-SOAFB	Select this code when using Cantilever armrests with Reclining backs ▲ The maximum user weight capacity is 300 lb				N/C

46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)

SOADPDH	Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPRH	Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPLH	Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168
SOADPRH	Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 1/4" to 15 1/4" with 1/2" increments)	\$168

47. RECLINING FLIP-UP, HEIGHT ADJUSTABLE ARMRESTS (CHOOSE 2)

SOANIDLH	Reclining armrest – Desk Length Left Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments) –	STD Length 12"	N/C		
SOANIDRH	Reclining armrest – Desk Length Right Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments) –	STD Length 12"	N/C		
SOANIFLH	Reclining armrest – Full Length Left Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments) –	STD Length 16"	N/C		
SOANIFRH	Reclining armrest – Full Length Right Side – (Height Adjustable from 10 1/2" to 14" with 1/2" increments) –	STD Length 16"	N/C		
CUST-48247	Custom Armrest Length – Select Length for the Left side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C
CUST-48247	Custom Armrest Length – Select Length for the Right side:	Custom Length 18"	Custom Length 20"	Custom Length 22"	N/C

48. ARMREST HEIGHT ADJUSTMENT AT FACTORY

- ▶ Set the Armrest Height at: _____ "
- ▲ The Armrest Height is taken from the seat pan to the top of the armpad. (Not required for Dual Post.)

49. ARMPADS

▶ **FOAM TYPE SELECTION**

SOWAPD	Desk Length Waterfall – (8 1/2" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	STD
	▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports.				
SOWAPF	Full Length Waterfall – (13" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	STD
	▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports.				

▶ **HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)**

SOGELD	Desk Length Hybrid Gel and Foam Flat armpad – (9 3/4" x 2")	Select side:	Left Hand Side	Right Hand Side	\$40
SOGELF	Full Length Hybrid Gel and Foam Flat armpad – (13 1/2" x 2")	Select side:	Left Hand Side	Right Hand Side	\$40
SOWGELD	Desk Length Hybrid Gel & Foam Waterfall armpad – (8 1/2" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports.				
SOWGELF	Full Length Hybrid Gel and Foam Waterfall armpad – (13" x 2 3/4")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ The Waterfall armpads may interfere with flipping back the armrests when combined with Lateral Trunk supports.				

▶ **SPECIALTY TYPE SELECTION**

SOCAR	Arm Trough with Hand Support (18" x 4 1/2")	Select side:	Left Hand Side	Right Hand Side	\$170
	▲ Arm Trough Armpads may interfere with Armrests Flipping back.				



50. FOOTRESTS

SOSAAB	NEW Signature Swing-Away Footrests	Set at:	60°	70°	STD
		Select Footrest Range:	SOFLA12 (12 ½" - 17 ½")	SOFLA14 (14 ½" - 19 ½")	
SOAAF	NEW Signature Heavy Duty Swing-Away Footrests	Set at:	60°	70°	\$175
		Select Footrest Range:	SOFLHD12 (12 ½" - 17 ½")	SOFLHD14 (14 ½" - 19 ½")	
ADDML	NEW Signature Manual Elevating & Articulating Legrests - Pair				\$595
		Select Footrest Range:	SOFLME12 (12 ½" - 17 ½")	SOFLME14 (14 ½" - 19 ½")	
	▲ Does not include calf pads. Please select calf pads in section 51.				
ADDPL	NEW Signature Power Elevating & Articulating Legrests				\$ 2,495
		Select Function:	Individual	Combined	
		Select Footrest Range:	SOFLPW12 (12 ½" - 17 ½")	SOFLPW14 (14 ½" - 19 ½")	
	▲ This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.				

51. FOOTPLATES & ACCESSORIES

SOSFAF	Flip-up Footplates				STD
SOAFP	Angle Adjustable & Flip-up Footplates – (Each)				\$65
SOAHL	Heel Loops – (Each)				\$20
SOACPF	Calf Pads for Adult Heavy-Duty Swing-Away footrests (SOAAF) with adjustable hardware – (Each)				\$185
SOAMA	Swing-away Amputee Support	Select :	Left leg	Right leg	\$455

52. FOOTREST LENGTH ADJUSTMENT AT FACTORY

- ▶ Set the Footrest Length at: _____ "
- ▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

53. CENTER MOUNTS

SOAAF90	90° Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°)				200\$
		Select Foot Platform Range:	SOFBA12 (12 ½" - 16 ½")	SOFBA16 (16 ½" - 19 ½")	
	▲ When selecting SOFBA16 (16 ½" - 19 ½") range, the lowest possible seat to floor height is 18.75"				
	▲ Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).				
	▲ With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.				
ADDPF	90° Power Center Mount Elevating & Articulating Foot Platform				\$ 2,495
		Select Foot Platform Range:	SOFBPW12 (10 ½"- 16")	SOFBPW16 (16"- 20")	
	▲ When selecting SOFBPW16 (16" - 20") range, the lowest possible seat to floor height is 18.75"				
	▲ This option includes calf pads. The power center mount counts as one power function.				
	▲ Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).				
	▲ For leg lengths between 10 ½" and 12 ½", there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.				
	▲ The starting angle is reduced from 90° to 85° for some seat configurations when used with the narrow base configuration (FS2424B)				
	▲ The E-BX (MAC) option is required with the combined selection of the 90° Power Center Mount Elevating & Articulating Foot Platform (ADDPF) and the narrow base configuration (FS2424B)				

54. FOOT PLATFORMS & ACCESSORIES

SOFUFP	Flip-up one piece Foot Platform				\$140
SOAFPFB	Angle Adjustable & Flip-up Split Footplates for center mount – (Each)				\$65
SOACPCM	Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each)				\$93

55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY

- ▶ Set the Foot Platform Height at: _____ "
- ▲ The Foot Platform Height is taken from the top of the seat pan to the top of the foot platform or split foot plates.

56. SEAT POSITIONING (NOT COMPATIBLE WITH THE WC19-ISO OPTION)

SOAUTO	55" Seat Positioning Strap				STD
SOBPHP2C	Bodypoint padded Hip Belt (Two Point – Center Pull)				\$170
SOPHB	Padded Hip Belt, 2-PT Magnetic Center Pull				\$194
SOBPCS	Bodypoint padded Chest Strap				\$250
SOPCS	Padded Chest Strap, 2-PT Magnetic Center Pull				\$262



57. HEADREST

▼ Select one of the headrest options with the mounting hardware and specify where the headrest will be installed.		
SONUAHR	Using headrest mounting hardware other than SOMAH located bellow. ◀ Specify: _____	N/C
SOFHRS	Small Curved Headrest Pad – 4" x 9"	\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"	\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"	\$155
SOMAH	Multi-axis Removable Mounting Hardware	\$220
SOHRBC	Bracket to Mount on the Contour Back Support (Without Recline)	N/C
SOHRPR	Bracket to Mount on the Seat Frame (For Reclining Backs)	N/C

58. LATERAL TRUNK SUPPORTS

▼ Select one of the lateral pads and the adjustable & removable hardware.		
SOSLS3	Adjustable Laterals with Small Curved Pads – (Each) – 3" x 5"	\$100
SOSLM4	Adjustable Laterals with Medium Curved Pads – (Each) – 4" x 6"	\$100
SOSLL5	Adjustable Laterals with Large Curved Pads – (Each) – 5" x 7"	\$100
SOSLHM	Adjustable Swing-Away Hardware Mounts – (Each)	\$220

59. HIP & THIGH SUPPORTS

▼ Select one of the hip & thigh support pads and the adjustable & removable hardware.		
SOHTSS10	Adjustable Hip & Thigh Support Pads – (Each) – Small 10" x 4"	\$130
SOHTSM12	Adjustable Hip & Thigh Support Pads – (Each) – Medium 12" x 5"	\$130
SOHTSL14	Adjustable Hip & Thigh Support Pads – (Each) – Large 14" x 6"	\$130
SOHTSARH	Adjustable and Removable Hip & Thigh Support Hardware – (Each)	\$136

60. ACCESSORIES

SOCH	Cup Holder – (Self-Leveling)	Install on:	Left Side	Right Side	\$115
SOPHP	Ram Mounts X-Grip Cell Phone Holder	Install on:	Left Side	Right Side	\$190
SOTHS	Ram Mounts X-Grip Small Tablet Holder	Install on:	Left Side	Right Side	\$205
SOTHL	Ram Mounts X-Grip Large Tablet Holder with Swivel	Install on:	Left Side	Right Side	\$360
▲ Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.					
▲ With multiple accessories installed on the same side, a custom upgrade will be required.					

61. REFLECTOR KIT

REFLECT	Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner's Manual)	N/C
----------------	---	-----

62. BACK PACK AND ARMREST SIDE POUCHES

SOBP	Back Pack – (Comes with durable straps that can be hooked around the back canes)	\$90
SOAP	Armrest Side Pouch – (Comes with durable straps that can be hooked around the armrest tubes)	Order: Left Side Right Side \$53
▲ If only choosing one (1) armrest pouch, please show what side it should be installed		

63. SWITCHES AND OPTIONS

ASCPWIFI	^{NEW} PGDT - WIFI Programmer Interface Kit	\$844
ASCPCP	PGDT R-net Dongle and PC programmer kit	\$600
ASCPVR2	PGDT VR2 with PPI Dongle and PC programmer kit	\$1,450
SSBS	Buddy Button Switch	\$230
SSES	Egg Switch-Black	\$116
SSML	Ultra Light Switch	\$295
ESS-RED	Emergency STOP Switch - Red	\$314



SPECIALTY DRIVER CONTROL OPTIONS

▼ Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)

64. THERAFIN SIP-N-PUFF KIT

SCSNPK	Therafin Sip-N-Puff Kit			\$1,975
SCSNPBTk	Therafin Sip-N-Puff Breath Tube Kit			N/C
	▲ Required with the Therafin Sip-N-Puff kit.			
SCSRMSNP	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	SWL (Swings on Left)	SWR (Swings on Right)	\$300
	▲ Required with the Therafin Sip-N-Puff kit.			
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.			
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation	\$314
	▲ Required with the Therafin Sip-N-Puff kit.			
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Therafin Sip-N-Puff kit.			
▶ CHOOSE ONE (1) HEADREST PAD				
SOFHRS	Small Curved Headrest Pad – 4" x 9"			\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"			\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"			\$155

65. PROXIMITY HEAD ARRAY

SCSIHAA	Adult – Switch-It 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse			\$5,995
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.			
ESS-RED	Emergency STOP Switch - Red	User Activation	Attendant Activation	\$314
	▲ Required with the Proximity Head Array.			
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Proximity Head Array.			

66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL

SCSIMP	Switch-It Light touch mini Proportional Joystick			\$4,495
SCSRMMP	Stealth Rock-n-Lock Swingaway Arm for Chin MP	SWL (Swings on Left)	SWR (Swings on Right)	\$870
	▲ Required with the Light touch mini Proportional Joystick.			
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.			
JCC	Chin Cup			\$55
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Light touch mini Proportional Joystick.			
▶ CHOOSE ONE (1) HEADREST PAD				
SOFHRS	Small Curved Headrest Pad – 4" x 9"			\$155
SOFHRM	Medium Curved Headrest Pad – 5" x 10"			\$155
SOFHRL	Large Curved Headrest Pad – 6" x 11"			\$155



SPECIALTY DRIVER CONTROL OPTIONS (CONTINUED)

▼ Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)

67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT

SCPGMJS	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$2,369
SCGMJSM	Gatlin Midline Swing-away Mount for JC200	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
SCSIMP	Switch-It Light touch mini Proportional Joystick			\$4,495
SCGMMPM	Gatlin Midline Swing-away Mount for Light touch mini proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT

SCPGMJS-M	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$2,369
SCARMJSM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ Must be selected for use with the PG Drives JC200 Mini Joystick.			
SCSIMP	Switch-It Light touch mini Proportional Joystick			\$4,495
SCARMPM	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$375
	▲ Must be selected for use with the Light touch mini Proportional Joystick.			
SCBTM	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ This option is only compatible with the Light touch Mini Proportional Joystick.			

69. PROPORTIONAL RIM HEAD CONTROL PACKAGE

SCPHC	Proportional RIM Head Control Package			\$2,795
SCSRM	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
SSES	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ Required with the Stealth Rock-n-Lock Swing Away mount.			
	▲ The Mode Switch Type can be changed from this selection in SWITCHES & OPTIONS section.			
SOMAH	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ Required with the Proportional RIM Head Control.			



CUSTOM SECTION

► Lead times will increase depending on the custom order pipeline at the time the purchase order is received. Please call for Quote and specifications for these options.

70. O² HOLDER

SUPOXY	Custom Non-Articulating Oxygen Tank Holder ▲ Not available with a power recline	MLH (Mount on Left)	MRH (Mount on Right)	\$295
SUPOXYA	Custom Articulating Oxygen Tank Holder ▲ Required with a power recline D-Size Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve E-Size Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve	MLH (Mount on Left)	MRH (Mount on Right)	\$1,195

71. VENT TRAY (PLEASE COMPLETE IN FULL)

▼ Amylor does not supply the battery, charger or the wiring for the Vent. These items are normally provided by the Vent Manufacturer.

SOVENT	Custom Vent Tray Articulating Vent Box – (Required with Power Tilt and/or Recline) Fixed Vent Box			\$1,360 \$400
	SPECIFY MODEL Trilogy 100/200 LTV LTV with Sprintpack Other, specify: _____			
SOVTBE	Extended Battery Compartment ▲ The Vent Tray option is not available on narrow options N2424 & N2221 (FS2424B & FS2221 B) The Vent is being powered by 1 U1 Battery and the chair by 2 Group 22NF batteries The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 22NF Batteries The Vent is being powered by 1 Group 24 Battery and the chair Powered by 2 Group 22NF Batteries The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries			\$495

72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)

▼ Please specify details and specifications for these options.

CUSTOM	CUSTOM Lithium ion - 54 Amps (represents one 24V Battery pack, equivalent to two Group 22 Batteries)			RFQ
CUSTOM	CUSTOM Lithium ion - 78 Amps (represents one 24V Battery pack, equivalent to two Group 24 Batteries)			RFQ
CUSTOM	CUSTOM 8A Offboard Battery Charger for Lithium ion batteries ▲ This option is required when choosing custom Lithium ion batteries			N/C
CUSTOM	CUSTOM Lithium ion Battery capacity LCD display gauge ▲ This option is required when choosing custom Lithium ion batteries			N/C
CUSTOM	Rear view Camera with 5" Monitor			RFQ
CUSTOM	USB Charger Port ◀ Specify the location for the installation: _____			RFQ
CUSTOM	Folding Backrest for Additional Storage possibilities and Travel accessibility			RFQ
CUSTOM	Large Ball Handle			RFQ
CUSTOM	Tilt + Elevate Function Through VR2 Non-Expandable Drive Control			RFQ
CUSTOM	Manual Wheel Locks			RFQ
CUSTOM	Power Anterior Tilt			RFQ
CUSTOM	Power Lateral Tilt 15° Raise the Left Side only Raise the Right Side only Raise the Left & Right Side			RFQ
CUSTOM	Power Latitude System – (Features 14" of forward travel and a 5" seat to floor height from the pan to the ground)			RFQ
CUSTOM	Power Elevating Stretcher – (Please call CSR for availability with configuration and user weight capacity)			RFQ
CUSTOM	Rearview mirror Left Side Right Side			RFQ
CUSTOM	Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity)			RFQ
CUSTOM	Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity)			RFQ
CUSTOM	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan			RFQ
CUSTOM	Endomorphic Back Canes – (Pair) 1" per side towards the inside = 2" 2" per side towards the inside = 4"			RFQ
CUSTOM	Transfer Handles – (Pair) Select a height range: 5" to 8" 6" to 9"			RFQ
CUSTOM	Foot Box Length: _____ Width: _____ Please confirm with CSR on sizing possibilities			RFQ
CUSTOM # 1	Other custom options. Specify in notes section below			RFQ
CUSTOM # 2	Other custom options. Specify in notes section below			RFQ



73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS
