



**1. ACCOUNT INFORMATION**

P.O. # or Quote: \_\_\_\_\_  
 Dealer Contact Name: \_\_\_\_\_  
 Email Address: \_\_\_\_\_  
 Phone Number: \_\_\_\_\_ Fax Number: \_\_\_\_\_  
 Bill to (Name): \_\_\_\_\_  
 Bill to (Location): \_\_\_\_\_  
 Ship to (Name): \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_  
 Province: \_\_\_\_\_ Postal Code: \_\_\_\_\_  
 User Name or Number: \_\_\_\_\_  
 User Weight: \_\_\_\_\_  
 User Height: \_\_\_\_\_  
 Seat Cushion Thickness: \_\_\_\_\_  
 Knee to Heel Length: Left: \_\_\_\_\_ " Right: \_\_\_\_\_ "



**THE ORDER CANNOT BE PROCESSED WITHOUT THE FOLLOWING:**

- A user name or number, the client's measurements and weight - Amylor is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Amylor summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

*\*IN RED SHOWS THE ADDITION OF NEW OPTIONS, UPGRADES, PRICING MODIFICATIONS OR ADDED NOTES*

**2. STANDARD OPTIONS**

Alltrack HD3 Power Base, 6-wheel independent wheel suspension, an additional seat suspension system, tie-down hardware for strap systems, 8A offboard charger, 8" front casters, and standard auto style buckle seat positioning strap, 5.4mph motors, Weight Capacity = 550 lb

**3. ALLTRACK HD3 POWER WHEELCHAIR**

<b>ATHD3</b>	ALLTRACK HD3 Power Wheelchair	\$11,295
--------------	-------------------------------	----------

**3.5 OWNER'S MANUAL**

<b>SOOMATENG</b>	English / Anglais	N/C	<b>SOOMATFR</b>	Français / French	N/C
------------------	-------------------	-----	-----------------	-------------------	-----

**4. USER WEIGHT PACKAGES**

<b>WCA55</b>	550 lb Maximum user weight capacity	STD
--------------	-------------------------------------	-----

**5. STANDARD COLORS**

<b>SCMOBM</b>	Onyx Black Matte	N/C	<b>SCMCAB</b>	Cabernet	N/C
<b>SCMCTB</b>	Cobalt Blue	N/C	<b>SCMBGM</b>	Bubble Gum Metallic	N/C
<b>SCMIOB</b>	Indigo Blue	N/C	<b>SCMHAB</b>	Habanero Metallic	N/C
<b>SCMCRE</b>	Candy Red	N/C	<b>SCMLIM</b>	SubLime	N/C

**▼ OPTIONAL COLORS**

<b>FCHDD-M</b>	Desert Digital Camo HydroGraphics	\$795
<b>FCHCF-M</b>	Carbon fiber HydroGraphics	\$795
<b>FCHWL-M</b>	Woodland Camo HydroGraphics	\$795



**5.5 WHEEL INSERT COLORS**

<b>SCCTB-WI</b>	Cobalt Blue	N/C	<b>SCHAB-WI</b>	Habanero Metallic	N/C
<b>SCIOB-WI</b>	Indigo Blue	N/C	<b>SCLIM-WI</b>	SubLime	N/C
<b>SCCRE-WI</b>	Candy Red	N/C	<b>SCBGM-WI</b>	Bubble Gum Metallic	N/C
<b>SCCAB-WI</b>	Cabernet	N/C			
▲ <i>Wheel insert colors are not compatible with optional Hydrographic colors (FCHDD, FCHWL and FCHCF)</i>					
▲ <i>Wheel insert colors are not compatible with the optional 14" X 4" Extra Wide Cast Wheels (DW144PT and DW144PI)</i>					
▲ <i>Wheel inserts are not available in Onyx Black Matte</i>					

**7. TRANSPORT OPTION**

<b>TTD4</b>	Tie-down attachment points compliant with ISO 7176-19. See owner's manual for additional details				STD
-------------	--	--	--	--	-----

**8. MOTOR PACKAGE**

<b>MP544B</b>	5.4 MPH 4-Pole Motor Package / with Brake Release				N/C
▲ <i>Not compatible with the Smart-Track Assistive Driving Package option (SMARTT)</i>					
<b>MP544-OPT</b>	5.4 MPH 4-Pole Motor Package / with COG Release				\$395

**9. ASSISTIVE DRIVING PACKAGE**

<b>SMARTT</b>	Smart-Track Assistive Driving Package				\$1,420
▲ <i>Smart-Track technology provides increased tracking capabilities with enhanced traction and steering precision. This is ideal for those with marginal hand function or switch control.</i>					
▲ <i>Requires the R-net PM 120 Amp Expandable Controller (ECR120)</i>					
▲ <i>This option requires the MP544-OPT motor package.</i>					
▲ <i>Recommended for non-proportional Drive Controls</i>					

**10. BATTERY CHARGER**

<b>BCOB</b>	8A Offboard Battery Charger				N/C
-------------	-----------------------------	--	--	--	-----

**11. BATTERY TRAY AND BASE WIDTH**

<b>FS24B</b>	Battery Tray for Group 24 Batteries – (Total base width = 25 3/8")				N/C
--------------	--	--	--	--	-----

**12. BATTERIES - REQUIRES A SELECTION**

<b>BAT24GEL</b>	Premium CG Power - GEL Battery – Group 24 – 70 Amps (Each)				\$465
<b>BATO</b>	Omission of the batteries				N/C

**14. DRIVE WHEEL**

<b>DW14PT</b>	14" Cast Wheel Pneumatic Tires				STD
<b>DW14PI</b>	14" Cast Wheel Pneumatic with Foam Filled Tires				N/C
<b>DW144PT</b>	14" x 4" Extra Wide Cast Wheels, Pneumatic Tires				\$200
<b>DW144PI</b>	14" x 4" Extra Wide Cast Wheels, Pneumatic with Foam Filled Tires				\$595
▲ <i>Increases the overall width of the base by 2 3/4"</i>					

**15. FRONT CASTERS**

<b>FC8PTM</b>	8" Pneumatic Casters				STD
<b>FC8PIM</b>	8" Pneumatic Casters with Foam Filled Tires				N/C



**16. SEAT TO FLOOR HEIGHT (MEASUREMENT TO THE TOP OF SEAT PAN ± 1/4")**

SEAT TO FLOOR HEIGHT TABLE	NO POWER OPTIONS			TILT OR TILT/RECLINE
	0°	3°	6°	
STFH1700 (17" STFH)				
STFH1800 (18" STFH)				
STFH1900 (19" STFH)				
STFH2000 (20" STFH)				

**17. SEAT TO FLOOR HEIGHT SPECIFICATION (IF NOT STANDARD)**

<b>SOMHSS</b>	Specific seat to floor height higher than standard. It's not possible to lower the seat to floor height below minimum ranges.	RFQ
▲ <b>Required at:</b> _____ " (± 1/4" Range)		

**18. PRE-SET SEAT FRAME ANGLE**

<b>SFE-0</b>	0° Pre-set seat frame seat angle (Must be selected when using power options)	STD
<b>SFE-3</b>	3° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C
<b>SFE-6</b>	6° Pre-set seat frame seat angle (Only available with a seat configuration without power options)	N/C

**19. LIGHTING PACKAGE**

<b>EOLP</b>	LED Light Package – (Includes front and rear lights with integrated flashers)	\$1,425
▲ <b>Requires the upgraded Joystick with Light Operation or the OMNI 2 (ECVJ6KL, ECRJCL, ECRJC35B or EDIB).</b>		

**20. NON-EXPANDABLE ELECTRONICS - PF DRIVES VR2 CONTROLLER AND JOYSTICKS**

▶ **VR2 Electronics can only accommodate specific power seat functions and combinations such as Tilt / Elevate / Recline / Power Center Mount / Combined Power Legs / Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs / Tilt + Elevate (Custom)**

<b>ECV902A</b>	VR2 90 Amp Drive Controller – (Requires the selection of VR2 joystick below)	STD
<b>ECVJ4K</b>	VR2 Standard 4-Key Joystick	STD
<b>ECVJ6K</b>	VR2 6-Key Joystick (Allows for single or select dual power options to operate through this joystick)	N/C
<b>ECVJ9KL</b>	VR2 9-Key Joystick with Lighting Operation (Allows for single or select dual power options to operate through this joystick)	\$150

**21. EXPANDABLE ELECTRONICS - PG DRIVES R-NET CONTROLLERS**

<b>ECEL90</b>	R-net EL 90 Amp Expandable Controller	\$175
<b>ECR120</b>	R-net PM 120 Amp Expandable Controller	\$275

**22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS**

▶ **SELECTION OF ONE JOYSTICK OPTION BELOW IS REQUIRED WITH THE EXPANDABLE CONTROLLER ECEL90 OR ECR120**

<b>ECO</b>	Omit Joystick – (For use with an Omni 2 Display)	N/C
<b>ECELJ</b>	R-net - (LED) – Joystick	STD
▲ <b>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</b>		
<b>ECRJC</b>	R-net Joystick with Color screen	\$155



**22. EXPANDABLE ELECTRONICS - PG DRIVES R-NET JOYSTICKS (CONTINUED)**

▶ **JOYSTICKS WITH LIGHTING OPERATION**

<b>ECELJL</b>	R-net – (LED) – Joystick with Lighting Operation	\$125
	▲ <i>This joystick is not compatible with the Memory Seat Program Option (SACMEMORY and MACMEMORY)</i>	
<b>ECRJCL</b>	R-net Joystick with Color screen and Lighting Operation	\$180
▶ <b>ADVANCED JOYSTICK</b>		
<b>ECRJ35B</b>	R-net Advanced Bluetooth Joystick (3.5" Screen). Lighting Operation, paddle switches, shortcut keys and stereo jacks	\$1,105

**23. JOYSTICK MOUNT**

<b>JSMF</b>	Fixed Joystick Mount	N/C
<b>JMSA</b>	Swing Away Joystick Mount – Height adjustable, Angle adjustable and Multi-Axis	\$235
<i>Select:</i>	<i>Set-up with standard installation to swing-away sideways</i>	
	<i>Set-up at 35° to swing-away sideways and downwards</i>	
	▲ <i>Set-up at 35° angle may cause interference with the user's limbs when the new Signature cantilever armrest (SOAFBxxx) is internally rotated.</i>	
<b>JSMCM</b>	Bodypoint Tri-Lock Midline Mounted Swing Away Joystick Mount	\$275

**24. JOYSTICK MOUNT SUPPORT**

<b>JSMFH</b>	Joystick mount support Bracket	STD
--------------	--------------------------------	-----

**25. JOYSTICK HAND MOUNT**

<b>HML</b>	Install Joystick on Left hand side	<b>HMR</b>	Install Joystick on Right hand side
------------	------------------------------------	------------	-------------------------------------

**26. JOYSTICK HANDLE**

<b>JSH</b>	Standard Joystick Gimble	STD	<b>JBAT</b>	Stick Handle 3" in Length	\$130
<b>JBH</b>	Ball Handle	\$85	<b>JFBH</b>	Foam Ball Handle	\$70
<b>JTH</b>	T-Handle	\$155	<b>JMUSH</b>	Bodypoint Mushroom Handle	\$210
<b>JGP</b>	Bodypoint 4" U-shaped Handle	\$210			

**27. ATTENDANT DRIVE CONTROL OPTION**

<b>EADC</b>	Attendant Hand Control	\$1,150
-------------	------------------------	---------

**28. ATTENDANT JOYSTICK HAND MOUNT**

<b>EADCL</b>	Install Attendant Joystick on Rear Left hand side	<b>EADCR</b>	Install Attendant Joystick on Rear Right hand side
--------------	---	--------------	--

**29. R-NET OMNI 2 - ENHANCED DISPLAY & DIGITAL INTERFACE BOX (REQUIRED WITH SPECIALTY CONTROLS)**

<b>EDIB</b>	Omni 2 – (Enhanced Display & Digital Interface Box with Infrared)	\$1,915
	▲ <i>Allows the switch operation of alternative devices through driver controls. Requires the R-net fully Expandable Controller (ECR120)</i>	
	▲ <i>The Omni 2 display mount must be selected from the option below.</i>	

**30. OMNI 2 ENHANCED DISPLAY MOUNT**

<b>EDSAM</b>	User Swing Away Mount	N/C
	<b>EDML – (Install Display on LH side)</b>	<b>EDMR – (Install Display on RH side)</b>



**31. ENVIRONMENTAL CONTROL MODULE**

<b>ECU</b>	Environment control module – (Output Module)	\$995
	▲ <i>Requires one of the fully expandable R-net controllers (ECEL90 or ECR120).</i>	

**33. SIGNATURE TAILORED POWER POSITIONING SOLUTIONS / REHAB SEAT**

<b>SAABAHD-550</b>	Bariatric <b>Signature</b> Rehab Seat – (Seat Frame only)	\$1,045
<b>SAPTXVHD</b>	Bariatric <b>Signature</b> CG Power Tilt	\$8,995
	▲ <i>This option includes the <b>Signature</b> Rehab Seat – (SAABAHD)</i>	
<b>SAPRXVHD</b>	Bariatric <b>Signature</b> 165° Power Recline with Mechanical Shear Reduction (MSR)	\$9,845
	▲ <i>This option includes the <b>Signature</b> Rehab Seat – (SAABAHD)</i>	
<b>SAPTPRXVHD</b>	Bariatric <b>Signature</b> CG Power Tilt and 165° Power Recline with Mechanical Shear Reduction (MSR)	\$14,495
	▲ <i>This option includes the <b>Signature</b> Rehab Seat – (SAABAHD)</i>	

**34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER**

▶ <b>SINGLE SEAT FUNCTION USING AN EXTERNAL TOGGLE OR OPTIONAL SWITCH ONLY (For both expandable and non-expandable systems)</b>									
<b>SOLO</b>	Single function controller (Includes an external toggle switch)								STD
▼ <b>OPTIONAL SWITCHES</b>									
<b>SWTOGGLE</b>	Single Toggle Switch							Attendant Activation	\$100
<b>SWPUSH1</b>	Single Push Button Switch – (Single Mode: SAC)					User Activation	Attendant Activation	\$90	
<b>SWPUSH2</b>	Dual Push Button Switch – (Dual Mode: SAC)					User Activation	Attendant Activation	\$180	
<b>SSBS</b>	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$230	
<b>SSES</b>	Egg Switch-Black					User Activation	Attendant Activation	\$116	
<b>SSML</b>	Ultra Light Switch					User Activation	Attendant Activation	\$295	
▶ <b>SINGLE SEAT FUNCTION THROUGH DRIVE CONTROL</b>									
<b>EPSV</b>	Single Function Through VR2 Non-Expandable Drive Control								N/C
	▲ <i>Allows for activation of a single power option through VR2 controller</i>								
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>								
	▲ <i>VR2 Electronics can only accommodate specific power seat functions such as Tilt / Recline / Power Center Mount / Combined Power Legs</i>								
	▲ <i>Choose EPSV2 to operate select power combinations</i>								
<b>EPSEL</b>	Single Function Through Drive Control with R-net 90-amp controller								N/C
	▲ <i>Allows for activation of single power option through R-net 90-amp controller</i>								
	▲ <i>Compatible with all R-net joysticks</i>								
	▲ <i>Does not allow for control of power option through drive control and separate switches</i>								
	▲ <i>Can not operate more than one power option</i>								
	▲ <i>If selecting memory seating, choose SAC below</i>								
	▲ <i>If selecting 120-amp R-net controller, choose SAC below</i>								
<b>SAC</b>	Single Function Through Drive Control with R-net 90-amp or 120-amp controller								N/C
	▲ <i>Allows for activation of single power option through R-net drive controls</i>								
	▲ <i>Compatible with 90-amp or 120-amp controller, all R-net joysticks, Omni 2 and specialty controls</i>								
▶ <b>SINGLE SEAT FUNCTION OPTIONS (OPTIONS BELOW ARE COMPATIBLE ONLY WHEN CHOOSING SAC ABOVE)</b>									
<b>SACMEMORY</b>	Memory Seat Program								\$1,495
	▲ <i>This option is not compatible with any VR2 or LED joysticks (ECELJ and ECELJL)</i>								
	▲ <i>Memory seat programming is not compatible with the Very Heavy Duty Tilt option (SAPTXVHD and SAPTPRXVHD)</i>								

MKT.F74 | OF | PD-A20602 | R07



**34. ELECTRONICS - SINGLE (1 ACTUATOR) POWER SEAT FUNCTION CONTROLLER (CONTINUED)**

<b>SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH</b>								
<b>SACTOGGLE</b>	Single Toggle Switch		Attendant Activation	\$410				
<b>SACPUSH1</b>	Single Push Button Switch – (Single Mode: SAC)	User Activation	Attendant Activation	\$400				
<b>SACPUSH2</b>	Dual Push Button Switch – (Dual Mode: SAC)	User Activation	Attendant Activation	\$490				
<b>SACSSBS</b>	Buddy Button Switch	Black	Red	Green	Yellow	User Activation	Attendant Activation	\$540
<b>SACSSES</b>	Egg Switch-Black	User Activation	Attendant Activation	\$426				
<b>SACSSML</b>	Ultra Light Switch	User Activation	Attendant Activation	\$605				

**35. ELECTRONICS - MULTIPLE SEAT FUNCTION CONTROL TYPE**

<b>DUAL SEAT FUNCTIONS THROUGH DRIVE CONTROL</b>				
<b>EPSV2</b>	Dual Function Through Drive Control with VR2 Non-Expandable Electronics			N/C
	▲ <i>Allows for the activation of 2 power options through the VR2 Controller. Combination of TILT and RECLINE is not possible with EPSV2</i>			
	▲ <i>Does not allow for control of power options through drive control and separate switches</i>			
	▲ <i>VR2 Electronics can only accommodate specific power seat combinations such as Tilt + Power Center Mount / Tilt + Combined Power Legs / Recline + Power Center Mount / Recline + Combined Power Legs</i>			

**MULTIPLE SEAT FUNCTIONS THROUGH R-NET**

<b>MAC</b>	E-BX Multiple Function through Drive Control with Expandable Electronics – Powers iQ Actuators	\$2,295
------------	--	---------

**MAC OPTIONS**

<b>MACMEMORY</b>	Memory Seat Program	\$175
	▲ <i>Memory seat programming is not compatible with the Very Heavy Duty Tilt option (SAPTXVHD and SAPTPRXVHD)</i>	

**SELECTION OF SWITCH BELOW GIVES ACCESS TO POWER OPTIONS THROUGH DRIVE CONTROLS AND THROUGH EXTERNAL SWITCH**

<b>MACKEYPAD-P</b>	Keypad w/ 5 Push Button Switches	User Activation	Attendant Activation	\$445
<b>MACKEYPAD-T</b>	Keypad w/ 5 Toggle Switches	User Activation	Attendant Activation	\$445
<b>MACMULTI</b>	4-Way Toggle Switch – (Up to 4 Activations)	User Activation	Attendant Activation	\$180
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
<b>MACMULTIT</b>	Triple Toggle Switch – (Up to 5 Activations)	User Activation	Attendant Activation	\$135
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
<b>MACPUSH4</b>	4-Way Push Button Switch	User Activation	Attendant Activation	\$225
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			
<b>MACPUSH5</b>	5-Way Push Button Switch	User Activation	Attendant Activation	\$240
	▲ <i>By selecting this switch, memory recall is only available through the Joystick.</i>			

**36. SWITCH ACTIVATION HAND MOUNT**

Switch Installation:	<b>TSML – (Install Switch on LH side)</b>	<b>TSMR – (Install Switch on RH side)</b>	N/C
----------------------	---	---	-----

**37. SIGNATURE TAILORED SOLUTIONS - PRE-SET BACK ANGLE (FOR NON-RECLINING BACKS)**

<b>SOBAA</b>	Select one:	82°	90°	98°	106°	N/C
	▲ <i>Back angle of 106° is not compatible with the OMNI 2.</i>					



**39. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT WIDTH RANGE BETWEEN 20" AND 28" (510 - 710 MM)**

<b>STS20-25-22</b>	Width adjustment from 20" to 25" (510-635 MM) and Depth range of 18" to 22" (455 - 560 MM)	N/C			
	▲ <i>The minimum seat depth available with these options is 18" (455 MM)</i>				
<b>STS23-28-22</b>	Width adjustment from 23" to 28" (585-710 MM) and Depth range of 18" to 22" (455 - 560 MM)	N/C			
	▲ <i>The minimum seat depth available with these options is 18" (455 MM)</i>				
<b>STS23-28-24</b>	Width adjustment from 23" to 28" (585-710 MM) and Depth range of 18" to 24" (455 - 610 MM)	N/C			
	▲ <i>The minimum seat depth available with these options is 20" (510 MM)</i>				
<b>20" (SOSW20)</b>	<b>21" (SOSW21)</b>	<b>22" (SOSW22)</b>	<b>23" (SOSW23)</b>	<b>24" (SOSW24)</b>	<b>25" (SOSW25)</b>
<b>26" (SOSW26)</b>	<b>27" (SOSW27)</b>	<b>28" (SOSW28)</b>			

**40. SIGNATURE TAILORED SOLUTIONS - ADJUSTABLE SEAT DEPTH RANGE BETWEEN 18" AND 24" (455 - 610 MM)**

<b>18" (SODN18)</b>	<b>19" (SODN19)</b>	<b>20" (SODN20)</b>	<b>21" (SODN21)</b>	<b>22" (SODN22)</b>	<b>23" (SODN23)</b>	N/C
<b>24" (SODN24)</b>						N/C

**41. SIGNATURE TAILORED SOLUTIONS - EXTRUDED BACK CANE HEIGHT (THE MEASUREMENT IS TAKEN FROM THE SEAT PAN TO THE TOP OF THE BACK CANE)**

► *Not compatible with a power or manual recline*

<b>16" (SOBCH16)</b>	<b>17" (SOBCH17)</b>	<b>18" (SOBCH18)</b>	<b>19" (SOBCH19)</b>
----------------------	----------------------	----------------------	----------------------

▼ *20" is the minimum back height required with a power or manual recline*

<b>20" (SOBCH20)</b>	<b>21" (SOBCH21)</b>	<b>22" (SOBCH22)</b>	<b>23" (SOBCH23)</b>	<b>24" (SOBCH24)</b>	<b>25" (SOBCH25)</b>
<b>26" (SOBCH26)</b>	<b>27" (SOBCH27)</b>	<b>28" (SOBCH28)</b>			

**41.5 SIGNATURE TAILORED SOLUTIONS - BACK CANE OPTIONS**

▼ **PUSH HANDLES**

<b>SOPH</b>	Push Handles (pair)	N/C
-------------	---------------------	-----

▼ **ENDOMORPH BACK CANES FOR NON-RECLINING BACKS**

<b>STSEND0-F1</b>	Endomorph Back Canes for non-reclining backs, 1 IN (25 MM) total	\$120
<b>STSEND0-F2</b>	Endomorph Back Canes for non-reclining backs, 2 IN (50 MM) total	\$120
<b>STSEND0-F3</b>	Endomorph Back Canes for non-reclining backs, 3 IN (75 MM) total	\$120
<b>STSEND0-F4</b>	Endomorph Back Canes for non-reclining backs, 4 IN (100 MM) total	\$120

▼ **ENDOMORPH BACK CANES FOR RECLINING BACKS**

<b>STSEND0-R1</b>	Endomorph Back Canes for reclining backs, 1 IN (25 MM) total	N/C
<b>STSEND0-R2</b>	Endomorph Back Canes for reclining backs, 2 IN (50 MM) total	N/C

▼ **FOLDING BACKREST (FOR NON-RECLINING BACKS)**

<b>STSFb</b>	Folding Backrest for Storage/Travel Access for non-reclining backs	\$240
--------------	--	-------

**42. BACK SUPPORT**

<b>DSTR-wwhh</b>	Adjustable Tension Back Upholstery, WW x HH	N/C
	▲ <i>Only available with Dual Post armrests (SOADP)</i>	
<b>STDSB-wwhh</b>	Contour back support up to 21" wide, with a maximum shell height of 22".	\$605
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
<b>STDSB-wwhh</b>	Contour back support up to 21" wide for non-standard sizes with shell height ≥ 23".	\$755
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	
<b>STDSB-wwhh</b>	Contour back support for Heavy-Duty from 22" wide and higher, with any height.	\$1,075
	▲ <i>The minimum space between the back and seat pan is 4" (± 1") with a fixed back and 6" (± 1") with a reclining back</i>	



**42. BACK SUPPORT (CONTINUED)**

▶ **MATRX BACKS**

<b>PBE-wwhh-HD</b>	MATRX HEAVY DUTY ELITE BACK – (21” to 24” Wide - Back height of 16” and 20”)	\$1,150
	▲ <i>Back shell height is available in 16” or 20”. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i>	
	▲ <i>Weight capacity up to 500 lb</i>	
<b>EDB-wwhh-HD</b>	MATRX ELITE DEEP HD BACK – (21” to 24” Wide - Back height of 16”, 18” and 20”)	\$1,195
	▲ <i>Back shell height is available in 16”, 18” or 20”. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Weight capacity up to 500 lb</i>	

▶ **NXT BACKS**

<b>NOACB-wwhh</b>	NXT OPTIMA FH AC CONTOUR BACK – (15” to 21” Wide)	\$690
	▲ <i>Weight capacity up to 300 lb</i>	
<b>NOACB-wwhh</b>	NXT OPTIMA FH AC CONTOUR BACK – (22” to 28” Wide)	\$720
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>Back shell heights only come in 13”, 16” or 19” high. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i>	
<b>NOPCB-wwhh</b>	NXT OPTIMA FH PC DEEP CONTOUR BACK – (15” to 21” Wide)	\$790
	▲ <i>Weight capacity up to 300 lb</i>	
<b>NOPCB-wwhh</b>	NXT OPTIMA FH PC DEEP CONTOUR BACK – (22” to 28” Wide)	\$820
	▲ <i>Weight capacity up to 550 lb</i>	
	▲ <i>This product is not compatible with Lateral Trunk Supports</i>	
	▲ <i>Back shell heights only come in 13”, 16” or 19” high. The minimum space between the back and seat pan is 4” (± 1”) with a fixed back and 6” (± 1”) with a reclining back</i>	

▶ **OMISSION OF THE BACK SUPPORT**

<b>STDSBO</b>	Omission of the Contour Back Support ◀ <b>Specify:</b> _____	N/C
<b>STDSBO-UNIVBKT</b>	Universal 4-Point Fixed Brackets for Aftermarket Backs	\$185

**43. BACK SUPPORT TOTAL HEIGHT (MEASURED FROM THE SEAT PAN TO THE TOP OF THE BACK)**

			▼ <i>20” is the minimum back height required with a reclining back.</i>				
<b>16” (SOBH16)</b>	<b>18” (SOBH18)</b>	<b>20” (SOBH20)</b>	<b>22” (SOBH22)</b>	<b>24” (SOBH24)</b>	<b>26” (SOBH26)</b>	<b>28” (SOBH28)</b>	
<i>This height is not available with the NXT Optima AC and PC options</i>							▲

**44. SEAT CUSHIONS**

▶ **PERFORMA CG AIR - THERAPEUTIC AIR CUSHION**

Custom Cushion (Please use the **Performa CG Air** order form)

Cell Height (Choose 1):  2”  3”  4”

▶ **SINGLE (1) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING**

▶ **PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER**

<b>CGA1VSCxxxx-x</b>	14” to 20” (355 mm to 510 mm) wide	\$696
<b>CGA1VSCxxxx-x</b>	22” + (560 mm +) wide	\$783
	▲ <i>Cushion sizes only come in even numbers for width and depth measurements.</i>	



**44. SEAT CUSHIONS**

<b>▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER</b>		
<b>CGA1V3Dxxxx-x</b>	14" to 20" (355 mm to 510 mm) wide	\$735
<b>CGA1V3Dxxxx-x</b>	22" + (560 mm +) wide	\$825
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER</b>		
<b>CGA1VICxxxx-x</b>	14" to 20" (355 mm to 510 mm) wide	\$785
<b>CGA1VICxxxx-x</b>	22" + (560 mm +) wide	\$830
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>▶ DUAL (2) VALVE - DESIRED COVER TYPE, DIMENSIONS AND PRICING</b>		
<b>▶ PERFORMA CG AIR CUSHION WITH STANDARD SPACER KNIT COVER</b>		
<b>CGA2VSCxxxx-x</b>	14" to 20" (355 mm to 510 mm) wide	\$696
<b>CGA2VSCxxxx-x</b>	22" + (560 mm +) wide	\$783
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>▶ PERFORMA CG AIR CUSHION WITH 3D SPACER KNIT DOUBLE LAYER COVER</b>		
<b>CGA2V3Dxxxx-x</b>	14" to 20" (355 mm to 510 mm) wide	\$735
<b>CGA2V3Dxxxx-x</b>	22" + (560 mm +) wide	\$825
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>▶ PERFORMA CG AIR CUSHION WITH EXTERNAL INCONTINENCE COVER</b>		
<b>CGA2VICxxxx-x</b>	14" to 20" (355 mm to 510 mm) wide	\$785
<b>CGA2VICxxxx-x</b>	22" + (560 mm +) wide	\$830
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>▶ DUAL VALVE CUSHION SEGMENTATION (PLEASE SELECT 1 BOX)</b>		
<b>CGA-2VLR</b>	CG Air Dual Valves with Left / Right Location	N/C
<b>CGA-2VF6</b>	CG Air Dual Valves with Front / Back Location, Front distance 6 in (150 mm) from the front	N/C
<b>CGA-2VF8</b>	CG Air Dual Valves with Front / Back Location, Front distance 8 in (205 mm) from the front	N/C
<b>CGA-2VF10</b>	CG Air Dual Valves with Front / Back Location, Front distance 10 in (255 mm) from the front	N/C
<b>▶ PERFORMA CG AIR OPTIONS</b>		
<b>CGA-GZ</b>	GripZone Anti-Slip Technology	\$75
<b>CGA-TT-WW</b>	Transfer Trellis	\$125
<b>CGA-AC</b>	Anatomic Contour Seat Insert	\$140
<b>CGA-SP</b>	Seat Pouch	\$44
<b>▶ CONTOUR SEAT CUSHION</b>		
<b>SCPSC-WWDD</b>	Contour Seat Cushion up to 19" Wide	\$620
<b>SCPSC-WWDD</b>	Contour Seat Cushion for Non-Standard Sizes of 20" to 24" Wide	\$735
<b>SCPSC-WWDD</b>	Contour Seat Cushion for Heavy-Duty of ≥ 25" Wide ▲ Specify: Width: _____" Depth: _____"	\$865
<b>▶ MATRX CUSHION</b>		
<b>PSP-WWDD</b>	MATRX PSP Seat Cushion – (11" to 20" Wide)	\$645
▲ <b>The maximum user weight capacity is 300 lb</b>		
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		
<b>PSP-WWDD</b>	MATRX PSP Seat Cushion – (22" to 26" Wide)	\$995
▲ <b>The maximum user weight capacity is 550 lb</b>		
▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>		



**44. SEAT CUSHIONS (CONTINUED)**

▶ **NXT BIO-FIT CUSHION**

<b>NXTBIOFIT-wwdd</b>	NXT BIO FIT Cushion – (15" to 21" Wide)	\$595
	▲ <b>The maximum user weight capacity is 300 lb</b>	
<b>NXTBIOFIT-wwdd</b>	NXT BIO FIT Cushion – (22" to 26" Wide)	\$625
	▲ <b>The maximum user weight capacity is 400 lb</b>	
	▲ <b>Cushion sizes only come in even numbers for width and depth measurements.</b>	

**45. SIGNATURE TAILORED SOLUTIONS - CANTILEVER ARMRESTS - ADJUSTABLE HEIGHT, INWARD/ OUTWARD ROTATION AND VERTICAL ANGLE. ALLOWS ARMPAD TO MOVE 1" FORWARD (CHOOSE 2)**

**FOR NON-RECLINING BACKS**

▼	<b>Height Adjustable from 7.25" (185 mm) to 15.25" (385 mm) in 1" (25 mm) increments</b>	
▼	<b>Not compatible with Waterfall Armpads (SOWAPx et SOWGELx)</b>	
<b>SOAFBDLH</b>	Signature Cantilever Armrest – Desk Length 12" (305 mm) - Left Side	N/C
<b>SOAFBDRH</b>	Signature Cantilever Armrest – Desk Length 12" (305 mm) - Right Side	N/C
<b>SOAFBFLH</b>	Signature Cantilever Armrest – Full Length 16" (405 mm) - Left Side	N/C
<b>SOAFBFRH</b>	Signature Cantilever Armrest – Full Length 16" (405 mm) - Right Side	N/C
<b>SOAFBELH</b>	Signature Cantilever Armrest – Extended Length 18" (455 mm) - Left Side	N/C
<b>SOAFBERH</b>	Signature Cantilever Armrest – Extended Length 18" (455 mm) - Right Side	N/C
▲	<b>The extended-length 18" (455 mm) armrests are only compatible with the extra-long armpads (16" x 3 ½") (SOXWFAP16 or SOXWGAP16)</b>	
▼	<b>CUSTOM LENGTH</b>	
<b>SOAFBCLH</b>	Signature Cantilever Armrest – Custom Length 20" (510 mm) - Left Side	\$95
<b>SOAFBCRH</b>	Signature Cantilever Armrest – Custom Length 20" (510 mm) - Right Side	\$95
<b>SOAFBCLH</b>	Signature Cantilever Armrest – Custom Length 22" (560 mm) - Left Side	\$95
<b>SOAFBCRH</b>	Signature Cantilever Armrest – Custom Length 22" (560 mm) - Right Side	\$95
▲	<b>The custom-length 20" (510 mm) and 22" (560 mm) armrests are only compatible with the extra-long armpads (16" x 3 ½") (SOXWFAP16 or SOXWGAP16)</b>	

**46. DUAL POST – HEIGHT ADJUSTABLE ARMRESTS – INCLUDES TRANSFER GRIPS (FOR NON-RECLINING BACKS - CHOOSE 2)**

▼	<b>NOT compatible with the extra-long armpads (16" x 3 ½") (SOXWFAP16 or SOXWGAP16) and the Arm Trough with Hand Support (SOCAR)</b>	
<b>SOADPDLH</b>	Dual Post Height Adjustable – Desk Length Left Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
<b>SOADPDRH</b>	Dual Post Height Adjustable – Desk Length Right Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
<b>SOADPFLH</b>	Dual Post Height Adjustable – Full Length Left Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168
<b>SOADPFRH</b>	Dual Post Height Adjustable – Full Length Right Side – (Height Adjustable from 10 ¼" to 15 ¼" with ½" increments)	\$168

**47. SIGNATURE TAILORED SOLUTIONS - RECLINING ARMRESTS WITH VERTICAL POST - FLIP-UP, HEIGHT ADJUSTABLE (CHOOSE 2)**

▼	<b>Height Adjustable from 10.25" (260 mm) to 16.25" (415 mm) in 1" (25 mm) increments</b>	
▼	<b>Choose a Height Range :</b>	<b>Low Range 10.25" to 13.25" (260 - 335 mm)      High Range 14.25" to 16.25" (360 - 415 mm)</b>
▼	<b>Not compatible with Waterfall Armpads (SOWAPx and SOWGELx)</b>	
<b>SOANIDLH</b>	Signature Reclining Armrest with Vertical Post – Desk Length 12" (305 mm) - Left Side	N/C
<b>SOANIDRH</b>	Signature Reclining Armrest with Vertical Post – Desk Length 12" (305 mm) - Right Side	N/C
<b>SOANIFLH</b>	Signature Reclining Armrest with Vertical Post – Full Length 16" (405 mm) - Left Side	N/C



**47. SIGNATURE TAILORED SOLUTIONS - RECLINING ARMRESTS WITH VERTICAL POST - FLIP-UP, HEIGHT ADJUSTABLE (CHOOSE 2) (CONTINUED)**

<b>SOANIFRH</b>	Signature Reclining Armrest with Vertical Post – Full Length 16" (405 mm) - Right Side	N/C
<b>SOANIELH</b>	Signature Reclining Armrest with Vertical Post – Extended Length 18" (455 mm) - Left Side	N/C
<b>SOANIERH</b>	Signature Reclining Armrest with Vertical Post – Extended Length 18" (455 mm) - Right Side	N/C
▲ <i>The extended-length 18" (455 mm) armrests are only compatible with the extra-long armpads (16" x 3 ½") (SOXWFAP16 or SOXWGAP16)</i>		
▼ <b>CUSTOM LENGTH</b>		
<b>SOANICLH</b>	Signature Reclining Armrest with Vertical Post– Custom Length 20" (510 mm) - Left Side	\$95
<b>SOANICRH</b>	Signature Reclining Armrest with Vertical Post– Custom Length 20" (510 mm) - Right Side	\$95
<b>SOANICLH</b>	Signature Reclining Armrest with Vertical Post– Custom Length 22" (560 mm) - Left Side	\$95
<b>SOANICRH</b>	Signature Reclining Armrest with Vertical Post– Custom Length 22" (560 mm) - Right Side	\$95
▲ <i>The custom-length 20" (510 mm) and 22" (560 mm) armrests are only compatible with the extra-long armpads (16" x 3 ½") (SOXWFAP16 or SOXWGAP16)</i>		

**48. ARMREST HEIGHT ADJUSTMENT AT FACTORY**

▶	Set the Armrest Height at: _____ "
▲	The Armrest Height is taken from the seat pan to the top of the armpad.
▲	Not required for Dual Post.

**49. SIGNATURE TAILORED SOLUTIONS - ARMPADS**

**FOAM TYPE SELECTION**

<b>SOFAPD</b>	Desk Length Flat Armpad – (10 ½" x 2 ½")	Select side:	Left Hand Side	Right Hand Side	STD
<b>SOFAPF</b>	Full Length Flat Armpad – (14" x 2 ½")	Select side:	Left Hand Side	Right Hand Side	STD
<b>SOXWFAP10</b>	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 10" (255 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C
<b>SOXWFAP14</b>	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 14" (355 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C
<b>SOXWFAP16</b>	X-Wide 3.5" (90 mm) Flat FOAM Armpad, 16" (405 mm) Long	Select side:	Left Hand Side	Right Hand Side	N/C
▲ <i>The X-Wide 16" x 3 ½" Armpads require Cantilever (SOAFBExx) or Reclining armrest (SOANIEEx) with an extended length of 18"(455 mm) or more</i>					
▲ <i>The X-Wide 16" x 3 ½" Armpads are NOT compatible with Dual Post Armrests (SOADPxx)</i>					
<b>SOWAPD</b>	Desk Length Waterfall – (8 ½" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	N/C
▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>					
▲ <i>NOT compatible with Cantilever (SOAFBExx) or Reclining armrests (SOANIEEx)</i>					
<b>SOWAPF</b>	Full Length Waterfall – (13" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	N/C
▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>					
▲ <i>NOT compatible with Cantilever (SOAFBxxx) or Reclining armrests (SOANIxxx)</i>					



**49. SIGNATURE TAILORED SOLUTIONS - ARMPADS (CONTINUED)**

**HYBRID GEL AND FOAM TYPE SELECTION (THE GEL ALWAYS STAYS SOFT)**

<b>SOGELD</b>	Desk Length Hybrid Gel and Foam Flat armpad – (9 ¾" x 2")	Select side:	Left Hand Side	Right Hand Side	\$36
<b>SOGELF</b>	Full Length Hybrid Gel and Foam Flat armpad – (13 ½" x 2")	Select side:	Left Hand Side	Right Hand Side	\$36
<b>SOXWGAP10</b>	X-Wide 3.5" (90 mm) Flat GEL Armpad, 10" (255 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$95
<b>SOXWGAP14</b>	X-Wide 3.5" (90 mm) Flat GEL Armpad, 14" (355 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$95
<b>SOXWGAP16</b>	X-Wide 3.5" (90 mm) Flat GEL Armpad, 16" (405 mm) Long	Select side:	Left Hand Side	Right Hand Side	\$105

- ▲ *The X-Wide 16" x 3 ½" Armpads require Cantilever (SOAFBExx) or Reclining armrest (SOANIExx) with an extended length of 18" (455 mm) or more*
- ▲ *The X-Wide 16" x 3 ½" Armpads are NOT compatible with Dual Post Armrests (SOADPxx)*

<b>SOWGELD</b>	Desk Length Hybrid Gel & Foam Waterfall armpad – (8 ½" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx) or Reclining armrests (SOANlxxx)</i>				
<b>SOWGELF</b>	Full Length Hybrid Gel and Foam Waterfall armpad – (13" x 2 ¾")	Select side:	Left Hand Side	Right Hand Side	\$95
	▲ <i>May interfere with flipping back the armrests when combined with Lateral Trunk supports.</i>				
	▲ <i>NOT compatible with Cantilever (SOAFBxxx) or Reclining armrests (SOANlxxx)</i>				

**► SPECIALTY TYPE SELECTION**

<b>SOCAR</b>	Arm Trough with Hand Support (18" x 4 ½")	Select side:	Left Hand Side	Right Hand Side	\$290
	▲ <i>Arm Trough Armpads may interfere with Armrests Flipping back.</i>				
	▲ <i>SOCAR armpads require Cantilever (SOAFBExx) or Reclining armrest (SOANIExx) with an extended length of 18" (455 mm) or more</i>				

**50. SIGNATURE TAILORED SOLUTIONS - FOOTRESTS**

<b>SOAAF</b>	Signature Heavy Duty Swing-Away Footrests				\$175
	Set at:	60°		70°	
	Select Footrest Range:	SOFLHD12 (12 ½" – 17 ½")		SOFLHD14 (14 ½" – 19 ½")	
<b>ADDML</b>	Signature Manual Elevating & Articulating Legrests - Pair				\$595
	Select Footrest Range:	SOFLME12 (12 ½" – 17 ½")		SOFLME14 (14 ½" – 19 ½")	
	▲ <i>Does not include calf pads. Please select calf pads in section 51.</i>				
<b>ADDPLHD</b>	Signature Heavy Duty Power Elevating & Articulating Legrests				\$2,745
	Select Function:	Individual		Combined	
	Select Footrest Range:	SOFLPW12 (12 ½" – 17 ½")		SOFLPW14 (14 ½" – 19 ½")	
	▲ <i>This option includes calf pads. Individual legs count as two power functions. Combined legs count as one power function.</i>				

**51. SIGNATURE TAILORED SOLUTIONS - FOOTPLATES & ACCESSORIES**

<b>SOSFAF</b>	Flip-up Footplates				STD
<b>SOAFP</b>	Angle Adjustable & Flip-up Footplates – (Each)				\$65
<b>SOAHL</b>	Heel Loops – (Each)				\$20
<b>SOACPFR</b>	Calf Pads for Adult Heavy-Duty Swing-Away footrests – (SOAAF) – with adjustable hardware – (Each)				\$93
<b>SOAMA</b>	Swing-away residual limb Support	Select:	Left leg	Right leg	\$455



**52. FOOTREST LENGTH ADJUSTMENT AT FACTORY**

- ▶ Set the Footrest Length at: \_\_\_\_\_ “
- ▲ The Footrest Length is taken from the top of the seat pan to the top of the footplate.

**53. SIGNATURE TAILORED SOLUTIONS - CENTER MOUNTS**

<b>SOAAF90</b>	Fixed Center Mount Interface Bracket – (Offers Manual adjustability from 94° to 70°)	\$200
	Select Foot Platform Range:	
	SOFBA12 (12 ½” – 16 ½”)	SOFBA16 (16 ½” – 19 ½”)
	▲ <i>When selecting SOFBA16 (16 ½” - 19 ½”) range, the lowest possible seat to floor height is 18.75”</i>	
	▲ <i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>	
	▲ <i>With the split footplates option (SOAFPFB), this bracket offers independent height adjustability for each footplate.</i>	
<b>ADDPFVHD</b>	90° Heavy Duty Power Center Mount Elevating & Articulating Foot Platform	\$2,745
	Select Footrest Range:	
	SOFBPW12 (10 ½” - 16”)	SOFBPW16 (16” - 20”)
	▲ <i>When selecting SOFBPW16 (16”- 20”) range, the lowest possible seat to floor height is 18.75”</i>	
	▲ <i>This option includes calf pads. The power center mount counts as one power function.</i>	
	▲ <i>Select the one piece center mount foot platform (SOFUFP) or the split footplates (SOAFPFB).</i>	
	▲ <i>For leg lengths between 10 ½” and 12 ½”, there is a partial conflict where the calf pads prevent the foot platform or foot plates from folding completely.</i>	

**54. SIGNATURE TAILORED SOLUTIONS - FOOT PLATFORMS & ACCESSORIES**

<b>SOFUFP</b>	Flip-up one piece Foot Platform	\$140
<b>SOAFPFB</b>	Angle Adjustable & Flip-up Split Footplates for center mount – (Each)	\$65
<b>SOACPCM</b>	Calf Pads for 90° Fixed Center Mount – (SOAAF90) with adjustable hardware – (Each)	\$93

**55. FOOT PLATFORM HEIGHT ADJUSTMENT AT FACTORY**

- ▶ Set the Foot Platform Height at: \_\_\_\_\_ “
- ▲ The Foot Platform Height is taken from the top of the seat pan to the top of the foot platform or split foot plates.

**56. SEAT POSITIONING**

<b>SOAUTO80-STD</b>	80” Seat Positioning Strap	STD
<b>SOBPHP2C</b>	Bodypoint padded Hip Belt – (Two Point – Center Pull)	\$170
<b>SOBPHP2C-L82</b>	XL 82” Bodypoint padded Hip Belt – (Two Point – Center Pull)	\$190
<b>SOPHB</b>	Padded Hip Belt 60”, 2-PT Magnetic Center Pull	\$194
<b>SOPCS</b>	Padded Chest Strap 60”, 2-PT Magnetic Center Pull	\$262

**57. HEADREST**

▼ **SELECT ONE OF THE HEADREST OPTIONS WITH THE MOUNTING HARDWARE AND SPECIFY WHERE THE HEADREST WILL BE INSTALLED.**

<b>SONUAHR-12</b>	Using headrest mounting hardware other than SOMAH below. <b>Available with both Fixed and Reclining backrests. 1/2 IN (12.7 MM) Mounts</b>	N/C
<b>SONUAHR-15</b>	Using headrest mounting hardware other than SOMAH below. <b>Available with Reclining backrests only. 19/32 IN (15 MM) Mounts</b>	N/C
<b>SOFHRS</b>	Small Curved Headrest Pad – 4” x 9”	\$155
<b>SOFHRM</b>	Medium Curved Headrest Pad – 5” x 10”	\$155
<b>SOFHRL</b>	Large Curved Headrest Pad – 6” x 11”	\$155
<b>SOMAH</b>	Multi-axis Removable Mounting Hardware	\$220
<b>SOHRBC</b>	Bracket to Mount on the Contour Back Support – (Without Recline)	N/C
<b>SOHRPR</b>	Bracket to Mount on the Seat Frame – (For Reclining Backs)	N/C



**58. LATERAL TRUNK SUPPORTS**

▼ **SELECT ONE OF THE LATERAL PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

<b>SOSLS3</b>	Adjustable Laterals with Small Curved Pads – (Each) – 3” x 5”	\$100
<b>SOSLM4</b>	Adjustable Laterals with Medium Curved Pads – (Each) – 4” x 6”	\$100
<b>SOSLL5</b>	Adjustable Laterals with Large Curved Pads – (Each) – 5” x 7”	\$100
<b>SOSLHM</b>	Adjustable Swing-Away Hardware Mounts – (Each)	\$220
▲	<i>The Lateral Trunk Supports are not compatible with the NXT Optima FH PC 6” Deep Contour Fixed Thoracic Back Support (NOPCB).</i>	
▲	<i>Waterfall armpads (SOWAPx &amp; SOWGELx) may interfere with flipping back the armrests when combined with Lateral Trunk supports</i>	

**59. HIP & THIGH SUPPORTS**

▼ **SELECT ONE OF THE HIP & THIGH SUPPORT PADS AND THE ADJUSTABLE & REMOVABLE HARDWARE.**

<b>SOHTSS10</b>	Adjustable Hip & Thigh Support Pads – (Each) – Small 10” x 4”	\$130
<b>SOHTSM12</b>	Adjustable Hip & Thigh Support Pads – (Each) – Medium 12” x 5”	\$130
<b>SOHTSL14</b>	Adjustable Hip & Thigh Support Pads – (Each) – Large 14” x 6”	\$130
<b>SOHTSARH</b>	Adjustable and Removable Hip & Thigh Support Hardware – (Each)	\$136

**60. ACCESSORIES**

<b>SOCH</b>	Cup Holder – (Self-Leveling)	Install on:	Left Side	Right Side	\$115
<b>SOCPH</b>	Ram Mounts X-Grip Cell Phone Holder	Install on:	Left Side	Right Side	\$190
<b>SOTHS</b>	Ram Mounts X-Grip Small Tablet Holder	Install on:	Left Side	Right Side	\$205
<b>SOTHL</b>	Ram Mounts X-Grip Large Tablet Holder with Swivel	Install on:	Left Side	Right Side	\$360
▲	<i>Standard installation is on the opposite side of joystick. Installing accessories on the same side as the joystick may require a custom upgrade.</i>				
▲	<i>With multiple accessories installed on the same side, a custom upgrade will be required.</i>				

**61. REFLECTOR KIT**

<b>REFLECT</b>	Red & White Reflector Kit Stickers (They will not be installed – Included in a bag with the Owner’s Manual)	N/C
----------------	---	-----

**62. BACK PACK AND ARMREST SIDE POUCHES**

<b>SOAP</b>	Armrest Side Pouch – (Comes with durable straps that can be hooked around the armrest tubes)	Order:	Left Side	Right Side	\$53
▲	<i>If only choosing one (1) armrest pouch, please show what side it should be installed</i>				
<b>S OBP</b>	Back Pack – (Comes with durable straps that can be hooked around the back canes)				\$90

**63. SWITCHES AND OPTIONS**

<b>ASCPWIFI</b>	PGDT - Wireless programmer, WIFI dongle	\$844
<b>ASCPCP</b>	PGDT R-NET Dongle and PC programmer kit	\$600
<b>ASCPVR2</b>	PGDT VR2 with PPI Dongle and PC programmer kit	\$1,450
<b>SSBS</b>	Buddy Button Switch	\$230
<b>SSES</b>	Egg Switch-Black	\$116
<b>SSML</b>	Ultra Light Switch	\$295
<b>ESS-RED</b>	Emergency STOP Switch - Red	\$314



**SPECIALTY DRIVER CONTROL OPTIONS**

► Specialty Drive Control options all require the Omni 2 Enhanced Display & Digital Interface Box (EDIB)

**64. THERAFIN SIP-N-PUFF KIT**

<b>SCSNPK</b>	Therafin Sip-N-Puff Kit				\$1,975
<b>SCSNPBTk</b>	Therafin Sip-N-Puff Breath Tube Kit				N/C
	▲ <i>Required with the Therafin Sip-N-Puff kit.</i>				
<b>SCSRMSNP</b>	Stealth Rock-n-Lock Swingaway Arm for Sip-N-Puff	SWL (Swings on Left)	SWR (Swings on Right)		\$300
	▲ <i>Required with the Therafin Sip-N-Puff kit.</i>				
<b>SCSRM</b>	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
<b>SSES</b>	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i>				
	▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES &amp; OPTIONS section</i>				
<b>ESS-RED</b>	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ <i>Required with the Therafin Sip-N-Puff kit.</i>				
<b>SOMAH</b>	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ <i>Required with the Therafin Sip-N-Puff kit.</i>				
► <b>CHOOSE ONE (1) HEADREST PAD</b>					
<b>SOFHRS</b>	Small Curved Headrest Pad – 4" x 9"				\$155
<b>SOFHRM</b>	Medium Curved Headrest Pad – 5" x 10"				\$155
<b>SOFHRL</b>	Large Curved Headrest Pad – 6" x 11"				\$155

**65. PROXIMITY HEAD ARRAY**

<b>SCSIHAA</b>	Adult - Switch-It 3 Switch Proximity Head Array with 1/8" Mono Jack for mode/reverse				\$5,995
<b>SCSRM</b>	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510
<b>SSES</b>	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)		\$116
	▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i>				
	▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES &amp; OPTIONS section.</i>				
<b>ESS-RED</b>	Emergency STOP Switch - Red	User Activation	Attendant Activation		\$314
	▲ <i>Required with the Proximity Head Array.</i>				
<b>SOMAH</b>	Multi-axis Headrest Removable Mounting Hardware				\$220
	▲ <i>Required with the Proximity Head Array.</i>				

**66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL**

<b>SCSIMP</b>	Switch-It Light touch mini Proportional Joystick				\$4,495
<b>SCSRMMP</b>	Stealth Rock-n-Lock Swingaway Arm for Chin MP	SWL (Swings on Left)	SWR (Swings on Right)		\$870
	▲ <i>Required with the Light touch mini Proportional Joystick.</i>				
<b>SCSRM</b>	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)		\$510



**66. LIGHT TOUCH PROPORTIONAL MINI JOYSTICK FOR CHIN / LIP CONTROL (CONTINUED)**

<b>SSES</b>	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i>			
	▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES &amp; OPTIONS section</i>			
<b>JCC</b>	Chin Cup			\$55
<b>SOMAH</b>	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ <i>Required with the Light touch mini Proportional Joystick.</i>			
▶ <b>CHOOSE ONE (1) HEADREST PAD</b>				
<b>SOFHRS</b>	Small Curved Headrest Pad – 4" x 9"			\$155
<b>SOFHRM</b>	Medium Curved Headrest Pad – 5" x 10"			\$155
<b>SOFHRL</b>	Large Curved Headrest Pad – 6" x 11"			\$155

**67. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR MIDLINE MOUNT**

<b>SCPGMJS</b>	PG Drives JC200 Proportional Compact Remote Joystick for Midline Mount			\$2,369
<b>SCGMJSM</b>	Gatlin Midline Swing-away Mount for JC200	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i>			
<b>SCSIMP</b>	Switch-It Light touch mini Proportional Joystick			\$4,495
<b>SCGMPPM</b>	Gatlin Midline Swing-away Mount for Light touch mini proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i>			
<b>SCBTM</b>	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i>			

**68. PROPORTIONAL COMPACT OR MINI JOYSTICK FOR ARMREST MOUNT**

<b>SCPGMJS-M</b>	PG Drives JC200 Proportional Compact Remote Joystick for Armrest Mount			\$2,369
<b>SCARMJSM</b>	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$325
	▲ <i>Must be selected for use with the PG Drives JC200 Mini Joystick.</i>			
<b>SCSIMP</b>	Switch-It Light touch mini Proportional Joystick			\$4,495
<b>SCARMPM</b>	Swing-away and Multi-Axis armrest mount	MLH (Mount on Left)	MRH (Mount on Right)	\$375
	▲ <i>Must be selected for use with the Light touch mini Proportional Joystick.</i>			
<b>SCBTM</b>	Bullet Hand Tray for Light Touch mini Proportional Joystick	MLH (Mount on Left)	MRH (Mount on Right)	\$750
	▲ <i>This option is only compatible with the Light touch Mini Proportional Joystick.</i>			

**69. PROPORTIONAL RIM HEAD CONTROL PACKAGE**

<b>SCPHC</b>	Proportional RIM Head Control Package			\$2,795
<b>SCSRM</b>	Stealth Rock-n-Lock Swingaway Mount for Mode Switch	SWL (Swings on Left)	SWR (Swings on Right)	\$510
<b>SSES</b>	Egg Switch for Changing Modes – Black	MLH (Mount on Left)	MRH (Mount on Right)	\$116
	▲ <i>Required with the Stealth Rock-n-Lock Swing Away mount.</i>			
	▲ <i>The Mode Switch Type can be changed from this selection in SWITCHES &amp; OPTIONS section</i>			
<b>SOMAH</b>	Multi-axis Headrest Removable Mounting Hardware			\$220
	▲ <i>Required with the Proportional RIM Head Control.</i>			



**CUSTOM SECTION**

▶ **Lead times will increase depending on the custom order pipeline at the time the purchase order is received.**  
 ▶ **Please call for Quote and specifications for these options.**

**70. O<sup>2</sup> HOLDER**

<b>SUPOXY</b>	Custom Non-Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$295
	▲ <b>Not available with a power recline</b>			
<b>SUPOXYA</b>	Custom Articulating Oxygen Tank Holder	MLH (Mount on Left)	MRH (Mount on Right)	\$1,195
	▲ <b>Required with a power recline</b>			
	D-Size	Oxygen Tank – (* 22" Tall x 4.3" Diameter) – * The height includes the valve		
	E-Size	Oxygen Tank – (* 31" Tall x 4.3" Diameter) – * The height includes the valve		

**71. VENT TRAY (PLEASE COMPLETE IN FULL)**

▼ **AMYLOR DOES NOT SUPPLY THE BATTERY, CHARGER OR THE WIRING FOR THE VENT. THESE ITEMS ARE NORMALLY PROVIDED BY THE VENT MANUFACTURER.**

<b>SOVENT</b>	Custom Vent Tray			
	Articulating Vent Box (Required with Power Tilt and/or Recline)			\$1,360
	Fixed Vent Box			\$400
	<b>SPECIFY MODEL</b>			
	Trilogy 100/200			
	LTV			
	LTV with Sprintpack			
	Other, specify: _____			
<b>SOVTBE</b>	Extended Battery Compartment			\$495
	The Vent is being powered by 1 U1 Battery and the chair by 2 Group 24 batteries			
	The Vent is being powered by 1 Group 22NF Battery and the chair Powered by 2 Group 24 batteries			

**72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE)**

<b>CUSTOM</b>	Rear view Camera with 5" Monitor				RFQ
<b>CUSTOM</b>	USB Charger Port	◀ Specify the location for the installation: _____			RFQ
<b>CUSTOM</b>	<b>Folding Backrest for Storage/Travel Access for Reclining backs</b>				RFQ
<b>CUSTOM</b>	Large Ball Handle				RFQ
<b>CUSTOM</b>	Manual Wheel Locks				RFQ
<b>CUSTOM</b>	Power Anterior Tilt				RFQ
<b>CUSTOM</b>	Power Flip-up Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				RFQ
<b>CUSTOM</b>	Power Vertical Lifting Foot Platform – (Please call CSR for availability with configuration and user weight capacity)				RFQ
<b>CUSTOM</b>	Rearview mirror	Left Side	Right Side		RFQ
<b>CUSTOM</b>	Independent Depth Adjustability built into the Footrest Receiver that is attached to the Seat Pan				RFQ
<b>CUSTOM</b>	Transfer Handles – (Pair)	Select a height range:	5" to 8"	6" to 9"	RFQ



**72. OTHER CUSTOM OPTIONS (REQUEST FOR QUOTE) (CONTINUED)**

<b>CUSTOM</b>	Foot Box	◀ Specify	Length: _____	Width: _____	Please confirm with CSR on sizing possibilities		RFQ
<b>CUSTOM</b>	Accessory Hooks located at the Rear of the seat structure			2 Hooks	4 Hooks		RFQ
<b>CUSTOM</b>	Heavy-Duty Wheelchair Cover, Large						RFQ
<b>CUSTOM #1</b>	Other custom options. Specify in notes section below						RFQ
<b>CUSTOM #2</b>	Other custom options. Specify in notes section below						RFQ

**73. NOTES / SPECIAL REQUESTS / CUSTOM SPECIFICATIONS**