



A	103641	925-L TRUNK SUPPORT
B	103637	924-L PELVIC SUPPORT
C	103633	930-L PELVIC HARNESS
D	103613	809-L ARM SUPPORT
E	103599	926-L HANDLEBAR FOR FRONT DRIVE
F	303235	FRAME FOR FRONT DRIVE
G	102810	890-C Large LEG DIVIDER+ABDUCTOR
N	103656	928 med/large REAR WHEELS NARROW BASE WITH BRAKE PEDALS
R	105605	810 M-L PAIR OF WEIGHT BARS
S	102813	890-SC Large LEG DIVIDER
T	102970	939 ERGONOMIC HANDLES
V	103616	943 M-L ARM STRAPS
POS.	PART NR.	DESCRIPTION

Spare Parts FRONT DRIVE GRILLO size large



F1	301797	PAIR OF KNOBS TO FIX HORIZONT.SUPPORT.
F2	302284	PAIR OF HORIZ.BARS FOR HANDLEBAR GRILLO large(5)
F3	302450	RGT-LEFT SUPPORTS OF 809 GRILLO
F4	301799	PAIR OF M8X35 HANDLES
F5	301800	PAIR OF M6X30 HANDLES
F6	301801	SET OF LEVERS TO LOCK THE HANDLEBAR VERTICALLY
POS.	PART NR.	DESCRIPTION



Spare Parts **HANDLEBAR SUPPORT + HORIZONTAL BARS**
FRONT-DRIVE GRILLO size Large

pag. 2

ORMESA®




GRILLO

THIS DRAWING BELONGS TO ORMESA: IT IS FORBIDDEN TO REPRODUCE, CHANGE OR INFORM ANY THIRD PARTY ABOUT ITS CONTENTS





F7	303236	UPPER CROSSPIECE SUPPORT+LOCKING LEVERS
F8	302287	LEFT UPPER CROSSPIECE SUPPORT+LOCKING LEVER
F9	302286	RIGHT UPPER CROSSPIECE SUPPORT+LOCKING LEVER
F10	302288	SET OF LEVERS TO LOCK TRUNK SUPPORT VERTICALLY
F11	301808	PAIR OF KNOBS FOR GAS SPRING GRILLO
F12	302290	COMMAND SUPPORT + LEVERS AND KNOBS
F13	303406	CABLES SET+ DOUBLER
POS.	PART NR.	DESCRIPTION

Spare Parts UPPER PART OF FRONT-DRIVE
GRILLO BASE size Large




F13	303406	CABLES SET+ DOUBLER FOR GRILLO large(5)
F14	302292	COMPLETE RGT SIDE GRILLO large(5)
F15	302293	COMPLETE LEFT SIDE GRILLO large(5)
F16	302033	UPPER GAS SPRING GRILLO large(5)
F17	302395	LOWER GAS SPRING GRILLO large(5)
F18	303237	FRONT-DRIVE GRILLO Large BASE + WHEELS
POS.	PART NR.	DESCRIPTION

Spare Parts FRONT-DRIVE GRILLO SIDES size Large



F19	303508	PAIR OF DIRECT.CAPS FOR WHEELS GRILLO
F20	301812	PAIR OF UPPER/LOWER SUPPORTS OF SIDE GRILLO
F21	303208	PAIR OF WHEEL SUPPORTS
F22	303211	KIT OF LOWER CROSSPIECE SUPPORT
F23	303238	PAIR OF BASE SIDES
F24	303239	LOWER CROSSPIECE+SCREWS
F26	303476	PAIR OF COMPLETE REAR WHEELS
F27	302546	PAIR OF DIRECTION LOCKS
F28	303387	PAIR OF COMPLETE FRONT WHEELS
POS.	PART NR.	DESCRIPTION


Spare Parts LOWER PART OF FRONT-DRIVE
GRILLO BASE size Large



A1	301832	SET OF CLAMPING LEVER TO ROTATE TRUNK SUPPORT
A2	160571P	CENTRAL PADDING
A3	303190	PAIR OF TRUNK HALF SUPPORT
A4	303240	EXTERNAL ROTATION SUPPORT + WIDTH ADJUSTMENT
A5	302298	PAIR OF PROFILE TO SLIDE TRUNK/PELVIC SUPPORT
A6	303241	PAIR OF TRUNK PADDING + STRAPS AND RIVETS
A7	301842	SET OF LEVERS FOR TRUNK SUPPORT WIDTH
A8	303223	PAIR OF TRUNK SUPPORT JOINTS
A9	301843	SET OF TRUNK SUPPORT JOINT LEVERS
A10	303242	FRAME+JOINTS TRUNK SUPPORT



POS. PART NR. DESCRIPTION

Spare Parts TRUNK SUPPORT size Large





B1	160571P	CENTRAL PADDING
B2	303189	PAIR OF PELVIC HALF SUPPORT
B3	303244	EXTERNAL SUPPORT + PELVIC WIDTH ADJUSTMENT
B4	303245	PAIR OF PELVIC PADDING + STRAPS AND RIVETS
B5	303246	PELVIC TRAVERSE SUPPORT
B6	302298	PAIR OF PROFILE TO SLIDE TRUNK SUPPORT
B7	301842	SET OF LEVERS FOR TRUNK SUPPORT WIDTH
POS.	PART NR.	DESCRIPTION

Spare Parts PELVIC SUPPORT size Large



D1	171254	809 medium/large ARM SUPPORT UPHOLSTERY
D2	301850	PAIR OF CURVED PLATE 809 medium/large
D3	301851	PAIR OF LOWER/UPPER TOOTHED SUPPORTS 809
D4	303108	PAIR OF M6X25 LEVERS
D5	301852	SET OF TOOTHED SUPP.LEVERS 809
D6	301853	PAIR OF HANDLEBAR BARS 809
D7	301855	PAIR OF VERT.BARS 809
POS.	CODICE	DESCRIZIONE

Spare Parts ARM SUPPORTS size Medium and Large



G1	171535	890 large(5) UPHOLSTERY for GRILLO
G2	302933	PAIR OF CURVED PLATE 890 large(5) GRILLO
G3	301858	PAIR OF ABDUCT.SLIDING PROFILE GRILLO
G4	301859	PAIR OF ABDUCT. 890 STOP GRILLO
G5	302304	HORIZONTAL BAR 890 large(5) + CAP GRILLO
G6	301862	890 FIXING BAR +LEVER/SCREWS
G7	303232	890 VERTICAL BAR + CAP
G8	303233	FIXING LOCK 890 + SCREWS
G9	301866	SET OF HEIGHT ADJUST.LEVER 890 GRILLO
POS.	PART NR.	DESCRIPTION

Spare Parts LEG DIVIDER+ABDUCTORS size Large

pag. 9

ORMESA®

GRILLO

THIS DRAWING BELONGS TO ORMESA: IT IS FORBIDDEN TO REPRODUCE, CHANGE OR INFORM ANY THIRD PARTY ABOUT ITS CONTENTS